




Products for DNA Research

2019 Catalog



Oligo synthesis success. The first time and every time.

TABLE OF CONTENTS

INTRODUCTION	5
ABOUT US	5
CATALOG	6
STERLING	7
QUALITY AND PERFORMANCE ASSURED	7
APPLIED BIOSYSTEMS INSTRUMENTS	8
STERLING CE PHOSPHORAMIDITES	8
STERLING SOLVENTS/REAGENTS	8
STERLING SUPPORTS	9
AB 3900 POLYSTYRENE MODIFIER COLUMNS	11
EXPEDITE™ INSTRUMENTS	12
STERLING CE PHOSPHORAMIDITES	12
STERLING SOLVENTS/REAGENTS	12
STERLING SUPPORTS	13
DNA PHOSPHORAMIDITES - SPECIAL PACKAGING	15
MERMADE INSTRUMENTS	16
STERLING CE PHOSPHORAMIDITES	16
STERLING SOLVENTS/REAGENTS	16
STERLING SUPPORTS	17
GE HEALTHCARE LIFE SCIENCES INSTRUMENTS	18
STERLING CE PHOSPHORAMIDITES	18
STERLING SOLVENTS/REAGENTS	19
DR. OLIGO INSTRUMENTS	20
STERLING CE PHOSPHORAMIDITES	20
STERLING SOLVENTS/REAGENTS	20
STERLING SUPPORTS	21
OLIGONUCLEOTIDE PURIFICATION	21
ALTERNATIVE PROTECTING GROUPS	22
DEPURINATION RESISTANT CE PHOSPHORAMIDITES	22
ULTRAMILD CE PHOSPHORAMIDITES	23
ULTRAMILD SUPPORTS	23
ULTRAMILD SOLVENTS/REAGENTS	23
ULTRAMILD DNA SYNTHESIS	23
SUPPORTS	24
GLEN UNYSUPPORT	24
GLEN UNYSUPPORT FC	25
UNIVERSAL SUPPORT III	26
Q-SUPPORTS	27
HIGH LOAD CPG	29
REAGENTS	30
ALTERNATIVE SOLVENTS/REAGENTS	30
CSO FOR NON-AQUEOUS OXIDATION	32
UNICAP PHOSPHORAMIDITE	32
BACKBONE MODIFICATION	33
SULFURIZING REAGENTS	33
5'-CE PHOSPHORAMIDITES	34
5'-SUPPORTS	35
METHYL PHOSPHONAMIDITES	36
PACE PHOSPHORAMIDITES	37
METHYL PHOSPHORAMIDITES	38

TABLE OF CONTENTS

ULTRAMILD SOLVENTS/REAGENTS	38
H-PHOSPHONATE MONOMERS	39
H-PHOSPHONATE REAGENTS	39
THIOPHOSPHORAMIDITES	40
LOCKED ANALOG PHOSPHORAMIDITES	41
OLIGONUCLEOTIDE-DIRECTED MUTAGENESIS	42
TRIMER PHOSPHORAMIDITES	42
DUPLEX STABILITY MODIFICATION	46
BASES AFFECTING DUPLEX STABILITY	46
ZIP NUCLEIC ACIDS (ZNA®)	48
CDPI, MGB™ LABELING	48
SELECTIVELY BINDING COMPLEMENTARY (SBC) OLIGOS	48
UNNATURAL BASE PAIRS	48
CAPS FOR INCREASED DUPLEX STABILITY AND BASE-PAIRING FIDELITY	49
EPIGENETICS	50
DNA METHYLATION	50
PCR/SEQUENCING APPLICATIONS	51
DUPLEX EFFECTS	51
T _m MODULATION	53
CLEANAMP™ MONOMERS	54
CHAIN TERMINATORS	55
STRUCTURAL STUDIES	57
STRUCTURE/ACTIVITY RELATIONSHIP	57
HALOGENATED NUCLEOSIDES	60
DNA DAMAGE/REPAIR	61
CLICK DNA AND RNA LIGATION	64
5'-LABELING OF MicroRNAs	64
2'-5' LINKED OLIGONUCLEOTIDES	65
MUTAGENESIS	66
<i>IN SITU</i> SYNTHESIS OF DNA ANALOGS	67
PROBING DNA STRUCTURE WITH FLUORESCENT NUCLEOSIDES	68
PHOTO-REGULATION OF DNA FUNCTION	70
INHIBITION OF DNA METHYLTRANSFERASES	71
LARGE SCALE SYNTHESIS	71
NON-CANONICAL STRUCTURES	72
G-QUADRUPLEX	72
TRIPLEX-FORMING OLIGONUCLEOTIDES	72
i-MOTIF DNA STRUCTURES	72
APTAMER DEVELOPMENT	73
MODIFIERS	74
TERMINUS MODIFIERS	74
SEQUENCE MODIFIERS	77
3'-MODIFIERS	79
CHEMICAL PHOSPHORYLATION	82
ALDEHYDE MODIFICATION	83
SPACER MODIFIERS	84
DENDRIMERS	85
BRANCHING PHOSPHORAMIDITE	85
PHOTOCLEAVABLE MONOMERS	86
CONJUGATION USING CLICK CHEMISTRY	87
OLIGO-CLICK KITS	90
COPPER-FREE CLICK CHEMISTRY	91
SERINOL REAGENTS FOR MODIFICATION AND LABELING	94
COT SERINOL PHOSPHORAMIDITE	97
DABCYL LABELING	98
BIOTIN LABELING	99
FLUORESCHEIN LABELING	102
FLUORESCHEIN LABELING (SIMA)	105

TABLE OF CONTENTS

CYANINE LABELING	106
ELITECHGROUP DYES AND QUENCHER	108
BLACK HOLE QUENCHER DYES	110
BLACKBERRY® QUENCHER (BBQ-650®)	112
RHODAMINE (TAMRA) LABELING	113
ACRIDINE LABELING	114
DNP LABELING	114
CHOLESTEROL LABELING	115
TOCOPHEROL LABELING	115
STEARYL LABELING	115
N-ACETYL GALACTOSAMINE (GalNAc) LABELING	116
CDPI, MGB™ LABELING	117
PSORALEN LABELING	118
EDTA LABELING	118
FERROCENE LABELING	119
METHYLENE BLUE LABELING	119
LABELING WITH METAL CHELATES	120
LABELING WITH POLYAROMATIC HYDROCARBONS	120
PUROMYCIN CPG	121
QUENCHED AUTOLIGATION (QUAL) PROBES	121
LABELING FOR PHOTO-REGULATION OF OLIGONUCLEOTIDES	122
LABELLING WITH ULTRAFAST PHOTO CROSS-LINKER	123
RNA SUPPORTS	125
RNA SUPPORTS FOR 3' MODIFICATION	125
RNA SYNTHESIS	126
TOM-PROTECTED RNA PHOSPHORAMIDITES	126
RNA SUPPORTS FOR TOM RNA SYNTHESIS	126
TBDMS-PROTECTED RNA PHOSPHORAMIDITES	128
RNA PHOSPHORAMIDITES - SPECIAL PACKAGING	128
ULTRAMILD TBDMS RNA PHOSPHORAMIDITES	129
TBDMS RNA SUPPORTS	129
ULTRAMILD SOLVENTS/REAGENTS	130
MINOR RNA BASES	131
MINOR RNA PHOSPHORAMIDITES (TOM PROTECTED)	131
RNA SEQUENCE MODIFIER (TOM PROTECTED)	132
MINOR RNA PHOSPHORAMIDITES (TBDMS PROTECTED)	133
MINOR RNA TRIPHOSPHATES	136
2'-OME-RNA SYNTHESIS	137
2'-OME-RNA PHOSPHORAMIDITES	137
ULTRAMILD 2'-OME-RNA	138
ULTRAMILD SOLVENTS/REAGENTS	138
2'-OME-RNA SUPPORTS	139
MINOR 2'-OME-RNA PHOSPHORAMIDITES	140
2'-OME-THIOPHOSPHORAMIDITES	142
2'-F RNA SYNTHESIS	143
2'-F-RNA PHOSPHORAMIDITES	143
2'-F ANA SYNTHESIS	144
2'-F-ARABINONUCLEIC ACID (2'-F-ANA)	144
2'-OME-RNA-PACE SYNTHESIS	145
2'-OME-RNA-PACE PHOSPHORAMIDITES	145
PURIFICATION	147
GLEN-PAK™ PURIFICATION	147
POLY-PAK™ PURIFICATION	149
GLEN GEL-PAK™ DESALTING	150
OLIGO-AFFINITY SUPPORT	151

PHYSICAL DATA	152
PHYSICAL DATA	152
INDEX	162
GENERAL INFORMATION	172
ORDERING	172
DISCOUNTS	172
TERMS AND CONDITIONS OF SALE	172
PATENTS	172

ABOUT US

Glen Research develops, manufactures and markets reagents for oligonucleotide synthesis, modification, labeling and purification. The company serves customers worldwide involved in basic research, diagnostics and therapeutics. Although Glen Research's original mission was to provide state-of-the-art reagents to researchers, the company also began offering standard reagents for oligonucleotide synthesis but with the innovation that every batch was accompanied by a Certificate of Analysis. The analytical techniques and quality criteria used for the evaluation and acceptance of these reagents were to become an industry standard years later. The company is headquartered in Sterling, Virginia. A privately held company, Glen Research was acquired by Maravai LifeSciences in December 2017.

OVER 30 YEARS OF ASSURED QUALITY FOR OLIGO SYNTHESIS

<p>1987</p> <p>Glen Research was incorporated in the Commonwealth of Virginia</p>	<p>1991</p> <p>Company awarded SBIR grant for the investigation of large scale oligonucleotide synthesis using H-phosphonate chemistry</p>
<p>1993</p> <p>Glen Research introduced the Sterling line of products, a new standard of quality for oligonucleotide synthesis</p>	<p>1995</p> <p>Glen Research negotiated an exclusive agreement to supply 5'-biotin phosphoramidite worldwide</p>
<p>1996</p> <p>Company negotiated an exclusive license with Gilead Sciences to supply C5-propynyl pyrimidine nucleosides and G-Clamp phosphoramidites</p>	<p>1997</p> <p>Glen Research moves into a custom built building in Sterling, Virginia</p>
<p>1999</p> <p>Company awarded patents for a chemical phosphorylation reagent compatible with DMT-ON purification</p>	<p>2002</p> <p>Company made an agreement with Epoch Biosciences, Inc. to supply their proprietary dyes and nucleosides to the research market</p>
<p>2003</p> <p>Glen Research negotiated an agreement with GE Healthcare Biosciences Corp. to supply Cyanine Dyes to the research market</p>	<p>2004</p> <p>Company awarded patents for a truly universal support for oligonucleotide synthesis - US III.</p>
<p>2006</p> <p>In collaboration with Berry & Associates, Inc., Glen Research awarded patents for pyrrolo-C analogues (fluorescent C analogues).</p>	<p>2008</p> <p>Glen Research obtained a license for the sale of Glen UnySupport from Ionis Pharmaceuticals</p>
<p>2013</p> <p>In collaboration with Nelson Biotechnologies, Inc., company awarded patent for serinol phosphoramidites and supports</p>	<p>2017</p> <p>Glen Research is acquired by Maravai LifeSciences</p>

CATALOG

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers For Instrument type	Add
Expedite	E
MerMade	M
Columns For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

Welcome to the Glen Research Catalog containing the most complete selection of products for DNA and RNA research. The Table of Contents at the beginning and the Index at the end of the Catalog are the most comprehensive we have produced. There are always limitations to printed catalogs in a fast-moving technology sector and a complete and up-to-date catalog is also maintained on our web site.

All minor bases, modifiers and RNA products are packaged for Applied Biosystems instruments. We can provide vials and columns for a wide variety of other instruments. As shown in the table to the left, we can accommodate catalog numbers for unusual products to fit all popular instruments. The table to the left is reproduced on all relevant spreads of this catalog.

We are unique in conducting a QC test for supports to show the length of oligo that can be prepared before a drop-off in coupling due to steric effects begins to occur. The drop-off point is recorded in the Certificate of Analysis or Analytical Report. Unless otherwise specified, our minor base and modification supports are 1000Å CPG, which results in improved performance and the ability to make much longer oligos. Polystyrene supports are also available for some of our most popular items.

For reasons of quality assurance, we do not transfer powders or oils from stock Applied Biosystems vials to vials for other instruments. Powders may be hygroscopic and electrostatic, making transfer difficult, and oils have to be dissolved and the solvent evaporated. For best performance, it is preferable for the customer to dissolve the product and immediately transfer the solution to the correct instrument vial. Consequently, the product will be delivered in an industry-standard septum-capped vial along with a clean dry vial for the appropriate instrument.

Glen Research's distributors cover a very significant percentage of countries where oligonucleotide synthesis is commonly practiced. Our vast selection of unusual products is really only comprehensively stocked here in Virginia and some of our web viewers have asked us to set up a direct shipping channel. For them, we offer the eGlen program which is described in the following web link: <http://www.glenresearch.com/Reference/eGlen.html>.

Authorized distributors for Glen Research products are listed below. Other countries not listed are covered by direct sales from our Sterling, USA office.

UK and Ireland

Cambio Ltd
Telephone Number: +44 (0) 1954 210200
Fax Number: +44 (0) 1954 210300
e-mail addresses: support@cambio.co.uk and
orders@cambio.co.uk
Website: <http://www.cambio.co.uk/>

China

Beijing LeBo Biotech Co., Ltd
Telephone Number: +86-10-52405563
Fax Number: +86-10-58850899
email address: info@lab-bio.com
Website: <http://www.lab-bio.com/>

Netherlands

Eurogentec b.v.
Telephone Number: +31 43 352 06 98
Fax Number: +31 43 354 19 65
e-mail address: info@eurogentec.com

Nordic and Baltic Countries

BioNordika AS
Telephone Number: +47 23 03 58 00
Fax Number: +47 23 03 58 01
e-mail address: info@bionordika.no
Website: <http://www.bionordika.no/>

Belgium

Eurogentec S.A.
Telephone Number: +32 4 372 74 00
Fax Number: +32 4 372 75 00
e-mail address: info@eurogentec.com
Website: <http://www.eurogentec.com/>

Germany

Eurogentec GmbH
Telephone Number: +49 221 258 94 55
Fax Number: +49 221 258 94 54
e-mail address: info@eurogentec.com

Republic of Korea

Bosung Scientific Co., Ltd.
Telephone Number: +82-02-6105-5630
Fax Number: +82-02-6105-5680
email address: info@bosungsci.com
Website: <https://bosungsci.com/>

Japan

Nihon Techno Service Co., Ltd.
Telephone Number: +81 29 886 6811
Fax Number: +81 29 870 0210
e-mail address: info@ntsbio.com
Website: <http://www.ntsbio.com/>

Israel

Eisenberg Bros. Ltd.
Telephone Number: 972-3-9777000
Fax Number: 972-3-9777001
e-mail address: nicoles@eb1.co.il
Website: <http://www.eisenbros.co.il/>

France

Eurogentec s.a.
Telephone Number: +33 2 41 73 33 73
Fax Number: +33 2 41 73 10 26
e-mail address: info@eurogentec.com

QUALITY AND PERFORMANCE ASSURED

Glen Research has developed and implemented a Quality Management System (QMS) designed to enhance customer satisfaction by focusing on processes for continual improvement and on assurance of conformity to customer needs, with full consideration of applicable regulatory requirements.

STERLING QUALITY


The benchmark for excellence in DNA and RNA synthesis. All Sterling materials must pass stringent purity and identity tests prior to acceptance. Sterling products are formulated, filtered, and packaged in optimal environments using specially cleaned and dried glassware and columns. Color-coded labeling and post-packaging analysis guarantee accuracy and Sterling Quality.




STERLING is a trademark of Glen Research Corporation.

Glen Research offers the highest level of Quality Assurance for reagents for DNA and RNA synthesis - Sterling Quality and Performance. We now apply the Sterling criteria of quality and performance to all of Glen Research's established products.


The common monomers and supports, whose structures are illustrated below, are available for the variety of synthesizers listed on the following pages.




dA-CE Phosphoramidite




dC-CE Phosphoramidite




Ac-dC-CE Phosphoramidite




dG-CE Phosphoramidite




dT-CE Phosphoramidite




dA-Icaa-CPG




dC-Icaa-CPG



Ac-dC-Icaa-CPG



dG-Icaa-CPG



dT-Icaa-CPG




STERLING

STERLING CE PHOSPHORAMIDITES

Glen Research CE (β-cyanoethyl) Phosphoramidites are produced and packaged to ensure the highest performance on DNA synthesizers. Every Glen Research product is accompanied by a Certificate of Analysis and HPLC trace, showing the results of our QC testing. Every Glen Research monomer vial is specially cleaned to eliminate particulate contamination and tested to ensure a tight fit on synthesizers.

Item	Catalog No.	Pack	Price (\$)
dA-CE Phosphoramidite	10-1000-02	0.25g	12.50
	10-1000-05	0.5g	25.00
	10-1000-10	1.0g	50.00
	10-1000-20	2.0g	100.00
	10-1000-40	4.0g	200.00
dC-CE Phosphoramidite	10-1010-02	0.25g	12.50
	10-1010-05	0.5g	25.00
	10-1010-10	1.0g	50.00
	10-1010-20	2.0g	100.00
	10-1010-40	4.0g	200.00
Ac-dC-CE Phosphoramidite	10-1015-02	0.25g	12.50
	10-1015-05	0.5g	25.00
	10-1015-10	1.0g	50.00
	10-1015-20	2.0g	100.00
	10-1015-40	4.0g	200.00
dG-CE Phosphoramidite	10-1020-02	0.25g	12.50
	10-1020-05	0.5g	25.00
	10-1020-10	1.0g	50.00
	10-1020-20	2.0g	100.00
	10-1020-40	4.0g	200.00
dmf-dG-CE Phosphoramidite	10-1029-02	0.25g	12.50
	10-1029-05	0.5g	25.00
	10-1029-10	1.0g	50.00
	10-1029-20	2.0g	100.00
	10-1029-40	4.0g	200.00
dT-CE Phosphoramidite	10-1030-02	0.25g	12.50
	10-1030-05	0.5g	25.00
	10-1030-10	1.0g	50.00
	10-1030-20	2.0g	100.00
	10-1030-40	4.0g	200.00



dmf-dG-CE Phosphoramidite

ABI INSTRUMENTS

- 60mL septum-capped vials used on oldest ABI 380, 381 and 391 instruments. 200mL oxidizer and 450mL deblock screw-capped bottles also used on ABI 380, 381 and 391 instruments.
- Small screw-capped vials used on ABI 392 and 394 instruments.
- Larger screw-capped vials used on ABI 392, 394 and 3400 instruments.
- Large bottles used on ABI 3900 instruments.

SEE ALSO

Depurination Resistant dA on page 22

STERLING CE PHOSPHORAMIDITES (CONT.)

Item	Catalog No.	Pack	Price (\$)
<i>Cap Mix A</i> THF/Pyridine/Ac2O	40-4110-45 ¹	45mL	16.00
	40-4110-52 ²	200mL	30.00
	40-4110-57 ³	450mL	72.00
	40-4110-62 ⁴	2000mL	325.00
<i>Cap Mix B</i> 16% 1-Melm in THF <i>(This Cap B solution is identical to the formulation produced by Applied Biosystems.)</i>	40-4220-45 ¹	45mL	20.00
	40-4220-52 ²	200mL	40.00
	40-4220-62 ⁴	2000mL	425.00
<i>Oxidizing Solution</i> 0.02M I2 in THF/Pyridine/H2O	40-4330-52 ^{1,2}	200mL	30.00
	40-4330-57 ³	450mL	72.00
	40-4330-62 ⁴	2000mL	325.00
<i>Deblocking Mix</i> 3% TCA/DCM	40-4140-57 ^{1,2}	450mL	36.00
	40-4140-62 ^{3,4}	2000mL	144.00

STERLING SUPPORTS

All Glen Research CPG supports use the standard long chain alkylamino (lcaa) linker but differ in the glass pore size, 500Å, 1000Å or 2000Å. The 500Å support is appropriate for shorter sequences, while the 1000Å supports perform better in the synthesis of longer (>30-mer) DNA sequences. The 2000Å support is best for very long (>150-mer) oligonucleotides. We have instituted an additional QC test for supports to show the length of oligo that can be prepared before a drop-off in coupling due to steric effects begins to occur. The drop-off point is recorded in the Certificate of Analysis. All Glen Research supports are fully end-capped to ensure that the CPG surface is totally inert, thereby avoiding the introduction of impurity sequences containing deletions at the 3'-terminus.


Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)				
dA	dC	dG	dT	dA,dC,dG,dT (1 column of each base)	Ac-dC	dmf-dG		

500Å Columns

20-2100-42	20-2110-42	20-2120-42	20-2130-42	20-2140-42	20-2113-42	4x0.2µm	40.00
20-2100-41	20-2110-41	20-2120-41	20-2130-41	20-2140-41	20-2113-41	4x1.0µm	60.00
20-2100-13	20-2110-13	20-2120-13	20-2130-13		20-2113-13	1x10µm	100.00

1000Å Columns

20-2101-45	20-2111-45	20-2121-45	20-2131-45	20-2141-45	20-2115-45	20-2129-45	4x40nm	40.00
20-2101-42	20-2111-42	20-2121-42	20-2131-42	20-2141-42	20-2115-42	20-2129-42	4x0.2µm	40.00
20-2101-41	20-2111-41	20-2121-41	20-2131-41	20-2141-41	20-2115-41	20-2129-41	4x1.0µm	60.00
20-2101-13	20-2111-13	20-2121-13	20-2131-13		20-2115-13	20-2129-13	1x10µm	100.00



dmf-dG-CPG

ABBREVIATIONS

- Ac₂O = Acetic Anhydride
- CE = Cyanoethyl
- CPG = Controlled Pore Glass
- DCM = Dichloromethane
- dmf = dimethylformamide
- I₂ = Iodine
- lcaa = long chain alkylamino
- Melm = 1-Methylimidazole
- µm = micromole(s)
- nm = nanomole(s)
- TCA = Trichloroacetic Acid
- THF = Tetrahydrofuran

SEE ALSO

Alternative Solvents on page 30

STERLING SUPPORTS (CONT.)

Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	dG	dT	dA,dC,dG,dT (1 column of each base)	Ac-dC	dmf-dG		
2000Å Columns								
20-2102-42	20-2112-42	20-2122-42	20-2132-42	20-2142-42			4x0.2µm	40.00
Low Volume (LV) Polystyrene Columns								
26-2100-45	26-2110-45	26-2120-45	26-2130-45	26-2140-45			4x40nm	48.00
26-2100-42	26-2110-42	26-2120-42	26-2130-42	26-2140-42			4x0.2µm	48.00
AB 3900 Polystyrene Columns								
26-2600-65	26-2610-65		26-2630-65		26-2629-65	200x40nm		825.00
26-2600-62	26-2610-62		26-2630-62		26-2629-62	200x200nm		825.00
AB 3900 1000Å CPG Columns								
20-2101-65			20-2131-65		20-2115-65	20-2129-65	200x40nm	600.00
20-2101-62			20-2131-62		20-2115-62	20-2129-62	200x200nm	650.00
20-2101-61			20-2131-61		20-2115-61	20-2129-61	200x1.0µm	875.00
500Å Bulk CPG								
20-2000-01	20-2010-01	20-2020-01	20-2030-01		20-2013-01		0.1g	9.00
20-2000-02	20-2010-02	20-2020-02	20-2030-02		20-2013-02		0.25g	20.00
20-2000-10	20-2010-10	20-2020-10	20-2030-10		20-2013-10		1.0g	75.00
1000Å Bulk CPG								
20-2001-01	20-2011-01	20-2021-01	20-2031-01		20-2015-01	20-2029-01	0.1g	9.00
20-2001-02	20-2011-02	20-2021-02	20-2031-02		20-2015-02	20-2029-02	0.25g	20.00
20-2001-10	20-2011-10	20-2021-10	20-2031-10		20-2015-10	20-2029-10	1.0g	75.00
2000Å Bulk CPG								
20-2002-01	20-2012-01	20-2022-01	20-2032-01				0.1g	15.00
20-2002-02	20-2012-02	20-2022-02	20-2032-02				0.25g	30.00
20-2002-10	20-2012-10	20-2022-10	20-2032-10				1.0g	105.00

Item	Catalog No.	Pack	Price (\$)
Empty Synthesis Columns-TWIST 40nm, 0.2µm or 1µm	20-0030-00	Pack of 10	60.00
Empty Synthesis Columns - TWIST 10µm/15µm	20-0040-00	Pack of 10	300.00
Replacement Frits - TWIST 10µm/15µm	20-0040-0F	Pack of 20	30.00

Product structures are shown in page 5. TWIST is a trademark of Glen Research Corporation.

AB 3900 POLYSTYRENE MODIFIER COLUMNS

Some of our more popular minor base and modifier supports are available on polystyrene in columns fully compatible with the Applied Biosystems 3900 synthesizer. These include our popular Universal Support III, which will allow DNA, RNA or LNA oligos to be produced on the 3900 with ANY base at the 3' terminus. At the same time, we are offering 1 µmole columns of Universal Support III for the 3900 instrument. Structures and more complete descriptions are found in the relevant catalog sections for each item. AB 3900 columns can be prepared with virtually any of the CPG supports in this catalog. It is no longer necessary to adjust the flow using our AB 3900 CPG columns, as noted in the box to the right. Modified CPG columns are only available in 200 nmole size - simple add 'A' to the regular catalog number to order.

Item	Catalog No.	Pack	Price (\$)
Universal Support III PS			
200 nmole columns	26-5110-52	Pack of 10	100.00
40 nmole columns (AB 3900 Format)	26-5110-55	Pack of 10	100.00
Glen UnySupport™ PS			
200 nmole columns	26-5140-52	Pack of 10	100.00
40 nmole columns	26-5140-55	Pack of 10	100.00
3'-Phosphate PS			
200 nmole columns	26-2900-52	Pack of 10	150.00
40 nmole columns	26-2900-55	Pack of 10	150.00
3'-PT-Amino-Modifier C6 PS			
200 nmole columns	26-2956-52	Pack of 10	220.00
40 nmole columns	26-2956-55	Pack of 10	220.00
3'-(6-FAM) PS			
200 nmole columns	26-2961-52	Pack of 10	300.00
40 nmole columns	26-2961-55	Pack of 10	300.00
3'-DabcyI PS			
200 nmole columns	26-5912-52	Pack of 10	300.00
40 nmole columns	26-5912-55	Pack of 10	300.00
3'-TAMRA PS			
200 nmole columns	26-5910-52	Pack of 10	300.00
40 nmole columns	26-5910-55	Pack of 10	300.00
3'-BiotinTEG PS			
200 nmole columns	26-2955-52	Pack of 10	300.00
40 nmole columns	26-2955-55	Pack of 10	300.00

SEE ALSO
Universal Supports on page 24

AB 3900 1000Å CPG COLUMNS
Glen Research's AB 3900 1000Å CPG columns bring the lower cost of CPG to this platform while maintaining the high synthesis efficiency of 1000Å CPG. Our columns offer the following key attributes:

- No need to change instrument settings
- No need to change software parameters
- Easier handling post-synthesis compared to PS
- High quality 1000Å CPG for optimal synthesis results

STERLING CE PHOSPHORAMIDITES

Glen Research CE (β-cyanoethyl) Phosphoramidites are produced and packaged to ensure the highest performance on DNA synthesizers. Every Glen Research product is accompanied by a Certificate of Analysis and HPLC trace, showing the results of our QC testing. Every Glen Research monomer vial is specially cleaned to eliminate particulate contamination.

Item	Catalog No.	Pack	Price (\$)
dA-CE Phosphoramidite	10-1000-C2	0.25g	12.50
	10-1000-C5	0.5g	25.00
	10-1000-1C	1.0g	50.00
	10-1000-2C	2.0g	100.00
dC-CE Phosphoramidite	10-1010-C2	0.25g	12.50
	10-1010-C5	0.5g	25.00
	10-1010-1C	1.0g	50.00
	10-1010-2C	2.0g	100.00
Ac-dC-CE Phosphoramidite	10-1015-C2	0.25g	12.50
	10-1015-C5	0.5g	25.00
	10-1015-1C	1.0g	50.00
	10-1015-2C	2.0g	100.00
dG-CE Phosphoramidite	10-1020-C2	0.25g	12.50
	10-1020-C5	0.5g	25.00
	10-1020-1C	1.0g	50.00
	10-1020-2C	2.0g	100.00
dmf-dG-CE Phosphoramidite	10-1029-C2	0.25g	12.50
	10-1029-C5	0.5g	25.00
	10-1029-1C	1.0g	50.00
	10-1029-2C	2.0g	100.00
dT-CE Phosphoramidite	10-1030-C2	0.25g	12.50
	10-1030-C5	0.5g	25.00
	10-1030-1C	1.0g	50.00
	10-1030-2C	2.0g	100.00

STERLING SOLVENTS/REAGENTS

All solvents and reagents are prepared to our exacting specifications to ensure the highest synthesis efficiency and are passed through a 0.2 micron filter during packaging to eliminate particulate contamination. Glen Research uses freshly sublimed 1H-tetrazole for premium performance on Expedite synthesizers.

Item	Catalog No.	Pack	Price (\$)
Activator Tetrazole in Acetonitrile	30-3102-66 ¹	60mL	50.00
	30-3102-52 ²	200mL	100.00
	30-3100-57 ²	450mL	200.00
Diluent Acetonitrile, anhydrous	40-4050-45	60mL	12.00
	40-4050-50	100mL	16.00

STERLING SOLVENTS/REAGENTS (CONT.)

Item	Catalog No.	Pack	Price (\$)
Anhydrous Wash Acetonitrile, anhydrous	40-4050-53 ¹	300mL	40.00
	40-4050-57 ²	450mL	50.00
Cap Mix A THF/Ac2O	40-4012-66 ¹	60mL	15.00
	40-4012-52 ²	200mL	30.00
	40-4012-57 ²	450mL	72.00
Cap Mix B 10% 1-Melm in THF/Pyridine	40-4122-66 ¹	60mL	20.00
	40-4122-52 ²	200mL	40.00
	40-4122-57 ²	450mL	96.00
Oxidizing Solution 0.02M I2 in THF/H2O/Pyridine	40-4132-66 ¹	60mL	20.00
	40-4132-52 ²	200mL	40.00
	40-4132-57 ²	450mL	96.00
Deblocking Mix 3% TCA/DCM	40-4140-68 ¹	180mL	18.00
	40-4140-71 ²	1L	80.00

STERLING SUPPORTS

All Glen Research supports use the standard long chain alkylamino (Icaa) linker but differ in the glass pore size, 500Å, 1000Å or 2000Å. The 500Å support is appropriate for shorter sequences, while the 1000Å supports perform better in the synthesis of longer (>30-mer) DNA sequences. The 2000Å support is best for very long (>150-mer) oligonucleotides. We have instituted an additional QC test for supports to show the length of oligo that can be prepared before a drop-off in coupling due to steric effects begins to occur. The drop-off point is recorded in the Certificate of Analysis. All Glen Research supports are fully end-capped to ensure that the CPG surface is totally inert, thereby avoiding the introduction of impurity sequences containing deletions at the 3'-terminus.

Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	dG	dT	dA,dC,dG,dT (1 column of each base)	Ac-dC	dmf-dG		
500Å Columns								
20-2200-42	20-2210-42	20-2220-42	20-2230-42	20-2240-42	20-2213-42		4x0.2µm	40.00
20-2200-41	20-2210-41	20-2220-41	20-2230-41	20-2240-41	20-2213-41		4x1.0µm	60.00
20-2200-14	20-2210-14	20-2220-14	20-2230-14		20-2213-14		1x15µm	150.00
1000Å Columns								
20-2201-45	20-2211-45	20-2221-45	20-2231-45	20-2241-45	20-2215-45	20-2229-45	4x40nm	40.00
20-2201-42	20-2211-42	20-2221-42	20-2231-42	20-2241-42	20-2215-42	20-2229-42	4x0.2µm	40.00
20-2201-41	20-2211-41	20-2221-41	20-2231-41	20-2241-41	20-2215-41	20-2229-41	4x1.0µm	60.00
20-2201-14	20-2211-14	20-2221-14	20-2231-14		20-2215-14	20-2229-14	1x15µm	150.00

ABBREVIATIONS

Ac₂O = Acetic Anhydride
 CE = Cyanoethyl
 CPG = Controlled Pore Glass
 DCM = Dichloromethane
 dmf = dimethylformamide
 I₂ = Iodine
 Icaa = long chain alkylamino
 Melm = 1-Methylimidazole
 µm = micromole(s)
 nm = nanomole(s)
 TCA = Trichloroacetic Acid
 THF = Tetrahydrofuran

SEE ALSO

Alternative Solvents on page 30

BULK CPG LOADING

500Å supports 35-50µmoles/g
 1000Å supports 25-40µmoles/g

QUALITY ASSURANCE

Every batch of these CE Phosphoramidites is tested as follows:

- HPLC**
a) Identity is confirmed by comparison with a reference sample.
b) Purity is determined by HPLC to be ≥98.0%.
- TLC**
Purity is verified by TLC.
- ³¹P NMR**
Purity is determined by ³¹P NMR to be ≥98%.
- Coupling Test**
Coupling efficiency is determined to be ≥99%.
- Solution Test**
A 0.1M solution is determined to be clear and free of particulate contamination.
- Loss on Drying**
Volatile contaminants are determined to be ≤2%.

SEE ALSO

Depurination Resistant dA on page 22

STERLING SUPPORTS (CONT.)

Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	dG	dT	dA,dC,dG,dT (1 column of each base)	Ac-dC	dmf-dG		
2000Å Columns								
20-2202-42	20-2212-42	20-2222-42	20-2232-42	20-2242-42			4x0.2µm	40.00
500Å Bulk CPG								
20-2000-01	20-2010-01	20-2020-01	20-2030-01		20-2013-01		0.1g	9.00
20-2000-02	20-2010-02	20-2020-02	20-2030-02		20-2013-02		0.25g	20.00
20-2000-10	20-2010-10	20-2020-10	20-2030-10		20-2013-10		1.0g	75.00
1000Å Bulk CPG								
20-2001-01	20-2011-01	20-2021-01	20-2031-01		20-2015-01	20-2029-01	0.1g	9.00
20-2001-02	20-2011-02	20-2021-02	20-2031-02		20-2015-02	20-2029-02	0.25g	20.00
20-2001-10	20-2011-10	20-2021-10	20-2031-10		20-2015-10	20-2029-10	1.0g	75.00
2000Å Bulk CPG								
20-2002-01	20-2012-01	20-2022-01	20-2032-01				0.1g	15.00
20-2002-02	20-2012-02	20-2022-02	20-2032-02				0.25g	30.00
20-2002-10	20-2012-10	20-2022-10	20-2032-10				1.0g	105.00

SEE ALSO

Universal Supports on page 24
 Q-Supports on page 27
 High Load Supports on page 29

Item	Catalog No.	Pack	Price (\$)
Empty Synthesis Columns, 40nm, 0.2µm Expedite Style	20-0021-02	Pack of 10	48.00
Empty Synthesis Columns, 1µm Expedite Style	20-0021-01	Pack of 10	48.00
Replacement Filters-Expedite	20-0021-0F	Pack of 20	20.00
Empty Synthesis Columns - TWIST 10µm/15µm	20-0040-00	Pack of 10	300.00
Replacement Frits - TWIST 10µm/15µm	20-0040-0F	Pack of 20	30.00

Product structures are shown in page 5. TWIST is a trademark of Glen Research Corporation.
 Expedite is a trademark of Applied Biosystems.

DNA PHOSPHORAMIDITES - SPECIAL PACKAGING

We offer our high quality DNA phosphoramidites specifically packaged for high throughput and large-scale synthesis customers. These customers normally require high quality materials produced under the guidelines of a validated quality management system while still being priced aggressively. These products include the usual Glen Research certification and guarantees and they are available in larger packs or in bulk. The core catalog numbers for regular DNA phosphoramidites are shown below. For these products, please request a quote.

Item	Catalog No.	Pack	Price (\$)
dA-CE Phosphoramidite	10-1000-SP		
dC-CE Phosphoramidite	10-1010-SP		
Ac-dC-CE Phosphoramidite	10-1015-SP		
dG-CE Phosphoramidite	10-1020-SP		
dmf-dG-CE Phosphoramidite	10-1029-SP		
dT-CE Phosphoramidite	10-1030-SP		

INSTRUMENT TYPES

Glen Research packages these monomers in a variety of industry-standard vials and bottles. Please provide the exact specification of the bottle required prior to receiving a quotation.

STERLING CE PHOSPHORAMIDITES

MerMade synthesizers belong to a family of synthesizers, including the column-based MerMade 4, MerMade 6 and 12 instruments and the parallel array synthesizers, MerMade 192 and MerMade 192E, manufactured by BioAutomation Corporation. Their web site can be found at: <http://www.BioAutomation.com>. Phosphoramidite monomers are packaged in 30mL and 240mL amber bottles for dissolving at a concentration of 1g/20mL and are connected directly to the instrument. Some instruments may also be configured to accept Applied Biosystems serum vials, as shown on page 6.

Item	Catalog No.	Pack	Price (\$)
dA-CE Phosphoramidite	10-1000-02M	0.25g	12.50
	10-1000-05M	0.5g	25.00
	10-1000-10M	1.0g	50.00
	10-1000-5S	5.0g	250.00
	10-1000-1S	10.0g	500.00
dC-CE Phosphoramidite	10-1010-02M	0.25g	12.50
	10-1010-05M	0.5g	25.00
	10-1010-10M	1.0g	50.00
	10-1010-5S	5.0g	250.00
	10-1010-1S	10.0g	500.00
Ac-dC-CE Phosphoramidite	10-1015-02M	0.25g	12.50
	10-1015-05M	0.5g	25.00
	10-1015-10M	1.0g	50.00
	10-1015-5S	5.0g	250.00
	10-1015-1S	10.0g	500.00
dG-CE Phosphoramidite	10-1020-02M	0.25g	12.50
	10-1020-05M	0.5g	25.00
	10-1020-10M	1.0g	50.00
	10-1020-5S	5.0g	250.00
	10-1020-1S	10.0g	500.00
dmf-dG-CE Phosphoramidite	10-1029-02M	0.25g	12.50
	10-1029-05M	0.5g	25.00
	10-1029-10M	1.0g	50.00
	10-1029-5S	5.0g	250.00
	10-1029-1S	10.0g	500.00
dT-CE Phosphoramidite	10-1030-02M	0.25g	12.50
	10-1030-05M	0.5g	25.00
	10-1030-10M	1.0g	50.00
	10-1030-5S	5.0g	250.00
	10-1030-1S	10.0g	500.00

STERLING SOLVENTS/REAGENTS

All solvents and reagents are prepared to our exacting specifications to ensure the highest synthesis efficiency and are passed through a 0.2 micron filter during packaging to eliminate particulate contamination. Parallel synthesizers typically use 5-ethylthio-1H-tetrazole (ETT) as activator to minimize the chance of crystallization. ETT is used at a concentration of 0.25M in acetonitrile, which is far below the level at which crystallization may occur.

Item	Catalog No.	Pack	Price (\$)
Activator 0.25M 5-Ethylthio-1H-Tetrazole in Acetonitrile	30-3140-57	450mL	200.00
	30-3140-61	960mL	365.00
	30-3140-62	2000mL	760.00

QUALITY ASSURANCE

- Every batch of these CE Phosphoramidites is tested as follows:
- HPLC**
 - Identity is confirmed by comparison with a reference sample.
 - Purity is determined by HPLC to be $\geq 98.0\%$.
 - TLC**

Purity is verified by TLC.
 - ^{31}P NMR**

Purity is determined by ^{31}P NMR to be $\geq 98\%$.
 - Coupling Test**

Coupling efficiency is determined to be $\geq 99\%$.
 - Solution Test**

A 0.1M solution is determined to be clear and free of particulate contamination.
 - Loss on Drying**

Volatile contaminants are determined to be $\leq 2\%$.

SEE ALSO

Depurination Resistant dA on page 22

SEE ALSO

Alternative Activators on page 30

STERLING SOLVENTS/REAGENTS (CONT.)

Item	Catalog No.	Pack	Price (\$)
Diluent Acetonitrile, anhydrous	40-4050-50	100mL	16.00
Cap Mix A THF/2,6-Lutidine/Ac2O	40-4010-57	450mL	72.00
	40-4010-61	960mL	154.00
	40-4010-62	2000mL	325.00
Cap Mix B 16% 1-Melm in THF	40-4220-57	450mL	96.00
	40-4220-61	960mL	204.00
	40-4220-62	2000mL	425.00
Oxidizing Solution 0.02M I ₂ in THF/Pyridine/H ₂ O	40-4330-57	450mL	72.00
	40-4330-61	960mL	154.00
	40-4330-62	2000mL	325.00
Deblocking Mix 3% Dichloroacetic acid in DCM	40-4040-57	450mL	36.00
	40-4040-61	960mL	75.00
	40-4040-62	2000mL	144.00
3% TCA/DCM	40-4140-57	450mL	36.00
	40-4140-61	960mL	75.00
	40-4140-62	2000mL	144.00

STERLING SUPPORTS

Columns containing 1000Å CPG are available in packs of 200 to fit MerMade plates. Regular 500Å or 1000Å supports, listed on page 8, may also be used to fill the wells of regular 96 well plates. However, this requires each plate to be prepared with each nucleoside accurately in all wells. A universal support clearly removes the need for four specific supports and makes preparing plates straightforward. Glen UnySupport™ 40 nmole frits, as described on page 22, can also be used.

Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	dG	dT	Ac-dC	dmf-dG		
Mermade 1000Å Columns							
20-2001-65		20-2021-65	20-2031-65	20-2015-65	20-2029-65	200x50nm	750.00
20-2001-62		20-2021-62	20-2031-62	20-2015-62	20-2029-62	200x200nm	750.00
20-2001-61		20-2021-61	20-2031-61	20-2015-61	20-2029-61	48x1.0µm	300.00

Item	Catalog No.	Pack	Price (\$)
Glen UnySupport™ 1000 1 µmole columns 200 nmole columns 40 nmole columns	20-5141-91	Pack of 96	375.00
	20-5141-92	Pack of 96	250.00
	20-5141-95	Pack of 96	250.00
Empty MerMade Columns Empty MerMade Columns (50nm) Empty MerMade Columns (200nm and 1µm)	20-0050-05	Pack of 48	200.00
	20-0050-02	Pack of 48	200.00

ABBREVIATIONS

- Ac₂O = Acetic Anhydride
- CE = Cyanoethyl
- CPG = Controlled Pore Glass
- DCM = Dichloromethane
- dmf = dimethylformamide
- I₂ = Iodine
- Melm = 1-Methylimidazole
- TCA = Trichloroacetic Acid
- THF = Tetrahydrofuran

SEE ALSO

Alternative Solvents on page 30

SEE ALSO

Universal Supports on page 24
Q-Supports on page 27
High Load Supports on page 29

STERLING CE PHOSPHORAMIDITES

Glen Research CE (β-cyanoethyl) Phosphoramidites are produced and packaged to ensure the highest performance on DNA synthesizers. Every Glen Research product is accompanied by a Certificate of Analysis and HPLC trace, showing the results of our QC testing. Every Glen Research monomer vial is specially cleaned to eliminate particulate contamination.

Item	Catalog No.	Pack	Price (\$)
ÅKTA oligopilot dA-CE Phosphoramidite	10-1000-20	2.0g	100.00
	10-1000-50	5.0g	250.00
dC-CE Phosphoramidite	10-1010-20	2.0g	100.00
	10-1010-50	5.0g	250.00
Ac-dC-CE Phosphoramidite	10-1015-20	2.0g	100.00
	10-1015-50	5.0g	250.00
dG-CE Phosphoramidite	10-1020-20	2.0g	100.00
	10-1020-50	5.0g	250.00
dmf-dG-CE Phosphoramidite	10-1029-20	2.0g	100.00
	10-1029-50	5.0g	250.00
dT-CE Phosphoramidite	10-1030-20	2.0g	100.00
	10-1030-50	5.0g	250.00

QUALITY ASSURANCE

Every batch of these CE Phosphoramidites is tested as follows:

- HPLC**
a) Identity is confirmed by comparison with a reference sample.
b) Purity is determined by HPLC to be ≥98.0%.
- TLC**
Purity is verified by TLC.
- ³¹P NMR**
Purity is determined by ³¹P NMR to be ≥98%.
- Coupling Test**
Coupling efficiency is determined to be ≥99%.
- Solution Test**
A 0.1M solution is determined to be clear and free of particulate contamination.
- Loss on Drying**
Volatile contaminants are determined to be ≤2%.

SEE ALSO

Depurination Resistant dA on page 22

STERLING SOLVENTS/REAGENTS

All solvents and reagents are prepared to our exacting specifications to ensure the highest synthesis efficiency and are passed through a 0.2 micron filter during packaging to eliminate particulate contamination.

Item	Catalog No.	Pack	Price (\$)
Diluent Acetonitrile, anhydrous	40-4050-45	60mL	12.00
	40-4050-50	100mL	16.00
ÅKTA oligopilot			
Activator 0.40M Tetrazole in Acetonitrile	30-3105-71	1L	380.00
	Cap Mix A Acetonitrile/Melm		
Cap Mix B* Acetonitrile/Ac2O/Lutidine			40-4028-71
Oxidizing Solution 0.05M I2 in Pyridine/H2O			40-4035-71
Deblocking Mix 3% Dichloroacetic acid in DCM			40-4040-71
3% TCA/DCM			40-4140-71
3% DCA in Toluene			40-4240-71

ABBREVIATIONS

Ac₂O = Acetic Anhydride
 CE = Cyanoethyl
 CPG = Controlled Pore Glass
 DCA = Dichloroacetic Acid
 DCM = Dichloromethane
 I₂ = Iodine
 Melm = 1-Methylimidazole
 μm = micromole(s)

SEE ALSO

Alternative Solvents on page 30

* Cap Mix B is a two part formulation that is combined immediately before shipment.

STERLING CE PHOSPHORAMIDITES

Dr. Oligo synthesizers belong to a family of synthesizers, including the parallel array synthesizers, Dr. Oligo 96, Dr. Oligo 192, Dr. Oligo 384 and Dr. Oligo 768, manufactured by Biolytic® Lab Performance, Inc. in Fremont, CA. Their web site can be found at: <http://www.biolytic.com>. Phosphoramidite monomers are packaged in 30mL and 240mL amber bottles for dissolving at a concentration of 1g/20mL and are connected directly to the instrument. Some instruments may also be configured to accept Applied Biosystems serum vials, as shown on page 6.

Item	Catalog No.	Pack	Price (\$)
dA-CE Phosphoramidite	10-1000-02M	0.25g	12.50
	10-1000-05M	0.5g	25.00
	10-1000-10M	1.0g	50.00
	10-1000-5S	5.0g	250.00
	10-1000-1S	10.0g	500.00
dC-CE Phosphoramidite	10-1010-02M	0.25g	12.50
	10-1010-05M	0.5g	25.00
	10-1010-10M	1.0g	50.00
	10-1010-5S	5.0g	250.00
	10-1010-1S	10.0g	500.00
Ac-dC-CE Phosphoramidite	10-1015-02M	0.25g	12.50
	10-1015-05M	0.5g	25.00
	10-1015-10M	1.0g	50.00
	10-1015-5S	5.0g	250.00
	10-1015-1S	10.0g	500.00
dG-CE Phosphoramidite	10-1020-02M	0.25g	12.50
	10-1020-05M	0.5g	25.00
	10-1020-10M	1.0g	50.00
	10-1020-5S	5.0g	250.00
	10-1020-1S	10.0g	500.00
dmf-dG-CE Phosphoramidite	10-1029-02M	0.25g	12.50
	10-1029-05M	0.5g	25.00
	10-1029-10M	1.0g	50.00
	10-1029-5S	5.0g	250.00
	10-1029-1S	10.0g	500.00
dT-CE Phosphoramidite	10-1030-02M	0.25g	12.50
	10-1030-05M	0.5g	25.00
	10-1030-10M	1.0g	50.00
	10-1030-5S	5.0g	250.00
	10-1030-1S	10.0g	500.00

STERLING SOLVENTS/REAGENTS

All solvents and reagents are prepared to our exacting specifications to ensure the highest synthesis efficiency and are passed through a 0.2 micron filter during packaging to eliminate particulate contamination. Parallel synthesizers typically use 5-ethylthio-1H-tetrazole (ETT) as activator to minimize the chance of crystallization. ETT is used at a concentration of 0.25M in acetonitrile, which is far below the level at which crystallization may occur.

Item	Catalog No.	Pack	Price (\$)
<i>Activator</i> 0.25M 5-Ethylthio-1H-Tetrazole in Acetonitrile	30-3140-57	450mL	200.00
	30-3140-62	2000mL	760.00

STERLING SOLVENTS/REAGENTS (CONT.)

Item	Catalog No.	Pack	Price (\$)
<i>Diluent</i> Acetonitrile, anhydrous	40-4050-50	100mL	16.00
<i>Cap Mix A</i> THF/2,6-Lutidine/Ac2O	40-4010-57	450mL	72.00
	40-4010-62	2000mL	325.00
<i>Cap Mix B</i> 16% 1-Melm in THF	40-4220-57	450mL	96.00
	40-4220-62	2000mL	425.00
<i>Oxidizing Solution</i> 0.02M I2 in THF/Pyridine/H2O	40-4330-57	450mL	72.00
	40-4330-62	2000mL	325.00
<i>Deblocking Mix</i> 3% Dichloroacetic acid in DCM	40-4040-57	450mL	36.00
	40-4040-62	2000mL	144.00
	40-4140-57	450mL	36.00
	40-4140-62	2000mL	144.00

STERLING SUPPORTS

Dr. Oligo instruments are designed for flexibility in the use of supports and columns. They can use fritted plates with loose CPG (page 8) and AB 3900 style polystyrene and CPG columns. Glen UnySupport™ 40 nmole frits can also be used.

Dr. Oligo instruments are designed for flexibility in the use of supports and columns. They can use fritted plates with loose CPG (page 8) and AB 3900 style polystyrene and CPG columns. Glen UnySupport™ 40 nmole frits can also be used.

Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	dG	dT	Ac-dC	dmf-dG		
<i>AB 3900 Polystyrene Columns</i>							
26-2600-65	26-2610-65		26-2630-65		26-2629-65	200x40nm	825.00
26-2600-62	26-2610-62		26-2630-62		26-2629-62	200x200nm	825.00
<i>AB 3900 1000Å CPG Columns</i>							
20-2101-65			20-2131-65	20-2115-65	20-2129-65	200x40nm	600.00
20-2101-62			20-2131-62	20-2115-62	20-2129-62	200x200nm	650.00
20-2101-61			20-2131-61	20-2115-61	20-2129-61	200x1.0µm	875.00

OLIGONUCLEOTIDE PURIFICATION

Biolytic Labs. also offers the innovative Dr. Oligo Processor for high throughput purification of oligonucleotides using Glen-Pak™ DNA Purification Cartridges: <https://www.biolytic.com/p-6814-dr-oligo-processor-fully-automated.aspx>.

QUALITY ASSURANCE

Every batch of these CE Phosphoramidites is tested as follows:

- HPLC**
 - Identity is confirmed by comparison with a reference sample.
 - Purity is determined by HPLC to be ≥98.0%.
- TLC**

Purity is verified by TLC.
- ³¹P NMR**

Purity is determined by ³¹P NMR to be ≥98%.
- Coupling Test**

Coupling efficiency is determined to be ≥99%.
- Solution Test**

A 0.1M solution is determined to be clear and free of particulate contamination.
- Loss on Drying**

Volatile contaminants are determined to be ≤2%.

SEE ALSO

Depurination Resistant dA on page 22

SEE ALSO

Alternative Activators on page 30

ABBREVIATIONS

Ac₂O = Acetic Anhydride
 CE = Cyanoethyl
 CPG = Controlled Pore Glass
 DCM = Dichloromethane
 dmf = dimethylformamide
 I₂ = Iodine
 Melm = 1-Methylimidazole
 TCA = Trichloroacetic Acid
 THF = Tetrahydrofuran

SEE ALSO

Alternative Solvents on page 30

SEE ALSO

Universal Supports on page 24
Q-Supports on page 27
High Load Supports on page 29
Glen-Pak™ DNA on page 147

DEPURINATION RESISTANT CE PHOSPHoramidites

Depurination is defined as the cleavage of the glycosidic bond attaching a purine base to the sugar moiety. Electron withdrawing acyl protecting groups like benzoyl and isobutyryl on the purine amino group(s) destabilize the glycosidic bond, whereas electron donating formamidine protecting groups stabilize the glycosidic bond. The consequence of depurination during oligonucleotide synthesis is the loss of the purine base to form an internucleotide linkage containing the abasic sugar at that position. This site is stable during further synthesis cycles but, upon deprotection with basic reagents, the oligonucleotide is cleaved at that position leading to two shorter fragments. The fragment towards the 5' terminus still contains the DMT group. If DMT-ON purification is being used, the depurinated fragments are co-purified along with the full length product as truncated oligonucleotides.




The most commonly used dA-CE Phosphoramidite containing benzoyl protecting groups suffers substantial degradation by depurination after excessive exposure to TCA. At the same time, two depurination resistant dA monomers, protected with diethylformamidine (def) and dimethylacetamidine (dma), are essentially stable to depurination during the same exposure to TCA.

Both new depurination resistant dA monomers (def and dma protected), were rapidly deprotected in ammonium hydroxide and are fully compatible with regular deprotection strategies. Def-protected-dA was rapidly deprotected with AMA at 65° in 20 minutes, which makes it fully compatible with regular AMA deprotection. In contrast, the dma-protected-dA required 80 minutes with AMA at 65° for complete deprotection.

Dmf-dG is also a depurination resistant CE Phosphoramidite with the isobutyryl group of the original monomer replaced with dimethylformamidine (dmf).

Although depurination does occur in regular oligonucleotide synthesis, the degradation is at an extremely low level. However in certain other circumstances, depurination may become more significant, such as synthesis of long oligos, chip-based synthesis, and large-scale synthesis

Item	Catalog No.	Pack	Price (\$)
def-dA-CE Phosphoramidite	10-1504-02	0.25g	15.00
	10-1504-05	0.5g	30.00
	10-1504-10	1.0g	60.00
dma-dA-CE Phosphoramidite <i>Please inquire.</i>	10-1505		
dmf-dG-CE Phosphoramidite	10-1029-02	0.25g	12.50
	10-1029-05	0.5g	25.00
	10-1029-10	1.0g	50.00
	10-1029-20	2.0g	100.00
	10-1029-40	4.0g	200.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

ULTRAMILD CE PHOSPHoramidites

An alternative protecting scheme for the normal CE phosphoramidites should allow UltraMILD deprotection and should not react with a wider variety of tags and labels. A set of monomers using phenoxyacetyl (Pac) protected dA and 4-isopropylphenoxyacetyl (iPr-Pac) protected dG, along with acetyl protected dC, met the desired criteria for UltraMILD deprotection.

We recommend the use of phenoxyacetic anhydride (Pac₂O) in Cap A. This modification removes the possibility of exchange of the iPr-Pac protecting group on the dG with acetate from the acetic anhydride capping mix. Cleavage and deprotection can be carried out in 2 hours at room temperature with ammonium hydroxide or 4 hours with 0.05M potassium carbonate in methanol.

Item	Catalog No.	Pack	Price (\$)
Pac-dA-CE Phosphoramidite	10-1601-02	0.25g	15.00
	10-1601-05	0.5g	30.00
	10-1601-10	1.0g	60.00
Ac-dC-CE Phosphoramidite	10-1015-02	0.25g	12.50
	10-1015-05	0.5g	25.00
	10-1015-10	1.0g	50.00
iPr-Pac-dG-CE Phosphoramidite	10-1621-02	0.25g	15.00
	10-1621-05	0.5g	30.00
	10-1621-10	1.0g	60.00

ULTRAMILD SUPPORTS

Item	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
	Pac-dA	Ac-dC	iPr-Pac-dG		
UltraMild CPG (Bulk)	20-2601-01	Listed	20-2621-01	0.1g	18.00
	20-2601-02	on	20-2621-02	0.25g	40.00
	20-2601-10	Page 8	20-2621-10	1.0g	150.00
ABI Columns	20-2701-45	20-2115-45	20-2721-45	4X40nm	40.00
	20-2701-42	20-2115-42	20-2721-42	4X0.2µm	40.00
	20-2701-41	20-2115-41	20-2721-41	4X1µm	60.00
	20-2701-13	20-2115-13	20-2721-13	10µm	100.00
	Expedite Columns	20-2801-45	20-2215-45	20-2821-45	4X40nm
20-2801-42		20-2215-42	20-2821-42	4X0.2µm	40.00
20-2801-41		20-2215-41	20-2821-41	4X1µm	60.00
20-2801-14		20-2215-14	20-2821-14	15µm	150.00

ULTRAMILD SOLVENTS/REAGENTS

Item	Catalog No.	Pack	Price (\$)
<i>Cap Mix A</i> THF/Pyridine/Pac ₂ O <i>(Applied Biosystems)</i>	40-4210-52	200mL	140.00
	40-4210-57	450mL	300.00
THF/Pac ₂ O <i>(Expedite)</i>	40-4212-52	200mL	140.00
	40-4212-57	450mL	300.00
<i>Deprotection Solution</i> 0.05M Potassium Carbonate in Methanol	60-4600-30	30mL	30.00

SEE ALSO

Universal Support III on page 26

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.



Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)



GLEN UNYSUPPORT

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M

Columns	
For Instrms type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

REFERENCES


- (1) A.P. Guzaev, and M. Manoharan, *J Am Chem Soc*, 2003, **125**, 2380-2381.
 (2) R.K. Kumar, A.P. Guzaev, C. Rentel, and V.T. Ravikumar, *Tetrahedron*, 2006, **62**, 4528.

ELIMINATION CONDITIONS

Reagent	Conditions
Ammonium hydroxide	80°C/2h 55°C/8h
Ammonium hydroxide/ 40% Methylamine (AMA)	80°C/0.5h 65°C/1h 55°C/8h
Methylamine Gas	65°C/0.5h/30psi
Potassium Carbonate in Methanol	RT/17h
t-Butylamine/Water (1:3 v/v)	60°C/4h

INTELLECTUAL PROPERTY

This product is covered by US Patent 7,202,264 owned by Ionis Pharmaceuticals, Inc..



Glen UnySupport

24

A recent development has been the use of a support based on a molecule which is “conformationally preorganized” to accelerate the dephosphorylation reaction.^{1,2} By using a rigid bicyclic molecule on the support, the rate of elimination is markedly faster than the original Universal Support. The structure of Glen UnySupport™ is shown below. The N-phenyl version, developed at Isis Pharmaceuticals as UnyLinker™, is available from several companies for large scale oligo synthesis. Glen UnySupport is the N-methyl version, which is preferred for high throughput oligonucleotide synthesis since methylamine rather than aniline is formed on deprotection. We are happy to offer Glen UnySupport in a variety of popular formats under license from Ionis Pharmaceuticals.

Item	Catalog No.	Pack	Price(\$)
Bulk Supports			
Glen UnySupport (500Å CPG)	20-5040-01	0.1g	11.00
	20-5040-02	0.25g	25.00
	20-5040-10	1.0g	95.00
Glen UnySupport (1000Å CPG)	20-5041-01	0.1g	11.00
	20-5041-02	0.25g	25.00
	20-5041-10	1.0g	95.00
High Load Glen UnySupport	25-5040-01	0.1g	15.00
	25-5040-02	0.25g	30.00
	25-5040-10	1.0g	115.00
Glen UnySupport PS	26-5040-01	0.1g	16.00
	26-5040-02	0.25g	35.00
	26-5040-10	1.0g	125.00
Columns			
The 1000Å columns and frits below are routinely stocked.			
ABI Format (not LV)			
1 µmole columns	20-5141-41	Pack of 4	60.00
0.2 µmole columns	20-5141-42	Pack of 4	40.00
40 nmole columns	20-5141-45	Pack of 4	40.00
10 µmole column (TWIST Format)	20-5141-13	Pack of 1	100.00
40 nmole frits	20-5441-95	Pack of 96	150.00
Female-Female Luer Adapter for 40 nmole frits	20-0060-00	Pack of 10	20.00
AB 3900 Format			
Glen UnySupport PS			
200 nmole columns	26-5140-52	Pack of 10	100.00
40 nmole columns	26-5140-55	Pack of 10	100.00
Expedite Format			
1 µmole columns	20-5241-41	Pack of 4	60.00
0.2 µmole columns	20-5241-42	Pack of 4	40.00
40 nmole columns	20-5241-45	Pack of 4	40.00
15 µmole column (TWIST Format)	20-5241-14	Pack of 1	150.00
96 Well Format (MerMade, etc.)			
1 µmole columns	20-5141-91	Pack of 96	375.00
200 nmole columns	20-5141-92	Pack of 96	250.00
40 nmole columns	20-5141-95	Pack of 96	250.00

GLEN UNYSUPPORT FC

The extended time required to cleave the succinate linkage of the original Glen UnySupport can be problematical, especially in high-throughput production of oligos, due to the outgassing of ammonia and/or methylamine. This reduction in concentration of gas can necessitate the evaporation of the cleavage solution and addition of fresh Ammonium Hydroxide:MethylAmine 1:1 (AMA) or ammonium hydroxide (NH₄OH) to ensure complete deprotection and dephosphorylation of the product oligos. Using a diglycolate linkage in Glen UnySupport FC instead of the succinate in Glen UnySupport, a significant increase in the rate of cleavage has been achieved. The minimum cleavage times for both versions are as follows:

	AMA	NH ₄ OH
Glen UnySupport	10 min.	40 min.
Glen UnySupport FC	2 min.	5 min.

With the cleavage time of Glen UnySupport FC reduced to less than 5 minutes, there is minimal loss of volatile gas and, therefore, no need to evaporate the cleavage solution and replenish with fresh AMA or ammonium hydroxide solutions.

We offer Glen UnySupport FC attached to 1000Å CPG in a variety of formats suited to high throughput synthesis, as well as in bulk for more routine use.


Item	Catalog No.	Pack	Price(\$)
Bulk Support			
Glen UnySupport FC (1000Å CPG)	22-5041-01	0.1g	11.00
	22-5041-02	0.25g	25.00
	22-5041-10	1.0g	95.00
ABI Format (not LV)			
1 µmole columns	22-5141-41	Pack of 4	60.00
0.2 µmole columns	22-5141-42	Pack of 4	40.00
40 nmole columns	22-5141-45	Pack of 4	40.00
10 µmole column (TWIST Format)	22-5141-13	Pack of 1	100.00
AB 3900 Format			
Glen UnySupport CPG			
200 nmole columns	22-5141-52	Pack of 10	100.00
40 nmole columns	22-5141-55	Pack of 10	100.00
Expedite Format			
1 µmole columns	22-5241-41	Pack of 4	60.00
0.2 µmole columns	22-5241-42	Pack of 4	40.00
40 nmole columns	22-5241-45	Pack of 4	40.00
15 µmole column (TWIST Format)	22-5241-14	Pack of 1	150.00
96 Well Format (MerMade, etc.)			
1 µmole columns	22-5141-91	Pack of 96	375.00
200 nmole columns	22-5141-92	Pack of 96	250.00
40 nmole columns	22-5141-95	Pack of 96	250.00

ELIMINATION CONDITIONS

Reagent	Conditions
Ammonium hydroxide	80°C/2h 55°C/8h
Ammonium hydroxide/ 40% Methylamine (AMA)	80°C/0.5h 65°C/1h 55°C/8h
Methylamine Gas	65°C/0.5h/30psi
Potassium Carbonate in Methanol	RT/17h
t-Butylamine/Water (1:3 v/v)	60°C/4h

INTELLECTUAL PROPERTY

This product is covered by US Patent 7,202,264 owned by Ionis Pharmaceuticals, Inc..



Glen UnySupport FC

25

REFERENCES

- (1) A.V. Azhayev, *Tetrahedron*, 1999, **55**, 787-800.
 (2) A.V. Azhayev and M. Antopolsky, *Tetrahedron*, 2001, **57**, 4977-4986.

INTELLECTUAL PROPERTY

This product is covered by US Patent No.: 6,770,754 and European Patent No.: 1404695.

CLEAVAGE AND DEPROTECTION


1. Cleavage
 For standard and UltraFast deprotection protocols, cleave the oligo from the support using 2M ammonia in methanol at room temperature for 30 minutes. (Only for oligonucleotides greater than 50 nucleotides in length, rinse the support with a further volume of water. Combine the two washes and evaporate to dryness.)

2. Deprotection
Standard
 Add 1 volume of 30% ammonium hydroxide, seal and deprotect using the conditions appropriate for removal of the protecting groups on the nucleobases.


UltraFast
 Add 1 volume of AMA (ammonium hydroxide/40% aqueous methylamine 1:1) seal and deprotect at 65°C for 10 minutes.

UltraMild Using Ammonium Hydroxide
 Add 1 volume of ammonium hydroxide, seal and leave at room temperature for 8 hours.

UltraMild Cleavage and Deprotection Using Potassium Carbonate in Methanol
 Cleave the oligo from the support using 50 mM potassium carbonate in methanol at room temperature for 30 minutes. Seal and leave overnight at room temperature.



Universal Support (1)



Universal Support III (2)

UNIVERSAL SUPPORT III

The key step in the use of any universal support in oligonucleotide synthesis is the dephosphorylation of the 3'-phosphate group to form the desired 3'-hydroxyl group. Azhayev^{1,2} has excelled in the investigation of neighboring group assistance in the dephosphorylation reaction. Amide groups may be considered to be weak N-H acids and can display basic properties in ammonium hydroxide or aqueous methylamine. In the original work^{1,2}, (±)-3-amino-1,2-propanediol was used to form a novel universal support (1). A succinate linker attaches the 3-amino group to the support and the 2-OH is protected with a base-labile group to set up an amide assisted elimination in mildly basic conditions. In this way, the dephosphorylation reaction would eliminate the desired 3'-OH oligonucleotide into solution and the product of any β-elimination competing side reaction would remain bound to the support. A further improvement has been achieved by using a carbamate group to connect the universal linker to the support, as in our product Universal Support III (2). Using Universal Support III, an oligo yield of >80% can be achieved on polymeric supports, with purity equivalent to the same oligo prepared normally.

Conditions for Cleavage and Deprotection are outlined in the table opposite. Universal Support III has been shown to generate oligonucleotides with the same efficacy in polymerase extension reactions as regular oligos. Despite the mild elimination reaction, oligonucleotides up to 75mer in length can be prepared routinely without loss of oligo during the synthesis cycles. This support is also used for the production of siRNA oligos.

Item	Catalog No.	Pack	Price(\$)
Bulk Support			
Universal Support III PS	26-5010-01	0.1g	16.00
	26-5010-02	0.25g	35.00
	26-5010-10	1.0g	125.00
ABI Format (not LV)			
Universal Support III PS			
1 μmole columns	26-5110-41	Pack of 4	60.00
0.2 μmole columns	26-5110-42	Pack of 4	40.00
40 nmole columns	26-5110-45	Pack of 4	40.00
10 μmole column (TWIST Format)	26-5110-13	Pack of 1	100.00
Expedite Format			
1 μmole columns	26-5210-41	Pack of 4	60.00
0.2 μmole columns	26-5210-42	Pack of 4	40.00
40 nmole columns	26-5210-45	Pack of 4	40.00
15 μmole column (TWIST Format)	26-5210-14	Pack of 1	150.00
96 Well Format (MerMade, etc.)			
Universal Support III PS			
1 μmole columns	26-5110-91	Pack of 96	375.00
200 nmole columns	26-5110-92	Pack of 96	250.00
40 nmole columns	26-5110-95	Pack of 96	250.00
AB 3900 Format			
Universal Support III PS			
200 nmole columns	26-5110-52	Pack of 10	100.00
40 nmole columns	26-5110-55	Pack of 10	100.00

Q-SUPPORTS

Oligonucleotides are routinely prepared on supports to which the first nucleoside is attached via a succinate linkage. Over the years, the succinate linkage has demonstrated stability during the synthesis process but has sufficient lability to be cleaved quickly in the deprotection step. However, if the cleavage step is carried out with ammonium hydroxide manually or on the synthesizer, it consumes one hour of precious time while releasing only about 80% of the oligonucleotide. This step is, therefore, a bottleneck in the productivity of many synthesis groups.

Is it possible to find a replacement to the succinate group which offers good stability to the synthesis reagents while offering a much faster cleavage step? The oxalate group has been shown to be very labile during cleavage but its stability to the normal synthesis reagents is not good, requiring changes for successful use. In a practical but elegant study¹ of various bifunctional carboxylic acids, Richard Pon's group identified hydroquinone-O,O'-diacetic acid as the most satisfactory alternative to the succinate group. Nucleosides with this linker arm (Q-linker) are attached to supports with the same ease as the succinyl linker arm.

The cleavage time in ammonium hydroxide at room temperature was found to be 2 minutes, but what about the stability during synthesis? How significant was premature cleavage of oligonucleotide on the synthesizer because of the basic reagents in the capping mixes and oxidizer? Pon showed that the Q-linker is stable to the capping reagents but very slightly labile to the oxidizer (8% cleavage in overnight exposure which would correspond to about 2,000 normal synthesis cycles).

We tested the significance of premature cleavage by preparing sixteen 20mer oligonucleotides on a 0.2 μmole scale, eight with succinate and eight with Q-linkers. The succinate supported oligos were cleaved for 1 hour at room temperature, while those on the Q-support were cleaved for 2 minutes. Both sets were then deprotected normally with ammonium hydroxide. The Q-supports actually gave 5% better yields of product than the succinate supports. Oligo purities were equivalent in both sets.

The Q-linker is absolutely compatible with all hydrolytic cleavage procedures, but especially mild procedures like potassium carbonate in methanol. Pon also showed that it is preferable for RNA supports, improving the cleavage time for 2'-silyl protected nucleoside supports from 2 hours (60-65% cleavage) to 5 minutes (95% cleavage).

We are offering Q-linkers of the four regular nucleosides on 500Å CPG in 0.2 and 1 μmole scales.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type [Add](#)

Expedite E
 MerMade M

Columns
For Instrument type [Add](#)

Expedite E
 Applied Biosystems 3900 A
 MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

REFERENCE

- (1) R.T. Pon and S.Y. Yu, *Tetrahedron Lett*, 1997, **38**, 3327-3330.





Q/SUCCINATE COMPARISON

Q-Support	Succinate
(2 minutes cleavage)	(60 minutes cleavage)
132 ODU*	125 ODU*

**Average crude yield from eight 1 μmole columns deprotected normally.*

Q-SUPPORTS (CONT.)





Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
dA	dC	Ac-dC	dmf-dG	dT		
500Å Bulk Support						
21-2000-01	21-2010-01	21-2013-01	21-2029-01	21-2030-01	0.1g	11.00
21-2000-02	21-2010-02	21-2013-02	21-2029-02	21-2030-02	0.25g	25.00
21-2000-10	21-2010-10	21-2013-10	21-2029-10	21-2030-10	1.0g	95.00
ABI Format (not LV)						
21-2100-41	21-2110-41	21-2113-41	21-2129-41	21-2130-41	4X1μm	60.00
21-2100-42	21-2110-42	21-2113-42	21-2129-42	21-2130-42	4X0.2μm	40.00
Expedite Format						
21-2200-41	21-2210-41	21-2213-41	21-2229-41	21-2230-41	4X1μm	60.00
21-2200-42	21-2210-42	21-2213-42	21-2229-42	21-2230-42	4X0.2μm	40.00



HIGH LOAD CPG

Our high loading support is based on controlled pore silica and it retains the usual 500Å pores. The spacer is also conventional. The only significant difference is the loading which is in the range 80 - 130μmoles/g or about 2.5 times the loading of normal 500Å CPG. Typical loadings for our high load CPG are in the 100 - 120μmoles/g range. As a consequence of the high loading, this support should not be used for sequences longer than 40mers. This high loading support is available in columns for most synthesizers. The 2.5μmole column is identical to our standard 1μmole column (with the exception of the loading). It should be used on occasions when greater than 1μmole is desired but when a 10 or 15μmole synthesis is too high. It should be run using the 1μmole cycle. The 25μmole column is identical to the 10μmole column used on Applied Biosystems synthesizers. It is run using the 10μmole cycle. The 35μmole column is used as an alternative to the 15μmole Expedite column. Again no changes to the standard cycle are recommended. The support is of course available in bulk for use on large-scale synthesizers. A word of caution is in order. When using a column with a higher load than recommended by the instrument manufacturer, there is a much smaller margin for error. All reagents must be fresh and anhydrous diluent and activator must be used. Should you decide to prepare higher-loading columns, ensure that the molar excess of monomer to support nucleoside is at least 5X and preferably 10X.

Item	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Pack	Price(\$)
	dA	dC	dG	dT		
Columns						
(ABI)	25-2100-46	25-2110-46	25-2120-46	25-2130-46	4X2.5μm	75.00
	25-2100-17	25-2110-17	25-2120-17	25-2130-17	1X25μm	125.00
(Expedite)	25-2200-46	25-2210-46	25-2220-46	25-2230-46	4X2.5μm	75.00
	25-2200-18	25-2210-18	25-2220-18	25-2230-18	1X35μm	185.00
Bulk						
	25-2000-02	25-2010-02	25-2020-02	25-2030-02	0.25g	25.00
	25-2000-10	25-2010-10	25-2020-10	25-2030-10	1.0g	90.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

SEE ALSO


Glen UnySupport on page 24

ALTERNATIVE SOLVENTS/REAGENTS


Glen Research offers alternative solvents and reagents in suitable bottles and formulations for use on various DNA synthesizers. All solvents and reagents are prepared to our exacting specifications to ensure the highest coupling efficiencies and are passed through a 0.2 micron filter during packaging to eliminate particulate contamination. Glen Research offers the activators below in powder form for later dissolution in anhydrous acetonitrile or as a prepared solution.

ABBREVIATIONS


Ac₂O = Acetic Anhydride
 DCA = Dichloroacetic Acid
 DCM = Dichloromethane
 DMAP = Dimethylaminopyridine
 I₂ = Iodine
 Melm = 1-Methylimidazole
 TCA = Trichloroacetic Acid
 THF = Tetrahydrofuran




5-Ethylthio-1H-tetrazole



DCI



5-Benzylthio-1H-tetrazole



Saccharin 1-Methylimidazole

INTELLECTUAL PROPERTY

SMI is sold under license from Avecia Biotechnology Inc.

Item	Catalog No.	Pack	Price (\$)
Activator			
5-Ethylthio-1H-tetrazole (ETT)	30-3040-10	1g	35.00
<i>(Dissolve 1g in 31mL anhydrous acetonitrile for a 0.25M solution)</i>	30-3040-20	2g	60.00
	30-3040-25	25g	500.00
0.25M 5-Ethylthio-1H-tetrazole in Acetonitrile <i>(Applied Biosystems)</i>	30-3140-45	45mL	40.00
	30-3140-52	200mL	100.00
	30-3140-57	450mL	200.00
	30-3140-62	2L	760.00
<i>(Expedite)</i>	30-3142-52	200mL	100.00
	30-3140-57	450mL	200.00
4,5-Dicyanoimidazole (DCI), crystalline <i>(Dissolve 1g in 34mL anhydrous acetonitrile for a 0.25M solution)</i>	30-3050-10	1g	35.00
	30-3050-25	25g	500.00
0.25M DCI in Acetonitrile <i>(Applied Biosystems)</i>	30-3150-45	45mL	40.00
	30-3150-52	200mL	100.00
	30-3150-57	450mL	200.00
	30-3150-62	2L	760.00
	30-3152-52	200mL	100.00
	30-3150-57	450mL	200.00
<i>(Expedite)</i>			
5-Benzylthio-1H-tetrazole (BTT) <i>(Dissolve 1g in 21.3mL anhydrous acetonitrile for a 0.25M solution)</i>	30-3070-10	1g	35.00
	30-3070-20	2g	60.00
	30-3070-25	25g	500.00
0.25M 5-Benzylthio-1H-tetrazole in Acetonitrile <i>(Applied Biosystems)</i>	30-3170-45	45mL	40.00
	30-3170-52	200mL	100.00
	30-3170-57	450mL	200.00
	30-3170-62	2L	760.00
<i>(Expedite)</i>	30-3172-52	200mL	100.00
	30-3170-57	450mL	200.00
Saccharin 1-Methylimidazole (SMI) <i>(Dissolve 1g in 31mL anhydrous acetonitrile for a 0.2M solution)</i>	30-3080-10	1g	35.00
	30-3080-20	2g	60.00
	30-3080-25	25g	500.00
0.2M Saccharin 1-Methylimidazole (SMI) in Acetonitrile <i>(Applied Biosystems)</i>	30-3180-45	45mL	40.00
	30-3180-52	200mL	100.00
	30-3180-57	450mL	200.00
	30-3180-62	2L	760.00
<i>(Expedite)</i>	30-3182-52	200mL	100.00
	30-3180-57	450mL	200.00

ALTERNATIVE SOLVENTS/REAGENTS (CONT.)

Item	Catalog No.	Pack	Price (\$)
Cap Mix A			
THF/Lutidine/Ac ₂ O	40-4010-52	200mL	30.00
	40-4010-57	450mL	72.00
	40-4010-62	2L	325.00
THF/Ac ₂ O (9:1)	40-4012-62	2L	275.00
Cap Mix B			
6.5% DMAP in THF <i>(Cap B solutions containing DMAP are preferred by some researchers for preparing long oligos.)</i>	40-4020-52	200mL	42.00
10% Melm in THF	40-4120-52	200mL	30.00
	40-4120-57	450mL	72.00
	40-4120-62	2L	325.00
10% Melm in THF/Pyridine (8:1)	40-4122-62	2L	325.00
Oxidizing Solution			
0.02M I ₂ in THF/Pyridine/H ₂ O	40-4132-62	2L	325.00
Deblocking Mix			
3% DCA/DCM <i>(DCA solutions are more mildly acidic than the TCA equivalents, possibly causing less depurination of dA sites.)</i>	40-4040-57	450mL	36.00
	40-4040-62	2L	144.00
2.5% DCA/DCM	40-4042-57	450mL	36.00
	40-4042-62	2L	144.00

5'-CE PHOSPHORAMIDITES

Glen Research 5'-CE (β -cyanoethyl) Phosphoramidites are designed for the production of 5'-5' or 3'-3' linkages, useful in antisense studies, or to synthesize oligonucleotide segments in the opposite sense from normal synthesis (Reverse Synthesis), for structural studies. These monomers are packaged in ABI-style vials (see note box).

Item	Catalog No.	Pack	Price (\$)
dA-5'-CE Phosphoramidite	10-0001-02	0.25g	75.00
	10-0001-05	0.5g	150.00
	10-0001-10	1.0g	300.00
dC-5'-CE Phosphoramidite	10-0101-02	0.25g	75.00
	10-0101-05	0.5g	150.00
	10-0101-10	1.0g	300.00
dmf-dG-5'-CE Phosphoramidite	10-9201-02	0.25g	75.00
	10-9201-05	0.5g	150.00
	10-9201-10	1.0g	300.00
dT-5'-CE Phosphoramidite	10-0301-02	0.25g	75.00
	10-0301-05	0.5g	150.00
	10-0301-10	1.0g	300.00

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type [Add](#)

Expedite E
MerMade M


Columns
For Instrument type [Add](#)

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)




dA-5'-CE Phosphoramidite



dC-5'-CE Phosphoramidite



Dmf-dG-5'-CE Phosphoramidite




dT-5'-CE Phosphoramidite


5'-SUPPORTS

The following supports are used to produce oligonucleotides with nuclease resistant 3'-3' linkages at the 3' terminus (by attaching regular 3'-CE phosphoramidites) or to produce oligonucleotide sections in the opposite sense (by attaching 5'-CE phosphoramidites). ABI-style columns are supplied unless otherwise requested (see note box).


Item	Catalog No.	Pack	Price (\$)
dA-5'-CPG	20-0002-01	0.1g	50.00
	20-0002-10	1.0g	375.00
	1 μ mole columns	Pack of 4	100.00
	0.2 μ mole columns	Pack of 4	75.00
	10 μ mole column (ABI)	Pack of 1	225.00
15 μ mole column (Expedite)	20-0012-13	Pack of 1	225.00
	20-0012-14	Pack of 1	300.00
dC-5'-CPG	20-0102-01	0.1g	50.00
	20-0102-10	1.0g	375.00
	1 μ mole columns	Pack of 4	100.00
	0.2 μ mole columns	Pack of 4	75.00
	10 μ mole column (ABI)	Pack of 1	225.00
15 μ mole column (Expedite)	20-0112-13	Pack of 1	225.00
	20-0112-14	Pack of 1	300.00
dG-5'-CPG	20-0202-01	0.1g	50.00
	20-0202-10	1.0g	375.00
	1 μ mole columns	Pack of 4	100.00
	0.2 μ mole columns	Pack of 4	75.00
	10 μ mole column (ABI)	Pack of 1	225.00
15 μ mole column (Expedite)	20-0212-13	Pack of 1	225.00
	20-0212-14	Pack of 1	300.00
dT-5'-CPG	20-0302-01	0.1g	50.00
	20-0302-10	1.0g	375.00
	1 μ mole columns	Pack of 4	100.00
	0.2 μ mole columns	Pack of 4	75.00
	10 μ mole column (ABI)	Pack of 1	225.00
15 μ mole column (Expedite)	20-0312-13	Pack of 1	225.00
	20-0312-14	Pack of 1	300.00




dA-5'-CPG



dC-5'-CPG



dG-5'-CPG




dT-5'-CPG


METHYL PHOSPHONAMIDITES

Methyl Phosphonamidites may be used in DNA synthesizers following conventional CE Phosphoramidite protocols to produce oligonucleotides containing one or more methyl phosphonate linkages. However, deprotection and purification techniques differ and a description of the procedures is included in the Technical Bulletin. We also offer the dC monomer with acetyl base protection.¹ This protecting group is removed with ammonium hydroxide during the cleavage step, eliminating modification at the dC sites during the deprotection step using ethylenediamine in ethanol.


Item	Catalog No.	Pack	Price (\$)
dA-Me Phosphonamidite	10-1100-02	0.25g	50.00
	10-1100-05	0.5g	100.00
Ac-dC-Me Phosphonamidite	10-1115-02	0.25g	50.00
	10-1115-05	0.5g	100.00
dG-Me Phosphonamidite	10-1120-02	0.25g	50.00
	10-1120-05	0.5g	100.00
dT-Me Phosphonamidite	10-1130-02	0.25g	50.00
	10-1130-05	0.5g	100.00




dA-Me Phosphonamidite



Ac-dC-Me Phosphonamidite



dG-Me Phosphonamidite




dT-Me Phosphonamidite

PACE PHOSPHORAMIDITES


Phosphonoacetate (PACE) modified oligonucleotides show great potential as biological modifiers in a wide variety of research applications. PACE monomers are part of a family of Phosphonocarboxylate monomers. The monomers can be easily incorporated into complex oligonucleotides and are compatible with a wide variety of other sugar or heterobase modifications. PACE DNA can be conjugated through the carboxylic acid functional group. They have been shown to be active in siRNA duplexes and accelerate the initial rate of cleavage by RNase H-1 when incorporated with phosphorothioates. However, the most interesting observation to date is that they exhibit an unprecedented enhancement in penetration of cultured cells.

PACE monomers are fully soluble in acetonitrile at a recommended concentration of 0.1M and are compatible with standard DNA synthesizers. As an optimal cycle, we recommend using DCI as an activator (30-3150-XX) and a 15 minute coupling time. Following coupling, cap using Unicap (10-4410-XX) with a regular coupling time and then oxidize using 0.5 M CSO for 3 minutes. Alternatively, a 33 minute coupling time using 0.45 M tetrazole, oxidation using low-water iodine (40-4032-XX) followed by capping with 6.5% DMAP as Cap B will give acceptable results. For deprotection, pre-treat the synthesis column with 1.5% DBU in anhydrous acetonitrile for 60 minutes at room temperature to remove 1,1-dimethyl-2-cyanoethyl protecting groups. Rinse the column with acetonitrile, dry under argon and complete the deprotection with 40% aqueous methylamine for 2 hours at room temperature.


Item	Catalog No.	Pack	Price(\$)
dA-PACE Phosphoramidite	10-1140-02	0.25g	100.00
	10-1140-05	0.5g	200.00
	10-1140-10	1.0g	400.00
Ac-dC-PACE Phosphoramidite	10-1150-02	0.25g	100.00
	10-1150-05	0.5g	200.00
	10-1150-10	1.0g	400.00
dG-PACE Phosphoramidite	10-1160-02	0.25g	100.00
	10-1160-05	0.5g	200.00
	10-1160-10	1.0g	400.00
dT-PACE Phosphoramidite	10-1170-02	0.25g	100.00
	10-1170-05	0.5g	200.00
	10-1170-10	1.0g	400.00




dA-PACE Phosphoramidite



Ac-dC-PACE Phosphoramidite



dG-PACE Phosphoramidite



dT-PACE Phosphoramidite

INTELLECTUAL PROPERTY

These products are covered by patents, US 6,693,187 and 7,067,641, and patents pending owned by Metasense Technologies. Purchase of all or any of these products includes a limited license to use the products solely for the manufacture of oligonucleotides for research use only. This license specifically excludes the use of the product or oligonucleotides containing the product for: (a) therapeutic or diagnostic applications (including kits, pools, libraries and other products or services that incorporate oligonucleotides containing the product), (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), or (c) resale (including sale of kits, pools, libraries and other products or services that incorporate the product or oligonucleotides containing the product). If such activities have commercial application, a separate license is required from Metasense Technologies. Neither the product nor any product created through its use may be used in human clinical trials.

A simple agreement must be signed before end-users and custom oligo services may purchase these products for use as defined above.
<http://www.glenresearch.com/Reference/PACE.pdf>

SEE ALSO

DCI on page 30
 UniCap on page 32
 0.5M CSO on page 32
 2'-OMe-PACE on page 145


METHYL PHOSPHORAMIDITES

For many years, Glen Research has supplied methyl phosphoramidites in addition to β-cyanoethyl (CE) phosphoramidites for the few situations where the more labile cyanoethyl group is not an advantage. Some of our customers, probably remembering that the methyl group was removed specifically with thiophenol, have tried to use these monomers to prepare the interesting, uncharged, and nuclease-resistant methyl phosphotriester linkage. Unfortunately, this linkage is labile to ammonium hydroxide and the regular phosphodiester linkage is formed (along with a small amount of chain scission). We offer UltraMild methyl phosphoramidites for this application. Oligos produced from these monomers can be deprotected with potassium carbonate in methanol to produce methyl phosphotriester linkages. Since these linkages are diastereomeric and uncharged, the oligos may be hard to handle. Consequently, it is likely that chimeras will be produced using these monomers along with the regular UltraMild CE phosphoramidites. If many dG residues are included in the oligonucleotide, we recommend the use of phenoxyacetic anhydride (Pac₂O) in Cap A. This modification removes the possibility of exchange of the isopropyl-phenoxyacetate (iPr-Pac) protecting group on the dG with acetate from the acetic anhydride capping mix.


Item	Catalog No.	Pack	Price(\$)
Pac-dA-Me Phosphoramidite	10-1301-02	0.25g	25.00
	10-1301-05	0.5g	50.00
	10-1301-10	1.0g	100.00
Ac-dC-Me Phosphoramidite	10-1315-02	0.25g	25.00
	10-1315-05	0.5g	50.00
	10-1315-10	1.0g	100.00
iPr-Pac-dG-Me Phosphoramidite	10-1321-02	0.25g	25.00
	10-1321-05	0.5g	50.00
	10-1321-10	1.0g	100.00
dT-Me Phosphoramidite	10-1330-02	0.25g	25.00
	10-1330-05	0.5g	50.00
	10-1330-10	1.0g	100.00

ULTRAMILD SOLVENTS/REAGENTS


Item	Catalog No.	Pack	Price (\$)
Cap Mix A THF/Pyridine/Pac ₂ O (Applied Biosystems)	40-4210-52	200mL	140.00
	40-4210-57	450mL	300.00
THF/Pac ₂ O (Expedite)	40-4212-52	200mL	140.00
	40-4212-57	450mL	300.00
Deprotection Solution 0.05M Potassium Carbonate in Methanol	60-4600-30	30mL	30.00




Pac-dA-Me Phosphoramidite



Ac-dC-Me Phosphoramidite



iPr-Pac-dG-Me Phosphoramidite



dT-Me Phosphoramidite

H-PHOSPHONATE MONOMERS


Glen Research H-Phosphonates are analyzed by HPLC and are synthesis-tested. H-Phosphonates are especially useful for the preparation of modified internucleotide linkages which are unattainable by phosphoramidite chemistry. The most popular application is the preparation of radiolabeled phosphorothioates, since the sulfurization reaction is carried out off the synthesizer. These monomers are packaged in ABI-style vials (see note box).

Item	Catalog No.	Pack	Price(\$)
dA-H-Phosphonate, TEA Salt	10-1200-02	0.25g	40.00
	10-1200-05	0.5g	80.00
dC-H-Phosphonate, DBU Salt	10-1210-02	0.25g	40.00
	10-1210-05	0.5g	80.00
dG-H-Phosphonate, TEA Salt	10-1220-02	0.25g	40.00
	10-1220-05	0.5g	80.00
dT-H-Phosphonate, TEA Salt	10-1230-02	0.25g	40.00
	10-1230-05	0.5g	80.00


H-PHOSPHONATE REAGENTS

Our H-Phosphonate solvents and reagents have been discontinued. H-Phosphonate reagents are easily prepared using high purity products and the formulations shown below.


Item	Catalog No.	Pack	Price (\$)
1-Adamantanecarbonyl chloride is available from Aldrich, Catalog No. 117722. Dilute to 0.1M. (Activator for monomers and capping reagent)			
Acetonitrile/Pyridine (50:50), anhydrous (Monomer Diluent)			
Acetonitrile/Pyridine (95:5), anhydrous (Activator Diluent)			
1% Isopropyl Phosphite in Acetonitrile/Pyridine (50:50) (Capping Reagent)			
Acetonitrile/Pyridine (50:50) (Neutralizer and Wash Solvent)			
4% I ₂ in Pyridine/H ₂ O/THF (10:10:80) THF/H ₂ O/TEA (80:10:10) (Both reagents are required for oxidation of H-phosphonate linkages)			




dA-H-Phosphonate



dC-H-Phosphonate



dG-H-Phosphonate



dT-H-Phosphonate

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

ABBREVIATIONS

I₂ = Iodine
TEA = Triethylamine
THF = Tetrahydrofuran

THIOPHOSPHORAMIDITES

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)





Replacing two non-bridging oxygen atoms with sulfur atoms in a DNA phosphodiester linkage creates a phosphorodithioate (PS2) linkage.¹ Like natural DNA, the phosphorodithioate linkage is achiral at phosphorus. This analog is completely resistant to nuclease degradation and forms complexes with DNA and RNA with somewhat reduced stabilities.² Moreover, it has been found that PS2-ODNs bind proteins with a higher affinity than their phosphodiester analogues²⁻⁶ suggesting that PS2-ODNs may have additional utility in the form of sulfur-modified phosphate ester aptamers (thioaptamers)^{3,6-8} for therapeutic and diagnostic applications. Thiophosphoramidites are now commercially available after recent work at AM Biotechnologies (<http://www.thioaptamer.com/>).

- 1) Thiophosphoramidites (ThioPAs) are not soluble in anhydrous acetonitrile diluent. Rather, 10% DCM (v/v) in acetonitrile is an ideal diluent for all four of the thioPAs for a final amidite concentration of 0.15 M.
- 2) ThioPAs are somewhat less stable than normal DNA phosphoramidites in anhydrous acetonitrile containing 10% DCM; however, the coupling efficiency of all four thioPAs is not reduced after two days in solution at room temperature.
- 3) After synthesis, the thioPA bottle on the synthesizer should be replaced with one containing acetonitrile diluent and the synthesizer line flushed with acetonitrile.

A typical cycle for the solid-phase synthesis of a PS2 linkage is different from a standard cycle for the synthesis of normal phosphate linkages. After coupling, the resulting thiophosphite triester is then sulfurized with DDTT. Capping is carried out AFTER sulfurization.

Upon completion of the automated synthesis, deprotection is carried out using a concentrated ammonia:ethanol (3:1, v:v) mix containing 20 mM DTT at 55 °C for 15-16 h.

Item	Catalog No.	Pack	Price(\$)
dA-Thiophosphoramidite	10-1700-90	100 μmole	150.00
	10-1700-02	0.25g	360.00
dC-Thiophosphoramidite	10-1710-90	100 μmole	150.00
	10-1710-02	0.25g	360.00
dG-Thiophosphoramidite	10-1720-90	100 μmole	150.00
	10-1720-02	0.25g	360.00
dT-Thiophosphoramidite	10-1730-90	100 μmole	150.00
	10-1730-02	0.25g	360.00







LOCKED ANALOG PHOSPHORAMIDITES

Locked Nucleic Acid (LNA) was first described by Wengel and co-workers in 1998¹ as a novel class of conformationally restricted oligonucleotide analogues. LNA is a bicyclic nucleic acid where a ribonucleoside is linked between the 2'-oxygen and the 4'-carbon atoms with a methylene unit. Oligonucleotides containing LNA exhibit unprecedented thermal stabilities towards complementary DNA and RNA², which allows excellent mismatch discrimination. In fact, the high binding affinity of LNA oligos allows for the use of short probes in, for example, SNP genotyping³, allele specific PCR and mRNA sample preparation. LNA is recommended for use in any hybridization assay that requires high specificity and/or reproducibility, e.g., dual labelled probes, *in situ* hybridization probes, molecular beacons and PCR primers. Furthermore, LNA offers the possibility to adjust T_m values of primers and probes in multiplex assays. LNA can be mixed with DNA and RNA, as well as other nucleic acid analogues, modifiers and labels. LNA oligonucleotides are water soluble, and can be separated by gel electrophoresis and precipitated by ethanol.

Glen Research is pleased to offer these highly useful reagents - Locked Analog (LA) Phosphoramidites - as tools for this technology.

Item	Catalog No.	Pack	Price(\$)
Bz-A-LA-CE Phosphoramidite	10-2000-05	0.5g	75.00
	10-2000-10	1.0g	150.00
5-Me-Bz-C-LA-CE Phosphoramidite	10-2011-05	0.5g	75.00
	10-2011-10	1.0g	150.00
dmf-G-LA-CE Phosphoramidite	10-2029-05	0.5g	75.00
	10-2029-10	1.0g	150.00
T-LA-CE Phosphoramidite	10-2030-05	0.5g	75.00
	10-2030-10	1.0g	150.00



REFERENCES

- (1) A.A. Koshkin, S.K. Singh, P. Nielsen, V.K. Rajwanshi, R. Kumar, M. Meldgaard, C.E. Olsen, and J. Wengel, *Tetrahedron*, 1998, **54**, 3607-3630.
- (1b) S.K. Singh, P. Nielsen, A.A. Koshkin, and J. Wengel, *Chem. Comm.*, 1998, (4), 455-456.
- (2) L. Kvaernø and J. Wengel, *Chem. Comm.*, 1999, (7), 657-658.
- (3) P. Mouritzen, A.T. Nielsen, H.M. Pfrundheller, Y. Choleva, L. Kongsbak, and S. Møller, *Expert Review of Molecular Diagnostics*, 2003, **3**(1), 27-38.

REFERENCES

- (1) J. Nielsen, W.K.D. Brill, and M.H. Caruthers, *Tetrahedron Letters*, 1988, **29**, 2911-2914.
- (2) L. Cummins, D. Graff, G. Beaton, W.S. Marshall, and M.H. Caruthers, *Biochemistry*, 1996, **35**, 8734-41.
- (3) X. Yang, and D.G. Gorenstein, *Curr Drug Targets*, 2004, **5**, 705-15.
- (4) W.S. Marshall, and M.H. Caruthers, *Science*, 1993, **259**, 1564-70.
- (5) J.L. Tonkinson, et al., *Antisense Research and Development*, 1994, **4**, 269-278.
- (6) X. Yang, et al., *Bioorg Med Chem Lett*, 1999, **9**, 3357-62.
- (7) X. Yang, et al., *Ann N Y Acad Sci*, 2006, **1082**, 116-9.
- (8) X. Yang, et al., *Nucleic Acids Res*, 2002, **30**, e132.

SEE ALSO

2'-OMe-RNA
Thiophosphoramidites on page 168

TRIMER PHOSPHoramidites

Trimer phosphoramidites¹⁻⁴ have proven to be extremely valuable because they allow codon-based mutagenesis, which circumvents the common problems of codon-bias, frame-shift mutations, and the introduction of nonsense or stop codons.⁵ This is accomplished by introducing a mixture of all 20 amino acid codons (or subset thereof) at any location within the sequenced to be mutated. This leads to the production of clonal libraries of exceptional diversity with order-of-magnitude increases in amino acid sequence variance while either maintaining a uniform amino acid distribution⁶ or one that is biased toward a desired set of amino acids.⁷

However, difficulties arise when trying to introduce mutations in multiple distal regions of a gene simultaneously. The synthesis of long oligonucleotides is required, which inevitably leads to lower sequence fidelity due to deletion mutants, depurination events and, to a lesser extent, mutations arising from deamination of cytidine, for example.

An elegant solution to this problem is the use of Antisense Trimer Phosphoramidites. These trimers are the reverse complement of the canonical 'sense' codons. When these antisense codons are put into the noncoding strand of a template DNA and amplified by PCR, they will code for the sense codon in the opposite strand of DNA. This allows the powerful technique of PCR Assembly⁸ to generate not only kilobase-sized genes from short 50mer oligonucleotides, but to simultaneously mutate multiple distal regions of that gene, as shown in Figure 1.

The sense and their corresponding antisense codons are listed in Table 1. Conveniently, many of our existing sense trimers can act as antisense codons. For example, AAC, which codes for asparagine, has the anticodon GTT, which is the sense codon for valine. However, some of the existing trimers, while they can act as an antisense codon, are not good choices for use. For example, TGG, which codes for tryptophan, could be used as an antisense codon for proline because CCA is one of proline's synonymous codons. However, CCA has a relatively low Codon Adaptation Index (CAI) value⁹ in *E. coli*, which could limit protein expression in that commonly used organism. For this reason, the anticodon CGG was chosen for optimal expression in *E. coli*, as were the other new antisense codons shown in bold in Table 1.


Included in Table 1 are the reaction factors (RFs) for each of the sense and antisense trimers. The reaction factor is critical since the trimers will likely be mixed and they exhibit different rates of reaction when coupling during oligonucleotide synthesis. An example where the RF is used to compensate for differing rates of coupling follows. The RF for AAC is 1.0 and for TAC is 1.6. Therefore, 1.6 equivalents of TAC are needed for every 1.0 equivalent of AAC for equal coupling rates. So to obtain 25 umoles of trimer mix that yields, on average, a 1:1 ratio of AAC/TAC at the mutation site, 9.6 umoles of AAC would be added to 15.4 umoles of TAC.

All of the trimers are available individually so the researchers can prepare custom trimer mixes. Two pre-made catalog trimer mixes are available: 13-1991-xx, for incorporating all 20 amino acid codons equally into a sequence and 13-1992-xx, for incorporating 19 amino acid codons (-Cys). For a custom trimer mix of a particular subset of codons or a trimer mix that represents a set of trimers that is biased toward a particular codon or codons, please contact support@glenresearch.com for a quotation and projected delivery date.

There is a concern that the sequence of the trimers has to be verified. For example, CAT coding for histidine, has to be differentiated from TAC, coding for tyrosine. These two trimers have virtually identical lipophilicity and their identity cannot be clearly confirmed by HPLC. This problem has been solved⁴ using HPLC electrospray mass spectrometric analysis of the trimers, which provides data confirming molecular weight and sequence.

REFERENCES

- (1) A.L. Kayushin, M.D. Korosteleva, A.I. Miroshnikov, W. Kosch, D. Zubov, and N. Piel, *Nucleic Acids Research*, 1996, **24**, 3748-3755.
- (2) A. Kayushin, et al., *Nucleos Nuclcot*, 1999, **18**, 1531-1533.
- (3) A. Kayushin, M. Korosteleva, and A. Miroshnikov, *Nucleos Nuclcot Nucleic Acids*, 2000, **19**, 1967-1976.
- (4) T. Mauriala, S. Auriola, A. Azhayev, A. Kayushin, M. Korosteleva, and A. Miroshnikov, *J Pharm Biomed Anal*, 2004, **34**, 199-206.
- (5) C. Neylon, *Nucleic Acids Res*, 2004, **32**, 1448-59.
- (6) L.R. Krumpe, K.M. Schumacher, J.B. McMahon, L. Makowski, and T. Mori, *BMC Biotechnol*, 2007, **7**, 65-72.
- (7) F.A. Fellouse, et al., *J Mol Biol*, 2007, **373**, 924-40.
- (8) W.P. Stemmer, A. Cramer, K.D. Ha, T.M. Brennan, and H.L. Heyneker, *Gene*, 1995, **164**, 49-53.
- (9) P.M. Sharp, and W.H. Li, *Nucleic Acids Res*, 1987, **15**, 1281-95.



General Structure of Trimer Phosphoramidites, where B=A⁺, C⁺, G⁺, T

Figure 1: Simultaneous Mutation of Multiple Distal Regions of Gene





TABLE 1: RF of Trimer Phosphoramidites

Sense codons (5'→3')	Reaction Factor (RF)	Antisense codons (5'→3')	Reaction Factor (RF)
AAA (Lys)	1.10	TTT	1.70
AAC (Asn)	1.00	GTT	1.90
ACT (Thr)	1.60	GGT	1.10
ATC (Ile)	1.50	GAT	1.40
ATG (Met)	1.30	CAT	1.30
CAG (Gln)	2.00	CTG	1.20
CAT (His)	1.30	ATG	1.30
CCG (Pro)	1.80	CGG	0.80
CGT (Arg)	1.40	GCG	0.60
CTG (Leu)	1.20	CAG	2.00
GAA (Glu)	1.40	TTC	1.30
GAC (Asp)	1.60	ATC	1.50
GCT (Ala)	1.50	TGC	1.50
GGT (Gly)	1.10	ACC	0.90
GTT (Val)	1.90	AAC	1.00
TAC (Tyr)	1.60	GTA	1.50
TCT (Ser)	1.30	AGA	1.40
TGC (Cys)	1.50	GCA	1.00
TGG (Trp)	1.10	CCA	1.10
TTC (Phe)	1.30	GAA	1.40



ATC Trimer

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type [Add](#)

Expedite E
MerMade M

Columns
For Instrument type [Add](#)

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

Item	Catalog No.	Pack	Price(\$)
Sense Trimers			
AAA Trimer Phosphoramidite (Lys)	13-1000-95 13-1000-90	50 µm 100 µm	385.00 770.00
AAC Trimer Phosphoramidite (Asn)	13-1001-95 13-1001-90	50 µm 100 µm	385.00 770.00
ACT Trimer Phosphoramidite (Thr)	13-1013-95 13-1013-90	50 µm 100 µm	385.00 770.00
ATC Trimer Phosphoramidite (Ile)	13-1031-95 13-1031-90	50 µm 100 µm	385.00 770.00
ATG Trimer Phosphoramidite (Met)	13-1032-95 13-1032-90	50 µm 100 µm	385.00 770.00
CAG Trimer Phosphoramidite (Gln)	13-1102-95 13-1102-90	50 µm 100 µm	385.00 770.00
CAT Trimer Phosphoramidite (His)	13-1103-95 13-1103-90	50 µm 100 µm	385.00 770.00
CCG Trimer Phosphoramidite (Pro)	13-1112-95 13-1112-90	50 µm 100 µm	385.00 770.00
CGT Trimer Phosphoramidite (Arg)	13-1123-95 13-1123-90	50 µm 100 µm	385.00 770.00
CTG Trimer Phosphoramidite (Leu)	13-1132-95 13-1132-90	50 µm 100 µm	385.00 770.00
GAA Trimer Phosphoramidite (Glu)	13-1200-95 13-1200-90	50 µm 100 µm	385.00 770.00
GAC Trimer Phosphoramidite (Asp)	13-1201-95 13-1201-90	50 µm 100 µm	385.00 770.00
GCT Trimer Phosphoramidite (Ala)	13-1213-95 13-1213-90	50 µm 100 µm	385.00 770.00
GGT Trimer Phosphoramidite (Gly)	13-1223-95 13-1223-90	50 µm 100 µm	385.00 770.00
GTT Trimer Phosphoramidite (Val)	13-1233-95 13-1233-90	50 µm 100 µm	385.00 770.00
TAC Trimer Phosphoramidite (Tyr)	13-1301-95 13-1301-90	50 µm 100 µm	385.00 770.00
TCT Trimer Phosphoramidite (Ser)	13-1313-95 13-1313-90	50 µm 100 µm	385.00 770.00
TGC Trimer Phosphoramidite (Cys)	13-1321-95 13-1321-90	50 µm 100 µm	385.00 770.00
TGG Trimer Phosphoramidite (Trp)	13-1322-95 13-1322-90	50 µm 100 µm	385.00 770.00
TTC Trimer Phosphoramidite (Phe)	13-1331-95 13-1331-90	50 µm 100 µm	385.00 770.00
Trimer Phosphoramidite Mix 1 (Mix of above 20 trimers)	13-1991-95 13-1991-90	50 µm 100 µm	570.00 1140.00
Trimer Phosphoramidite Mix 2 (Mix of above 20 trimers less TGC-Cys)	13-1992-95 13-1992-90	50 µm 100 µm	570.00 1140.00

Item	Catalog No.	Pack	Price(\$)
Antisense Trimers			
AAC Trimer Phosphoramidite (Anti Val)	13-1001-95 13-1001-90	50 µm 100 µm	385.00 770.00
ACC Trimer Phosphoramidite (Anti Gly)	13-1011-95 13-1011-90	50 µm 100 µm	385.00 770.00
AGA Trimer Phosphoramidite (Anti Ser)	13-1020-95 13-1020-90	50 µm 100 µm	385.00 770.00
ATC Trimer Phosphoramidite (Anti Asp)	13-1031-95 13-1031-90	50 µm 100 µm	385.00 770.00
ATG Trimer Phosphoramidite (Anti His)	13-1032-95 13-1032-90	50 µm 100 µm	385.00 770.00
CAG Trimer Phosphoramidite (Anti Leu)	13-1102-95 13-1102-90	50 µm 100 µm	385.00 770.00
CAT Trimer Phosphoramidite (Anti Met)	13-1103-95 13-1103-90	50 µm 100 µm	385.00 770.00
CCA Trimer Phosphoramidite (Anti Trp)	13-1110-95 13-1110-90	50 µm 100 µm	385.00 770.00
CGG Trimer Phosphoramidite (Anti Pro)	13-1122-95 13-1122-90	50 µm 100 µm	385.00 770.00
GAA Trimer Phosphoramidite (Anti Phe)	13-1200-95 13-1200-90	50 µm 100 µm	385.00 770.00
GAT Trimer Phosphoramidite (Anti Ile)	13-1203-95 13-1203-90	50 µm 100 µm	385.00 770.00
GCA Trimer Phosphoramidite (Anti Cys)	13-1210-95 13-1210-90	50 µm 100 µm	385.00 770.00
GCG Trimer Phosphoramidite (Anti Arg)	13-1212-95 13-1212-90	50 µm 100 µm	385.00 770.00
GGT Trimer Phosphoramidite (Anti Thr)	13-1223-95 13-1223-90	50 µm 100 µm	385.00 770.00
GTA Trimer Phosphoramidite (Anti Tyr)	13-1230-95 13-1230-90	50 µm 100 µm	385.00 770.00
TGC Trimer Phosphoramidite (Anti Ala)	13-1321-95 13-1321-90	50 µm 100 µm	385.00 770.00
TTC Trimer Phosphoramidite (Anti Glu)	13-1331-95 13-1331-90	50 µm 100 µm	385.00 770.00
TTT Trimer Phosphoramidite (Anti Lys)	13-1333-95 13-1333-90	50 µm 100 µm	385.00 770.00


BASES AFFECTING DUPLEX STABILITY

Substitution of C-5 propynyl-dC (pdC) for dC and C-5 propynyl-dU (pdU) for dT are effective strategies to enhance base pairing. Using these base substitutions, duplex stability and melting temperatures are raised by the following amounts: C-5 propynyl-C 2.8° per substitution; C-5 propynyl-U 1.7° per substitution. AP-dC (G-clamp) substitutes for dC and is another very important modified nucleoside that enhances hybridization by 7-21° per substitution depending upon the sequence and location of the AP-dC. The ability of these modified bases to enhance binding while maintaining specificity has proven useful in antisense research and in the synthesis of high affinity probes. AP-dC is also a fluorescent nucleoside and should find uses in DNA structural research.

Item	Catalog No.	Pack	Price(\$)
pdC-CE Phosphoramidite	10-1014-90	100 μmole	85.00
	10-1014-02	0.25g	245.00
	10-1014-05	0.5g	490.00
pdU-CE Phosphoramidite	10-1054-90	100 μmole	65.00
	10-1054-02	0.25g	195.00
	10-1054-05	0.5g	390.00
AP-dC-CE Phosphoramidite (G-Clamp)	10-1097-95	50 μmole	230.00
	10-1097-90	100 μmole	460.00
	10-1097-02	0.25g	1175.00

C-5 methyl pyrimidine nucleosides are known to stabilize duplexes relative to the non-methylated bases. Therefore, enhanced binding can be achieved using 5-methyl-dC in place of dC, duplex melting temperature being increased by 1.3°. Ac-5-Me-dC-CE Phosphoramidite is fully compatible with AMA deprotection and none of the N4-Me transamination mutant is observed on deprotection.

Item	Catalog No.	Pack	Price(\$)
5-Me-dC-CE Phosphoramidite	10-1060-90	100 μmole	50.00
	10-1060-02	0.25g	120.00
Ac-5-Me-dC-CE Phosphoramidite	10-1560-90	100 μmole	50.00
	10-1560-02	0.25g	120.00



pdC

pdU

AP-dC

5-Me-dC

Ac-5-Me-dC

BASES AFFECTING DUPLEX STABILITY (CONT.)


The simplest approach to the design of high affinity primers and probes is to substitute A sites with 2-amino-A, since the 2-amino-A-T base pair is equivalent in strength to the G-T base pair. 2-Amino-A also destabilizes A-G wobble mismatches, thus increasing specificity. In 1998, we introduced a 2-amino-dA monomer which exhibits fast and effective deprotection in ammonium hydroxide and it is stabilized to depurination during synthesis. We now recommend the use of 0.5 M CSO in anhydrous acetonitrile (40-4632-xx) for best results with multiple additions of 2-amino-dA. This is because the bis formamidinium protected 2-amino-dA leads to significant strand scission when standard iodine oxidation is used during synthesis. For this reason, we have also added Pac-2-Amino-dA, a monomer with optimized protection to meet the following criteria: stable during oligonucleotide synthesis, oxidation, and detritylation; labile towards common deprotection conditions (NH₃, AMA, MeNH₂); and the nucleobase protecting groups are cleaved under fairly mild conditions.

Item	Catalog No.	Pack	Price(\$)
2-Amino-dA-CE Phosphoramidite (2,6-diaminopurine)	10-1085-95	50 μmole	70.00
	10-1085-90	100 μmole	125.00
	10-1085-02	0.25g	250.00
Pac-2-Amino-dA-CE Phosphoramidite (2,6-diaminopurine)	10-1585-95	50 μmole	70.00
	10-1585-90	100 μmole	125.00
	10-1585-02	0.25g	250.00

Sequences with high GC content may contain mismatches and still hybridize because of the high stability of the G-C base pair. The N4-ethyl analogue of dC (N4-Et-dC) hybridizes specifically to natural dG but the stability of the base pair is reduced to about the level of an AT base pair.

Coupling N6-Me-dA (10-1003) and N4-Et-dC (10-1068) with 1H-tetrazole leads to a trace of branching at the secondary amine positions, while DCI leads to around 15% branching. In collaboration with Berry and Associates, the acetyl protected monomers were prepared. Acetyl protection was chosen since it would block branching reactions. Oligonucleotides synthesized using these monomers proved to be compatible with all popular deprotection strategies from UltraMild to UltraFast. When the acetyl protected monomers were compared with the unprotected monomers using DCI as activator, branching was reduced from 15% to zero.

Item	Catalog No.	Pack	Price(\$)
N4-Et-dC-CE Phosphoramidite	10-1068-95	50 μmole	125.00
	10-1068-90	100 μmole	225.00
	10-1068-02	0.25g	675.00
N4-Ac-N4-Et-dC-CE Phosphoramidite	10-1513-95	50 μmole	125.00
	10-1513-90	100 μmole	225.00
	10-1513-02	0.25g	675.00
N6-Me-dA-CE Phosphoramidite	10-1003-90	100 μmole	162.50
	10-1003-02	0.25g	495.00
N6-Ac-N6-Me-dA-CE Phosphoramidite	10-1503-90	100 μmole	162.50
	10-1503-02	0.25g	495.00



2-Amino-dA

Pac-2-Amino-dA

N4-Et-dC

N4-Ac-N4-Et-dC

N6-Me-dA

N6-Ac-N6-Me-dA

SEE ALSO

0.5M CSO on page 32
N6-Me-dA on page 62

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

ZIP NUCLEIC ACIDS (ZNA®)

INTELLECTUAL PROPERTY

"Spermine phosphoramidite" synthon is the subject matter of U.S. Divisional Patent Application No. 14/745,871, European Patent No. 1973927 and foreign equivalents for which Polyplus-transfection is the co-owner. Product is sold for research purposes only. Product shall not be used to manufacture oligospermine-oligonucleotide conjugates for use in diagnostics, clinical or commercial applications including use in humans. There is no implied license to manufacture oligospermine-oligonucleotide conjugates for diagnostic, clinical, or commercial applications, including but not limited to contract research. Please contact Polyplus-transfection at licensing@polyplus-transfection.com to obtain a license for such use.

ZNA® is a registered trademark of Polyplus-transfection SA.

Spermine phosphoramidite is used to produce oligospermine-oligonucleotide conjugates - Zip Nucleic Acids (ZNA®) Oligos. The name reflects the presumed mode of action. The conjugates are believed to use the oligospermine to seek out and move along (scan) oligonucleotide strands until the probe complementary sequence is located. The oligospermine then performs the function of stabilizing the formed duplex by reducing electrostatic repulsion, thereby leading to significantly increased binding affinities. ZNA® Oligos have found use in the following applications: Multiplex PCR; PCR of AT-rich Regions; RT qPCR; Detection of MicroRNA; Improved SNP Discrimination; and Antisense and Antigene Effects. Spermine phosphoramidite is simple to use in oligonucleotide synthesis and can be added multiple times at the 3' or 5' terminus. Deprotection and isolation are also straightforward. HPLC analysis of the conjugates requires high pH to suppress the ionization of the spermine residues.

Item	Catalog No.	Pack	Price(\$)
Spermine Phosphoramidite	10-1939-95	50 µmole	145.00
	10-1939-90	100 µmole	270.00
	10-1939-02	0.25g	525.00

CDPI₃ MGB™ LABELING


Synthetic oligonucleotides with covalently-attached CDPI₃ have enhanced DNA affinity and improved the hybridization properties of sequence-specific DNA probes. Short CDPI₃-oligonucleotides hybridize with single-stranded DNA to give more stable DNA duplexes than unmodified ODNs of similar length. The simplest approach to MGB probe design is to use an MGB support, add a quencher molecule as the first addition and complete the synthesis with a 5'-fluorophore. Alternatively, a fluorophore support could be used with the 5' terminus containing a quencher molecule followed by a final MGB addition at the 5' terminus. Glen Research offers 5'-CDPI₃ MGB™ Phosphoramidite and 3'-CDPI₃ MGB™ CPG.

SELECTIVELY BINDING COMPLEMENTARY (SBC) OLIGOS

SBC oligos exhibit high affinity for natural oligonucleotides but they show little affinity for other SBC oligos even of a complementary sequence. Oligos in which A has been replaced with 2-amino-A and T with 2-thio-T represent an excellent example of SBC oligos. While 2-amino-A forms a very stable base pair with T containing three hydrogen bonds, the stability of the base pair with 2-thio-T is greatly diminished. However, 2-thio-T base pairs perfectly well with A. As an example, SBC 20mers annealed against a DNA 20mer target exhibited Tm values 10 °C higher than the corresponding DNA-DNA hybrid, whereas the SBC-SBC hybrid yielded Tm values 30 °C lower.

UNNATURAL BASE PAIRS

Unnatural base pairs display unique abilities in duplex DNA and in nucleic acid and protein biosyntheses. A standard Watson and Crick base pair is formed between iso-C and iso-G, but the hydrogen bonding pattern is quite different from the natural base pairs A-T and C-G. Iso-bases can, therefore, increase specificity of nucleic acid hybridization when introduced as a third base pair. It has also been demonstrated that iso-bases 5-Me-iso-dC and iso-dG can function as degenerate pyrimidine and purine bases, respectively. Iso-dG further functioned as a degenerate base opposite B (C, T, and G) ambiguous sites.



CAPS FOR INCREASED DUPLEX STABILITY AND BASE-PAIRING FIDELITY

New cap structures allow for the preparation of hybridization probes with increased affinity for complementary sequences. The monomers used to prepare capped oligonucleotides are phosphoramidites that can be readily introduced via automated DNA synthesis at the end of solid phase syntheses. The caps favor the formation of stable Watson-Crick duplexes by stacking on the terminal base pair (Figures 1 and 2).




FIGURE 1: STACKING OF CAP ON 5' TERMINAL BASE PAIR







FIGURE 2: STACKING OF Uaq CAP ON 3' TERMINAL BASE PAIR

Melting point increases of over 10 °C per modification can be realized for short duplexes.^{1,2} The caps fit canonical Watson-Crick base pairs and do not stack well on mismatched base pairs. This leads to increased base pairing selectivity at the terminal and the penultimate position of oligonucleotides featuring the caps. Base pairing fidelity is usually low at the termini, where fraying occurs frequently in the absence of caps. The beneficial effects of the caps are also realized when longer target strands are bound, so there is no need for blunt ends for the duplexes formed.^{1,2} The caps, when attached to the 5' terminus of an oligonucleotide, also facilitate purification as their lipophilicity leads to prolonged retention on reversed phase columns or cartridges. Finally, capping of termini may discourage the degradation of oligonucleotides by exonucleases.

3'-Uaq Cap CPG, a Uridine support modified with a 2'- anthraquinone residue, is the most effective oligonucleotide cap known to date.^{3,4} For short hybrid duplexes between DNA probes and RNA target strands, the increase in Tm is up to 18 °C and the modification is effective in increasing the Tm of DNA:DNA, RNA:RNA, and DNA:RNA hybrid duplexes. 3'-Uaq Cap also increases probe specificity by depressing the melting point of terminal mismatches.

Item	Catalog No.	Pack	Price(\$)
5'-Trimethoxystilbene Cap Phosphoramidite	10-1986-90	100 µmole	195.00
	10-1986-02	0.25g	495.00
5'-Pyrene Cap Phosphoramidite	10-1987-90	100 µmole	195.00
	10-1987-02	0.25g	495.00
3'-Uaq Cap CPG	20-2980-01	0.1g	180.00
	20-2980-10	1.0g	1500.00
	20-2980-41	Pack of 4	300.00
	20-2980-42	Pack of 4	150.00
	20-2980-13	Pack of 1	750.00
	20-2980-14	Pack of 1	1125.00



REFERENCES

- (1) Dogan, Z.; Paulini, R.; Rojas Stütz, J. A.; Narayanan, S.; Richert, C. J. *Amer. Chem. Soc.* **2004**, *126*, 4762-4763.
- (2) Narayanan, S.; Gall, J.; Richert, C. *Nucleic Acids Res.* **2004**, *32*, 2901-2911.
- (3) A. Patra, C. Richert, *J. Amer. Chem. Soc.*, 2009, **131**, 12671-12681.
- (4) C. Ahlborn, K. Siegmund, C. Richert, *J. Amer. Chem. Soc.*, 2007, **129**, 15218-15232.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)


DNA METHYLATION

One of the fastest growing fields in biology and cancer research is epigenetics. While the underlying genetic code defines which proteins and gene products are synthesized, it is epigenetic control that defines when and where they are expressed. This dynamic control of gene expression is essential for X chromosome inactivation, embryogenesis, cellular differentiation and appears integral to memory formation and synaptic plasticity.

In 2009, two reports^{1,2} described the discovery of 5-hydroxymethyl-2'-deoxycytidine (hmdC), a novel dC modification in Purkinje neurons and embryonic stem cells. Later, a third report found this modification to be strongly enriched in brain tissues associated with higher cognitive functions.³ This dC modification is generated by the action of α -ketoglutarate dependent ten eleven translocation (TET) enzymes, which oxidizes 5-Me-dC to hmdC. This finding stimulated discussion about active demethylation pathways that could occur, e.g., *via* base excision repair (BER), with the help of specialized DNA glycosylases. Alternatively, one could envision a process in which the hydroxymethyl group of hmdC is further oxidized to 5-formyl-dC (fdC) or 5-carboxy-dC (cadC) followed by elimination of either formic acid or carbon dioxide^{4,5}.

Glen Research has supported this research since its inception by providing the building blocks for the synthesis of oligonucleotides containing all the new dC derivatives - hmdC, fdC and cadC. The first generation hmdC phosphoramidite was fairly very well accepted but requires fairly harsh deprotection conditions. Therefore, a second generation building block (5-Hydroxymethyl-dC II) developed by Carell and co-workers that is compatible with UltraMild deprotection was introduced.⁶ 5-Formyl-dC III has been designed to meet all of the requirements to prepare an oligo containing all of the methylated variants.⁷

Item	Catalog No.	Pack	Price(\$)
5-Hydroxymethyl-dC-CE Phosphoramidite	10-1062-95	50 μ mole	335.00
	10-1062-90	100 μ mole	650.00
	10-1062-02	0.25g	1675.00
5-Carboxy-dC-CE Phosphoramidite	10-1066-95	50 μ mole	230.00
	10-1066-90	100 μ mole	450.00
	10-1066-02	0.25g	1200.00
5-Formyl-dC-CE Phosphoramidite	10-1514-95	50 μ mole	610.00
	10-1514-90	100 μ mole	1200.00
	10-1514-02	0.25g	3225.00
5-Hydroxymethyl-dC II-CE Phosphoramidite	10-1510-95	50 μ mole	345.00
	10-1510-90	100 μ mole	670.00
	10-1510-02	0.25g	2100.00
5-Formyl-dC III-CE Phosphoramidite	10-1564-95	50 μ mole	360.00
	10-1564-90	100 μ mole	700.00
	10-1564-02	0.25g	1800.00



5-Hydroxymethyl-dC

5-Hydroxymethyl-dC II

5-Carboxy-dC

5-Formyl-dC


5-Formyl-dC III

DUPLIX EFFECTS


The design of primers is frequently complicated by the degeneracy of the genetic code. Three strategies are now available to confront this problem. In the first, a mixed base addition (N) is used to form the degenerate site. This approach is best if the number of degenerate sites is small. A second option is the use of 2'-deoxyinosine or 2'-deoxyNebularine which exhibit low, but unequal, hydrogen bonding to the other four bases. The third option is the use of a universal nucleoside. In this strategy, the base analog does not hybridize significantly to the other four bases and makes up some of the duplex destabilization by acting as an intercalating agent. 3-Nitropyrrole 2'-deoxynucleoside (M) is the first example of a set of universal bases. Subsequently, 5-nitroindole was determined to be an effective universal base and to be superior to 3-nitropyrrole, based on duplex melting experiments.

The modified bases designated P and K show considerable promise as degenerate bases. The pyrimidine derivative P, when introduced into oligonucleotides, base pairs with either A or G, while the purine derivative K base pairs with either C or T. A dP+dK mix also can serve as a mixed base with much less degeneracy than dA+dC+dG+dT (N).

Item	Catalog No.	Pack	Price(\$)
dA+dG-CE Phosphoramidites	10-1002-02	0.25g	40.00
dC+dT-CE Phosphoramidites	10-1013-02	0.25g	40.00
dA+dC+dG+dT-CE Phosphoramidites	10-1023-02	0.25g	40.00
Other pack sizes, mixed base combinations and custom doping of individual monomers are available on request. Also, mixed base columns are available in 0.2 and 1.0 μ mole sizes on request.			
dI-CE Phosphoramidite	10-1040-90	100 μ mole	50.00
	10-1040-02	0.25g	120.00
dI-CPG 500	20-2040-01	0.1g	30.00
	20-2190-41	Pack of 4	120.00
	20-2190-42	Pack of 4	72.00
dI-CPG 1000	20-2041-01	0.1g	30.00
	20-2191-41	Pack of 4	120.00
	20-2191-42	Pack of 4	72.00
dU-CE Phosphoramidite	10-1050-90	100 μ mole	35.00
	10-1050-02	0.25g	100.00
dU-CPG 500	20-2050-01	0.1g	30.00
	20-2150-41	Pack of 4	120.00
	20-2150-42	Pack of 4	72.00
dU-CPG 1000	20-2051-01	0.1g	50.00
	20-2151-41	Pack of 4	200.00
	20-2151-42	Pack of 4	120.00



2'-deoxyinosine



2'-deoxyuridine

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

DUPLEX EFFECTS (CONT.)

Item	Catalog No.	Pack	Price(\$)
2'-DeoxyNebularine-CE Phosphoramidite (Purine)	10-1041-90	100 μmole	105.00
	10-1041-02	0.25g	255.00
5-Nitroindole-CE Phosphoramidite	10-1044-90	100 μmole	125.00
	10-1044-02	0.25g	325.00
dP-CE Phosphoramidite	10-1047-90	100 μmole	195.00
	10-1047-02	0.25g	595.00
dK-CE Phosphoramidite	10-1048-90	100 μmole	195.00
	10-1048-02	0.25g	595.00
dP+dK-CE Phosphoramidite	10-1049-90	100 μmole	195.00
	10-1049-02	0.25g	595.00

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)



DUPLEX EFFECTS (CONT.)

Unnatural base pairs display unique abilities in duplex DNA and in nucleic acid and protein biosyntheses. A standard Watson and Crick base pair is formed between iso-C and iso-G, but the hydrogen bonding pattern is quite different from the natural base pairs A-T and C-G. (The 5-methyl analogue was chosen as the synthetic target due to the reported instability of 2'-deoxyisocytidine caused by deamination during oligonucleotide synthesis or deprotection.)

Item	Catalog No.	Pack	Price(\$)
dmf-5-Me-isodC-CE Phosphoramidite	10-1065-90	100 μmole	100.00
	10-1065-02	0.25g	255.00
dmf-isodG-CE Phosphoramidite	10-1078-90	100 μmole	165.00
	10-1078-02	0.25g	355.00


Tm MODULATION

Any technique that involves hybridization of multiple sequences simultaneously, as in DNA chip and reverse hybridization technologies, is subject to inaccuracies due to differences in GC content. Sequences with high GC content may contain mismatches and still hybridize, whereas a low GC content probe may match perfectly and yet disassociate from the target, leading to false positives and negatives, respectively.

An elegant way of circumventing this problem would be to use a modified base that normalized the stability of the GC and AT base pairs. The N4-ethyl analogue (N4-Et-dC) hybridizes specifically to natural dG but the stability of the base pair is reduced to about the level of an AT base pair. In a series of probes whose GC content ranged from 0 to 100%, the range in Tm values when N4-Et-dC was used was only 7 °C; when dC was used, that range was 39 °C.

SEE ALSO

N4-Et-dC on page 47
N4-Ac-N4-Et-dC on page 47



CLEANAMP™ MONOMERS

CleanAmp™ Primers offer an alternative to other Hot Start technologies and allow greater control of primer hybridization and extension during PCR. It has been demonstrated that CleanAmp Primers outperform other technologies in multiple applications. Indeed, over a broad range of applications, CleanAmp Primers reduce or eliminate off-target amplification. Greater amplicon yield is also achieved, due to improvement in specificity and sensitivity. By using either the slow-releasing Precision primers with two CleanAmp phosphotriester linkages or the faster-releasing Turbo Primers with a single CleanAmp phosphotriester linkage, the rate of formation of unmodified primer can be controlled to suit reaction needs. A table to aid in the selection of Turbo and Precision Primers for specific applications is shown below.

Turbo Primers

Fast cycling
 Multiplex PCR
 Improves amplicon yield
 Reduces mis-priming/ primer dimer formation

Precision Primers

Standard cycling
 One-step reverse-transcription PCR
 Improved specificity and limit of detection
 Greatest reduction in mis-priming/primer dimer formation

Synthesis of CleanAmp Primers requires the use of UltraMild Chemistry.





CleanAmp™ Primers and monomers are available from TriLink BioTechnologies.

SEE ALSO

UltraMild DNA Synthesis on page 23

INTELLECTUAL PROPERTY







CleanAmp™ is a trademark of TriLink BioTechnologies, Inc. CleanAmp™ products or portions thereof are covered by TriLink pending Patent Applications, US 2007281308 and WO2007139723, US Provisional Patent Application Serial # 61/056,324 and US Patent 6762298 licensed from the Department of Health and Human Services. CleanAmp™ products are sold exclusively for R & D use by the purchaser. They may not be resold, distributed or re-packaged. No license is granted or implied with the purchase of this product. Amplification methods used in connection with Polymerase Chain Reaction ("PCR") Process are covered by many patents. It may be necessary to obtain a separate license for certain patented applications in which the product is used. CleanAmp™ Licenses are available directly from TriLink BioTechnologies. www.trilinkbiotech.com



CHAIN TERMINATORS

In situations where ligation must be blocked at the 5' terminus, 5'-OMe-dT may be used. 5'-OMe modification of a strand of siRNA using 5'-OMe-T can control guide strand selection and targeting specificity.¹ 5'-Amino-dT terminates an oligonucleotide with a 5'-amino group which may be used for attaching a peptide or a PNA sequence. To avoid polymerase extension at the 3' terminus, 2',3'-dideoxynucleoside and 3'-deoxynucleoside CPGs have proved to be effective. 2',3'-Phosphoramidites are designed to be used with the 5'-phosphoramidites and supports. Since these phosphoramidites have no DMT group, they are not compatible with purification by the DMT-on technique. Ion exchange HPLC or PAGE should be used to purify these dideoxy terminated oligos to ensure that shorter sequences (containing 3'-OH) groups are removed. (3'-Termination can also be effected using a 3'-3' linkage formed using 5'-supports, or 3'-spacer C3 CPG.)

Item	Catalog No.	Pack	Price(\$)
5'-OMe-dT-CE Phosphoramidite	10-1031-90	100 μmole	135.00
	10-1031-02	0.25g	355.00
5'-Amino-dT-CE Phosphoramidite	10-1932-90	100 μmole	150.00
	10-1932-02	0.25g	400.00
3'-dA-CPG	20-2004-01	0.1g	400.00
	20-2104-41	Pack of 4	675.00
	20-2104-42	Pack of 4	225.00
3'-dC-CPG	20-2064-01	0.1g	300.00
	20-2164-41	Pack of 4	600.00
	20-2164-42	Pack of 4	200.00
3'-dG-CPG	20-2074-01	0.1g	300.00
	20-2174-41	Pack of 4	600.00
	20-2174-42	Pack of 4	200.00
3'-dT-CPG	20-2084-01	0.1g	300.00
	20-2184-41	Pack of 4	600.00
	20-2184-42	Pack of 4	200.00



SEE ALSO

5'-Phosphoramidites on page 34
 5'-Supports on page 35
 3'-Spacer C3 CPG on page 84

REFERENCE

(1) PY. Chen, et al., *RNA*, 2008, **14**, 263-274..

OTHER INSTRUMENT TYPES






All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
<i>For Instrument type</i>	<i>Add</i>
Expedite	E
MerMade	M
Columns	
<i>For Instrument type</i>	<i>Add</i>
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

CHAIN TERMINATORS (CONT.)







Item	Catalog No.	Pack	Price(\$)
2',3'-ddC-CPG	20-2017-01	0.1g	300.00
1 μmole columns	20-2117-41	Pack of 4	600.00
0.2 μmole columns	20-2117-42	Pack of 4	200.00
2',3'-ddA-CE Phosphoramidite	10-7001-90	100 μmole	130.00
	10-7001-02	0.25g	545.00
2',3'-ddC-CE Phosphoramidite	10-7101-90	100 μmole	130.00
	10-7101-02	0.25g	545.00
2',3'-ddG-CE Phosphoramidite	10-7201-90	100 μmole	145.00
	10-7201-02	0.25g	675.00
2',3'-ddT-CE Phosphoramidite	10-7301-90	100 μmole	130.00
	10-7301-02	0.25g	545.00



STRUCTURE/ACTIVITY RELATIONSHIP

The following products are used to investigate the effect on the activity of an oligonucleotide when key structural elements are changed. The 7-deaza purine monomers lack groups critical for hydrogen bonding. 7-Deaza-8-aza-A and 7-deaza-8-aza-G (PPG) monomers are isomers of A and G and have similar electron density. Their presence in oligos is slightly stabilizing relative to A and G. Unlike G, PPG does not lead to aggregation and G-rich oligos can be easily prepared and isolated. 5'-Fluorescein oligos with PPG at the 5'-terminus are much less quenched than the equivalent G oligos. As a purine analogue of Thymidine, 7-deaza-2'-deoxyXanthosine (7-deaza-dX) promises to have interesting effects on DNA structure of triplexes. 7-Deaza-dX also forms a non-standard base pair with a 2,4-diaminopyrimidine nucleoside analogue. Standard nucleobases have an unshared pair of electrons that project into the minor groove of duplex DNA. Enzymes that interact with DNA, polymerases, reverse transcriptases, restriction enzymes, etc., may use a hydrogen bond donating group to contact the hydrogen bond acceptor in the minor groove. 3-Deaza-2'-deoxyadenosine is very interesting in that it maintains the ability for regular Watson-Crick hydrogen bonding to T but is lacking the electron pair at the 3-position normally provided by N3.

Item	Catalog No.	Pack	Price(\$)
7-Deaza-dA-CE Phosphoramidite	10-1001-95	50 μmole	177.50
	10-1001-90	100 μmole	355.00
	10-1001-02	0.25g	975.00
7-Deaza-8-aza-dA-CE Phosphoramidite	10-1083-95	50 μmole	177.50
	10-1083-90	100 μmole	355.00
	10-1083-02	0.25g	975.00
7-Deaza-dG-CE Phosphoramidite	10-1021-95	50 μmole	177.50
	10-1021-90	100 μmole	355.00
	10-1021-02	0.25g	975.00
7-Deaza-8-aza-dG-CE Phosphoramidite (PPG)	10-1073-95	50 μmole	207.50
	10-1073-90	100 μmole	395.00
	10-1073-02	0.25g	1150.00
7-deaza-dX-CE Phosphoramidite	10-1076-95	50 μmole	177.50
	10-1076-90	100 μmole	355.00
	10-1076-02	0.25g	975.00
3-Deaza-dA-CE Phosphoramidite	10-1088-95	50 μmole	177.50
	10-1088-90	100 μmole	355.00
	10-1088-02	0.25g	975.00



STABILITY NOTES

7-Deaza-dG is unstable to iodine oxidation. Add a maximum of 2 times when using iodine oxidation or use 0.5M (10-camphorsulfonyl)-oxaziridine (CSO) in anhydrous acetonitrile and 3 min. oxidation time. (See Glen Report-Vol.9, No.1, 1996, page 8.)

INTELLECTUAL PROPERTY

The use of PPG is subject to proprietary rights of ELITechGroup and it is sold under license from ELITechGroup.


(1) I.V. Kutyavin, et al., *Nucleic Acids Res.*, 2002, **30**, 4952-4959.

STRUCTURE/ACTIVITY RELATIONSHIP (CONT.)

The C-nucleoside 2'-deoxypseudouridine, in contrast to dU, forms stable C:pseudoU-A triplets. 2-Aminopurine lacks groups critical for hydrogen bonding and is a mildly fluorescent base.


Demand for sulfur modified bases continues to expand for investigations of oligonucleotide structure, but primarily for cross-linking purposes. 6-Thio-dG, 4-Thio-dT and 4-thio-dU are very useful modifications for photo cross-linking and photoaffinity labeling experiments. Oligos containing 2-thio-dT are useful in examining protein-DNA interaction by acting as photosensitizing probes. The thiocarbonyl group in 2-thio-dT is especially interesting in that it is available to react with compounds associating with the minor groove of DNA. 2-Amino-A forms a very stable base pair with T containing three hydrogen bonds but the stability of the base pair with 2-thio-T is greatly diminished. Due to steric interactions between the 2-thio group of thymidine and the 2-amino group of 2-amino-A, the base pair contains only a single hydrogen bond. Oligos containing 2-amino-dA and 2-thio-dT exhibit high affinity for natural oligonucleotides but show little affinity for other similar oligos even of a complementary sequence.

Item	Catalog No.	Pack	Price(\$)
2'-deoxypseudoU-CE Phosphoramidite	10-1055-95	50 μ mole	177.50
	10-1055-90	100 μ mole	355.00
	10-1055-02	0.25g	975.00
2-Aminopurine-CE Phosphoramidite	10-1046-90	100 μ mole	135.00
	10-1046-02	0.25g	355.00
6-Thio-dG-CE Phosphoramidite	10-1072-95	50 μ mole	177.50
	10-1072-90	100 μ mole	355.00
	10-1072-02	0.25g	975.00
4-Thio-dT-CE Phosphoramidite	10-1034-95	50 μ mole	165.00
	10-1034-90	100 μ mole	295.00
	10-1034-02	0.25g	675.00
4-Thio-dU-CE Phosphoramidite	10-1052-95	50 μ mole	165.00
	10-1052-90	100 μ mole	295.00
	10-1052-02	0.25g	675.00
2-Thio-dT-CE Phosphoramidite	10-1036-95	50 μ mole	165.00
	10-1036-90	100 μ mole	295.00
	10-1036-02	0.25g	675.00




2'-dpseudoU


58




2-Aminopurine




6-Thio-dG



4-Thio-dT



4-Thio-dU



2-Thio-dT

STRUCTURE/ACTIVITY RELATIONSHIP (CONT.)


8-Amino-dA and 8-amino-dG are useful in triplex formation due to the presence of the additional amino groups.

2'-DeoxyXanthosine (dX) is a naturally occurring nucleoside that may be derived from oxidative deamination of 2'-deoxyGuanosine (dG). dX has a similar bonding pattern to thymidine and it may base pair with dA, with such purine-purine interactions causing duplex distortion. dX also featured in attempts to extend the genetic alphabet with a new base pair of dX and pyrimidine-2,4-diamine nucleoside. dX has also interested researchers in the field of DNA damage and repair since it is a product of nitric oxide-induced mutagenesis.


Item	Catalog No.	Pack	Price(\$)
8-Amino-dA-CE Phosphoramidite	10-1086-95	50 μ mole	177.50
	10-1086-90	100 μ mole	355.00
	10-1086-02	0.25g	975.00
8-Amino-dG-CE Phosphoramidite	10-1079-95	50 μ mole	177.50
	10-1079-90	100 μ mole	355.00
	10-1079-02	0.25g	975.00
2'-dX-CE Phosphoramidite	10-1537-95	50 μ mole	105.00
	10-1537-90	100 μ mole	200.00
	10-1537-02	0.25g	420.00

STABILITY NOTE


Synthetic oligonucleotides containing 8-amino-dG must be cleaved and deprotected with ammonium hydroxide containing 0.25M 2-mercaptoethanol to avoid oxidative degradation of 8-amino-dG sites.



8-Amino-dA



8-Amino-dG




dX

59


HALOGENATED NUCLEOSIDES

Brominated and iodinated nucleosides are used in X-ray crystallography studies of oligonucleotide structure. They are also photolabile and are used for cross-linking studies to probe the structure of protein-DNA complexes. Antibodies exist to Br-dU and oligonucleotides containing Br-dU can be used as probes. 5-Fluoro-dU can be used as a non-photoreactive alternative to 5-Br-dU with similar electron density. 5-F-dU base pairs more strongly than T to both dA and the dG mismatch. It is also useful for probing DNA structure using ¹⁹F NMR spectroscopy.


Item	Catalog No.	Pack	Price(\$)
8-Br-dA-CE Phosphoramidite	10-1007-90	100 μmole	115.00
	10-1007-02	0.25g	295.00
8-Br-dG-CE Phosphoramidite	10-1027-90	100 μmole	105.00
	10-1027-02	0.25g	255.00
5-Br-dC-CE Phosphoramidite	10-1080-90	100 μmole	60.00
	10-1080-02	0.25g	160.00
5-I-dC-CE Phosphoramidite	10-1081-90	100 μmole	135.00
	10-1081-02	0.25g	355.00
5-Br-dU-CE Phosphoramidite	10-1090-90	100 μmole	60.00
	10-1090-02	0.25g	160.00
5-I-dU-CE Phosphoramidite	10-1091-90	100 μmole	60.00
	10-1091-02	0.25g	160.00
5-F-dU-CE Phosphoramidite	10-1092-90	100 μmole	135.00
	10-1092-02	0.25g	355.00
5-Br-dU-CPG	20-2090-01	0.1g	50.00
1 μmole columns	20-2090-41	Pack of 4	200.00
0.2 μmole columns	20-2090-42	Pack of 4	120.00




8-Bromo-2'-deoxyAdenosine




8-Bromo-2'-deoxyGuanosine




5-Bromo-2'-deoxyCytidine




5-Iodo-2'-deoxyCytidine



5-Bromo-2'-deoxyUridine



5-Iodo-2'-deoxyUridine




5-Fluoro-2'-deoxyUridine


DNA DAMAGE/REPAIR

Cellular DNA is constantly being damaged by oxidation and alkylation, by free radicals, and by ultraviolet and ionizing radiation. The body has therefore evolved a number of repair enzyme systems to excise and repair these lesions. The 8-oxo purine monomers allow investigation of the structure and activity of oligonucleotides containing an 8-oxo mutation which is formed naturally when DNA is subjected to oxidative conditions or ionizing radiation. 5,6-Dihydro pyrimidines are naturally occurring compounds that are structural components of alanine transfer RNA. Dihydrouracil and the hydroxy pyrimidines are major base damage products formed by exposure of DNA to ionizing radiation.


Item	Catalog No.	Pack	Price(\$)
8-Oxo-dA-CE Phosphoramidite	10-1008-90	100 μmole	135.00
	10-1008-02	0.25g	355.00
8-Oxo-dG-CE Phosphoramidite	10-1028-95	50 μmole	177.50
	10-1028-90	100 μmole	355.00
	10-1028-02	0.25g	975.00
	10-1530-90	100 μmole	195.00
5,6-Dihydro-dT-CE Phosphoramidite	10-1530-02	0.25g	595.00
	10-1550-90	100 μmole	195.00
5,6-Dihydro-dU-CE Phosphoramidite	10-1550-02	0.25g	595.00
	10-1063-90	100 μmole	275.00
5-OH-dC-CE Phosphoramidite	10-1063-02	0.25g	775.00
	10-1053-90	100 μmole	225.00
5-OH-dU-CE Phosphoramidite	10-1053-02	0.25g	675.00
	10-1093-90	100 μmole	225.00
5-Hydroxymethyl-dU-CE Phosphoramidite	10-1093-02	0.25g	675.00




8-oxo-2'-deoxyAdenosine




8-oxo-2'-deoxyGuanosine




5,6-Dihydro-dT




5,6-Dihydro-dU



5-OH-dC



5-OH-dU



5-Hydroxymethyl-dU

STABILITY NOTES

Synthetic oligonucleotides containing 8-oxo-dG must be cleaved and deprotected with ammonium hydroxide containing 0.25M 2-mercaptoethanol to avoid oxidative degradation of 8-oxo-dG sites.

Oligonucleotides synthesized using 5,6-dihydro-dU or 5,6-dihydro-dT and UltraMILD monomers can be cleaved using either concentrated ammonium hydroxide or 50 mM potassium carbonate in anhydrous methanol. Complete cleavage and deprotection can be accomplished at room temperature in 2-4 hours without damaging either the dihydro-dU or dihydro-dT bases.

SEE ALSO

5-Hydroxymethyl-dC on page 50
dX on page 59

DNA DAMAGE/REPAIR (CONT.)

STABILITY NOTES

Synthetic oligonucleotides containing 8-amino-dG must be cleaved and deprotected with ammonium hydroxide containing 0.25M 2-mercaptoethanol to avoid oxidative degradation of 8-amino-dG sites.

Oligonucleotides synthesized using Thymidine Glycol and UltraMILD monomers can be cleaved using either concentrated ammonium hydroxide or 50 mM potassium carbonate in anhydrous methanol. Complete cleavage and deprotection can be accomplished at room temperature in 2-4 hours without damaging Thymidine Glycol base. The best method to remove the TBDMS groups was achieved using TEA.3HF at 40°C overnight.

REFERENCE

(1) K. Groebke, and C.J. Leumann, *Helv Chim Acta*, 1990, **73**, 608-617.

SEE ALSO

dSpacer on page 84
Pyrrolidine on page 63

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add


Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)


8-Amino-G is formed along with 8-oxo-G as the major mutagenic lesions formed in DNA damage caused by 2-nitropropane. 2-Nitropropane is an industrial solvent and a component of paints, dyes and varnishes, and is also present in cigarette smoke. Thymine glycol (5,6-dihydroxy-5,6-dihydrothymine) is formed when thymine is subjected to oxidative stress, including ionizing radiation. Oxidation of the 5,6 double bond of Thymidine generates two chiral centers at C5 and C6. The cis-5R,6S form is generated as the predominant product along with the other diastereomer, the cis-5S,6R form. The presence of thymidine glycol in DNA has significant biological consequences and many organisms possess specific repair enzymes for the excision of this lesion.

Hydrolysis of nucleoside residues in DNA occurs to generate abasic sites. Most commonly, dA sites are hydrolyzed causing depurination and leading to abasic residues. For researchers trying to determine if their source of depurination in chemical synthesis of DNA is reagent, fluidics or protocol-based, we offer a depurination-resistant dA monomer. A new chemical method allows the generation of abasic sites in double and single stranded oligonucleotides using very mild specific conditions and with very low probability of side reactions. Abasic II Phosphoramidite¹ has the advantage of simplicity in that the silyl group is removed post-synthesis using aqueous acetic acid. dSpacer has also been used successfully as a mimic of the highly base-labile abasic site.


Item	Catalog No.	Pack	Price(\$)
8-Amino-dG-CE Phosphoramidite	10-1079-95	50 µmole	177.50
	10-1079-90	100 µmole	355.00
	10-1079-02	0.25g	975.00
Thymidine Glycol CE Phosphoramidite	10-1096-95	50 µmole	180.00
	10-1096-90	100 µmole	360.00
	10-1096-02	0.25g	975.00
Abasic II Phosphoramidite (dR Precursor)	10-1927-95	50 µmole	80.00
	10-1927-90	100 µmole	150.00
	10-1927-02	0.25g	475.00



8-Amino-dG



Thymidine Glycol



Abasic II Phosphoramidite



DNA DAMAGE/REPAIR (CONT.)

One of the major sources of DNA damage in all organisms is the UV component of sunlight. The predominant reaction induced by UV light on DNA is dimerization of adjacent pyrimidine bases leading to cyclobutane dimers (CPDs). The dimers formed in the most significant quantity are the cis-syn cyclobutane dimer of two thymine bases. Although formed routinely, these dimer products are efficiently excised and repaired enzymatically by nucleotide excision repair (NER) or the dimerization is reversed by photolase enzymes. A further mode of oxidative damage is radiation-induced damage of DNA, which has been shown to lead to bridged cyclonucleosides. The purines, cyclo-dA and cyclo-dG, are predominantly formed, although the cyclo pyrimidines have also been detected. Cyclo-dA is doubly intriguing since it contains both damaged base and damaged sugar residues and, as such, should have a considerable biological impact. In a manner analogous to thymine dimer, cyclo purines cause significant distortion of the regular DNA helix and these lesions are repaired not by base excision repair (BER) but by NER.


Item	Catalog No.	Pack	Price(\$)
<i>Cis-syn</i> Thymine Dimer Phosphoramidite	11-1330-95	50 µmole	2100.00
	11-1330-90	100 µmole	4200.00
	11-1330-02	0.25g	10200.00
5',8-Cyclo-dA CE Phosphoramidite	10-1098-95	50 µmole	950.00
	10-1098-90	100 µmole	1850.00
	10-1098-02	0.25g	5350.00
5',8-Cyclo-dG CE Phosphoramidite	10-1598-95	50 µmole	1250.00
	10-1598-90	100 µmole	2450.00

Base excision repair (BER) is one of the most studied repair mechanisms. In this pathway, DNA glycosylases recognize the damaged bases and catalyze their excision through hydrolysis of the N-glycosidic bond. Attempts to understand the structural basis for DNA damage recognition by DNA glycosylases have been hampered by the short-lived association of these enzymes with their DNA substrates. To overcome this problem, the Verdine group at Harvard synthesized a pyrrolidine analog that mimics the charged transition state of the enzyme-substrate complex. When incorporated into double-stranded DNA, they found the pyrrolidine analog (PYR), introduced as the Pyrrolidine-CE Phosphoramidite, forms an extremely stable complex with the DNA glycosylase AlkA, exhibiting a dissociation constant in the pM range and potently inhibited the reaction catalyzed by the enzyme.


Item	Catalog No.	Pack	Price(\$)
Pyrrolidine-CE Phosphoramidite (PYR)	10-1915-95	50 µmole	190.00
	10-1915-90	100 µmole	380.00
	10-1915-02	0.25g	1085.00

*Cis-syn* Thymine Dimer

5',8-Cyclo-dA



5',8-Cyclo-dG




Pyrrolidine

INSTRUMENT TYPES

For these very expensive phosphoramidites, an ABI septum vial is the standard vial. Add E to the catalog no. for an Expedite vial or V to the catalog no. for an Expedite V vial.

BIOCOMPATIBLE TRIAZOLE LINKAGE



SEE ALSO

5'-I-dT in *Click Chemistry* on page 90

Click Chemistry on page 88

REFERENCES - CLICK LIGATION

- (1) A.H. El-Sagheer, A.P. Sanzone, R. Gao, A. Tavassoli, and T. Brown, *Proc Natl Acad Sci U S A*, 2011, **108**, 11338-43.
- (2) A.H. el-Sagheer, and T. Brown, *Chem Commun (Camb)*, 2011, **47**, 12057-8.
- (3) A.P. Sanzone, A.H. El-Sagheer, T. Brown, and A. Tavassoli, *Nucleic Acids Res*, 2012.
- (4) A. Dallmann, et al., *Chemistry*, 2011, **17**, 14714-7.
- (5) A.H. El-Sagheer, and T. Brown, *Proc Natl Acad Sci U S A*, 2010, **107**, 15329-34.

REFERENCES - MicroRNA Labeling

- (1) H. Vogel, and C. Richert, *ChemBioChem*, 2012, **13**, 1474-82.
- (2) R. Eisenhuth, and C. Richert, *Journal of Organic Chemistry*, 2008, **74**, 26-37.
- (3) E. Kervio, A. Hochgesand, U.E. Steiner, and C. Richert, *Proc Natl Acad Sci U S A*, 2010, **107**, 12074-9.

CLICK DNA AND RNA LIGATION


Ligation of an oligo containing a 5'-azide with an oligo containing a 3'-propargyl group using Click Chemistry leads to a triazole linkage that has been shown to have *in vivo* biocompatibility. This technique has been used to synthesize DNA constructs up to 300 bases in length. When the resultant triazole linkage was placed in a PCR template, various polymerases were able to copy the sequence correctly. The linkage has also been shown to be compatible with transcription and rolling circle amplification, as well as gene expression in *E. coli*. In the RNA world, a hammerhead ribozyme containing the triazole linkage at the substrate cleavage site has been shown to retain its activity. A large variety of applications is envisaged for this biocompatible chemical ligation. Support for this technology is offered with the help of Tom Brown's group at the University of Southampton.

Item	Catalog No.	Pack	Price(\$)	
3'-Propargyl-5-Me-dC CPG	20-2982-01	0.1g	180.00	
	20-2982-10	1.0g	1500.00	
	1 μmole columns	20-2982-41	Pack of 4	300.00
	0.2 μmole columns	20-2982-42	Pack of 4	150.00
	10 μmole column (ABI)	20-2982-13	Pack of 1	750.00
15 μmole column (Expedite)	20-2982-14	Pack of 1	1125.00	


5'-LABELING OF MicroRNAs

Several methods have been developed for the detection of miRNAs, however, few allow the simultaneous detection of multiple miRNAs. To overcome this analytical deficiency, the Richert group at the University of Stuttgart has recently developed an ingenious method to selectively detect miRNAs on microarrays without interference from endogenous pre-mRNAs, mRNAs and other RNA species. In this method, a short oligonucleotide containing 3'-amino-dT and a 5' reporter molecule is chemically ligated to the microRNA in a one-step procedure by *in situ* activation of the microRNA. This is specifically achieved by taking advantage of the fact that miRNAs, unlike other RNAs, are 5'-phosphorylated. The reaction is template-directed (and thus sequence specific) and can be performed together with enzymatic 3'-extension/labeling, either in solution or on a support. The short DNA labeling strand may feature one of a variety of different labels, such as a biotin group or a fluorophore.

Item	Catalog No.	Pack	Price(\$)	
3'-Amino-dT CPG	20-2981-01	0.1g	120.00	
	20-2981-10	1.0g	995.00	
	1 μmole columns	20-2981-41	Pack of 4	200.00
	0.2 μmole columns	20-2981-42	Pack of 4	120.00
	10 μmole column (ABI)	20-2981-13	Pack of 1	500.00
15 μmole column (Expedite)	20-2981-14	Pack of 1	750.00	



3'-Propargyl-5-Me-dC CPG




3'-Amino-dT CPG


2'-5' LINKED OLIGONUCLEOTIDES

Cellular DNA and RNA are made up of ribo- and 2'-deoxyribonucleic acids linked together via 3'-5' phosphodiester linkages and by far comprise the bulk of polynucleic acids found in cells. Much less common are oligonucleotides which have 2'-5' linkages. However, a unique feature of 2'-5' linked oligonucleotides is their ability to bind selectively to complementary RNA. These features suggest a number of interesting uses for 2'-5' linked oligos such as their use as RNA specific probes or in antisense oligos. Recently, oligos have been synthesized using 3'-deoxy-2'-phosphoramidites and 2'-deoxy-3'-phosphoramidites to produce chimeras with 2'-5' linked ends and 3'-5' linked central regions. It was found that 2'-5' phosphorothioate oligos: 1) bind selectively to complementary RNA with the same affinity as phosphodiester oligos; 2) exhibit much less nonspecific binding to cellular proteins; 3) do not activate RNase H. A 3'-deoxynucleoside at the 3'-terminus of an otherwise normal oligonucleotide effectively blocks polymerase extension.


Item	Catalog No.	Pack	Price(\$)	
3'-dA-CE Phosphoramidite	10-1004-95	50 μmole	177.50	
	10-1004-90	100 μmole	355.00	
	10-1004-02	0.25g	975.00	
3'-dC-CE Phosphoramidite	10-1064-95	50 μmole	177.50	
	10-1064-90	100 μmole	355.00	
	10-1064-02	0.25g	975.00	
3'-dG-CE Phosphoramidite	10-1074-95	50 μmole	177.50	
	10-1074-90	100 μmole	355.00	
	10-1074-02	0.25g	975.00	
3'-dT-CE Phosphoramidite	10-1084-95	50 μmole	177.50	
	10-1084-90	100 μmole	355.00	
	10-1084-02	0.25g	975.00	
3'-dA-CPG	20-2004-01	0.1g	300.00	
	1 μmole columns	20-2104-41	Pack of 4	600.00
	0.2 μmole columns	20-2104-42	Pack of 4	200.00
3'-dC-CPG	20-2064-01	0.1g	300.00	
	1 μmole columns	20-2164-41	Pack of 4	600.00
	0.2 μmole columns	20-2164-42	Pack of 4	200.00




3'-dA




3'-dC




3'-dG



3'-dT



3'-dA-CPG



3'-dC-CPG

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

SEE ALSO

3'-deoxynucleoside CPG on page 51


2'-5' LINKED OLIGONUCLEOTIDES (CONT.)

Item	Catalog No.	Pack	Price(\$)
3'-dG-CPG	20-2074-01	0.1g	300.00
1 μmole columns	20-2174-41	Pack of 4	600.00
0.2 μmole columns	20-2174-42	Pack of 4	200.00
3'-dT-CPG	20-2084-01	0.1g	300.00
1 μmole columns	20-2184-41	Pack of 4	600.00
0.2 μmole columns	20-2184-42	Pack of 4	200.00

MUTAGENESIS

Cellular polynucleotides are alkylated by endogenous components, such as S-adenosylmethionine, or after reacting with two general classes of environmental and laboratory chemicals. SN1 chemical agents include alkylnitrosourea and N-alkyl-N-nitro-N-nitrosoguanidine that react with the N7 position of guanine, N3 of adenine, O6 of guanine, O2 or O4 of pyrimidines, and the non-phosphodiester oxygen atoms of the phosphate backbone. In contrast, SN2 chemical agents such as methyl methanesulfonate and dimethyl sulfate react primarily with the N1 position of adenine (1-Methyl-2'-deoxyadenosine) and N3 of cytosine. To avoid chain branching during synthesis when using DCI as activator, N6-Me-dA is offered with acetyl protection.


Item	Catalog No.	Pack	Price(\$)
06-Me-dG-CE Phosphoramidite	10-1070-90	100 μmole	105.00
	10-1070-02	0.25g	255.00
N6-Me-dA-CE Phosphoramidite	10-1003-90	100 μmole	162.50
	10-1003-02	0.25g	495.00
N6-Ac-N6-Me-dA-CE Phosphoramidite	10-1503-90	100 μmole	162.50
	10-1503-02	0.25g	495.00
O4-Me-dT-CE Phosphoramidite	10-1032-90	100 μmole	135.00
	10-1032-02	0.25g	355.00
1-Me-dA-CE Phosphoramidite	10-1501-95	50 μmole	125.00
	10-1501-90	100 μmole	250.00
	10-1501-02	0.25g	750.00



IN SITU SYNTHESIS OF DNA ANALOGS

The convertible nucleoside strategy is one of the most versatile methods for producing modifications in bases to examine their effects on DNA structure and activity. In some cases, with versatility comes difficulty in that the convertible base is modified after oligonucleotide synthesis. The chemistry is sometimes complex and base composition analysis of the final oligonucleotide is required to verify structure. The convertible dU monomer can be used to introduce a variety of modifications at the convertible position, including N, O and S modifications. Convertible F-dC is by far the simplest approach to the preparation of oligonucleotides containing F-dC - normal ammonium hydroxide treatment effects the conversion to F-dC. Convertible dA has been used to prepare oligonucleotides containing multiple points for attachment to solid supports. In this way, high capacity affinity supports for the purification of DNA binding proteins have been prepared. 2-F-dI is a convertible nucleoside for the preparation of 2'-dG derivatives following the displacement of the 2-fluorine by primary amines.

Item	Catalog No.	Pack	Price(\$)
TMP-F-dU-CE Phosphoramidite (Convertible F-dC)	10-1016-90	100 μmole	195.00
	10-1016-02	0.25g	495.00
O6-Phenyl-dI-CE Phosphoramidite (Convertible dA)	10-1042-90	100 μmole	135.00
	10-1042-02	0.25g	355.00
O4-Triazolyl-dU-CE Phosphoramidite (Convertible dU)	10-1051-90	100 μmole	135.00
	10-1051-02	0.25g	355.00
2-F-dI-CE Phosphoramidite (Convertible dG)	10-1082-95	50 μmole	180.00
	10-1082-90	100 μmole	360.00
	10-1082-02	0.25g	975.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

ABBREVIATION

TMP = 2,4,6-trimethylphenyl

SEE ALSO

2-Aminopurine on page 58
 AP-dC (G-Clamp) on page 46
 UltraMild Chemistry on page 23
 Pyrrolo-C on page 132
 Pyrrolo-CTP on page 136

INTELLECTUAL PROPERTY

Pyrrolo-dC is a joint development project of Berry & Associates, Inc. and Glen Research Corporation. Pyrrolo-dC is covered by US Patent No.: 7,144,995.

SPECTRAL PROPERTIES

The spectral properties of pyrrolo-dC, coupled with its unique base-pairing ability, make this fluorescent analog extremely valuable in probing DNA structure. When the pyrrolo-dC is base-paired, its fluorescence is significantly quenched through what is most likely base stacking or dG interactions. The quantum yield of fluorescence for pyrrolo-dC is quite sensitive to its hybridization state, making it ideally suited for probing the dynamic structure of DNA.

	QY	λ	ϵ (L/mol.cm)
single-stranded	0.07	260nm 347nm	4000 3700
double-stranded	0.02		

(QY determined relative to quinine sulfate in 0.5M H2SO4)

REFERENCES


1. D.A. Berry, et al., *Tetrahedron Lett*, 2004, **45**, 2457-2461.
2. The Glen Report, 2007, **19**, 8-9.
3. P. Sandin, et al., *Nucleic Acids Res.*, 2008, **36**, 157-167.
4. P. Sandin, et al., *Nucleic Acids Res.*, 2005, **33**, 5019-5025.
5. K.C. Engman, et al., *Nucleic Acids Res.*, 2004, **32**, 5087-5095.

PROBING DNA STRUCTURE WITH FLUORESCENT NUCLEOSIDES


Etheno-dA is a fluorescent nucleoside which is especially useful in observing the transition between DNA structural types. It is quite base labile and should be deprotected with ammonium hydroxide at room temperature for 24 hours. Alternatively, UltraMild chemistry can be used. 2-Aminopurine and AP-dC (G-Clamp) are also useful fluorescent nucleosides.

Pyrrolo-dC is a fluorescent deoxycytidine analog that is an ideal probe of DNA structure and dynamics.^{1,2} It base-pairs as a normal dC nucleotide. An oligo fully substituted with pyrrolo-dC has the same T_m as the control dC oligo with the same specificity for dG. Its small size does not perturb the structure of the DNA helix and it is well tolerated by a number of DNA and RNA polymerases. It is highly fluorescent and its excitation and emission are well to the red of most fluorescent nucleotide analogs, which eliminates or reduces background fluorescence from proteins. Pyrrolo-dCTP has potential uses in biological assay development.


Item	Catalog No.	Pack	Price(\$)
Etheno-dA-CE Phosphoramidite	10-1006-90	100 μ mole	105.00
	10-1006-02	0.25g	255.00
Pyrrolo-dC-CE Phosphoramidite	10-1017-95	50 μ mole	110.00
	10-1017-90	100 μ mole	220.00
	10-1017-02	0.25g	675.00
Pyrrolo-dCTP (10 mM)	81-1017-01	100 μ L	150.00



Etheno-2'-deoxyAdenosine



Pyrrolo-dC



Pyrrolo-dCTP

PROBING DNA STRUCTURE WITH FLUORESCENT NUCLEOSIDES (CONT.)

By attaching pyrene or perylene to the 5 position of deoxyuridine through a triple bond, the fluorophore is electronically coupled to the deoxyuridine base. This electronic coupling of the base and the fluorophore makes the fluorescence sensitive to the base pairing of the dU portion of the molecule, allowing the discrimination between perfect and one base mismatched targets.

Item	Catalog No.	Pack	Price(\$)
Pyrene-dU-CE Phosphoramidite	10-1590-95	50 μ mole	105.00
	10-1590-90	100 μ mole	210.00
	10-1590-02	0.25g	550.00
Perylene-dU-CE Phosphoramidite	10-1591-95	50 μ mole	150.00
	10-1591-90	100 μ mole	300.00
	10-1591-02	0.25g	720.00

SPECTRAL PROPERTIES


	Absorbance Maximum	Emission Maximum
Pyrene-dU	402nm	472nm
Perylene-dU	473nm	490nm

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)



Pyrene-dU

Perylene-dU

PROBING DNA STRUCTURE WITH FLUORESCENT NUCLEOSIDES (CONT.)

The tricyclic fluorescent nucleoside analogues, 1,3-diaza-2-oxophenothiazine, tC, and 1,3-diaza-2-oxophenoxazine, tC^o, are deoxycytidine analogs that have been shown to base pair faithfully with dG with virtually no disruption of the normal duplex structure.³⁻⁵ This means that the stability of the DNA duplex is not compromised as compared to the control regardless of DNA sequence. The fluorescence quantum yield of tC is essentially unchanged between single stranded and double stranded DNA - 0.21 for single stranded DNA and 0.19 for duplex DNA. Also, the fluorescence characteristics of tC are not sensitive to neighboring base combinations. tC^o has been shown to be the brightest fluorescent nucleoside analogue in duplex context reported so far and even retains the majority of its fluorescence when surrounded by guanine residues. Indeed, tC^o has been reported to be 25-50 times brighter than 2-aminopurine.


The base analogue tC_{nitro} is a FRET-acceptor together with tC^o (or tC) as the donor molecule. This constitutes the first ever description of a nucleobase FRET-pair. This novel FRET-pair provides a unique tool for investigations of nucleic acid containing systems. tC_{nitro} is virtually non-fluorescent under normal conditions.

Item	Catalog No.	Pack	Price(\$)
tC-CE Phosphoramidite	10-1516-95	50 μmole	250.00
	10-1516-90	100 μmole	490.00
	10-1516-02	0.25g	1460.00
tC ^o -CE Phosphoramidite	10-1517-95	50 μmole	250.00
	10-1517-90	100 μmole	490.00
	10-1517-02	0.25g	1460.00
tC _{nitro} -CE Phosphoramidite	10-1518-95	50 μmole	265.00
	10-1518-90	100 μmole	520.00
	10-1518-02	0.25g	1460.00


PHOTO-REGULATION OF DNA FUNCTION

Glen Research's interest lies in the preparation of caged oligonucleotides whose function is restored after uncaging by UV light at a wavelength that causes no DNA damage. The Deiters group at North Carolina State University has described NPOM-Caged-dT, where the nucleobase is caged with the photolabile group, 6-nitropiperonyloxymethyl (NPOM), which can be removed using UV light at 365nm. Oligonucleotides containing NPOM-Caged-dT every five or six bases do not hybridize to their complementary strand. Photo-uncaging of the caged oligonucleotide is then easily carried out with UV light at 365 nm for seconds to minutes to restore the activity of the oligonucleotide.


Item	Catalog No.	Pack	Price(\$)
NPOM-Caged-dT-CE Phosphoramidite	10-1534-95	50 μmole	185.00
	10-1534-90	100 μmole	355.00
	10-1534-02	0.25g	895.00




tC



tC^o



tC_{nitro}



NPOM-Caged-dT

INHIBITION OF DNA METHYLTRANSFERASES

Zebularine (pyrimidin-2-one ribonucleoside) is a cytidine analogue that acts as a DNA demethylase inhibitor, as well as a cytidine deaminase inhibitor. This structure is very active biologically and Zebularine is now used as a potent anti-cancer drug. A 2'-deoxynucleoside analogue of Zebularine, 5-methyl-pyrimidin-2-one, 2'-deoxynucleoside, has been used to probe the initiation of the cellular DNA repair process by making use of its mildly fluorescent properties. This combination of biological activity and fluorescence properties would make 5-Me-2'-deoxyZebularine a strong addition to our array of nucleoside analogues.


Cytosine-5-methyltransferases are found in everything from archaeobacteria to mammals and when the regulation of cytosine-5-methyltransferases goes awry, cancer can result. The mechanism of action for this family of enzymes involves attack of a cysteine thiol group on the C6 position of cytosine, leading to a transient dihydrocytosine intermediate, which then facilitates the nucleophilic attack by C5 on the activated methyl group of the S-adenosyl-L-methionine cofactor. As with many enzymes, the intermediate can be trapped using a suicide substrate and 5-fluoro-cytosine has been used extensively in this role. An alternate strategy is to use a transition-state mimic that binds to the active site with high affinity. An excellent candidate was found in 5-aza-5,6-dihydrocytosine. Despite not being covalently bound to the enzyme, it was found^{1,2} to be a more potent inhibitor of cytosine-5-methyltransferases than 5-fluoro-cytosine. 5-Aza-5,6-dihydro-dC is compatible with standard oligonucleotide synthesis and deprotection conditions and is an excellent tool for use in methyltransferase research.

Item	Catalog No.	Pack	Price(\$)
5-Me-2'-deoxyZebularine-CE Phosphoramidite	10-1061-95	50 μmole	200.00
	10-1061-90	100 μmole	400.00
	10-1061-02	0.25g	975.00
5-Aza-5,6-dihydro-dC-CE Phosphoramidite	10-1511-95	50 μmole	180.00
	10-1511-90	100 μmole	360.00
	10-1511-02	0.25g	1120.00

LARGE SCALE SYNTHESIS

The most common side reaction during deprotection of oligonucleotides on a large scale is the alkylation of dT residues by acrylonitrile, formed by β-elimination of the cyanoethyl phosphate protecting groups, to generate N3-cyanoethyl-dT.

Item	Catalog No.	Pack	Price(\$)
N3-Cyanoethyl-dT	10-1531-90	100 μmole	200.00
	10-1531-02	0.25g	600.00



SEE ALSO

Ribo-tC^o on page 136

SPECTRAL PROPERTIES

Absorption and emission data for tC and tC^o are collected below:

tC	QY	λ	ε
		(nm)	(L/mol.cm)
single-stranded	0.21	385nm	4000
double-stranded	0.19		

tC ^o	QY	λ	ε
		(nm)	(L/mol.cm)
single-stranded	0.30	360nm	9000
double-stranded	0.21		

(QY determined relative to quinine sulfate in 0.5M H₂SO₄)

INTELLECTUAL PROPERTY

These products are offered in collaboration with ModyBase HB.

REFERENCES

- (1) G. Sheikhnjad, et al., *J Mol Biol*, 1999, **285**, 2021-2034.
- (2) V.E. Marquez, et al., *Antisense Nucleic Acid Drug D*, 1999, **9**, 415-421.

SEE ALSO

Convertible F-dC on page 65
5-Fluoro-2'-deoxyUridine on page 60
Pyrrolidine on page 63

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M

Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

NON-CANONICAL STRUCTURES

DNA and RNA structures are defined by Watson-Crick rules of hybridization. However, a variety of DNA and RNA structures have been defined which do not rely on simple A-T/U and G-C binding. Since these structures disobey the Watson-Crick canon, they are described as non-canonical. Non-canonical DNA and RNA segments are formed as a result of secondary structures. These include G-quadruplexes, triplex forming oligos, hairpins, cruciforms, and i-Motif structures.

G-QUADRUPLEX

Oligonucleotide structural analysis has demonstrated that DNA and RNA nucleic acid sequences containing G-tracts separated by other bases spontaneously fold into G-quadruplex structures. G-quadruplexes are formed when four adjacent guanine residues stack in a cyclic Hoogsteen hydrogen-bonding arrangement leading to four-stranded helical structures. The study of G-quadruplexes in basic genetic processes is an active area of research in telomerase activity, gene regulation, and functional genomics. Guanine analogues that have different hydrogen bonding characteristics - 7-deaza-8-aza-dG and 7-deaza-dG - have proved useful in analyzing G-quadruplex structures. Similarly, common DNA lesions - 8-oxo-dG and abasic sites - have been used to investigate their effect on G-quadruplex structure and activity.

TRIPLEX-FORMING OLIGONUCLEOTIDES

Triplex-forming oligonucleotides (TFO) bind in the major groove of duplex DNA in a sequence-specific manner through the formation of non Watson-Crick (Hoogsteen) hydrogen bonds. The formation of a triplex along the major groove competes with the binding of transcription factors and other proteins that are necessary for transcription, thereby inhibiting the expression of particular genes. A variety of nucleoside analogues have been used in TFO - 8-amino-dG, 8-amino-dA, 6-thio-dG and deoxypseudouridine.

i-MOTIF DNA STRUCTURES

Intercalated Motif (i-Motif) DNA structures may be formed in regions rich in 2'-deoxyCytidine. Especially at acidic pH, these structures could be described as C-Quadruplexes with two parallel stranded sequences also held together in an antiparallel orientation by cytosine-cytosine base pairs. Since these structures are stable at acidic pH, they can act as nanoswitches by change in pH. As they were not considered to be stable at physiological pH, they were not initially considered to be relevant to biological systems. However, the stability of the cytosine-cytosine base pair is enhanced by intercalating ligands and so a variety of i-Motif structures are now considered to be biologically significant. Since i-Motif structures have now been observed forming and dissolving in living cells, these structures are now the subject of active investigation of the meaning of their activity in human cells. Research is also being directed to the effect of common DNA lesions, like depurinated sites, 8-oxo-dG and 5-hydroxymethyl-dC, on these transient structures.

SEE ALSO

7-Deaza-8-Aza-2'-deoxyGuanosine on page 57
8-oxo-2'-deoxyGuanosine on page 61
7-Deaza-2'-deoxyGuanosine on page 57
Abasic II Phosphoramidite on page 62
dSpacer on page 84
8-Amino-dG on page 62
8-Amino-dA on page 59
6-Thio-dG on page 58
2'-deoxypseudoU-CE on page 58
5-Hydroxymethyl-dC on page 50

APTAMER DEVELOPMENT

Aptamers, generated through repetitive selection using SELEX or an equivalent *in vivo* procedure, are chosen for their ability to bind desired target molecules, which are frequently small molecules useful in therapeutics. In some ways, they may be described as chemically engineered versions of antibodies. Of course, nucleic acid aptamers have advantages over antibodies in that they can be developed rapidly by *in vitro* methods, with the reproducibility of chemical synthesis and inherent stability of modified oligonucleotides. A full battery of base, sugar and internucleotide modifications is available for aptamer development.

2'-F-RNA has been used extensively in aptamer development, as well as 2'-F-ANA more recently. An article in The Glen Report by Jeff Carter, Director, Process Chemistry, SomaLogic, Inc. described¹ the use of a DNA backbone with 5-substituted dU analogues as low off-rate modified aptamer (SOMAmer[®]) reagents to enable multiplexed screening of thousands of serum or plasma proteins. These aptamers also include PC Biotin along with a fluorophore, in this case Cyanine 3, for subsequent detection.

REFERENCE

(1) J. Carter, *The Glen Report*, 2015, **27.1**, 6-8.

SEE ALSO

2'-F-RNA Phosphoramidites on page 143
2'-F-Arabinonucleic Acid (2'-F-ANA) on page 144
PC Biotin Phosphoramidite on page 100
Cyanine 3 Phosphoramidite on page 106

TERMINUS MODIFIERS

INTELLECTUAL PROPERTY

5'-Carboxy-Modifier C10 is offered for sale under license from TriLink BioTechnologies, Inc. It is intended for research and development purposes only, and may not be used for commercial, clinical, diagnostic or any other use. It is covered under US Patent No. 6,320,041.

SEE ALSO

PC modifiers on page 86

ABBREVIATIONS








CNEt = Cyanoethyl
CPG = Controlled Pore Glass
DMT = 4,4'-Dimethoxytrityl
Fmoc = Fluorenylmethoxycarbonyl
iPr = Isopropyl
MMT = 4-Monomethoxytrityl
T = Trityl
TFA = Trifluoroacetyl

Glen Research 5'-Modifiers are designed for use in DNA synthesizers to functionalize the 5'-terminus of the target oligonucleotide. The 5'-Amino-Modifiers are available with a variety of chain lengths to fit exactly the desired application.

The DMS(O)MT-protected amino group is easier to deprotect compared to the MMT-protected one. The sulfoxy derivative survives conditions of oligonucleotide synthesis and can either be cleaved with standard deblock solution, or left intact for HPLC purification. At the same time, the DMS(O)MT group is fully compatible with cartridge purification. When detritylation on a cartridge is carried out, the DMS(O)MT+, which is more stable than MMT+, does not reattach itself to an amine. We now offer 5'-DMS(O)MT-Amino-Modifier C6 utilizing this new trityl based protecting group.

5'-Amino-Modifier TEG, a hydrophilic triethylene glycol ethylamine derivative, is 12 atoms in length and fully soluble in aqueous media.

Item	Catalog No.	Pack	Price (\$)
5'-Amino-Modifier C3-TFA	10-1923-90	100 μ mole	50.00
	10-1923-02	0.25g	175.00
5'-Amino-Modifier C6	10-1906-90	100 μ mole	60.00
	10-1906-02	0.25g	200.00
5'-Amino-Modifier C6-TFA	10-1916-90	100 μ mole	30.00
	10-1916-02	0.25g	100.00
5'-Amino-Modifier C12	10-1912-90	100 μ mole	90.00
	10-1912-02	0.25g	300.00
5'-Amino-Modifier 5	10-1905-90	100 μ mole	60.00
	10-1905-02	0.25g	200.00
5'-DMS(O)MT-Amino-Modifier C6	10-1907-90	100 μ mole	60.00
	10-1907-02	0.25g	200.00
5'-Amino-Modifier TEG	10-1917-90	100 μ mole	115.00
	10-1917-02	0.25g	500.00



TERMINUS MODIFIERS (CONT.)




Our more recent 5'-amino modifiers, protected by a novel phthalic acid diamide (PDA) protecting group, are stable solids. In contrast to the TFA protected amino modifiers, which are viscous oils, the analogous PDA protected compounds are granular powders. This important property of these compounds allows straightforward handling, storage and aliquoting and leads to a significant increase in stability.

Deprotection with methylamine in gas phase or aqueous solution or AMA leads to fast and complete removal of the PDA protecting group. However, ammonium hydroxide will not drive the equilibrium reaction to completion and only partial deprotection occurs - overnight deprotection with ammonium hydroxide will yield around 80% active amine.

We are offering three PDA Amino-Modifiers:

- 5'-Amino-Modifier C6-PDA
- Hydrophobic 5'-Amino-Modifier C12-PDA
- Hydrophilic 5'-Amino-Modifier-TEG-PDA

Item	Catalog No.	Pack	Price (\$)
5'-Amino-Modifier C6-PDA	10-1947-90	100 μ mole	30.00
	10-1947-02	0.25g	100.00
5'-Amino-Modifier C12-PDA	10-1948-90	100 μ mole	65.00
	10-1948-02	0.25g	240.00
5'-Amino-Modifier-TEG-PDA	10-1949-90	100 μ mole	105.00
	10-1949-02	0.25g	420.00



INTELLECTUAL PROPERTY

PDA amino-modifiers were developed by Stefan Pitsch and ReseaChem GmbH (S. Berger), Patent pending.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add
Expedite E
MerMade M

Columns
For Instrument type Add
Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

TERMINUS MODIFIERS (CONT.)

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M









Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

The disulfide thiol modifier may be used for introducing 3'- or 5'-thiol linkages. Dithiol Serinol, produced from lipoic acid and our patented serinol backbone, allows easy connection of multiply dithiol-labeled oligos to gold surfaces. 5'-Carboxy-Modifier C10 is a unique linker designed to be added at the terminus of an oligonucleotide synthesis. It generates an activated carboxylic acid N-hydroxysuccinimide (NHS) ester suitable for immediate conjugation on the synthesis column with molecules containing a primary amine, resulting in a stable amide linkage. An alternative carboxylate protecting group on an otherwise fully protected oligonucleotide. The 2-chlorotrityl group is also removed using the standard deblock cycle to generate a free carboxyl group on an otherwise fully protected oligonucleotide. The 2-chlorotrityl group is also removed during oligo deprotection with ammonium hydroxide or AMA and is incompatible with RP purification techniques. PC Amino-Modifier is a photocleavable C6 amino-modifier, part of our line of photocleavable (PC) modifiers. 5'-AminoOxy-Modifier 11 is based on a tetraethylene glycol linkage for improved solubility and for reducing the potential negative impact on hybridization of the oligo. The oxime formed from the conjugation of primary amines with aldehydes is not stable to acidic or basic conditions and requires subsequent reduction with borohydride to form stable amine conjugates. 5'-Maleimide Modifier Phosphoramidite, developed at the University of Barcelona, incorporates a maleimide cycloadduct that is stable to ammonium hydroxide at room temperature. This phosphoramidite can be incorporated into DNA and RNA with both phosphate and phosphorothioate linkages. A retro-Diels-Alder reaction deprotects the maleimide immediately prior to conjugation.

Item	Catalog No.	Pack	Price (\$)
5'-Thiol-Modifier C6	10-1926-90	100 µmole	60.00
	10-1926-02	0.25g	200.00
Thiol-Modifier C6 S-S	10-1936-90	100 µmole	150.00
	10-1936-02	0.25g	360.00
Dithiol Serinol Phosphoramidite	10-1991-95	50 µmole	120.00
	10-1991-90	100 µmole	215.00
	10-1991-02	0.25g	585.00
PC Amino-Modifier Phosphoramidite	10-4906-90	100 µmole	135.00
	10-4906-02	0.25g	395.00
5'-Carboxy-Modifier C10	10-1935-90	100 µmole	65.00
	10-1935-02	0.25g	265.00
5'-Carboxy-Modifier C5	10-1945-90	100 µmole	95.00
	10-1945-02	0.25g	330.00
5'-AminoOxy-Modifier 11	10-1919-95	50 µmole	140.00
	10-1919-90	100 µmole	265.00
	10-1919-02	0.25g	895.00
5'-Maleimide-Modifier Phosphoramidite	10-1938-90	100 µmole	70.00
	10-1938-02	0.25g	335.00




SEQUENCE MODIFIERS

Sequence Modifiers are designed for use in automated synthesis. The carboxy-dT is hydrolyzed during deprotection and can be coupled directly to a molecule containing a primary amino group by a standard peptide coupling or via the intermediate N-hydroxysuccinimide (NHS) ester. Amino-Modifier dA, Amino-Modifier dC, N2-Amino-Modifier dG and both Amino-Modifier dT products can be added in place of a dA, dC, dG and dT residue, respectively, during oligonucleotide synthesis. Corresponding Amino-Modifier supports can replace their respective deoxynucleoside supports. After deprotection, the primary amine on the C6 analogues is separated from the oligonucleotide by a spacer arm with a total of 7-10 atoms and can be labeled or attached to an enzyme. The C2 analogue is more suitable for the attachment of molecules designed to react with the oligonucleotide.

Item	Catalog No.	Pack	Price (\$)
Amino-Modifier C6 dA	10-1089-90	100 µmole	205.00
	10-1089-02	0.25g	455.00
Amino-Modifier C6 dC	10-1019-90	100 µmole	225.00
	10-1019-02	0.25g	450.00
N2-Amino-Modifier C6 dG	10-1529-95	50 µmole	240.00
	10-1529-90	100 µmole	480.00
	10-1529-02	0.25g	1100.00
Carboxy-dT	10-1035-90	100 µmole	180.00
	10-1035-02	0.25g	360.00
Amino-Modifier C2 dT	10-1037-90	100 µmole	180.00
	10-1037-02	0.25g	360.00
	10-1037-05	0.5g	720.00
Amino-Modifier C6 dT	10-1039-90	100 µmole	180.00
	10-1039-02	0.25g	360.00
	10-1039-05	0.5g	720.00

SEE ALSO
Amino-Modifier supports on page 79




SEQUENCE MODIFIERS (CONT.)

Our repertoire of NHS ester derivatives has been expanded to include the NHS-Carboxy-dT-CE Phosphoramidite. By making a dT analog of the Carboxy-Modifier C10, it is possible to label one or multiple sites within an oligonucleotide. This opens up the possibility to label any number of different dyes or molecules within an oligonucleotide when the phosphoramidite is unavailable. Doing so is straightforward and may be done manually off the synthesizer or even in a fully-automated manner on the DNA synthesizer.


We have never found conditions which allow the TFA group to be removed from an amino-modifier while the oligonucleotide remains attached to the support. We are able to solve this problem by using a 9-fluorenylmethoxycarbonyl (Fmoc) protecting group. The Fmoc group is removed using a two step procedure, the first to remove the cyanoethyl protection groups and flush the formed acrylonitrile from the synthesis column using 1% diisopropylamine in acetonitrile, and the second to remove the Fmoc group using 10% piperidine in DMF. The amino group so formed on the column can be reacted with a variety of activated esters. We offer Fmoc-Amino-Modifier C6 dT Phosphoramidite as a nucleosidic option and Amino-Modifier Serinol Phosphoramidite as a non-nucleosidic alternative. We also offer S-Bz-Thiol-Modifier C6-dT to join the ranks of thiol-modifiers for oligonucleotide synthesis. Thiol-Modifier C6-dT can be added as usual at the desired locations within a sequence.

Item	Catalog No.	Pack	Price (\$)
NHS-Carboxy-dT	10-1535-90	100 μ mole	210.00
	10-1535-02	0.25g	550.00
Fmoc-Amino-Modifier C6 dT	10-1536-90	100 μ mole	180.00
	10-1536-02	0.25g	360.00
S-Bz-Thiol-Modifier C6-dT	10-1538-95	50 μ mole	130.00
	10-1538-90	100 μ mole	245.00
	10-1538-02	0.25g	550.00
Amino-Modifier Serinol Phosphoramidite	10-1997-95	50 μ mole	125.00
	10-1997-90	100 μ mole	225.00
	10-1997-02	0.25g	595.00



NHS-Carboxy-dT

Fmoc-Amino-Modifier C6 dT




S-Bz-Thiol-Modifier C6-dT

Amino-Modifier Serinol Phosphoramidite


3'-MODIFIERS

3'-Amino-Modifier CPGs, containing amino groups protected with the base-labile Fmoc group, are designed to functionalize the 3'-terminus of the target oligonucleotide by the introduction of a primary amine. In an alternative approach, the nitrogen destined to become the 3'-amino group is included in a phthalimide (PT) group which is attached to the support through an amide group attached to the aromatic ring. This simple linkage is very stable to all conditions of oligonucleotide synthesis and contains no chiral center. Using an extended ammonium hydroxide treatment (55°C for 17 hours), the cleavage of the amine from the phthalimide is accomplished along with the deprotection of the oligonucleotide. ABI-style columns are supplied unless otherwise requested.

Item	Cat. No.	Pack	Price (\$)
3'-Amino-Modifier C7 CPG 1000	20-2958-01	0.1g	95.00
	20-2958-10	1.0g	675.00
	20-2958-41	Pack of 4	140.00
	20-2958-42	Pack of 4	85.00
	20-2958-13	Pack of 1	250.00
15 μ mole column (Expedite)	20-2958-14	Pack of 1	375.00
3'-Amino-Modifier Serinol CPG	20-2997-01	0.1g	95.00
	20-2997-10	1.0g	675.00
	20-2997-42	Pack of 4	85.00
	20-2997-41	Pack of 4	140.00
	20-2997-13	Pack of 1	250.00
15 μ mole column (Expedite)	20-2997-14	Pack of 1	375.00
3'-PT-Amino-Modifier C3 CPG	20-2954-01	0.1g	95.00
	20-2954-10	1.0g	675.00
	20-2954-41	Pack of 4	140.00
	20-2954-42	Pack of 4	85.00
	20-2954-13	Pack of 1	250.00
15 μ mole column (Expedite)	20-2954-14	Pack of 1	375.00
3'-PT-Amino-Modifier C6 CPG	20-2956-01	0.1g	95.00
	20-2956-10	1.0g	675.00
	20-2956-41	Pack of 4	140.00
	20-2956-42	Pack of 4	85.00
	20-2956-13	Pack of 1	250.00
15 μ mole column (Expedite)	20-2956-14	Pack of 1	375.00
3'-PT-Amino-Modifier C6 PS	26-2956-01	0.1g	125.00
	26-2956-10	1.0g	1025.00
	26-2956-52	Pack of 10	220.00
	26-2956-55	Pack of 10	220.00
200 nmole columns (AB 3900)			
40 nmole columns (AB 3900)			




3'-Amino-Modifier C7 CPG



3'-PT Amino-Modifier C3 CPG

3'-PT Amino-Modifier C6 CPG




3'-Amino-Modifier Serinol CPG


3'-MODIFIERS (CONT.)

The 3'-Thiol-Modifier S-S CPG supports are used to introduce 3'-thiol linkages with three and six atom spacers into oligonucleotides. 3'-Dithiol Serinol CPG is used to introduce a dithiol group at the 3'-terminus. In conjunction with Dithiol Serinol Phosphoramidite, it is simple to produce oligonucleotides with multiple thiol groups at the 3' terminus, which is ideal for conjugation to gold surfaces. With Glyceryl CPG the 3'-terminus of an oligonucleotide is readily oxidized by sodium periodate to form a 3'-phosphoglycaldehyde. The aldehyde may be further oxidized to the corresponding carboxylic acid. Either the aldehyde or the carboxylate may be used for subsequent conjugation to amine-containing products.


Item	Cat. No.	Pack	Price (\$)	
3'-Thiol-Modifier C3 S-S CPG	20-2933-01	0.1g	85.00	
	20-2933-10	1.0g	600.00	
	1 µmole columns	20-2933-41	Pack of 4	125.00
	0.2 µmole columns	20-2933-42	Pack of 4	75.00
	10 µmole column (ABI)	20-2933-13	Pack of 1	225.00
15 µmole column (Expedite)	20-2933-14	Pack of 1	350.00	
3'-Thiol-Modifier 6 S-S CPG	20-2938-01	0.1g	85.00	
	20-2938-10	1.0g	600.00	
	0.2 µmole columns	20-2938-42	Pack of 4	75.00
	1 µmole columns	20-2938-41	Pack of 4	125.00
	10 µmole column (ABI)	20-2938-13	Pack of 1	225.00
15 µmole column (Expedite)	20-2938-14	Pack of 1	350.00	
3'-Dithiol Serinol CPG	20-2991-01	0.1g	120.00	
	20-2991-10	1.0g	995.00	
	0.2 µmole columns	20-2991-42	Pack of 4	120.00
	1 µmole columns	20-2991-41	Pack of 4	200.00
	10 µmole column (ABI)	20-2991-13	Pack of 1	300.00
15 µmole column (Expedite)	20-2991-14	Pack of 1	450.00	
3'-Glyceryl CPG	20-2902-01	0.1g	85.00	
	20-2902-10	1.0g	600.00	
	1 µmole columns	20-2902-41	Pack of 4	125.00
	0.2 µmole columns	20-2902-42	Pack of 4	75.00
	10 µmole column (ABI)	20-2902-13	Pack of 1	225.00
15 µmole column (Expedite)	20-2902-14	Pack of 1	350.00	




3'-Thiol-Modifier C3 S-S CPG



3'-Thiol-Modifier 6 S-S CPG



3'-Dithiol Serinol CPG




3'-Glyceryl CPG


3'-MODIFIERS (CONT.)

3'-Amino-Modifier C6 dC CPG and 3'-Amino-Modifier C6 dT CPG replace a dC and T, respectively, at the 3'-terminus. These products allow convenient labeling at the 3' without blocking the terminus from desired enzymatic activity.

Item	Cat. No.	Pack	Price (\$)	
3'-Amino-Modifier C6 dC CPG	20-2019-01	0.1g	120.00	
	20-2019-10	1.0g	995.00	
	1 µmole columns	20-2019-41	Pack of 4	200.00
	0.2 µmole columns	20-2019-42	Pack of 4	120.00
	10 µmole column (ABI)	20-2019-13	Pack of 1	300.00
15 µmole column (Expedite)	20-2019-14	Pack of 1	450.00	
3'-Amino-Modifier C6 dT CPG	20-2039-01	0.1g	96.00	
	20-2039-10	1.0g	800.00	
	1 µmole columns	20-2039-41	Pack of 4	160.00
	0.2 µmole columns	20-2039-42	Pack of 4	96.00
	10 µmole column (ABI)	20-2039-13	Pack of 1	240.00
15 µmole column (Expedite)	20-2039-14	Pack of 1	360.00	



Amino-Modifier C6 dC CPG




Amino-Modifier C6 dT CPG

CHEMICAL PHOSPHORYLATION


Chemical Phosphorylation Reagent is most commonly used to phosphorylate the 5'-terminus of an oligonucleotide. Although this product is also successful in 3'-phosphorylation, 3'-Phosphate CPG allows direct preparation of oligonucleotides with a 3'-phosphate group. Chemical Phosphorylation Reagent II contains a DMT group on a side chain which is stable to base cleavage and can be left on the oligonucleotide for use in RP purification. The DMT group is later removed with aqueous acid and the side chain is eliminated after brief treatment with aqueous ammonium hydroxide to yield the 5'-phosphate.¹ Solid CPR II is similar in performance to CPR II but it is easier to prepare aliquots since it is a powder. Many researchers treat synthesis supports with a hindered base (e.g., diethylamine, diisopropylethylamine, or DBU) post-synthesis to eliminate and remove the cyanoethyl phosphate groups. In this way, the acrylonitrile formed in situ is removed from the support and is not available to alkylate dT residues at the N3 position in the oligos. Since the sulfonyl ethyl group in 3'-Phosphate CPG is also susceptible to β -elimination leading to oligo cleavage, this technique is not compatible with 3'-phosphate CPG. Using CPR II CPG, which is base labile but does not support β -elimination, the cyanoethyl groups can be removed from the oligo prior to cleavage and base deprotection. ABI-style vials and columns are supplied unless otherwise requested.

Item	Cat. No.	Pack	Price (\$)
Chemical Phosphorylation Reagent	10-1900-90	100 μ mole	50.00
	10-1900-02	0.25g	160.00
3'-Phosphate CPG	20-2900-01	0.1g	70.00
	20-2900-10	1.0g	480.00
	1 μ mole columns	Pack of 4	100.00
	0.2 μ mole columns	Pack of 4	60.00
	10 μ mole column (ABI)	Pack of 1	180.00
	15 μ mole column (Expedite)	Pack of 1	280.00
3'-Phosphate PS	26-2900-01	0.1g	75.00
	26-2900-10	1.0g	510.00
	200 nmole columns (AB 3900)	Pack of 10	150.00
	40 nmole columns (AB 3900)	Pack of 10	150.00
3'-Phosphate CPG (High Load)	25-2900-01	0.1g	85.00
	25-2900-10	1.0g	600.00
	2.5 μ mole columns	Pack of 4	120.00
Chemical Phosphorylation Reagent II (CPR II)	10-1901-90	100 μ mole	60.00
	10-1901-02	0.25g	200.00
Solid Chemical Phosphorylation Reagent II (Solid CPR II)	10-1902-90	100 μ mole	60.00
	10-1902-02	0.25g	200.00
3'-CPR II CPG	20-2903-01	0.1g	70.00
	20-2903-10	1.0g	480.00
	0.2 μ mole columns	Pack of 4	60.00
	1 μ mole columns	Pack of 4	100.00
	10 μ mole column (ABI)	Pack of 1	180.00
	15 μ mole column (Expedite)	Pack of 1	280.00




Chemical Phosphorylation Reagent


82




3'-Phosphate CPG



Chemical Phosphorylation Reagent II



Solid Chemical Phosphorylation Reagent II



3'-CPR II CPG


ALDEHYDE MODIFICATION

Aldehyde modifiers would be attractive electrophilic substitutions in oligonucleotides since they are able to react with amino groups to form a Schiff's base, with hydrazino groups to form hydrazones, and with semicarbazides to form semicarbazones. The Schiff's base is unstable and must be reduced with sodium borohydride to form a stable linkage but hydrazones and semicarbazides are very stable linkages.


Our collaboration with ELITechGroup, formerly Epoch Biosciences, has allowed us to offer 5'-Aldehyde-Modifier C2 Phosphoramidite. The acetal protecting group is sufficiently hydrophobic for use in RP HPLC and cartridge purification and is readily removed after oligonucleotide synthesis under standard oligonucleotide detritylation conditions with 80% acetic acid / 20% water or 2% aqueous trifluoroacetic acid during cartridge purification.

A formylindole nucleoside analogue has been used to introduce aldehyde groups within an oligonucleotide or at the 5' terminus. This product has no protecting group on the aldehyde, which means that deprotection of the modified oligonucleotide can be done without changing preferred conditions.

Item	Cat. No.	Pack	Price (\$)
5'-Aldehyde-Modifier C2 Phosphoramidite	10-1933-90	100 μ mole	85.00
	10-1933-02	0.25g	325.00
Formylindole CE Phosphoramidite	10-1934-90	100 μ mole	85.00
	10-1934-02	0.25g	325.00



5'-Aldehyde-Modifier C2



Formylindole

INTELLECTUAL PROPERTY

These Products are for research purposes only, and may not be used for commercial, clinical, diagnostic or any other use. These Products are subject to proprietary rights of ELITechGroup and are made and sold under license from ELITechGroup. There is no implied license for commercial use with respect to these Products and a license must be obtained directly from ELITechGroup with respect to any proposed commercial use of these Products. "Commercial use" includes but is not limited to the sale, lease, license or other transfer of the Product or any material derived or produced from it, the sale, lease, license or other grant of rights to use the Product or any material derived or produced from it, or the use of the Product to perform services for a fee for third parties (including contract research).

A simple agreement must be signed before end-users and custom oligo services may purchase these products for use as defined above.
<http://www.glenresearch.com/Reference/ELITechGroupProducts.pdf>

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

83

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers For Instrument type	Add
Expedite	E
MerMade	M

Columns For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

INTELLECTUAL PROPERTY

Glen Research offers PC Biotin, PC Amino-Modifier and PC Spacer products in association with AmberGen, Inc. and Link Technologies, Ltd. For a commercial application license, please contact AmberGen, Inc., +617-923-9990, (sales@ambergenc.com), http://www.ambergenc.com/.

PC Linker phosphoramidite is available from Glen Research in association with Link Technologies Ltd (Scotland).

SEE ALSO

5'-Biotin on page 99

REFERENCES


- P. Ordoukhanian and J.-S. Taylor, *J. Am. Chem. Soc.*, **117**, 9570-9571, 1995.
- (2a) F. Hausch and A. Jäschke, *Nucleic Acids Research*, 2000, **28**, e35.
- (2b) F. Hausch and A. Jäschke, *Tetrahedron*, 2001, **57**, 1261-1268.
- T. Wenzel, T. Elssner, K. Fahr, J. Bimmler, S. Richter, I. Thomas, and M. Kostrzewa, *Nucleosides, Nucleotides & Nucleic Acids*, 2003, **22**, 1579-1581.

PHOTOCLEAVABLE MONOMERS

PC Biotin Phosphoramidite can be used to prepare 5'-biotinylated oligonucleotides suitable for capture by streptavidin in a mode similar to our popular 5' Biotin Phosphoramidite. Amino- and thiol-modified oligonucleotides have proven to be very useful for the attachment of a variety of haptens and fluorophores, as well as for the tethering of the oligonucleotides to a diversity of beads and surfaces. PC Amino-Modifier Phosphoramidite is used to prepare 5'-amino-modified oligonucleotides suitable for subsequent photocleavage. PC Spacer Phosphoramidite can be used as an intermediary to attach any modification reagent, available as a phosphoramidite, to the terminus of oligonucleotides. After photocleavage, a 5'-phosphate is generated on the DNA, rendering it suitable for further biological transformations, such as gene construction and cloning after ligation.

A versatile photocleavable DNA building block has been described by researchers in Washington University, Missouri and used in phototriggered hybridization.¹ This reagent has also been used in the design of multifunctional DNA and RNA conjugates² for the in vitro selection of new molecules catalyzing biomolecular reactions. Researchers at Bruker Daltonik in Germany have also developed genoSNIP, a method for single-nucleotide polymorphism (SNP) genotyping by MALDI-TOF mass spectrometry.³ This method uses size reduction of primer extension products by incorporation of the photocleavable linker for phototriggering strand breaks near to the 3' end of the extension primer. PC Linker can be incorporated into oligonucleotides at any position by standard automated DNA synthesis methodology. PC Linker Phosphoramidite has the added advantage in that photocleavage results in monophosphate fragments at both the 3'- and 5'-termini of the oligonucleotide fragments.

Item	Catalog No.	Pack	Price(\$)
PC Biotin Phosphoramidite	10-4950-95	50 µmole	145.00
	10-4950-90	100 µmole	280.00
	10-4950-02	0.25g	675.00
PC Amino-Modifier Phosphoramidite	10-4906-90	100 µmole	135.00
	10-4906-02	0.25g	395.00
PC Spacer Phosphoramidite	10-4913-90	100 µmole	135.00
	10-4913-02	0.25g	395.00
PC Linker Phosphoramidite	10-4920-90	100 µmole	255.00
	10-4920-02	0.25g	795.00




CONJUGATION USING CLICK CHEMISTRY

The copper(I)-catalyzed azide-alkyne cycloaddition (CuAAC) reaction between azides and alkynes to form 1,2,3-triazoles, as reported¹ by Sharpless, was found to be so exquisitely regioselective and efficient at even the most mild conditions that Sharpless coined the term 'Click Chemistry' to describe it. The use of this method for DNA modification has been somewhat delayed by the fact that copper ions damage DNA, typically yielding strand breaks.² As these problems have now been overcome by the use of copper(I)-stabilizing ligands (e.g., tris(benzyltriazolylmethyl)amine, TBTA³), Carell et al. and Seela et al. discovered that the CuAAC reaction can be used to functionalize alkyne-modified DNA nucleobases with extremely high efficiency.⁴

Oligonucleotides bearing a single nucleosidic alkyne group can be prepared using a C8-Alkyne-dC or dT-CE Phosphoramidite. Purified oligonucleotides are usually modified with 2-5 equivalents of the corresponding marker-azide (e.g., fluorescent-dye azides). After the addition of precomplexed Cu(I), complete conversion to the labeled oligo is observed in a time span between 30 min and 4 hours. After a simple precipitation step, labeled oligonucleotides can be recovered in near quantitative yields. Using a combination of C8-Alkyne, C8-TIPS-Alkyne and C8-TMS-Alkyne, it is possible to label oligonucleotides in up to three separate click reactions. The alkyne groups on the last two monomers are protected, respectively, with triisopropylsilyl (TIPS) and trimethylsilyl (TMS) protecting groups.^{5,6} The first click reaction on solid phase on a C8-Alkyne yields the singly modified oligonucleotide with full retention of the TIPS and/or TMS protecting group. For double click, a C8-TIPS-Alkyne is used as the second nucleoside and the TIPS protecting group is cleaved with tetrabutylammonium fluoride (TBAF) without causing any damage to the DNA. The second click reaction in solution yields the doubly modified oligonucleotide in excellent yield. For the introduction of three different labels, all three nucleosides are introduced into oligonucleotides. The first reaction is performed directly on the resin. The singly modified oligonucleotide is subsequently cleaved from the support with concomitant cleavage of the TMS group and retention of the TIPS protecting group. The second click reaction is performed in solution. Precipitation of the doubly modified oligonucleotide, cleavage of the TIPS group with TBAF, and a subsequent third click reaction in solution furnishes the desired triply modified oligonucleotide in excellent overall yield.

Item	Catalog No.	Pack	Price(\$)
C8-Alkyne-dT-CE Phosphoramidite	10-1540-95	50 µmole	165.00
	10-1540-90	100 µmole	315.00
	10-1540-02	0.25g	900.00
C8-TIPS-Alkyne-dC-CE Phosphoramidite	10-1541-95	50 µmole	295.00
	10-1541-90	100 µmole	575.00
	10-1541-02	0.25g	1275.00
C8-TMS-Alkyne-dC-CE Phosphoramidite	10-1542-95	50 µmole	270.00
	10-1542-90	100 µmole	525.00
	10-1542-02	0.25g	1275.00
C8-Alkyne-dC-CE Phosphoramidite	10-1543-95	50 µmole	225.00
	10-1543-90	100 µmole	435.00
	10-1543-02	0.25g	1125.00



REFERENCES

- C.W. Tornøe, C. Christensen, M. Meldal, *J. Org. Chem.* 2002, **67**, 3057-3064; V. Rostovtsev, L. G. Green, V. V. Fokin, K. B. Sharpless, *Angew. Chem.* 2002, **114**, 2708-2711; *Angew. Chem. Int. Ed.* 2002, **41**, 2596-2599.
- C. J. Burrows, J. G. Muller, *Chem. Rev.* 1998, **98**, 1109 - 1151.
- T. R. Chan, R. Hilgraf, K. B. Sharpless, V. V. Fokin, *Org. Lett.* 2004, **6**, 2853 - 2855.
- J. Gierlich, G. A. Burley, P. M. E. Gramlich, D. M. Hammond, T. Carell, *Org. Lett.* 2006, **8**, 3639-3642. F. Seela, V. R. Sirivolu, *Chem. Biodiversity* 2006, **3**, 509-514.
- P. M. E. Gramlich, S. Warncke, J. Gierlich, T. Carell, *Angew. Chem.* 2008, **120**, 3491-3493; *Angew. Chem. Int. Ed.* 2008, **47**, 3442-3444.
- P. M. E. Gramlich, C. T. Wirges, A. Manetto, T. Carell, *Angew. Chem. Int. Ed.* 2008, **47**, 8350-8358.

INTELLECTUAL PROPERTY

baseclick GmbH has been granted the following patents (1-3) besides its further patent applications (4-5).

- WO 2006/117161 (New labeling strategies for the sensitive detection of analytes)
- WO 2008/952775 (Click chemistry for the production of reporter molecules)
- WO 2010/115957 (Click Chemistry on heterogeneous catalysts)
- PCT/EP 2013/064610 (Anandamide-modified nucleic molecules)
- PCT/EP 2015/056007 (Self-assembly of DNA Origami: a diagnostic tool)

baseclick GmbH holds a worldwide exclusive license for granted patent application WO 03/101972 (Copper-catalysed ligation of azides and acetylenes for the nucleic acid field) in the area of diagnostics and research.

As Glen Research and baseclick are partners, Glen Research is now able to help in sublicensing this outstanding technology.

CONJUGATION USING CLICK CHEMISTRY (CONT.)

5-Ethynyl-dU offers convenient click conjugation with an azide to generate a label rigidly attached to one of the oligonucleotide bases. 5-Ethynyl-dU is subject to base-catalyzed hydration during cleavage and deprotection, especially when using a strong base or heat. Hydration of an ethynyl group forms a methyl ketone which subsequently blocks potential click reactions. Mild deprotection conditions are necessary when using 5-Ethynyl-dU-CE Phosphoramidite to prevent this side reaction. TIPS-5-Ethynyl-dU-CE Phosphoramidite, containing a protected alkyne, offers broader compatibility with oligonucleotide synthesis and deprotection. Protecting the 5-ethynyl group with a triisopropylsilyl (TIPS) protecting group prevents acid or base catalyzed hydration during oligonucleotide synthesis and workup. A quick treatment with TBAF removes the TIPS protecting group.

SEE ALSO

3'-Propargyl-5-Me-dC CPG on page 64

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers

For Instrument type Add

Expedite E
MerMade M


Columns

For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)


Item	Catalog No.	Pack	Price(\$)
C8-TIPS-Alkyne-dT-CE Phosphoramidite	10-1544-95	50 μmole	220.00
	10-1544-90	100 μmole	425.00
	10-1544-02	0.25g	1020.00
C8-TMS-Alkyne-dT-CE Phosphoramidite	10-1545-95	50 μmole	205.00
	10-1545-90	100 μmole	395.00
	10-1545-02	0.25g	1050.00
5-Ethynyl-dU-CE Phosphoramidite	10-1554-95	50 μmole	130.00
	10-1554-90	100 μmole	245.00
	10-1554-02	0.25g	775.00
TIPS-5-Ethynyl-dU-CE Phosphoramidite	10-1555-95	50 μmole	195.00
	10-1555-90	100 μmole	370.00
	10-1555-02	0.25g	975.00
THPTA Ligand (Water soluble)	50-1004-92	25 μmole	50.00
	50-1004-90	100 μmole	180.00
Click-Solution (DMSO/t-BuOH)	50-1002-11	10 x 1.0mL	185.00



CONJUGATION USING CLICK CHEMISTRY (CONT.)

Oligonucleotides prepared using 5'-Hexynyl Phosphoramidite are stable to standard deprotection conditions and exhibit a slightly increased retention time on RP HPLC. Azides are not compatible with oligonucleotide synthesis using phosphoramidites so a post-synthesis reaction is required. Azidobutyrate NHS Ester is used¹ for azido-modification of amines at either the 3'-end or the 5'-end of an oligo and it can even be used for internal modification on an Amino-Modifier-C6 dX residue within the sequence. Specific to the 5'-terminus, 5'-Bromohexyl Phosphoramidite is added in the last cycle. This modifier can then be easily transformed into a 5'-azido group by displacement of bromide using sodium azide.² Alkyne NHS ester allows the functionalization of an amino moiety in a variety of molecules, including DNA and RNA oligonucleotides as well as peptides or proteins. We also offer two products for use in Click Chemistry based upon our 1,3-diol product portfolio with the serinol backbone - a phosphoramidite for adding an alkyne group at the 5' terminus or within the sequence, and a synthesis support for labeling the 3' terminus of oligonucleotides with an alkyne group.

Item	Catalog No.	Pack	Price(\$)
5'-Hexynyl Phosphoramidite	10-1908-90	100 μmole	60.00
	10-1908-02	0.25g	200.00
Azidobutyrate NHS Ester (Dissolve 2.3mg in 60μL of DMSO)	50-1904-23	2.3mg	60.00
	50-1904-24	23mg	300.00
5'-Bromohexyl Phosphoramidite	10-1946-90	100 μmole	60.00
	10-1946-02	0.25g	200.00
Alkyne-NHS Ester (Dissolve 2.3mg in 60μL of DMSO)	50-1905-23	2.3mg	60.00
	50-1905-24	23mg	300.00
Alkyne-Modifier Serinol Phosphoramidite	10-1992-95	50 μmole	100.00
	10-1992-90	100 μmole	185.00
	10-1992-02	0.25g	575.00
3'-Alkyne-Modifier Serinol CPG	20-2992-01	0.1g	105.00
	20-2992-10	1.0g	800.00
	0.2 μmole columns		
	20-2992-42	Pack of 4	100.00
	1 μmole columns		
	20-2992-41	Pack of 4	175.00
10 μmole column (ABI)			
20-2992-13	Pack of 1	260.00	
15 μmole column (Expedite)			
20-2992-14	Pack of 1	390.00	



REFERENCES

- (1) R. Kumar, et al., *Journal of the American Chemical Society*, 2007, **129**, 6859-6864.
- (2) J. Lietard, A. Meyer, J.J. Vasseur, and F. Morvan, *Tetrahedron Letters*, 2007, **48**, 8795-8798.

SEE ALSO

Serinol Products on page 94

CONJUGATION USING CLICK CHEMISTRY (CONT.)

1-Ethylnyl-dSpacer CE Phosphoramidite can be used in any position within an oligonucleotide while still retaining the high efficiency of click chemistry. The modifier is efficiently incorporated into oligonucleotides using standard phosphoramidite chemistry, is stable to common deprotection conditions, and is compatible with Glen-Pak™ purification. 1-Ethylnyl-dSpacer generates a substituted 1,2,3-triazole pseudo-nucleobase after click chemistry conjugation with an azide. The 1-ethylnyl-dSpacer modification exhibits similar duplex stability to the standard dSpacer (10-1914) and destabilizes the duplex when internally incorporated. Upon cycloaddition, the duplex stability is moderated by the resulting structure of the modification. Simple 1,2,3-triazoles were destabilizing, as were modifications that incorporated TEG linkers (6-FAM-TEG and Amino-TEG). Modifications that incorporated aromatic functional groups restored duplex stability to varying degrees with coumarin and psoralen significantly restoring stability. A 5'-iodo-modified oligonucleotide (prepared using 5'-Iodo-dT) can be quantitatively converted to the corresponding 5'-azide.

SEE ALSO

dSpacer on page 84

STABILITY NOTES


Oligonucleotides containing a 5'-iodo group are prepared conventionally with the exception that deprotection is carried out in ammonium hydroxide at room temperature for 24 hours. Under these conditions, degradation of the iodo group was less than 2%.

OLIGO-CLICK KITS


Oligo-Click Kits contain an air-stable, insoluble Cu(I) source in pellet form in a pre-loaded and ready-to-use vial. Within the kit, the TBTA ligand is replaced by an activator which is compatible with both aqueous and organic solvents. This innovative combination of catalyst and ligand/activator results in a much easier labeling work-flow of only three simple steps. The preparation of the oligonucleotide labeling via CuAAC now requires only a minimal hands-on time of a few minutes or even less and can be carried out in air without any additional precautions. Glen Research is offering the following kits in collaboration with baseclick GmbH.

- Oligo-Kit M Reload: This kit has sufficient reagents for conjugating up to nine alkyne-containing oligonucleotides on a 100 nmole scale or a single oligonucleotide on a 1 μmole scale. *The user must supply the azide and a solvent such as DMSO for dissolving the azide.*
- Oligo-Kit M Biotin, Oligo-Kit M Fluorescein and Oligo-Kit M TAMRA: Each kit has sufficient reagents for conjugating up to seven alkyne-containing oligonucleotides on a 100 nmole scale or a single oligonucleotide on a 1 μmole scale. *Each kit contains all of the ingredients necessary, including the azide and DMSO solvent.*

Item	Catalog No.	Pack	Price(\$)
baseclick Oligo-Click-M-Reload	50-2100-01	each	120.00
baseclick Oligo-Click-M-Biotin	50-2101-01	each	200.00
baseclick Oligo-Click-M-Fluorescein	50-2102-01	each	240.00
baseclick Oligo-Click-M-TAMRA	50-2103-01	each	270.00



1-Ethylnyl-dSpacer



5'-I-dT

COPPER-FREE CLICK CHEMISTRY


At Glen Research, our goal was to offer a copper-free click phosphoramidite reagent with the following properties:

- Simple to use
- Stable in solution on the synthesizer
- Stable to ammonium hydroxide and AMA
- Excellent click performance in 17 hours or less at room temperature


From the variety of cyclooctyne-based copper-free click reagents so far described, we have chosen to offer compounds based on a dibenzo-cyclooctyne (DBCO) structure. We are offering 5'-DBCO-TEG Phosphoramidite for preparing oligos with a 5'-DBCO modification and DBCO-dT-CE Phosphoramidite for inserting a DBCO group at any position within the oligonucleotide. In addition, we offer a further DBCO phosphoramidite – DBCO-Serinol Phosphoramidite. Using our proprietary serinol backbone as a non-nucleosidic spacer allows the DBCO group to be placed at any location within a sequence with multiple additions clearly possible. DBCO-sulfo-NHS Ester is also offered for post-synthesis conjugation reactions. DBCO-modified oligos may be conjugated with azides in organic solvents, such as DMSO, or aqueous buffers. Depending on the azide used, the reaction will go to completion in 4-17 hours at room temperature. Simple desalting on a Glen Gel-Pak™ leads to a product with virtually quantitative conjugation efficiency.

Note: We now recommend that synthesis of oligos containing DBCO-dT be completed using 0.5 M CSO in anhydrous acetonitrile (40-4632-xx). Acceptable results can be achieved with iodine oxidation if DBCO-dT is subjected to no more than 8-10 cycles.


Item	Catalog No.	Pack	Price(\$)
5'-DBCO-TEG Phosphoramidite	10-1941-95	50 μmole	125.00
	10-1941-90	100 μmole	230.00
	10-1941-02	0.25g	775.00
DBCO-dT-CE Phosphoramidite	10-1539-95	50 μmole	250.00
	10-1539-90	100 μmole	485.00
	10-1539-02	0.25g	975.00
DBCO-sulfo-NHS Ester (Dissolve 5.2mg in 60μL water or DMSO)	50-1941-23	5.2mg	60.00
	50-1941-24	52mg	300.00
DBCO-Serinol Phosphoramidite	10-1998-95	50 μmole	180.00
	10-1998-90	100 μmole	340.00
	10-1998-02	0.25g	895.00




5'-DBCO-TEG



DBCO-dT



DBCO-Serinol



DBCO-sulfo-NHS Ester

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers For Instrument type	Add
Expedite	E
MerMade	M
Columns For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

SEE ALSO

0.5M CSO on page 32
Serinol Products on page 94

CONJUGATION USING CLICK CHEMISTRY (CONT.)


Glen Research is offering first our most popular labels for general interest and, subsequently, we will add azide products that are not compatible with phosphoramidite chemistry.

Biotin is still our most commonly used label and biotinTEG, with its hydrophilic triethylene glycol spacer, is the most popular biotin product. Desthiobiotin is a biotin analogue that is well captured by streptavidin but the captured product can be easily released by applying a biotin solution to the streptavidin beads. 6-FAM is our most popular fluorescein derivative and we offer azides of both 6-FAM and pivaloyl-protected 6-FAM for situations where subsequent reactions require the 6-FAM to be protected. In both 6-FAM products, the hydrophilic TEG spacer is again used. The azides are offered in 25 and 100 μ mole packs for convenient oligonucleotide labeling.

7-Hydroxycoumarin, also known as umbelliferone, is a highly fluorescent, pH-sensitive fluorophore that emits in the blue region of the spectrum. However, its fluorescence is strongly quenched if the hydroxyl is alkylated or phosphorylated, making it useful in high-throughput screening for phosphatases and lipases. Interestingly, it was found that the 3-azido derivative is also highly quenched but, upon reaction with an alkyne in the presence of copper to form the triazole, the fluorescence is restored.¹ The clicked coumarin emits at a lambda max of 480 nm and absorbs at 358 nm.

HEX and TET are two of our most popular fluorescein-based dyes for labeling oligonucleotides. We are happy to offer 6-HEX and 6-TET Azides for use in click conjugations.

Item	Catalog No.	Pack	Price(\$)
BiotinTEG Azide	50-2000-92	25 μ mole	150.00
	50-2000-90	100 μ mole	450.00
DesthiobiotinTEG Azide	50-2001-92	25 μ mole	135.00
	50-2001-90	100 μ mole	400.00
Dipivaloyl 6-FAM-TEG Azide	50-2002-92	25 μ mole	230.00
	50-2002-90	100 μ mole	690.00
6-FAM-TEG Azide	50-2003-92	25 μ mole	180.00
	50-2003-90	100 μ mole	540.00
Coumarin Azide	50-2004-92	25 μ mole	115.00
	50-2004-90	100 μ mole	350.00
6-HEX Azide	50-2005-92	25 μ mole	150.00
	50-2005-90	100 μ mole	450.00
6-TET Azide	50-2006-92	25 μ mole	150.00
	50-2006-90	100 μ mole	450.00




CONJUGATION USING CLICK CHEMISTRY (CONT.)

Two nitroxide spin labels, TEMPO Azide and TEMPO-TEG Azide, for site directed spin labeling (SDSL) are now offered.

Click Chemistry with psoralen azide and one of our many nucleosidic and non-nucleosidic alkyne derivatives has the potential to generate a variety of practical cross-linkers. The well known reversible cross-linking behavior of psoralen with an adjacent thymidine residue could be very useful.

To better address applications in near-infrared (NIR) imaging, Glen Research is offering a water soluble Disulfo-Cyanine 7 azide that can be easily conjugated to DNA and RNA through standard click chemistry. This long wavelength dye offers the benefits of improved solubility, reduced aggregation, and improved stability in the near-infrared spectrum along with the convenience of click chemistry.

Item	Catalog No.	Pack	Price(\$)
TEMPO Azide	50-2007-92	25 μ mole	115.00
	50-2007-90	100 μ mole	350.00
TEMPO-TEG Azide	50-2008-92	25 μ mole	135.00
	50-2008-90	100 μ mole	400.00
Psoralen Azide	50-2009-92	25 μ mole	115.00
	50-2009-90	100 μ mole	350.00
Disulfo-Cyanine 7 Azide	50-2010-92	25 μ mole	325.00
	50-2010-90	100 μ mole	975.00



REFERENCE

(1) J. Gierlich, G.A. Burley, P.M. Gramlich, D.M. Hammond, and T. Carell, *Org Lett*, 2006, **8**, 3639-42.

SERINOL REAGENTS FOR MODIFICATION AND LABELING

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

INTELLECTUAL PROPERTY

Serinol Reagents for Modification and Labeling are covered by US Patent No.: 8,394,948.


Most popular non-nucleosidic phosphoramidites for modification and labeling are based on two structural types: 1,2-diols and 1,3-diols. Products based on a 1,2-diol backbone were first described to allow amino-modification and biotin labeling. Technically, the 1,2-diol backbone has some drawbacks relative to the 1,3-diol backbone. The 1,2-diol backbone can participate in a dephosphorylation reaction since the 1,2-diol can form a favored 5-membered cyclic phosphate intermediate. This reaction is competitive with simple hydrolysis of the protecting groups and leads to some loss of label. However, the degree of loss at the 3' terminus can be limited by the removal of the cyanoethyl protecting group using DBU or diethylamine prior to the cleavage and deprotection steps. Similarly, loss at the 5' terminus can be eliminated by retaining the DMT group until the oligo is fully deprotected. Fortunately, the elimination reaction is virtually non-existent in the 1,3-diol backbone since the cyclic intermediate would be a 6-membered ring which is not favored for a cyclic phosphate intermediate.

IVD customers have requested a new backbone based on a 1,3-diol that would overcome any technical or IP issues surrounding our current products. We now offer a line of products based on the serinol backbone, which have been developed in close collaboration between Glen Research and Nelson Biotechnologies. Protected Biotin Serinol Phosphoramidite and CPG are protected with a *t*-butylbenzoyl group on the biotin ring. This group is designed to stop any phosphoramidite reactions at this active position in biotin. This protection avoids branching when using nucleophilic activators like DCI. The protecting group is easily removed during oligonucleotide cleavage and deprotection. The BiotinLC versions are similarly protected and should be useful for the synthesis of highly sensitive biotinylated probes. 6-Fluorescein Serinol Phosphoramidite and CPG are designed to prepare oligonucleotides containing one or several 6-Fluorescein (6-FAM) residues. Amino-Modifier Serinol Phosphoramidite and CPG are used to add amino groups into one or several positions in oligonucleotides. The amino group is protected with Fmoc, which may be removed on the synthesis column prior to solid-phase conjugation to the amino groups, or which may be removed during deprotection for subsequent solution phase conjugation to the amino groups.


Combining lipoic acid and our patented serinol backbone, we now offer Dithiol Serinol Phosphoramidite and the related 3'-Dithiol Serinol CPG. This unique architecture moves the bulky dithiol away from the phosphate backbone, making it suitable for conjugation to gold surfaces. The long spacer arm of Dithiol Serinol also allows multiple consecutive incorporations of the modifier without the need for intermediate spacer phosphoramidite additions to achieve optimal stepwise coupling efficiency.

We offer three products for use in Click Chemistry based upon our 1,3-diol product portfolio with the serinol backbone - a phosphoramidite for adding an alkyne group at the 5' terminus or within the sequence, a synthesis support for labeling the 3' terminus of oligonucleotides with an alkyne group, and DBCO-Serinol phosphoramidite as a copper-free click reagent.

Item	Catalog No.	Pack	Price(\$)
Protected Biotin Serinol Phosphoramidite	10-1993-95	50 µmole	165.00
	10-1993-90	100 µmole	295.00
	10-1993-02	0.25g	675.00
6-Fluorescein Serinol Phosphoramidite	10-1994-95	50 µmole	165.00
	10-1994-90	100 µmole	295.00
	10-1994-02	0.25g	595.00



Protected Biotin Serinol Phosphoramidite




6-Fluorescein Serinol Phosphoramidite

SERINOL REAGENTS FOR MODIFICATION AND LABELING (CONT.)


Item	Catalog No.	Pack	Price(\$)
Protected BiotinLC Serinol Phosphoramidite	10-1995-95	50 µmole	205.00
	10-1995-90	100 µmole	365.00
	10-1995-02	0.25g	675.00
Amino-Modifier Serinol Phosphoramidite	10-1997-95	50 µmole	125.00
	10-1997-90	100 µmole	225.00
	10-1997-02	0.25g	595.00
Dithiol Serinol Phosphoramidite	10-1991-95	50 µmole	120.00
	10-1991-90	100 µmole	215.00
	10-1991-02	0.25g	585.00
Alkyne-Modifier Serinol Phosphoramidite	10-1992-95	50 µmole	100.00
	10-1992-90	100 µmole	185.00
	10-1992-02	0.25g	575.00
DBCO-Serinol Phosphoramidite	10-1998-95	50 µmole	180.00
	10-1998-90	100 µmole	340.00
	10-1998-02	0.25g	895.00

SEE ALSO


DBCO on page 91




Protected BiotinLC Serinol Phosphoramidite




Amino-Modifier Serinol Phosphoramidite



Dithiol Serinol



Alkyne-Modifier Serinol Phosphoramidite



DBCO-Serinol

SERINOL REAGENTS FOR MODIFICATION AND LABELING (CONT.)

Item	Catalog No.	Pack	Price(\$)
3'-Protected Biotin Serinol CPG	20-2993-01	0.1g	120.00
	20-2993-10	1.0g	995.00
	0.2 μmole columns	Pack of 4	120.00
	1 μmole columns	Pack of 4	200.00
	10 μmole column (ABI)	Pack of 1	300.00
15 μmole column (Expedite)	20-2993-14	Pack of 1	450.00
3'-6-Fluorescein Serinol CPG	20-2994-01	0.1g	120.00
	20-2994-10	1.0g	995.00
	0.2 μmole columns	Pack of 4	120.00
	1 μmole columns	Pack of 4	200.00
	10 μmole column (ABI)	Pack of 1	300.00
15 μmole column (Expedite)	20-2994-14	Pack of 1	450.00
3'-Protected BiotinLC Serinol CPG	20-2995-01	0.1g	120.00
	20-2995-10	1.0g	995.00
	0.2 μmole columns	Pack of 4	120.00
	1 μmole columns	Pack of 4	200.00
	10 μmole column (ABI)	Pack of 1	300.00
15 μmole column (Expedite)	20-2995-14	Pack of 1	450.00
3'-Amino-Modifier Serinol CPG	20-2997-01	0.1g	95.00
	20-2997-10	1.0g	675.00
	0.2 μmole columns	Pack of 4	85.00
	1 μmole columns	Pack of 4	140.00
	10 μmole column (ABI)	Pack of 1	250.00
15 μmole column (Expedite)	20-2997-14	Pack of 1	375.00

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)




SERINOL REAGENTS FOR MODIFICATION AND LABELING (CONT.)

Item	Catalog No.	Pack	Price(\$)
3'-Dithiol Serinol CPG	20-2991-01	0.1g	120.00
	20-2991-10	1.0g	995.00
	0.2 μmole columns	Pack of 4	120.00
	1 μmole columns	Pack of 4	200.00
	10 μmole column (ABI)	Pack of 1	300.00
15 μmole column (Expedite)	20-2991-14	Pack of 1	450.00
3'-Alkyne-Modifier Serinol CPG	20-2992-01	0.1g	105.00
	20-2992-10	1.0g	800.00
	0.2 μmole columns	Pack of 4	100.00
	1 μmole columns	Pack of 4	175.00
	10 μmole column (ABI)	Pack of 1	260.00
15 μmole column (Expedite)	20-2992-14	Pack of 1	390.00

COT SERINOL PHOSPHoramidite

Bright, long-lasting and non-phototoxic organic fluorophores are essential for the continued optimization of a diverse range of imaging applications. However, all currently available technologies remain susceptible to undesirable transitions to dark states. Dark states arise from non-fluorescent triplet electronic configurations from which the rate of return to the ground state is slow. When in the triplet state, the fluorophore is susceptible to photobleaching and fluorescence applications are compromised by unpredictably reducing the signal-to-noise ratio (SNR), as well as limiting the total duration of time over which information can be gathered. The direct conjugation of small-molecule protective agents (PAs) has enabled significant improvements through intra-molecular triplet quenching. Through a partnership with Lumidyne Technologies, Glen Research has created a novel PA-linked phosphoramidite using cyclooctatetraene (COT). COT Serinol Phosphoramidite provides a means to improve the photostability of virtually any fluorophore in a modular fashion. Our spectrofluorometric studies show that the presence of COT limited the amount of photobleaching of an oligo containing the cyanine 5 dye.

Item	Catalog No.	Pack	Price(\$)
COT Serinol Phosphoramidite	10-1996-95	50 μmole	310.00
	10-1996-90	100 μmole	600.00
	10-1996-02	0.25g	1800.00



INTELLECTUAL PROPERTY

This product is covered under US Patent 8,945,515 B2.

DABCYL LABELING


A molecular beacon probe¹ has its natural fluorescence quenched in solution unless it is hybridized to the target sequence. Consequently, the design of a molecular beacon requires a fluorophore to be in one part of the sequence and the quencher molecule to be in another, with both molecules being separated from the oligonucleotide by a hydrocarbon spacer. The Dabcyl group has been found to be a universal quencher. 3'-Dabsyl CPG and 3'-Dabcyl CPG are used to prepare probes with the quencher blocking the 3'-terminus. 5'-Dabcyl Phosphoramidite locates the quencher at the 5'-terminus and Dabcyl-dT places it within the sequence, leaving the 3'-terminus available for polymerase extension.

Item	Catalog No.	Pack	Price(\$)	
3'-Dabsyl CPG	20-5911-01	0.1g	120.00	
	20-5911-10	1.0g	975.00	
	1 μmole columns	20-5911-41	Pack of 4	200.00
	0.2 μmole columns	20-5911-42	Pack of 4	120.00
	10 μmole column (ABI)	20-5911-13	Pack of 1	350.00
15 μmole column (Expedite)	20-5911-14	Pack of 1	500.00	
3'-Dabcyl CPG	20-5912-01	0.1g	120.00	
	20-5912-10	1.0g	975.00	
	1 μmole columns	20-5912-41	Pack of 4	200.00
	0.2 μmole columns	20-5912-42	Pack of 4	120.00
	10 μmole column (ABI)	20-5912-13	Pack of 1	350.00
15 μmole column (Expedite)	20-5912-14	Pack of 1	500.00	
3'-Dabcyl PS	26-5912-01	0.1g	125.00	
	26-5912-10	1.0g	1025.00	
	200 nmole columns (AB 3900)	26-5912-52	Pack of 10	300.00
40 nmole columns (AB 3900)	26-5912-55	Pack of 10	300.00	
Dabcyl-dT	10-1058-95	50 μmole	180.00	
	10-1058-90	100 μmole	325.00	
	10-1058-02	0.25g	675.00	
5'-Dabcyl Phosphoramidite	10-5912-95	50 μmole	125.00	
	10-5912-90	100 μmole	225.00	
	10-5912-02	0.25g	650.00	


REFERENCE

(1) S. Tyagi and F.R. Kramer, *Nature Biotechnology*, 1996, **4**, 303-308.


DYE QUENCHER PLOT




http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf




Dabsyl CPG



Dabcyl CPG



Dabcyl-dT




5'-Dabcyl Phosphoramidite

BIOTIN LABELING


Glen Research biotin phosphoramidites for direct labeling of synthetic oligonucleotides exhibit the following features:

- All are soluble in acetonitrile at concentrations useful for DNA synthesis.
- All include a DMT group for cartridge purifications which is essential for the preparation of biotinylated PCR primers because of the potential for cross contamination in HPLC purifications.
- For the development of diagnostic probes, biotin phosphoramidite is capable of branching to allow multiple biotins to be introduced at the 3'- or 5'-terminus. BiotinTEG Phosphoramidite contains a 15 atom mixed polarity spacer arm based on a triethylene glycol.
- Protected Biotin Serinol Phosphoramidite and CPG are protected with a *t*-butylbenzoyl group on the biotin ring. This group is designed to stop any phosphoramidite reactions at this active position in biotin. This protection avoids branching when using nucleophilic activators like DCI. The protecting group is easily removed during oligonucleotide cleavage and deprotection. The BiotinLC versions are similarly protected and should be useful for the synthesis of highly sensitive biotinylated probes.


Item	Catalog No.	Pack	Price (\$)
Biotin Phosphoramidite	10-1953-95	50 μmole	165.00
	10-1953-90	100 μmole	295.00
	10-1953-02	0.25g	675.00
BiotinTEG Phosphoramidite	10-1955-95	50 μmole	165.00
	10-1955-90	100 μmole	295.00
	10-1955-02	0.25g	675.00
Protected Biotin Serinol Phosphoramidite	10-1993-95	50 μmole	165.00
	10-1993-90	100 μmole	295.00
	10-1993-02	0.25g	675.00
Protected BiotinLC Serinol Phosphoramidite	10-1995-95	50 μmole	205.00
	10-1995-90	100 μmole	365.00
	10-1995-02	0.25g	675.00




Biotin Phosphoramidite



BiotinTEG Phosphoramidite



Protected Biotin Serinol Phosphoramidite



Protected BiotinLC Serinol Phosphoramidite

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

BIOTIN LABELING (CONT.)





Biotin-dT can replace dT residues within the oligonucleotide sequence. 5'-Biotin phosphoramidite can be added ONLY ONCE to the 5'-terminus of an oligonucleotide. However, the DMT group on the biotin can be used in RP cartridge and HPLC purification techniques. PC Biotin is a photocleavable 5'-biotin phosphoramidite. BiotinTEG CPG and Protected BiotinLC Serinol CPG are designed for the direct synthesis of oligonucleotides containing biotin at the 3' terminus.

Desthiobiotin is a biotin analogue that exhibits lower binding to biotin-binding proteins such as streptavidin. This biotin analogue is lacking the sulfur group from the molecule and has a dissociation constant (Kd) several orders of magnitude less than biotin/streptavidin. As a result, biomolecules containing desthiobiotin are dissociated from streptavidin simply in the presence of buffered solutions of biotin. We offer desthiobiotinTEG phosphoramidite and the corresponding CPG.

ABI-style vials and columns are supplied unless otherwise requested (see note box).





SEE ALSO
PC Biotin on page 86

Item	Catalog No.	Pack	Price (\$)
5'-Biotin Phosphoramidite	10-5950-95	50 μmole	125.00
	10-5950-90	100 μmole	225.00
	10-5950-02	0.25g	650.00
Biotin-dT	10-1038-95	50 μmole	167.50
	10-1038-90	100 μmole	325.00
	10-1038-02	0.25g	625.00
PC Biotin Phosphoramidite	10-4950-95	50 μmole	145.00
	10-4950-90	100 μmole	280.00
	10-4950-02	0.25g	675.00
DesthiobiotinTEG Phosphoramidite	10-1952-95	50 μmole	185.00
	10-1952-90	100 μmole	335.00
	10-1952-02	0.25g	775.00



BIOTIN LABELING (CONT.)

Item	Catalog No.	Pack	Price (\$)
3'-BiotinTEG CPG	20-2955-01	0.1g	120.00
	20-2955-10	1.0g	995.00
	20-2955-42	Pack of 4	120.00
	20-2955-41	Pack of 4	200.00
	20-2955-13	Pack of 1	300.00
15 μmole column (Expedite)	20-2955-14	Pack of 1	450.00
3'-BiotinTEG PS	26-2955-01	0.1g	125.00
	26-2955-10	1.0g	1025.00
	26-2955-52	Pack of 10	300.00
40 nmole columns (AB 3900)	26-2955-55	Pack of 10	300.00
3'-Protected Biotin Serinol CPG	20-2993-01	0.1g	120.00
	20-2993-10	1.0g	995.00
	20-2993-42	Pack of 4	120.00
	20-2993-41	Pack of 4	200.00
	20-2993-13	Pack of 1	300.00
15 μmole column (Expedite)	20-2993-14	Pack of 1	450.00
3'-Protected BiotinLC Serinol CPG	20-2995-01	0.1g	120.00
	20-2995-10	1.0g	995.00
	20-2995-42	Pack of 4	120.00
1 μmole columns	20-2995-41	Pack of 4	200.00
	20-2995-13	Pack of 1	300.00
15 μmole column (Expedite)	20-2995-14	Pack of 1	450.00
DesthiobiotinTEG CPG	20-2952-01	0.1g	140.00
	20-2952-10	1.0g	1150.00
	20-2952-42	Pack of 4	140.00
	20-2952-41	Pack of 4	230.00
	20-2952-13	Pack of 1	345.00
15 μmole column (Expedite)	20-2952-14	Pack of 1	520.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)



FLUORESCIN LABELING

5'-Fluorescein phosphoramidite contains no 4,4'-dimethoxytrityl (DMT) group and can be added only once at the 5'-terminus, thereby terminating synthesis. This product is prepared using the 6-carboxyfluorescein derivative. The tetrachloro-, hexachloro- and dichloro-dimethoxy-fluorescein (TET, HEX and JOE, respectively) phosphoramidites are designed to take advantage of the multicolor detection capability of modern DNA sequencers and genetic analyzers. Fluorescein phosphoramidite is designed to produce the same fluorescein-type structure as had been previously prepared using fluorescein isothiocyanate (FITC). Our fluorescein phosphoramidite also contains a DMT group to allow quantification of coupling. The analogous structure, 6-Fluorescein Phosphoramidite, prepared using 6-FAM, is also available, along with 6-Fluorescein Serinol Phosphoramidite. Fluorescein-dT can be inserted into the desired sequence as a replacement for a dT residue.


We offer five fluorescein supports. Fluorescein CPG has traditionally been used to add the fluorescein label at the 3'-terminus. The analogous structure, 3'-(6-Fluorescein) CPG, prepared using 6-FAM, is now also available, along with 6-Fluorescein Serinol CPG. We also offer 3'-(6-FAM) CPG and Fluorescein-dT CPG, both derivatives of 6-carboxyfluorescein (6-FAM). Both are single isomers and use an amide linkage which is stable during cleavage and deprotection and does not allow isomer formation. 3'-(6-FAM) CPG allows effective blockage of the 3'-terminus from polymerase extension as well as exonuclease digestion. Fluorescein-dT CPG allows both of these enzymatic activities to proceed. Normal cleavage and deprotection with ammonium hydroxide readily generates the fluorescein labeled oligos.

The spectral characteristics of these dyes are detailed on the following page.


Item	Cat. No.	Pack	Price (\$)
5'-Fluorescein Phosphoramidite (6-FAM)	10-5901-95	50 µmole	110.00
	10-5901-90	100 µmole	215.00
	10-5901-02	0.25g	575.00
5'-Hexachloro-Fluorescein Phosphoramidite (HEX)	10-5902-95	50 µmole	190.00
	10-5902-90	100 µmole	375.00
	10-5902-02	0.25g	875.00
5'-Tetrachloro-Fluorescein Phosphoramidite (TET)	10-5903-95	50 µmole	180.00
	10-5903-90	100 µmole	350.00
	10-5903-02	0.25g	875.00
5'-Dichloro-dimethoxy-Fluorescein Phosphoramidite II (JOE)	10-5906-95	50 µmole	105.00
	10-5906-90	100 µmole	198.00
	10-5906-02	0.25g	495.00




5'-Fluorescein Phosphoramidite



5'-Hexachloro-Fluorescein Phosphoramidite



5'-Tetrachloro-Fluorescein Phosphoramidite



5'-Dichloro-dimethoxy-Fluorescein Phosphoramidite II

FLUORESCIN LABELING (CONT.)

Item	Cat. No.	Pack	Price (\$)
Fluorescein Phosphoramidite	10-1963-95	50 µmole	165.00
	10-1963-90	100 µmole	295.00
	10-1963-02	0.25g	595.00
6-Fluorescein Phosphoramidite	10-1964-95	50 µmole	165.00
	10-1964-90	100 µmole	295.00
	10-1964-02	0.25g	595.00
6-Fluorescein Serinol Phosphoramidite	10-1994-95	50 µmole	165.00
	10-1994-90	100 µmole	295.00
	10-1994-02	0.25g	595.00
Fluorescein-dT Phosphoramidite	10-1056-95	50 µmole	180.00
	10-1056-90	100 µmole	325.00
	10-1056-02	0.25g	675.00

FLUORESCENT DYES


	Absorbance Maximum	Emission Maximum	Color
Fluorescein	494nm	525nm	Green
Tetrachloro-Fluorescein	521nm	536nm	Orange
Hexachloro-Fluorescein	535nm	556nm	Pink
SIMA (HEX)	538nm	551nm	Pink
Dichloro-dimethoxy-Fluorescein	525nm	548nm	Orange/Pink
TAMRA	565nm	580nm	Rose
Cy3	546nm	563nm	Red
Cy3.5	588nm	604nm	Purple
Cy5	648nm	662nm	Violet
Cy5.5	683nm	707nm	Dark Blue
Yakima Yellow	530nm	549nm	Yellow
Redmond Red	579nm	595nm	Red

OTHER INSTRUMENT TYPES


All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M


(Please inquire for availability of vials and columns for other instrument types.)




Fluorescein Phosphoramidite



Fluorescein dT





6-Fluorescein Phosphoramidite



6-Fluorescein Serinol Phosphoramidite

FLUORESCIN LABELING (CONT.)

Item	Cat. No.	Pack	Price (\$)
3'-Fluorescein CPG	20-2963-01	0.1g	120.00
	20-2963-10	1.0g	995.00
	20-2963-41	Pack of 4	200.00
	20-2963-42	Pack of 4	120.00
	20-2963-13	Pack of 1	300.00
20-2963-14	Pack of 1	450.00	
3'-(6-Fluorescein) CPG	20-2964-01	0.1g	120.00
	20-2964-10	1.0g	995.00
	20-2964-41	Pack of 4	200.00
	20-2964-42	Pack of 4	120.00
	20-2964-13	Pack of 1	300.00
20-2964-14	Pack of 1	450.00	
3'-(6-FAM) CPG	20-2961-01	0.1g	120.00
	20-2961-10	1.0g	995.00
	20-2961-41	Pack of 4	200.00
	20-2961-42	Pack of 4	120.00
	20-2961-13	Pack of 1	300.00
20-2961-14	Pack of 1	450.00	
3'-(6-FAM) PS	26-2961-01	0.1g	130.00
	26-2961-10	1.0g	1045.00
	26-2961-52	Pack of 10	300.00
26-2961-55	Pack of 10	300.00	
3'-6-Fluorescein Serinol CPG	20-2994-01	0.1g	120.00
	20-2994-10	1.0g	995.00
	20-2994-42	Pack of 4	120.00
	20-2994-41	Pack of 4	200.00
	20-2994-13	Pack of 1	300.00
20-2994-14	Pack of 1	450.00	



http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf


FLUORESCIN LABELING (CONT.)

Item	Cat. No.	Pack	Price (\$)
3'-Fluorescein-dT CPG	20-2056-01	0.1g	120.00
	20-2056-10	1.0g	995.00
	20-2056-41	Pack of 4	200.00
	20-2056-42	Pack of 4	120.00
	20-2056-13	Pack of 1	300.00
20-2056-14	Pack of 1	450.00	

FLUORESCIN LABELING (SIMA)

Dichloro-diphenyl-fluorescein, SIMA (HEX) exhibits virtually identical absorbance and emission spectra to HEX. SIMA (HEX) is much more stable to basic deprotection conditions than HEX and oligonucleotides can be deprotected using ammonium hydroxide at elevated temperatures and even ammonium hydroxide/methylamine (AMA) at room temperature or 65°C for 10 minutes. SIMA absorption maximum was 3 nm blue-shifted compared to HEX at pH 7. The absorbance is broader, so the extinction coefficient is smaller than that of HEX, but when exciting at 500 nm where the absorbance was normalized, the emission was still 90% of HEX and the emission was red-shifted by 5 nm. A second SIMA (HEX) product, SIMA (HEX)-dT, can be used to introduce SIMA (HEX) in the synthetic oligonucleotide sequence, usually as a replacement for the native dT linkage. Again, this product is fully compatible with deprotection schemes using ammonium hydroxide at elevated temperatures or AMA at room temperature and 65°C.

Item	Cat. No.	Pack	Price (\$)
SIMA (HEX) Phosphoramidite	10-5905-95	50 μmole	90.00
	10-5905-90	100 μmole	175.00
	10-5905-02	0.25g	400.00
SIMA (HEX)-dT Phosphoramidite	10-5945-95	50 μmole	345.00
	10-5945-90	100 μmole	675.00
	10-5945-02	0.25g	995.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)



http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf

CYANINE LABELING

Two cyanine derivatives, Cyanine 3 and Cyanine 5, which differ in structure simply by the number of carbons in the conjugated poly-ene linkage, are joined by the closely related analogues, Cyanine 3.5 and Cyanine 5.5, and are available as phosphoramidites. Cyanine dyes are normally added once at the 5'-terminus and the MMT group should be removed on the synthesizer. The absorbance of the MMT cation (yellow) is noticeably different from the DMT cation (orange), and so, absorbance-based trityl monitors will detect it incorrectly as a low coupling. On the other hand, conductivity detectors will interpret the release more correctly. Cyanine dye phosphoramidites have also been used successfully adjacent to the 3'-terminus. Cyanine 3 and Cyanine 5 supports are also offered to allow simpler production of 3' cyanine dye-labeled oligonucleotides.

Deprotection of oligos containing Cyanine dyes may be carried out with ammonium hydroxide at room temperature, regardless of the base protecting groups on the monomers used. If there is a need to use ammonium hydroxide at elevated temperature, Cyanine 3 and Cyanine 3.5 are more stable than Cyanine 5 and Cyanine 5.5. However, it is always prudent to use monomers with base labile protecting groups to limit the exposure time to 2 hours or less at 65°C during deprotection.


To better address applications in near-infrared (NIR) imaging, Glen Research is offering a water soluble Disulfo-Cyanine 7 Azide. This long wavelength dye offers the benefits of improved solubility, reduced aggregation, and improved stability in the near-infrared spectrum along with the convenience of click chemistry.

SPECTRAL DATA FOR CYANINE DYES

	Absorbance Maximum	Emission Maximum	Color
Cyanine 3	546nm	563nm	Red
Cyanine 3.5	588nm	604nm	Purple
Cyanine 5	646nm	662nm	Violet
Cyanine 5.5	683nm	707nm	Dark Blue
Cyanine 7	750nm	773nm	Dark Green


(Measured in an oligo in 0.1M TEAA buffer, pH7.)

DYE QUENCHER PLOT




http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf


Item	Cat. No.	Pack	Price (\$)
Cyanine 3 Phosphoramidite	10-5913-95	50 µmole	205.00
	10-5913-90	100 µmole	375.00
	10-5913-02	0.25g	925.00
Cyanine 3.5 Phosphoramidite	10-5914-95	50 µmole	220.00
	10-5914-90	100 µmole	400.00
	10-5914-02	0.25g	925.00
Cyanine 5 Phosphoramidite	10-5915-95	50 µmole	205.00
	10-5915-90	100 µmole	375.00
	10-5915-02	0.25g	925.00
Cyanine 5.5 Phosphoramidite	10-5916-95	50 µmole	245.00
	10-5916-90	100 µmole	450.00
	10-5916-02	0.25g	925.00




Cyanine 3 Phosphoramidite



Cyanine 5 Phosphoramidite



Cyanine 3.5 Phosphoramidite



Cyanine 5.5 Phosphoramidite

CYANINE LABELING (CONT.)

Item	Cat. No.	Pack	Price (\$)
Cyanine 3 CPG	20-5913-01	0.1g	160.00
	20-5913-10	1.0g	1250.00
	1 µmole columns (TWIST format only)	20-5913-41	Pack of 4
0.2 µmole columns	20-5913-42	Pack of 4	70.00
Cyanine 5 CPG	20-5915-01	0.1g	160.00
	20-5915-10	1.0g	1250.00
	1 µmole columns (TWIST format only)	20-5915-41	Pack of 4
0.2 µmole columns	20-5915-42	Pack of 4	70.00
Disulfo-Cyanine 7 Azide	50-2010-92	25 µmole	325.00
	50-2010-90	100 µmole	975.00


OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)


DYE QUENCHER PLOT




http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf



Cyanine 3 CPG



Cyanine 5 CPG







Disulfo-Cyanine 7 Azide

ELITECHGROUP DYES AND QUENCHER

Glen Research's agreement with ELITechGroup, formerly Epoch Biosciences, allows us to offer several of their proprietary products designed for the synthesis of novel DNA probes. We are pleased to offer products based on ELITechGroup's Redmond Red®, Yakima Yellow® and AquaPhluor® 593 fluorophores and Eclipse® non-fluorescent quencher. Under our agreement we also supply PPG, a modified nucleoside, and 5'-Aldehyde-Modifier C2 Phosphoramidite. The fluorescent dyes, Yakima Yellow, Redmond Red and AquaPhluor 593, are available as phosphoramidites and supports. Yakima Yellow has an absorbance maximum at 530 nm and emission maximum at 549 nm, Redmond Red's absorbance and emission maxima are at 579 nm and 595 nm, respectively, and AquaPhluor 593 has an absorbance maximum at 593 nm and emission maximum at 613 nm.

The Eclipse quencher from ELITechGroup solves most of the problems inherent in the synthesis of molecular beacon and FRET probes. The Eclipse molecule is highly stable and can be used safely in all common oligo deprotection schemes. The absorbance maximum for Eclipse Quencher is at 522 nm, compared to 479 nm for dabcyl. In addition, the structure of the Eclipse Quencher is substantially more electron deficient than that of dabcyl and this leads to better quenching over a wider range of dyes, especially those with emission maxima at longer wavelengths (red shifted) such as Redmond Red and Cyanine 5. In addition, with an absorption range from 390 nm to 625 nm, the Eclipse Quencher is capable of effective performance in a wide range of colored FRET probes.

Item	Cat. No.	Pack	Price (\$)
Redmond Red® Phosphoramidite	10-5920-95	50 µmole	220.00
	10-5920-90	100 µmole	420.00
	10-5920-02	0.25g	1045.00
Yakima Yellow® Phosphoramidite	10-5921-95	50 µmole	230.00
	10-5921-90	100 µmole	440.00
	10-5921-02	0.25g	1045.00
5'-AquaPhluor® 593 Phosphoramidite	10-5923-95	50 µmole	405.00
	10-5923-90	100 µmole	795.00
	10-5923-02	0.25g	1575.00
Eclipse® Quencher Phosphoramidite	10-5925-95	50 µmole	250.00
	10-5925-90	100 µmole	480.00
	10-5925-02	0.25g	1185.00



SEE ALSO

PPG on page 57
5'-Aldehyde-Modifier C2 on page 83

FLUORESCENT DYES

	Absorbance Maximum	Emission Maximum	Color
Yakima Yellow	530nm	549nm	Yellow
Redmond Red	579nm	595nm	Red
AquaPhluor 593	593nm	613nm	Red

INTELLECTUAL PROPERTY





These Products are for research purposes only, and may not be used for commercial, clinical, diagnostic or any other use. These Products are subject to proprietary rights of ELITechGroup and are made and sold under license from ELITechGroup. There is no implied license for commercial use with respect to these Products and a license must be obtained directly from ELITechGroup with respect to any proposed commercial use of these Products. "Commercial use" includes but is not limited to the sale, lease, license or other transfer of the Product or any material derived or produced from it, the sale, lease, license or other grant of rights to use the Product or any material derived or produced from it, or the use of the Product to perform services for a fee for third parties (including contract research).

A simple agreement must be signed before end-users and custom oligo services may purchase these products for use as defined above.
<http://www.glenresearch.com/Reference/ELITechGroupProducts.pdf>

AquaPhluor®, Yakima Yellow®, Redmond Red® and Eclipse®, are registered Trademarks of ELITechGroup.

ELITECHGROUP DYES AND QUENCHER (CONT.)

Item	Cat. No.	Pack	Price (\$)	
Redmond Red® CPG	20-5920-01	0.1g	180.00	
	20-5920-10	1.0g	1500.00	
	1 µmole columns	20-5920-41	Pack of 4	300.00
	0.2 µmole columns	20-5920-42	Pack of 4	150.00
	10 µmole column (ABI)	20-5920-13	Pack of 1	750.00
15 µmole column (Expedite)	20-5920-14	Pack of 1	1125.00	
Yakima Yellow® CPG	20-5921-01	0.1g	180.00	
	20-5921-10	1.0g	1500.00	
	1 µmole columns	20-5921-41	Pack of 4	300.00
	0.2 µmole columns	20-5921-42	Pack of 4	150.00
	10 µmole column (ABI)	20-5921-13	Pack of 1	750.00
15 µmole column (Expedite)	20-5921-14	Pack of 1	1125.00	
AquaPhluor® 593 CPG	20-5923-01	0.1g	215.00	
	20-5923-10	1.0g	1800.00	
	1 µmole columns	20-5923-41	Pack of 4	325.00
	0.2 µmole columns	20-5923-42	Pack of 4	165.00
	10 µmole column (ABI)	20-5923-13	Pack of 1	925.00
15 µmole column (Expedite)	20-5923-14	Pack of 1	1395.00	
Eclipse® Quencher CPG	20-5925-01	0.1g	230.00	
	20-5925-10	1.0g	1925.00	
	1 µmole columns	20-5925-41	Pack of 4	350.00
	0.2 µmole columns	20-5925-42	Pack of 4	175.00
	10 µmole column (ABI)	20-5925-13	Pack of 1	995.00
15 µmole column (Expedite)	20-5925-14	Pack of 1	1495.00	




OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

DYE QUENCHER PLOT



http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf

BLACK HOLE QUENCHER DYES

With the growing popularity of red and near-infrared dyes, we are offering the Black Hole Quencher™ dyes (BHQs), whose physical properties are detailed in Table 1. BHQ dyes are robust dark quenchers that very nicely complement our existing product line. They are compatible with ammonium hydroxide deprotection, exhibit excellent coupling efficiencies, have large extinction coefficients and are completely non-fluorescent. Their absorbances are well-tuned to quench a variety of popular fluorophores – even those far into the red, such as Cy3 and Cy5. The dark quencher most typically used in a Molecular Beacon is Dabcyl. Because the quenching does not involve FRET, there is little, if any, dependence upon donor-acceptor spectral overlap. In a comprehensive paper by Marras, Kramer and Tyagi,¹ the ability of BHQ-1 and BHQ-2 to quench 22 different fluorophores was evaluated. For shorter wavelength fluorophores such as fluorescein, the quenching efficiency was roughly the same as Dabcyl (91% – 93%). However, for dyes emitting in the far red, such as Cy5, the BHQ dyes were far superior – quenching the Cy5 with 96% efficiency, compared to 84% with Dabcyl. This may reflect the BHQ's ability to form stable, non-fluorescent complexes which can be a plus even in FRET probes. Indeed, recent work suggests that these non-fluorescent complexes will form even in the absence of a hairpin stem structure used by Molecular Beacons.²

TABLE 1: BLACK HOLE QUENCHERS

Quencher	λ _{max} (nm)	E260 (L/mol.cm)	E _{max} (L/mol.cm)
BHQ-1	534	8,000	34,000
BHQ-2	579	8,000	38,000
BHQ-3	672	13,000	42,700

REFERENCES

- (1) S.A.E. Marras, F.R. Kramer, and S. Tyagi, *Nucleic Acids Res.*, 2002, **30**, E122.
- (2) M.K. Johansson, H. Fidder, D. Dick, and R.M. Cook, *J Am Chem Soc*, 2002, **124**, 6950-6956.


SEE OTHER QUENCHERS

Dabcyl on page 98
Eclipse™ on page 109
BBQ-650® on page 112

INTELLECTUAL PROPERTY

"Black Hole Quencher", "BHQ-0", "BHQ-1", "BHQ-2" and "BHQ-3" are trademarks of Biosearch Technologies, Inc., Novato, CA. The BHQ dye technology is the subject of pending patents and is licensed and sold under agreement with Biosearch Technologies, Inc.. Products incorporating the BHQ dye moiety are sold exclusively for R&D use by the end-user. They may not be used for clinical or diagnostic purposes and they may not be re-sold, distributed or re-packaged.

Item	Cat. No.	Pack	Price (\$)
5'-BHQ-1 Phosphoramidite	10-5931-95	50 μmole	100.00
	10-5931-90	100 μmole	200.00
	10-5931-02	0.25g	700.00
5'-BHQ-2 Phosphoramidite	10-5932-95	50 μmole	100.00
	10-5932-90	100 μmole	200.00
	10-5932-02	0.25g	700.00
BHQ-1-dT	10-5941-95	50 μmole	265.00
	10-5941-90	100 μmole	525.00
	10-5941-02	0.25g	925.00
BHQ-2-dT	10-5942-95	50 μmole	265.00
	10-5942-90	100 μmole	525.00
	10-5942-02	0.25g	925.00



BLACK HOLE QUENCHER DYES (CONT.)

Item	Cat. No.	Pack	Price (\$)
3'-BHQ-1 CPG	20-5931-01	0.1g	190.00
	20-5931-10	1.0g	1500.00
	20-5931-41	Pack of 4	300.00
	20-5931-42	Pack of 4	80.00
	20-5931-13	Pack of 1	575.00
3'-BHQ-2 CPG	20-5932-01	0.1g	190.00
	20-5932-10	1.0g	1500.00
	20-5932-41	Pack of 4	300.00
	20-5932-14	Pack of 1	825.00
3'-BHQ-3 CPG	20-5933-01	0.1g	190.00
	20-5933-10	1.0g	1500.00
	20-5933-41	Pack of 4	300.00
	20-5933-42	Pack of 4	80.00
	20-5933-14	Pack of 1	825.00


OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)

Dye Quencher Plot




http://www.glenresearch.com/ProductFiles/Dye_Quencher_Plot.pdf

BLACKBERRY® QUENCHER (BBQ-650®)


We are happy to offer several products containing the BlackBerry® Quencher (BBQ-650®), which exhibits a broad absorption profile from 550nm to 750nm, centered at 650nm. This range offers more effective quenching of some of our popular long wavelength dyes like TAMRA, Redmond Red, Cy dyes and DyLight dyes. We offer BBQ-650 products for the 3' and 5' termini, as well as BBQ-650-dT for inclusion within the oligonucleotide sequence, with the following properties:

- Quenches the fluorescence of long wavelength dyes
- Quenches in FRET and contact mode
- Absorbance maximum at ~650nm
- Quenching range – 550-750nm
- Compatible with standard oligo synthesis chemistry
- Compatible with regular deprotection but requires mild deprotection with AMA at room temperature
- Available for 3', 5', and internal substitution
- More stable than BHQ-3


Item	Cat. No.	Pack	Price (\$)
5'-BBQ-650® Phosphoramidite	10-5934-95	50 µmole	160.00
	10-5934-90	100 µmole	305.00
	10-5934-02	0.25g	925.00
BBQ-650®-dT	10-5944-95	50 µmole	280.00
	10-5944-90	100 µmole	545.00
	10-5944-02	0.25g	925.00
3'-BBQ-650® CPG	20-5934-01	0.1g	190.00
	20-5934-10	1.0g	1500.00
	20-5934-41	Pack of 4	300.00
	20-5934-42	Pack of 4	80.00
	20-5934-13	Pack of 1	575.00
	20-5934-14	Pack of 1	825.00



5'-BBQ-650™



3'-BBQ-650™ CPG



BBQ-650™-dT


INTELLECTUAL PROPERTY

BlackBerry® Quencher technology: US Patent 7,879,986. The purchase of BlackBerry® Quencher reagents includes a limited license to use these reagents exclusively for research and development purposes. They may not be used for clinical or diagnostic purposes and they may not be re-sold, distributed, or re-packaged without prior agreement and consent of Berry & Associates, Inc. Subsequent sale of products that are derived from BlackBerry® Quencher reagents is permitted so long as the following written disclaimer is included in written and electronic catalogs, in commercial advertisement, and in packages with containers of such derivative products: "BlackBerry is a trademark of Berry & Associates, Inc. Products derived from BlackBerry® Quencher reagents are sold exclusively for research and development use by the purchaser. They may not be used for clinical or diagnostic purposes without prior agreement and consent of Berry & Associates, Inc."


RHODAMINE (TAMRA) LABELING

Rhodamine derivatives are not sufficiently stable to survive conventional deprotection and these must be attached to amino-modified oligonucleotides using post-synthesis labeling techniques. Because Tetramethyl Rhodamine (TAMRA) is not base stable, the procedure to cleave and deprotect the labeled oligonucleotide must be carefully considered. Using the UltraMILD monomers and deprotection with potassium carbonate in methanol, TAMRA oligonucleotides can be fairly conveniently isolated. To streamline the preparation of TAMRA oligos, we offer 3'-TAMRA CPG for 3' labeling and TAMRA-dT for labeling within the sequence. We also offer TAMRA NHS ester for labeling amino-modified oligonucleotides.


Item	Cat. No.	Pack	Price (\$)
3'-TAMRA CPG	20-5910-01	0.1g	120.00
	20-5910-10	1.0g	995.00
	20-5910-41	Pack of 4	200.00
0.2 µmole columns	20-5910-42	Pack of 4	120.00
3'-TAMRA PS	26-5910-01	0.1g	130.00
	26-5910-10	1.0g	1045.00
	26-5910-52	Pack of 10	300.00
	26-5910-55	Pack of 10	300.00
TAMRA-dT	10-1057-95	50 µmole	250.00
	10-1057-90	100 µmole	495.00
	10-1057-02	0.25g	975.00
TAMRA NHS Ester (Solution in anhydrous DMSO)	50-5910-66	60 µL	240.00



TAMRA CPG



TAMRA-dT



TAMRA NHS Ester

SEE ALSO

UltraMILD monomers on page 23

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

ACRIDINE LABELING


Acridine phosphoramidite is designed to produce an oligonucleotide containing acridine at any position in the molecule. Acridine CPG is used to label the 3'-terminus. Acridine is an effective intercalating agent.

Item	Cat. No.	Pack	Price (\$)
Acridine Phosphoramidite	10-1973-95	50 μmole	165.00
	10-1973-90	100 μmole	295.00
	10-1973-02	0.25g	675.00
3'-Acridine CPG	20-2973-01	0.1g	120.00
	20-2973-10	1.0g	995.00
	20-2973-41	Pack of 4	200.00
	20-2973-42	Pack of 4	120.00
	20-2973-13	Pack of 1	300.00
	20-2973-14	Pack of 1	450.00

DNP LABELING

An analytical test based on detection of 2,4-dinitrophenyl (DNP) labeled oligonucleotides with anti-DNP antibodies has been proposed. We have chosen the branched triethylene glycol (TEG) spacer in our version of DNP phosphoramidite since it can be added once or several times to the 3' or 5' terminus.

Item	Catalog No.	Pack	Price(\$)
DNP-TEG Phosphoramidite	10-1985-95	50 μmole	165.00
	10-1985-90	100 μmole	295.00
	10-1985-02	0.25g	675.00



Acridine

Acridine CPG

DNP-TEG

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

CHOLESTEROL LABELING

Potential therapeutic oligonucleotides must permeate the cell membrane for optimal activity. The addition of lipophilic groups to an oligonucleotide would be expected to enhance cellular uptake/membrane permeation. The use of cholesteryl oligos and the consequent improvement in activity has been described. We have designed our Cholesteryl products with triethyleneglycol (TEG) spacers for maximum solubility.

Item	Catalog No.	Pack	Price(\$)
Cholesteryl-TEG Phosphoramidite	10-1975-95	50 μmole	140.00
	10-1975-90	100 μmole	265.00
	10-1975-02	0.25g	545.00
5'-Cholesteryl-TEG Phosphoramidite	10-1976-95	50 μmole	95.00
	10-1976-90	100 μmole	175.00
	10-1976-02	0.25g	525.00
3'-Cholesteryl-TEG CPG	20-2975-01	0.1g	85.00
	20-2975-10	1.0g	700.00
	20-2975-41	Pack of 4	140.00
	20-2975-42	Pack of 4	84.00
	20-2975-13	Pack of 1	210.00
	20-2975-14	Pack of 1	315.00

SEE ALSO

Spermine on page 48

TOCOPHEROL LABELING


Vitamin E is both lipophilic and non-toxic even at high doses so would be an excellent candidate as a lipophilic carrier for oligonucleotides. Therefore, as an addition to our cholesteryl product line, we offer simple α-tocopheryl (vitamin E) labeling. Totally synthetic α-tocopherol is racemic at its three chiral centers and is used to prepare this product.

Item	Catalog No.	Pack	Price(\$)
α-Tocopherol-TEG Phosphoramidite	10-1977-95	50 μmole	160.00
	10-1977-90	100 μmole	300.00
	10-1977-02	0.25g	575.00

STEARYL LABELING

We now offer a simple C18 lipid as an economical and effective carrier molecule. We envisage that the 5'-stearyl group will become a favored lipophilic carrier for experimentation with synthetic oligonucleotides.

Item	Catalog No.	Pack	Price(\$)
5'- Stearyl Phosphoramidite	10-1979-90	100 μmole	45.00
	10-1979-02	0.25g	180.00



Cholesteryl-TEG

5'-Cholesteryl-TEG

3'-Cholesteryl-TEG CPG

α-Tocopherol-TEG


5'- Stearyl

N-ACETYL GALACTOSAMINE (GalNAc) LABELING

A directed approach to the delivery of therapeutic oligonucleotides specifically to the liver has been to target the asialoglycoprotein receptor (ASGPR) using a suitable glycoconjugate. Indeed, ASGPR is the ideal target for delivery of therapeutic oligonucleotides to the liver since it combines tissue specificity, high expression levels and rapid internalization and turnover. The use of oligonucleotide glycoconjugates has led to significant advances in therapeutic delivery as evidenced by the work of Alnylam Pharmaceuticals and Ionis Pharmaceuticals using multivalent N-acetylgalactosamine (GalNAc) oligonucleotide conjugates.

Glen Research is delighted to introduce a GalNAc modification strategy using a monomeric GalNAc support and the equivalent GalNAc phosphoramidite. Our experimental work has shown that these products are fully compatible with regular oligonucleotide synthesis and deprotection. Oligonucleotides containing GalNAc can be deprotected using standard procedures during which the acetyl protecting groups on GalNAc are removed. We have demonstrated that 5'-GalNAc C3 phosphoramidite can be used to prepare oligonucleotides with multiple consecutive GalNAc additions at the 5' terminus. Glen Research offers these GalNAc C3 products under an agreement with AM Chemicals LLC.

Item	Catalog No.	Pack	Price(\$)	
5'-GalNAc C3 Phosphoramidite	10-1974-95	50 μmole	137.50	
	10-1974-90	100 μmole	255.00	
	10-1974-02	0.25g	500.00	
GalNAc C3 CPG	20-2974-01	0.1g	40.00	
	20-2974-10	1.0g	320.00	
	1 μmole columns	20-2974-41	Pack of 4	100.00
	0.2 μmole columns	20-2974-42	Pack of 4	60.00
	10 μmole column (ABI)	20-2974-13	Pack of 1	180.00
	15 μmole column (Expedite)	20-2974-14	Pack of 1	280.00



OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

CDPI₃ MGB™ LABELING

The tripeptide of dihydropyrroloindole-carboxylate (CDPI₃) is a minor groove binding (MGB) moiety derived from the natural product CC-1065 with strong DNA binding properties. Synthetic oligonucleotides with covalently-attached CDPI₃ have enhanced DNA affinity and have improved the hybridization properties of sequence-specific DNA probes. Short CDPI₃-oligonucleotides hybridize with single-stranded DNA to give more stable DNA duplexes than unmodified ODNs of similar length. CDPI₃ MGB-oligonucleotide conjugates have been found to be useful in the following applications:



- Arrest of primer extension and PCR blockers
- Short and fluorogenic PCR primers
- Real-time PCR probes
- miRNA Inhibitors

The simplest approach to MGB probe design is to use an MGB support, add a quencher molecule as the first addition and complete the synthesis with a 5'-fluorophore. Alternatively, a fluorophore support could be used with the 5' terminus containing a quencher molecule followed by a final MGB addition at the 5' terminus. Glen Research offers 5'-CDPI₃ MGB™ Phosphoramidite and 3'-CDPI₃ MGB™ CPG.

5'-CDPI₃ MGB phosphoramidite was found to be hydrophobic enough that it required 10% THF in ACN to go completely into solution at a 0.1 M concentration and required a 3 minute coupling time. Deprotection can be carried out in EtOH/NH₄OH 1:3 (v/v) 17 hr at 55 °C and CDPI₃ MGB is compatible with GlenPak™ purification.

With the CDPI₃ MGB CPG, optimal results are obtained if UltraMild monomers and Cap A are used during synthesis along with 0.5 M CSO oxidizer. However, the use of standard monomers with iodine oxidation followed by deprotection with EtOH/NH₄OH 1:3 (v/v) for 17 hr at 55 °C will give acceptable results.

Item	Catalog No.	Pack	Price(\$)	
5'-CDPI ₃ MGB™ Phosphoramidite	10-5924-95	50 μmole	705.00	
	10-5924-90	100 μmole	1390.00	
	10-5924-02	0.25g	2600.00	
CDPI ₃ MGB™ CPG	20-5924-01	0.1g	215.00	
	20-5924-10	1.0g	1800.00	
	1 μmole columns	20-5924-41	Pack of 4	325.00
	0.2 μmole columns	20-5924-42	Pack of 4	165.00
	10 μmole column (ABI)	20-5924-13	Pack of 1	925.00
	15 μmole column (Expedite)	20-5924-14	Pack of 1	1395.00



LEGAL NOTICE

NOTICE TO PURCHASER: LIMITED LICENSE

This product is sold under licensing arrangements between ELITechGroup Inc. and Glen Research. The purchase price of this product includes limited, nontransferable rights to use the product solely for activities of the purchaser which are directly related to human diagnostics. Other uses, including incorporation of the product into another commercial product, are prohibited without additional license rights. For information on purchasing a license to this product for purposes other than those stated above, contact:

ELITech Group Molecular Diagnostics,
21720 23rd Drive SE, Suite 150,
Bothell, WA 98021.
Phone (425) 482-5555.
Fax (425) 482-5550.
Email: mdx@elitechgroup.com.

This limited license permits the person or legal entity to which this product has been provided to use the product, and the data generated by use of the product, only for human diagnostics. Neither ELITechGroup Inc. nor its licensors grants any other licenses, expressed or implied for any other purposes.

Some components of nucleic acid analysis, such as specific methods and compositions for manipulating or visualizing nucleic acids for analysis, may be covered by one or more patents owned by other parties. Similarly, nucleic acids containing specific nucleotide sequences may be patented. Making, using, offering for sale, or selling such components or nucleic acids may require one or more licenses. Nothing in this document should be construed as an authorization or implicit license to make, use or sell any so covered component or nucleic acid under any such patents.

PSORALEN LABELING


Psoralen C2 at the 5'-terminus of an oligonucleotide serves effectively as a cross-linking reagent in double-stranded oligonucleotides. The 6 atom spacer arm of Psoralen C6 allows cross-linking with a triplex oligonucleotide strand. Click Chemistry with psoralen azide and one of our many nucleosidic and non-nucleosidic alkyne derivatives has the potential to generate a variety of practical cross-linkers. The well known reversible cross-linking behavior of psoralen with an adjacent thymidine residue could be very useful.

Item	Cat. No.	Pack	Price (\$)
Psoralen C2 Phosphoramidite	10-1982-90	100 μmole	195.00
	10-1982-02	0.25g	495.00
Psoralen C6 Phosphoramidite	10-1983-90	100 μmole	195.00
	10-1983-02	0.25g	495.00
Psoralen Azide	50-2009-92	25 μmole	115.00
	50-2009-90	100 μmole	350.00


EDTA LABELING

EDTA-C2-dT phosphoramidite contains the triethyl ester of EDTA which allows sequence-specific cleavage of single- and double-stranded DNA and RNA. The cleavage reaction is only initiated once Fe(II) and dithiothreitol are added and so is readily controlled. Coupling of EDTA-dT is normal but cleavage and deprotection should be carried out with sodium hydroxide in aqueous methanol (0.4M NaOH in methanol/water 4:1) overnight at room temperature.


Item	Cat. No.	Pack	Price (\$)
EDTA-C2-dT-CE Phosphoramidite	10-1059-95	50 μmole	250.00
	10-1059-90	100 μmole	495.00
	10-1059-02	0.25g	975.00




Psoralen C2



Psoralen C6



Psoralen Azide



EDTA-C2-dT

FERROCENE LABELING

With an excellent stability profile, ferrocene has always attracted considerable interest for DNA labeling to generate probes for electrochemical detection. Based on our Amino-Modifier C6-dT structure, Ferrocene-dT is easily added to oligonucleotides with no disruption of regular hybridization behavior. Multiple incorporations into an oligonucleotide probe are also simply achieved. Oligonucleotides are deprotected using standard techniques. Ferrocene oligonucleotides should be stored under Argon and aqueous solutions should be degassed immediately.


Item	Cat. No.	Pack	Price (\$)
Ferrocene-dT-CE Phosphoramidite	10-1576-95	50 μmole	170.00
	10-1576-90	100 μmole	330.00
	10-1576-02	0.25g	670.00

METHYLENE BLUE LABELING


Methylene Blue, which belongs to the phenothiazine family of dyes, is a unique dye with a variety of useful properties. Despite its high extinction coefficient in the visible region (81,000 L/mol.cm), it is weakly fluorescent due to its high rate of intersystem crossing from the S₁ excited state to the T₁ triplet state. This property makes it an excellent photosensitizer, and it has been used extensively to produce highly reactive singlet oxygen. Methylene blue has the ability to both intercalate in duplex DNA, preferring G:C over T:A base pairs, and can act as an electrochemical redox probe. Methylene blue has also been shown to be unmatched in performance as a redox-active reporter for electrochemical biosensors.

Earlier, we introduced Methylene Blue C3 Phosphoramidite but this product proved to have quite limited stability and has been discontinued. As an alternative option, we introduced Methylene Blue NHS Ester to allow researchers to label amino-modified oligonucleotides with this interesting dye. With the encouragement and technical expertise of Carole Chaix and her colleagues at the University of Lyon, we decided to prepare an alternative structure that seemed to have a much superior stability profile - Methylene Blue II Phosphoramidite. Fortunately, this structure did indeed prove more stable and we are now able to offer again a Methylene Blue Phosphoramidite.


Item	Cat. No.	Pack	Price (\$)
Methylene Blue NHS Ester (Dissolve 5.4mg in 60μL of DMSO)	50-1960-23	5.4mg	540.00
Methylene Blue II Phosphoramidite	10-5961-95	50 μmole	310.00
	10-5961-90	100 μmole	595.00
	10-5961-02	0.25g	1500.00



Ferrocene-dT



Methylene Blue II



Methylene Blue NHS Ester

INTELLECTUAL PROPERTY

Methylene Blue II is covered under patent applications FR12 51739 and PCT/FR2013/050356 and is sold under license from the University of Lyon.

LABELING WITH METAL CHELATES

2,2'-Dipicolylamine Phosphoramidite has been discontinued. This product was manufactured and developed by Syntrix Biosystems Inc. For further information, please contact:

Dean Y. Maeda, Ph.D., M.B.A.
 Director, Chemistry and Preclinical Development
 Syntrix Biosystems
 215 Clay St NW Ste B5
 Auburn, WA 98001
 tel: 253-833-8009 ext. 23
 fax: 253-833-8127
 Dmaeda@syntrixbio.com

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
 For Instrument type Add

Expedite E
 MerMade M


Columns
 For Instrument type Add

Expedite E
 Applied Biosystems 3900 A
 MerMade M


(Please inquire for availability of vials and columns for other instrument types.)

FLUORESCENT DYES


	Absorbance Maximum	Emission Maximum	Excimer
Pyrene-dU	402nm	472nm	486nm
Perylene-dU	473nm	490nm	Not Determined



Pyrene-dU



2,2'-Dipicolylamine



Perylene-dU

LABELING WITH POLYAROMATIC HYDROCARBONS

Pyrene and perylene are fluorescent polycyclic aromatic hydrocarbons that have the ability to form 'excited state dimers' known as excimers. This unstructured, long-wavelength emission arises from the formation of a charge-transfer complex between the excited state and the ground state of two fluorescent molecules. In Pyrene-dU and perylene-dU, the hydrocarbon is attached at the 5 position of deoxyuridine through a triple bond and is electronically coupled to the deoxyuridine base. This electronic coupling of the base and the hydrocarbon makes the fluorescence sensitive to the base pairing of the dU portion of the molecule, allowing the discrimination between perfect and one base mismatched targets.

Item	Cat. No.	Pack	Price (\$)
Pyrene-dU-CE Phosphoramidite	10-1590-95	50 μmole	105.00
	10-1590-90	100 μmole	210.00
	10-1590-02	0.25g	550.00
Perylene-dU-CE Phosphoramidite	10-1591-95	50 μmole	150.00
	10-1591-90	100 μmole	300.00
	10-1591-02	0.25g	720.00

PUROMYCIN CPG


One of the most challenging requirements associated with combinatorial chemistry is the recovery of sequence information of the oligonucleotide or peptide selected by the screening assay. A method¹ has been developed to generate a fusion product between mRNA and the polypeptide it encodes using *in vitro* translation of synthetic RNAs 3'-labeled with puromycin, an antibiotic that mimics transfer RNA. Puromycin binds in the ribosome's A site, forms a peptide bond with the growing peptide chain, and blocks further peptide elongation. By linking puromycin to mRNA, a peptide-RNA fusion product results from the translation of the message linking the encoding mRNA with its peptide product.

Item	Catalog No.	Pack	Price(\$)
Puromycin CPG	20-4040-01	0.1g	120.00
	20-4040-10	1.0g	995.00
	1 μmole columns	Pack of 4	200.00
	0.2 μmole columns	Pack of 4	120.00
	10 μmole column (ABI)	Pack of 1	360.00
15 μmole columns (Expedite)	20-4140-14	Pack of 1	540.00


QUENCHED AUTOLIGATION (QUAL) PROBES

QUAL probes¹ consist of two oligonucleotides, the first containing a nucleophilic group at the 3'-terminus, while the second has an electrophilic group at the 5'-terminus. When the probe pair finds the target, the oligos line up with the 3'-terminus of the first directly adjacent to the 5'-terminus of the second. An autoligation reaction then takes place to combine the two oligos into a single probe. As usual, the 3' nucleophilic group is the 3-thiophosphate, easily prepared using 3'-phosphate CPG with a sulfurizing step in the first cycle. In this case, the electrophilic group is a 5'-dabsyl group, which is an excellent leaving group as well as a fine quencher of fluorescence. The second oligo, therefore, contains a fluorophore which is quenched by the dabsyl group. A popular choice for fluorophore is fluorescein-dT but it is easy to imagine that a variety of fluorophores could be attached to any of the commercially available amino-modified nucleoside phosphoramidites.

Item	Catalog No.	Pack	Price(\$)
5'-Dabsyl-dT-CE Phosphoramidite	10-1532-90	100 μmole	250.00
	10-1532-02	0.25g	775.00



Puromycin CPG



5'-Dabsyl-dT

REFERENCE

(1) R.W. Roberts and J.W. Szostak, *Proc. Natl. Acad. Sci. USA*, 1997, **94**, 12297-302.

REFERENCE

(1) S. Sando and E.T. Kool, *J Amer Chem Soc*, 2002, **124**, 2096-2097.

SEE ALSO

3'-Phosphate CPG on page 82
 Sulfurizing Reagent on page 41
 Fluorescein-dT on page 103

LABELING FOR PHOTO-REGULATION OF OLIGONUCLEOTIDES


Photo-control, the use of ultraviolet or visible light to control a reaction, has a number of advantages over other external stimuli:

- Light does not introduce contaminants into the reaction system,
- Excitation wavelength can be controlled through the design of the photo-responsive molecule, and
- It is now straightforward to control irradiation time and/or local excitation.

When a photo-responsive molecule is directly attached to DNA as a receptor, photo-regulation of the bioprocess regulated by that DNA molecule could, in principle, be achieved. Such photo-responsive DNA could also be used as a switch in a DNA-based nano-machine. Professor Hiroyuki Asanuma and his group at the department of Molecular Design and Engineering of the Graduate School of Engineering of the Nagoya University (Japan) have developed an efficient method to achieve this goal. They have attached azobenzene to DNA and made it photo-responsive^{1,2}. Azobenzene is a typical photo-responsive molecule that isomerizes from its planar *trans*-form to the non-planar *cis*-form after UV-light irradiation with a wavelength between 300 nm and 400 nm (λ_{max} is around 330 nm). Interestingly, the system reverts from the *cis*-form to the *trans*-form after further irradiation with visible light (wavelength over 400 nm). This process is completely reversible, and the azobenzene group does not decompose or induce undesirable side reactions even on repeated *trans-cis* isomerization. By introducing azobenzenes into DNA through D-threoninol as a linker, Asanuma and co-workers succeeded in achieving photo-regulation of:

- Formation and dissociation of a DNA duplex^{3,4} and
- Transcription by T7-RNA polymerase reaction^{5,6,7}.

Item	Catalog No.	Pack	Price(\$)
Azobenzene Phosphoramidite	10-5800-95	50 μ mole	105.00
	10-5800-90	100 μ mole	200.00
	10-5800-02	0.25g	550.00




Azobenzene Phosphoramidite

LABELING WITH ULTRAFAST PHOTO CROSS-LINKER

When 3-cyanovinylcarbazole nucleoside (^{CNVK}) is incorporated into an oligonucleotide, very rapid photo cross-linking to the complementary strand can be induced at one wavelength and rapid reversal of the cross-link is possible at a second wavelength. Neither wavelength has the potential to cause significant DNA damage. Irradiation of a duplex containing a single incorporation of ^{CNVK} at 366nm led to 100% cross-linking to thymine base in 1 second, although complete cross-linking to cytosine takes 25 seconds.¹ A 30 second irradiation time should cover all situations. In addition, it was demonstrated that the purine bases were unreactive to cross-linking, allowing differentiation between pyrimidines and purines at the target site. The authors also determined the effect of sequence contexts around the ^{CNVK} site and demonstrated that the identity of bases on either side of the cross-linking site has little effect on the reaction. Once cross-linked, the UV melting temperature of the duplex was raised by around 30 °C relative to the duplex before irradiation. Complete reversal of the cross-link takes place at 312nm in 3 minutes. This facile reversal reaction is, therefore, accomplished with no damage to normal DNA.

In a later publication, a further application of this cross-linking technique was investigated.² When ^{CNVK} was cross-linked with a dC residue in duplex DNA, heating at 90°C for 3.5 hours led to deamination of the cytosine base to form uracil in the complementary strand. Reversal of the cross-link at 312nm led to a DNA strand in which dC had been converted to dU. The authors showed that this transformation is specific for the dC residue opposite the ^{CNVK} and any further adjacent dC residues are unaffected. Similarly, the authors have shown that ^{CNVK} can be cross-linked to an adjacent RNA strand.³

Item	Cat. No.	Pack	Price (\$)
3-Cyanovinylcarbazole Phosphoramidite (^{CNVK})	10-4960-95	50 μ mole	200.00
	10-4960-90	100 μ mole	390.00
	10-4960-02	0.25g	1125.00



3-Cyanovinylcarbazole

REFERENCES

(1) Y. Yoshimura, and K. Fujimoto, *Org Lett*, 2008, **10**, 3227-30.
 (2) K. Fujimoto, K. Konishi-Hiratsuka, T. Sakamoto, and Y. Yoshimura, *ChemBioChem*, 2010, **11**, 1661-4.
 (3) Y. Yoshimura, T. Ohtake, H. Okada, and K. Fujimoto, *ChemBioChem*, 2009, **10**, 1473-6.

REFERENCES

(1) H. Asanuma, et al., *Angew Chem Int Ed*, 2001, **40**, 2671-2673.
 (2) T. Takarada, et al., *Chem Lett.*, 2001, **30**, 732.
 (3) H. Asanuma, X.G. Liang, T. Yoshida, and M. Komiyama, *ChemBiochem*, 2001, **2**, 39-44.
 (4) H. Asanuma, D. Matsunaga, and M. Komiyama, *NUCLEIC ACIDS SYMP SER (OXF)*, 2005, **49**, 35.
 (5) H. Asanuma, et al., *ChemBiochem*, 2002, **3**, 786.
 (6) M. Liu, H. Asanuma, and M. Komiyama, *J. Amer. Chem. Soc.*, 2006, **128**, 1009.
 (7) H. Asanuma, et al., *Nature Protocols*, 2007, **2**, 203-212.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers	
For Instrument type	Add
Expedite	E
MerMade	M
Columns	
For Instrument type	Add
Expedite	E
Applied Biosystems 3900	A
MerMade	M

(Please inquire for availability of vials and columns for other instrument types.)





RNA SUPPORTS FOR 3' MODIFICATION

Glen Research offers RNA supports in which protected ribonucleosides are attached to CPG. With 5'-DMT protection, and all other protecting groups base-labile, the use of these supports is identical to DNA supports. These supports are suitable for use in producing oligodeoxynucleotides modified at the 3'-terminus or oligoribonucleotides. ABI-style columns are supplied unless otherwise requested (see note box).

Item	Catalog No.	Pack	Price (\$)	
Bz-A-RNA-CPG	20-3303-01	0.1g	40.00	
	20-3303-02	0.25g	95.00	
	20-3303-10	1.0g	355.00	
	1 μmole columns	20-3403-41	Pack of 4	100.00
	0.2 μmole columns	20-3403-42	Pack of 4	75.00
	10 μmole columns (ABI)	20-3403-13	Pack of 1	225.00
	15 μmole column (Expedite)	20-3403-14	Pack of 1	300.00
Ac-C-RNA-CPG	20-3315-01	0.1g	40.00	
	20-3315-02	0.25g	95.00	
	20-3315-10	1.0g	355.00	
	1 μmole columns	20-3415-41	Pack of 4	100.00
	0.2 μmole columns	20-3415-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3415-13	Pack of 1	225.00
	15 μmole column (Expedite)	20-3415-14	Pack of 1	300.00
Ac-G-RNA-CPG	20-3324-01	0.1g	40.00	
	20-3324-02	0.25g	95.00	
	20-3324-10	1.0g	355.00	
	1 μmole columns	20-3424-41	Pack of 4	100.00
	0.2 μmole columns	20-3424-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3424-13	Pack of 1	225.00
	15 μmole column (Expedite)	20-3424-14	Pack of 1	300.00
U-RNA-CPG	20-3330-01	0.1g	40.00	
	20-3330-02	0.25g	95.00	
	20-3330-10	1.0g	355.00	
	1 μmole columns	20-3430-41	Pack of 4	100.00
	0.2 μmole columns	20-3430-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3430-13	Pack of 1	225.00
	15 μmole column (Expedite)	20-3430-14	Pack of 1	300.00

ABBREVIATIONS

Ac = Acetyl
 Bz = Benzoyl
 CNEt = Cyanoethyl
 CPG = Controlled Pore Glass
 DMT = 4,4'-Dimethoxytrityl




TOM-PROTECTED RNA PHOSPHoramidites

RNA synthesis using monomers containing the 2'-O-TriisopropylsilylOxyMethyl (TOM) group (TOM-Protecting-Group™) is characterized by very high coupling efficiency along with fast, simple deprotection. High coupling efficiency is achieved because the TOM-Protecting-Group exhibits lower steric hindrance than the 2'-O-t-butylidimethylsilyl (TBDMS) group used in our alternative RNA monomers. Fast and reliable deprotection is achieved using methylamine in ethanol/water at room temperature. A further feature of the TOM-Protecting-Group is that during basic steps it can not undergo 2' to 3' migration. This migration under basic conditions leads to non-biologically active 2'-5' linkages when using the TBDMS group. These features allow the TOM-Protected monomers to produce longer oligonucleotides. TOM-Protected RNA monomers are also fully compatible with minor bases with 2'-O-TBDMS protection.


Item	Catalog No.	Pack	Price (\$)
A-TOM-CE Phosphoramidite	10-3004-02	0.25g	75.00
	10-3004-05	0.5g	150.00
	10-3004-10	1.0g	275.00
C-TOM-CE Phosphoramidite	10-3014-02	0.25g	75.00
	10-3014-05	0.5g	150.00
	10-3014-10	1.0g	275.00
G-TOM-CE Phosphoramidite	10-3024-02	0.25g	75.00
	10-3024-05	0.5g	150.00
	10-3024-10	1.0g	275.00
U-TOM-CE Phosphoramidite	10-3034-02	0.25g	75.00
	10-3034-05	0.5g	150.00
	10-3034-10	1.0g	275.00

RNA SUPPORTS FOR TOM RNA SYNTHESIS


Item	Catalog No.	Pack	Price (\$)	
Ac-A-RNA-CPG	20-3304-01	0.1g	40.00	
	20-3304-02	0.25g	95.00	
	20-3304-10	1.0g	355.00	
	1 µmole columns	20-3404-41	Pack of 4	100.00
	0.2 µmole columns	20-3404-42	Pack of 4	75.00
10 µmole column (ABI)	20-3404-13	Pack of 1	225.00	
15 µmole column (Expedite)	20-3404-14	Pack of 1	300.00	




A-TOM



C-TOM




G-TOM



U-TOM

INTELLECTUAL PROPERTY



TOM-Protecting-Group™

TOM-RNA Phosphoramidites are supplied under agreement with QIAGEN. RNA synthesis using the TOM-Protecting-Group is covered by US Patent No. 5,986,084.

TOM-Protecting-Group is a trademark of QIAGEN.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

RNA SUPPORTS FOR TOM RNA SYNTHESIS (CONT.)

Item	Catalog No.	Pack	Price (\$)	
Ac-C-RNA-CPG	20-3315-01	0.1g	40.00	
	20-3315-02	0.25g	95.00	
	20-3315-10	1.0g	355.00	
	1 µmole columns	20-3415-41	Pack of 4	100.00
	0.2 µmole columns	20-3415-42	Pack of 4	75.00
10 µmole column (ABI)	20-3415-13	Pack of 1	225.00	
15 µmole column (Expedite)	20-3415-14	Pack of 1	300.00	
Ac-G-RNA-CPG	20-3324-01	0.1g	40.00	
	20-3324-02	0.25g	95.00	
	20-3324-10	1.0g	355.00	
	1 µmole columns	20-3424-41	Pack of 4	100.00
	0.2 µmole columns	20-3424-42	Pack of 4	75.00
10 µmole column (ABI)	20-3424-13	Pack of 1	225.00	
15 µmole column (Expedite)	20-3424-14	Pack of 1	300.00	
U-RNA-CPG	20-3330-01	0.1g	40.00	
	20-3330-02	0.25g	95.00	
	20-3330-10	1.0g	355.00	
	1 µmole columns	20-3430-41	Pack of 4	100.00
	0.2 µmole columns	20-3430-42	Pack of 4	75.00
10 µmole column (ABI)	20-3430-13	Pack of 1	225.00	
15 µmole column (Expedite)	20-3430-14	Pack of 1	300.00	

TBDMS-PROTECTED RNA PHOSPHoramidites

Glen Research CE (β -cyanoethyl) Phosphoramidites for RNA synthesis are produced and packaged to ensure the highest performance on commercial synthesizers. Every batch is accompanied by a Certificate of Analysis and an HPLC trace, showing the results of our QC testing. RNA Phosphoramidites are synthesis-tested with a minimum coupling efficiency of 97%. Glen Research RNA monomers are packaged in industry standard vials which are specially cleaned to eliminate particulate contamination. These monomers are available in a variety of packs, including high throughput (HT) and low cost (LC). An UltraMild set is also available for situations where sensitive bases are in use. Dmf-G (10-3029) has been discontinued and may be substituted with Ac-G (10-3025).

Item	Catalog No.	Pack	Price (\$)
Bz-A-CE Phosphoramidite	10-3003-02	0.25g	40.00
	10-3003-05	0.5g	80.00
	10-3003-10	1.0g	160.00
Ac-C-CE Phosphoramidite	10-3015-02	0.25g	40.00
	10-3015-05	0.5g	80.00
	10-3015-10	1.0g	160.00
Ac-G-CE Phosphoramidite	10-3025-02	0.25g	40.00
	10-3025-05	0.5g	80.00
	10-3025-10	1.0g	160.00
U-CE Phosphoramidite	10-3030-02	0.25g	40.00
	10-3030-05	0.5g	80.00
	10-3030-10	1.0g	160.00


ABBREVIATIONS

Bz = Benzoyl
CNEt = Cyanoethyl
CPG = Controlled Pore Glass
dmf = Dimethylformamide
DMT = 4,4'-Dimethoxytrityl
iPr = Isopropyl
lcaa = long chain alkylamino
Pac = Phenoxyacetyl
PhOAc = Phenoxyacetyl
TBDMS = t-Butyl-dimethylsilyl


RNA PHOSPHoramidites - SPECIAL PACKAGING

We offer our high quality DNA phosphoramidites specifically packaged for high throughput and large-scale synthesis customers. These customers normally require high quality materials produced under the guidelines of a validated quality management system while still being priced aggressively. These products include the usual Glen Research certification and guarantees and they are available in larger packs or in bulk. The core catalog numbers for regular DNA phosphoramidites are shown below. For these products, please request a quote.


Item	Catalog No.	Pack	Price (\$)
Bz-A-CE Phosphoramidite	10-3003-SP		
Ac-C-CE Phosphoramidite	10-3015-SP		
Ac-G-CE Phosphoramidite	10-3025-SP		
U-CE Phosphoramidite	10-3030-SP		




Bz-A-CE Phosphoramidite



Ac-C-CE Phosphoramidite



Ac-G-CE Phosphoramidite



U-CE Phosphoramidite


ULTRAMILD TBDMS RNA PHOSPHoramidites

Item	Catalog No.	Pack	Price (\$)
Pac-A-CE Phosphoramidite	10-3000-02	0.25g	75.00
	10-3000-05	0.5g	150.00
	10-3000-10	1.0g	275.00
Ac-C-CE Phosphoramidite	10-3015-02	0.25g	40.00
	10-3015-05	0.5g	80.00
	10-3015-10	1.0g	160.00
iPr-Pac-G-CE Phosphoramidite	10-3021-02	0.25g	75.00
	10-3021-05	0.5g	150.00
	10-3021-10	1.0g	275.00
U-CE Phosphoramidite	10-3030-02	0.25g	40.00
	10-3030-05	0.5g	80.00
	10-3030-10	1.0g	160.00


TBDMS RNA SUPPORTS

ABI-style columns are supplied for 1 μ mole and 0.2 μ mole scales unless otherwise requested (see note box).


Item	Catalog No.	Pack	Price (\$)	
Pac-A-RNA-CPG	20-3300-01	0.1g	40.00	
	20-3300-02	0.25g	95.00	
	20-3300-10	1.0g	355.00	
	1 μ mole columns	20-3400-41	Pack of 4	100.00
	0.2 μ mole columns	20-3400-42	Pack of 4	75.00
	10 μ mole column (ABI)	20-3400-13	Pack of 1	225.00
	15 μ mole column (Expedite)	20-3400-14	Pack of 1	300.00
Bz-A-RNA-CPG	20-3303-01	0.1g	40.00	
	20-3303-02	0.25g	95.00	
	20-3303-10	1.0g	355.00	
	1 μ mole columns	20-3403-41	Pack of 4	100.00
	0.2 μ mole columns	20-3403-42	Pack of 4	75.00
	10 μ mole column (ABI)	20-3403-13	Pack of 1	225.00
	15 μ mole column (Expedite)	20-3403-14	Pack of 1	300.00




Pac-A-CE Phosphoramidite



Ac-C-CE Phosphoramidite



iPr-Pac-G-CE Phosphoramidite



U-CE Phosphoramidite

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

TBDMS RNA SUPPORTS (CONT.)

Item	Catalog No.	Pack	Price (\$)	
Ac-C-RNA-CPG	20-3315-01	0.1g	40.00	
	20-3315-02	0.25g	95.00	
	20-3315-10	1.0g	355.00	
	1 μmole columns	20-3415-41	Pack of 4	100.00
	0.2 μmole columns	20-3415-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3415-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3415-14	Pack of 1	300.00	
iPr-Pac-G-RNA-CPG	20-3321-01	0.1g	40.00	
	20-3321-02	0.25g	95.00	
	20-3321-10	1.0g	355.00	
	1 μmole columns	20-3421-41	Pack of 4	100.00
	0.2 μmole columns	20-3421-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3421-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3421-14	Pack of 1	300.00	
Ac-G-RNA-CPG	20-3324-01	0.1g	40.00	
	20-3324-02	0.25g	95.00	
	20-3324-10	1.0g	355.00	
	1 μmole columns	20-3424-41	Pack of 4	100.00
	0.2 μmole columns	20-3424-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3424-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3424-14	Pack of 1	300.00	
U-RNA-CPG	20-3330-01	0.1g	40.00	
	20-3330-02	0.25g	95.00	
	20-3330-10	1.0g	355.00	
	1 μmole columns	20-3430-41	Pack of 4	100.00
	0.2 μmole columns	20-3430-42	Pack of 4	75.00
	10 μmole column (ABI)	20-3430-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3430-14	Pack of 1	300.00	

ULTRAMILD SOLVENTS/REAGENTS

Item	Catalog No.	Pack	Price (\$)
<i>Cap Mix A</i> THF/Pyridine/Pac ₂ O (Applied Biosystems)	40-4210-52	200mL	140.00
	40-4210-57	450mL	300.00
THF/Pac ₂ O (Expedite)	40-4212-52	200mL	140.00
	40-4212-57	450mL	300.00
<i>Deprotection Solution</i> 0.05M Potassium Carbonate in Methanol	60-4600-30	30mL	30.00

MINOR RNA PHOSPHoramidites (TOM PROTECTED)

Glenn Research offers minor RNA phosphoramidites with either TOM or TBDMS protecting groups. 4-Thio-U, 5-Methyl-Cytidine, and 2-Amino-Adenosine are useful for analyzing RNA structure and activity relationships, for example, in ribozyme studies.

Pyrrolo-C is a fluorescent nucleoside whose fluorescence is sensitive to its environment and is ideal for probing RNA structure. It base-pairs as a normal C nucleotide. It is highly fluorescent and its excitation and emission are well to the red of most fluorescent nucleotide analogs, which eliminates or reduces background fluorescence from proteins. Pyrrolo-CTP has potential uses in biological assay development.

rSpacer is used to introduce an abasic site to an RNA sequence. The TOM protected version has been discontinued and is replaced with the TBDMS version.

The protecting scheme for 2,6-Diaminopurine has been changed and the original product (10-3084) has been replaced with the optimized product (10-3085) below.

Item	Catalog No.	Pack	Price(\$)
4-Thio-U-TOM-CE Phosphoramidite	10-3052-95	50 μmole	212.50
	10-3052-90	100 μmole	425.00
	10-3052-02	0.25g	975.00
5-Me-C-TOM-CE Phosphoramidite	10-3064-95	50 μmole	95.00
	10-3064-90	100 μmole	190.00
	10-3064-02	0.25g	475.00
2,6-Diaminopurine-TOM-CE Phosphoramidite (2-amino-A)	10-3085-95	50 μmole	212.50
	10-3085-90	100 μmole	425.00
	10-3085-02	0.25g	975.00

SEE ALSO

Minor TBDMS monomers on page 133
Pyrrolo-CTP on page 136
rSpacer TBDMS on page 134

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type *Add*

Expedite E
MerMade M


Columns
For Instrument type *Add*

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)



4-Thio-U-TOM



5-Me-C-TOM



2,6-diaminopurine-TOM


MINOR RNA PHOSPHoramidites (TOM PROTECTED) (CONT.)

Item	Catalog No.	Pack	Price(\$)
Pyrrolo-C-TOM-CE Phosphoramidite	10-3017-95	50 μmole	212.50
	10-3017-90	100 μmole	425.00
	10-3017-02	0.25g	975.00
rSpacer CE Phosphoramidite (See rSpacer TBDMS CE Phosphoramidite)	10-3914-95	50 μmole	DISCONTINUED
	10-3914-90	100 μmole	
	10-3914-02	0.25g	


RNA SEQUENCE MODIFIER (TOM PROTECTED)

Amino-Modifier C6-U has been added to the growing family of sequence modifiers and we envisage applications in RNA structural studies as well as for labeling siRNA to probe uptake and cellular distribution.


Item	Catalog No.	Pack	Price(\$)
Amino-Modifier C6-U Phosphoramidite	10-3039-95	50 μmole	360.00
	10-3039-90	100 μmole	720.00
	10-3039-02	0.25g	1475.00



Pyrrolo-C-TOM



rSpacer



Amino-Modifier C6-U

SEE ALSO

Pyrrolo-dC on page 68
Pyrrolo-CTP on page 136

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add


Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)


MINOR RNA PHOSPHoramidites (TBDMS PROTECTED)

Inosine and 5-Methyl-Uridine are useful for analyzing RNA structure and activity relationships. 5-Bromo-Uridine and 5-Iodo-Uridine have been used for crystallography studies and cross-linking experiments. 6-Thioguanosine (6-thio-G) has applications in ribozyme and siRNA research, as well as in RNA-protein interactions. The removal of the silyl protecting group without interfering with the sulfur is critical. This is removed¹ cleanly by triethylamine trihydrofluoride in DMSO but t-butylammonium fluoride (TBAF) leads to degradation of the thio-nucleotide analogue and should not be used. 2-Aminopurine riboside is useful for analyzing RNA structure and activity relationships, for example, in ribozyme studies.


Item	Catalog No.	Pack	Price(\$)
I-CE Phosphoramidite	10-3040-95	50 μmole	95.00
	10-3040-90	100 μmole	190.00
	10-3040-02	0.25g	475.00
5-Me-U-CE Phosphoramidite (T)	10-3050-95	50 μmole	95.00
	10-3050-90	100 μmole	190.00
	10-3050-02	0.25g	475.00
Br-U-CE Phosphoramidite	10-3090-95	50 μmole	98.00
	10-3090-90	100 μmole	195.00
	10-3090-02	0.25g	475.00
I-U-CE Phosphoramidite	10-3091-95	50 μmole	98.00
	10-3091-90	100 μmole	195.00
	10-3091-02	0.25g	475.00
6-Thio-G-CE Phosphoramidite	10-3072-95	50 μmole	250.00
	10-3072-90	100 μmole	500.00
	10-3072-02	0.25g	1200.00
2-Aminopurine-CE Phosphoramidite	10-3070-95	50 μmole	212.50
	10-3070-90	100 μmole	425.00
	10-3070-02	0.25g	975.00




Inosine




5-Me-U




Br-U



I-U



6-Thio-G



2-Aminopurine

SEE ALSO

Minor TOM monomers on page 131

REFERENCES

- (1) C.J. Adams, J.B. Murray, M.A. Farrow, J.R.P. Arnold, and P.G. Stockley, *Tetrahedron Lett.*, 1995, **36**, 5421-5424.
- (2) D.A. Berry, et al., *Tetrahedron Lett.*, 2004, **45**, 2457-2461.

MINOR RNA (TBDMS PROTECTED) (CONT.)


8-Aza-7-deaza-Adenosine is an isomer of Adenosine with virtually identical electron density. The N7 nitrogen is not available for hydrogen bonding.

Ribozyme activity is substantially affected by the substitution of modified pyrimidine bases. Zebularine (pyrimidin-2-one ribonucleoside) may be regarded as a Cytidine derivative lacking the exocyclic amino group. Zebularine and Pyridin-2-one Ribonucleoside, the 3-deaza analogue of Zebularine, are prime candidates for use in evaluating ribozyme activity and function. It should be noted that Zebularine is mildly fluorescent, absorbing at 298nm and emitting at 367nm.


PseudoUridine is one of the most common modified nucleosides found in RNA. The availability of a phosphoramidite will allow detailed research into the effects of this modified base on RNA structure and activity.

rSpacer is used to introduce an abasic site to an RNA sequence.


Item	Catalog No.	Pack	Price(\$)
Zebularine-CE Phosphoramidite	10-3011-95	50 μ mole	125.00
	10-3011-90	100 μ mole	250.00
	10-3011-02	0.25g	650.00
Pyridin-2-one-CE Phosphoramidite	10-3012-95	50 μ mole	210.00
	10-3012-90	100 μ mole	420.00
	10-3012-02	0.25g	1200.00
PseudoUridine-CE Phosphoramidite	10-3055-95	50 μ mole	175.00
	10-3055-90	100 μ mole	350.00
	10-3055-02	0.25g	995.00
8-Aza-7-deaza-A-CE Phosphoramidite	10-3083-95	50 μ mole	300.00
	10-3083-90	100 μ mole	600.00
	10-3083-02	0.25g	1500.00
rSpacer TBDMS CE Phosphoramidite	10-3915-95	50 μ mole	80.00
	10-3915-90	100 μ mole	135.00
	10-3915-02	0.25g	395.00




Zebularine




Pyridin-2-one




PseudoUridine




8-Aza-7-deaza-A




rSpacer TBDMS



1-Me-A



N6-Me-A



1-Me-Pseudouridine

MINOR RNA (TBDMS PROTECTED) (CONT.)

Methylation of adenosine at position 1 produces a drastic functional change in the nucleobase. 1-Methyladenosine (pK_a 8.25) is a much stronger base than adenosine (pK_a 3.5). N-1 methylation excludes participation of the adenine base in canonical Watson-Crick base pairing and provides a positive charge to the nucleobase. This modification also alters the hydrophobicity of the base, the stacking properties, the ordering of water molecules and the chelation properties. The base may become involved in non-canonical hydrogen bonding, in electrostatic interactions and, in general, it may contribute to the conformational dynamics of the tRNA.

In the central dogma of molecular biology, genetic information flows from DNA to RNA and then to protein. Reversible epigenetic modifications on genomic DNA and histone have been known to substantially regulate gene expression. On the other hand, there exists more than 100 naturally occurring chemical modifications in RNA; however, the functions of these RNA modifications are largely unknown. Whether some of these modifications in RNA can be reversed and could impact gene expression in the central dogma was unknown until the recent discovery of N6-methyladenosine (N6-Me-A) as the first example of reversible RNA methylation.¹ We offer the N6-Me-A RNA monomer with a phenoxyacetyl protecting group to minimize potential branching. We have shown N6-Me-A-CE Phosphoramidite to be completely compatible with all popular RNA synthesis and deprotection methods, from UltraMild to the most popular procedure using AMA for deprotection.

Item	Catalog No.	Pack	Price(\$)
1-Me-A-CE Phosphoramidite	10-3501-95	50 μ mole	190.00
	10-3501-90	100 μ mole	380.00
	10-3501-02	0.25g	975.00
N6-Me-A-CE Phosphoramidite	10-3005-95	50 μ mole	285.00
	10-3005-90	100 μ mole	550.00
	10-3005-02	0.25g	1295.00

RNA methylation occurs in a large selection of RNA nucleosides and this post transcriptional modification of RNA, carried out by a variety of RNA methyltransferases, appears in a wide variety of RNA species - including tRNA, mRNA, miRNA and RNA viruses. Over 90 methylated nucleosides have been found in tRNA and these play many significant roles in tRNA structure. In addition, methylation appears to mark the tRNA as mature, preventing its degradation as well as directing localization within the cell. mRNA, modified with 1-methylpseudouridine (1-Me- Ψ) alone or in combination with 5-methylcytidine (5-Me-C), significantly increases protein expression in cells and mouse models. 1-Me- Ψ is also a modified nucleobase that can greatly enhance the properties of mRNA by reducing immunogenicity and increasing stability.

Item	Catalog No.	Pack	Price(\$)
1-Me-Pseudouridine Phosphoramidite	10-3056-95	50 μ mole	420.00
	10-3056-90	100 μ mole	820.00
	10-3056-02	0.25g	2300.00

REFERENCE

(1) Y. Fu, D. Dominissini, G. Rechavi, and C. He, *Nat Rev Genet*, 2014, **15**, 293-306.

SEE ALSO

[5-Me-C on page 131](#)
[Pseudouridine on page 134](#)

MINOR RNA (TBDMS PROTECTED) (CONT.)



The bright fluorescent tricyclic cytosine analogues tC and tC^o stand out among fluorescent bases due to their virtually unquenched fluorescence inside single- or double-stranded DNA. Until recently, this family of tricyclic cytosines had only been studied and used in DNA contexts and, importantly, introduced as possible donors of the first DNA base analogue FRET-pair with tC_{nitro}. Fluorescent base analogues for RNA are limited in number compared to their DNA counterparts. To facilitate the application of such analogues, characterization of their structural and dynamics behavior in RNA compared to the corresponding natural nucleoside is important. We now introduce the tC^o ribonucleoside, which has been incorporated into a range of RNA sequences, where it was shown to be a very potent and useful fluorophore in this context.¹ Glen Research offers this useful fluorescent ribonucleoside analogue in cooperation with ModyBase HB.

Item	Catalog No.	Pack	Price(\$)
Ribo-tC ^o -CE Phosphoramidite	10-3517-95	50 μmole	245.00
	10-3517-90	100 μmole	470.00
	10-3517-02	0.25g	1195.00

MINOR RNA TRIPHOSPHATES

Pyrrolo-dC is a fluorescent nucleoside that codes as dC and base pairs efficiently with dG. Preliminary evidence indicates that pyrrolo-dC triphosphate is an excellent substrate for Taq, Pfu and Vent polymerases and is incorporated specifically opposite dG. Pyrrolo-dCTP has been available for some time and is in use in biological assays. Pyrrolo-CTP is a fluorescent ribonucleotide with fluorescence exquisitely sensitive to its environment and is of great interest for RNA structural research. The pyrrolo-C project is a joint development by Berry and Associates, Inc. and Glen Research Corporation.

Item	Catalog No.	Pack	Price(\$)
Pyrrolo-CTP 10mM	81-3017-01	100 μL	270.00

Ribo-tC^o

Pyrrolo-CTP

REFERENCE

(1) Füchtbauer, A.F., Preus, S., Börjesson, K., McPhee, S.A., Lilley D.M.J., Wilhelmsson, L.M., *Sci. Rep.*, 2017, **7**, 2393.

SEE ALSO

tC^o on page 70

INTELLECTUAL PROPERTY

These products are offered in collaboration with ModyBase HB.

SEE ALSO

Pyrrolo-dC on page 68
Pyrrolo-C on page 132

2'-OME-RNA PHOSPHoramidITES

Glen Research 2'-Ome-RNA CE (β-cyanoethyl) Phosphoramidites are designed to produce synthetic oligonucleotides containing nuclease resistant 2'-O-methyl ribonucleotide linkages. Deprotection, isolation and handling of 2'-O-methyl oligonucleotides are identical to the procedures for oligodeoxynucleotides.

Item	Catalog No.	Pack	Price(\$)
2'-Ome-A-CE Phosphoramidite	10-3100-90	100 μmole	20.00
	10-3100-02	0.25g	50.00
	10-3100-05	0.5g	100.00
	10-3100-10	1.0g	200.00
2'-Ome-C-CE Phosphoramidite	10-3110-90	100 μmole	20.00
	10-3110-02	0.25g	50.00
	10-3110-05	0.5g	100.00
	10-3110-10	1.0g	200.00
2'-Ome-Ac-C-CE Phosphoramidite	10-3115-90	100 μmole	20.00
	10-3115-02	0.25g	50.00
	10-3115-05	0.5g	100.00
	10-3115-10	1.0g	200.00
2'-Ome-iBu-G-CE Phosphoramidite	10-3120-90	100 μmole	20.00
	10-3120-02	0.25g	50.00
	10-3120-05	0.5g	100.00
	10-3120-10	1.0g	200.00
2'-Ome-G-CE Phosphoramidite	10-3121-90	100 μmole	20.00
	10-3121-02	0.25g	50.00
	10-3121-05	0.5g	100.00
	10-3121-10	1.0g	200.00
2'-Ome-U-CE Phosphoramidite	10-3130-90	100 μmole	20.00
	10-3130-02	0.25g	50.00
	10-3130-05	0.5g	100.00
	10-3130-10	1.0g	200.00

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers

For Instrument type [Add](#)


Expedite E
MerMade M

Columns


For Instrument type [Add](#)

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)




2'-Ome-A




2'-Ome-C




2'-Ome-Ac-C



2'-Ome-iBu-G



2'-Ome-G



2'-Ome-U

ULTRAMILD 2'-OME-RNA


The use of UltraMild monomers in oligonucleotide synthesis has allowed very sensitive dyes like TAMRA, HEX and Cy5 to be used virtually routinely. The DNA and RNA monomers are currently available and we also provide this set of 2'-OMe-RNA monomers. In our version of this chemistry, we use as protecting groups phenoxyacetyl (Pac) for A, acetyl (Ac) for C, and isopropyl-phenoxyacetyl (iPr-Pac) for G.

It has become clear that acetic anhydride in the conventional capping mix can cause transamidation in situations where an amine protecting group is quite labile. This leads to acetyl protection on the amino group that may be slow to be removed. Consequently, if many dG residues are included in the oligonucleotide, we recommend the use of phenoxyacetic anhydride (Pac₂O) in Cap A. This modification removes the possibility of exchange of the iPr-Pac protecting group on the dG with acetate from the acetic anhydride capping mix.


Item	Catalog No.	Pack	Price (\$)
2'-OMe-Pac-A-CE Phosphoramidite	10-3601-02	0.25g	62.50
	10-3601-05	0.5g	125.00
	10-3601-10	1.0g	250.00
2'-OMe-Ac-C-CE Phosphoramidite	10-3115-02	0.25g	50.00
	10-3115-05	0.5g	100.00
	10-3115-10	1.0g	200.00
2'-OMe-iPr-Pac-G-CE Phosphoramidite	10-3621-02	0.25g	62.50
	10-3621-05	0.5g	125.00
	10-3621-10	1.0g	250.00

ULTRAMILD SOLVENTS/REAGENTS


Cap Mix A THF/Pyridine/Pac ₂ O (Applied Biosystems)	40-4210-52	200mL	140.00
	40-4210-57	450mL	300.00
THF/Pac ₂ O (Expedite)	40-4212-52	200mL	140.00
	40-4212-57	450mL	300.00
Deprotection Solution 0.05M Potassium Carbonate in Methanol	60-4600-30	30mL	30.00



2'-OMe-Pac-A



2'-OMe-Ac-C



2'-OMe-iPr-Pac-G

2'-OME-RNA SUPPORTS

ABI-style columns are supplied for 1 μmole and 0.2 μmole scales unless otherwise requested (see note box).

Item	Catalog No.	Pack	Price (\$)
2'-OMe-A-RNA-CPG	20-3600-01	0.1g	40.00
	20-3600-02	0.25g	95.00
	20-3600-10	1.0g	355.00
	20-3700-41	Pack of 4	100.00
	20-3700-42	Pack of 4	75.00
10 μmole column (ABI)	20-3700-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3700-14	Pack of 1	300.00
2'-OMe-C-RNA-CPG	20-3610-01	0.1g	40.00
	20-3610-02	0.25g	95.00
	20-3610-10	1.0g	355.00
	20-3710-41	Pack of 4	100.00
	20-3710-42	Pack of 4	75.00
10 μmole column (ABI)	20-3710-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3710-14	Pack of 1	300.00
2'-OMe-Ac-C-RNA-CPG	20-3615-01	0.1g	40.00
	20-3615-02	0.25g	95.00
	20-3615-10	1.0g	355.00
	20-3715-41	Pack of 4	100.00
	20-3715-42	Pack of 4	75.00
10 μmole column (ABI)	20-3715-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3715-14	Pack of 1	300.00
2'-OMe-G-RNA-CPG	20-3621-01	0.1g	40.00
	20-3621-02	0.25g	95.00
	20-3621-10	1.0g	355.00
	20-3721-41	Pack of 4	100.00
	20-3721-42	Pack of 4	75.00
10 μmole column (ABI)	20-3721-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3721-14	Pack of 1	300.00
2'-OMe-U-RNA-CPG	20-3630-01	0.1g	40.00
	20-3630-02	0.25g	95.00
	20-3630-10	1.0g	355.00
	20-3730-41	Pack of 4	100.00
	20-3730-42	Pack of 4	75.00
10 μmole column (ABI)	20-3730-13	Pack of 1	225.00
15 μmole column (Expedite)	20-3730-14	Pack of 1	300.00

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)

MINOR 2'-OME-RNA PHOSPHoramidites


To aid in the evaluation of the structures of 2'-OME-RNA complexes, we offer the CE phosphoramidites listed below. 2'-OMe-T is useful in triplex studies while the 2-aminopurine derivative may be tested in ribozyme studies. By supporting an additional hydrogen bond, 2,6-diaminopurine (2-amino-adenosine) binds more strongly with uridine than does adenosine. Oligonucleotides containing 2'-OMe-5-Me-C and 2'-OMe-I would be of interest to researchers involved in triplex and antisense studies using 2'-OMe-RNA. The uses of 2'-OMe-5-bromo-U phosphoramidite range from crystallographic studies due to the heavy atom to cross-linking because of its photolability. 5-Fluoro-pyrimidine nucleosides have been useful as therapeutic agents and their effect on the structure and activity of oligonucleotides may be examined using the 2'-OMe-RNA derivatives. The 2,4,6-trimethylphenyl (TMP) protected 2'-OMe-U derivative is a convertible nucleoside and reaction with ammonia leads to the 5-fluoro-dC analogue. 2'-OMe-3-deaza-5-aza-C (Reverse C) derivative has the potential to mimic in oligonucleotides 5-azacytidine, a DNA methylase inhibitor. Its ability to bind as a C will likely be diminished.

ABI-style vials are supplied unless otherwise requested (see note box).


Item	Catalog No.	Pack	Price(\$)
2'-OMe-2-Aminopurine- CE Phosphoramidite (N-dmf-AP)	10-3123-95	50 μ mole	177.50
	10-3123-90	100 μ mole	355.00
	10-3123-02	0.25g	975.00
2'-OMe-2,6-Diaminopurine- CE Phosphoramidite (2-amino-A)	10-3124-95	50 μ mole	177.50
	10-3124-90	100 μ mole	355.00
	10-3124-02	0.25g	975.00
2'-OMe-5-Me-U-CE Phosphoramidite (2'-OMe-T)	10-3131-90	100 μ mole	150.00
	10-3131-02	0.25g	360.00



2'-OMe-2-AP




2'-OMe-2-amino-A




2'-OMe-5-Me-U

MINOR 2'-OME-RNA PHOSPHoramidites (CONT.)


Item	Catalog No.	Pack	Price(\$)
2'-OMe-I-CE Phosphoramidite	10-3140-90	100 μ mole	150.00
	10-3140-02	0.25g	360.00
2'-OMe-5-Me-C-CE Phosphoramidite	10-3160-90	100 μ mole	240.00
	10-3160-02	0.25g	675.00
2'-OMe-5-Br-U-CE Phosphoramidite	10-3190-90	100 μ mole	240.00
	10-3190-02	0.25g	675.00
2'-OMe-TMP-5-F-U-CE Phosphoramidite (2'-OMe-5-F-C Precursor)	10-3111-95	50 μ mole	177.50
	10-3111-90	100 μ mole	355.00
	10-3111-02	0.25g	975.00
2'-OMe-5-F-U-CE Phosphoramidite	10-3132-95	50 μ mole	177.50
	10-3132-90	100 μ mole	355.00
	10-3132-02	0.25g	975.00
2'-OMe-3-deaza-5-aza-C-CE Phosphoramidite has been discontinued.	10-3116		




2'-OMe-I




2'-OMe-5-Me-C




2'-OMe-5-Br-U



2'-OMe-TMP-5-F-U



2'-OMe-5-F-U



2'-OMe-3-deaza-5-aza-C

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.

Monomers
For Instrument type Add

Expedite E
MerMade M

Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M

(Please inquire for availability of vials and columns for other instrument types.)

2'-OME-THIOPHOSPHORAMIDITES

The phosphorodithioate linkage (PS2) is both achiral and essentially resistant to nucleases. Previous studies have shown very interesting results which include observations that DNA with PS2 linkages activates RNase H *in vitro*, strongly inhibits human immunodeficiency virus (HIV) reverse transcriptase, induces B-cell proliferation and differentiation, and is completely resistant to hydrolysis by various nucleases. 2'-OMe- RNA Thiophosphoramidites are RNA monomers designed to produce oligos combining the PS2 linkage with the 2'-O-methyl ribose modification. These PS2-modified RNA oligos have potential for use in siRNAs and dithiophosphate aptamers (thioaptamers).

Item	Catalog No.	Pack	Price(\$)
2'-OMe-A-Thiophosphoramidite	10-3170-90	100 μ mole	300.00
	10-3170-02	0.25g	720.00
2'-OMe-C-Thiophosphoramidite	10-3171-90	100 μ mole	300.00
	10-3171-02	0.25g	720.00
2'-OMe-G-Thiophosphoramidite	10-3172-90	100 μ mole	300.00
	10-3172-02	0.25g	720.00
2'-OMe-U-Thiophosphoramidite	10-3173-90	100 μ mole	300.00
	10-3173-02	0.25g	720.00

SEE ALSO

DNA Thiophosphoramidites on page 40

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type Add

Expedite E
MerMade M


Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)




2'-OMe-A-Thiophosphoramidite



2'-OMe-C-Thiophosphoramidite



2'-OMe-G-Thiophosphoramidite



2'-OMe-U-Thiophosphoramidite

2'-F-RNA PHOSPHORAMIDITES

2'-Deoxy-2'-fluoro-nucleosides adopt an RNA-type sugar conformation, presumably due to the high electronegativity of fluorine. Because of this sugar conformation, RNA duplexes (A-form) are generally more thermodynamically stable than DNA duplexes (B-form). As expected, the addition of 2'-F-RNA residues to oligodeoxynucleotides progressively increases the thermal stability of their duplexes with RNA. The stabilization is additive at approximately 2° per residue. This compares favorably with 2'-OMe-RNA at around 1.5° and RNA at 1.1° per residue. In the meantime, base pair specificity remains intact.

2'-F-RNA phosphodiester linkages are not nuclease resistant, although the corresponding phosphorothioate linkages are highly resistant. Researchers usually design antisense oligonucleotides to form duplexes with RNA, which are then substrates for RNase H. Uniformly modified 2'-F-RNA/RNA duplexes are not substrates for RNase H. However, it is straightforward to prepare chimeric 2'-F-RNA/DNA phosphorothioate oligonucleotides which exhibit enhanced binding to the RNA target, are substrates for RNase H, and are highly nuclease resistant.

Item	Catalog No.	Pack	Price(\$)
2'-F-A-CE Phosphoramidite	10-3400-02	0.25g	100.00
	10-3400-05	0.5g	200.00
2'-F-Ac-C-CE Phosphoramidite	10-3415-02	0.25g	50.00
	10-3415-05	0.5g	100.00
2'-F-G-CE Phosphoramidite	10-3420-02	0.25g	100.00
	10-3420-05	0.5g	200.00
2'-F-U-CE Phosphoramidite	10-3430-02	0.25g	50.00
	10-3430-05	0.5g	100.00

STABILITY NOTE

Synthetic oligonucleotides containing 2'-F-RNA linkages may be deprotected with ammonium hydroxide as normal. Deprotection using AMA at 65°C leads to some degradation and so we recommend the use of AMA at room temperature for 2 hours.

OTHER INSTRUMENT TYPES

All minor bases, RNA products and modifiers are packaged in septum-capped vials suitable for ABI and other instruments. If you would like another type of vial/column add the following to the end of the catalog number.


Monomers
For Instrument type Add

Expedite E
MerMade M


Columns
For Instrument type Add

Expedite E
Applied Biosystems 3900 A
MerMade M


(Please inquire for availability of vials and columns for other instrument types.)




2'-F-A



2'-F-Ac-C



2'-F-G



2'-F-U

2'-F-ARABINONUCLEIC ACID (2'-F-ANA)

REFERENCES

1. E. Viazovkina, M.M. Mangos, M.I. Elzagheid, and M.J. Damha, *Curr Protoc Nucleic Acid Chem*, 2002, **Chapter 4**, Unit 4 15.
2. J.K. Watts, and M.J. Damha, *Can. J. Chem.*, 2008, **86**, 641-656.
3. J.K. Watts, A. Katolik, J. Viladoms, and M.J. Damha, *Org Biomol Chem*, 2009, **7**, 1904-10.
4. A. Kalota, *et al.*, *Nucleic Acids Res.*, 2006, **34**, 451.
5. G.F. Deleavey, *et al.*, *Nucleic Acids Res.*, 2010, **38**, 4547-4557, J.K. Watts, *et al.*, *Nucleic Acids Res.*, 2007, **35**, 1441-1451, T. Dowler, *et al.*, *Nucleic Acids Res.*, 2006, **34**, 1669-1675.

INTELLECTUAL PROPERTY

2'-F-ANA is covered by intellectual property. Key patents covering siRNA and antisense applications are as follows:

WO/2009/146556 (siRNA); WO 03064441 and WO 0220773 (antisense).

STABILITY NOTE


Synthetic oligonucleotides containing 2'-F-RNA linkages may be deprotected with ammonium hydroxide as normal. Deprotection using AMA at 65°C leads to some degradation and so we recommend the use of AMA at room temperature for 2 hours.

Arabinonucleosides are epimers of ribonucleosides with the chiral switch being at the 2' position of the sugar residue. 2'-F-ANA adopts a more DNA-like B-type helix conformation, not through the typical C2'-endo conformation but, rather, through an unusual O4'-endo (east) pucker. However, the presence of the electronegative fluorine leads to a still significant increase (ΔT_m 1.2°C/mod) in melting temperature per modification.¹ 2'-F-ANA-containing oligonucleotides exhibit very high binding specificity to their targets. Indeed, a single mismatch in a 2'-F-ANA – RNA duplex leads to a ΔT_m of -7.2 °C and in a 2'-F-ANA - DNA duplex a ΔT_m of -3.9 °C.²


The presence of fluorine at the 2' position in 2'-F-ANA leads to increased stability to hydrolysis under basic conditions relative to RNA and even 2'-F-RNA.^{1,3} The stability of 2'-F-ANA to nucleases also makes this a useful modification for enhancing the stability of oligonucleotides in biological environments.² 2'-F-ANA hybridizes strongly to target RNA and, unlike most 2' modifications, induces cleavage of the target by RNase H. Phosphorothioate (PS) 2'-F-ANA is routinely used in these applications due to its increased nuclease resistance. Alternating 2'-F-ANA and DNA units provide among the highest potency RNase H-activating oligomers. Both the "altimer" and "gapmer" strand architectures consistently outperform PS-DNA and DNA/RNA gapmers.⁴

siRNA oligos were found to tolerate the presence of 2'-F-ANA linkages very well. High potency gene silencing was demonstrated⁵ with siRNA chimeras containing 2'-F-RNA and/or LNA and 2'-F-ANA. The high efficacy of these chimeras was attributed to the combination of the rigid RNA-like properties of 2'-F-RNA and LNA with the DNA-like properties of 2'-F-ANA.


Item	Catalog No.	Pack	Price(\$)
2'-F-A-ANA CE Phosphoramidite	10-3800-90	100 μ mole	150.00
	10-3800-02	0.25g	375.00
2'-F-Bz-C-ANA CE Phosphoramidite	10-3810-02	0.25g	200.00
	10-3810-05	0.5g	400.00
2'-F-Ac-C-ANA CE Phosphoramidite	10-3815-02	0.25g	200.00
	10-3815-05	0.5g	400.00
2'-F-G-ANA CE Phosphoramidite	10-3820-90	100 μ mole	165.00
	10-3820-02	0.25g	425.00
2'-F-U-ANA CE Phosphoramidite	10-3830-02	0.25g	125.00
	10-3830-05	0.5g	250.00




2'-F-A-ANA




2'-F-Bz-C-ANA



2'-F-Ac-C-ANA



2'-F-G-ANA




2'-F-U-ANA

2'-OME-RNA-PACE PHOSPHORAMIDITES


PACE modifications have enjoyed a resurgence in interest as applied to the field of CRISPR gene editing. In an initial publication, it was shown that single guide RNAs (sgRNA) provided significantly higher activity in cells when 2'-O-methylthiophosphonoacetates were incorporated on the ends of the guide RNA to protect against cellular nucleases.¹ In subsequent studies, 2'-OMe PACE modified sgRNAs were also shown to significantly increase on-target specificity of the CRISPR-Cas9 DNA cleavage in eukaryotic cells. In a recent paper, the incorporation of 2'-OMe PACE modified nucleotides in the 20-nucleotide guide region of the sgRNA was shown to decrease off-target cutting by over an order of magnitude while in most cases increasing the overall on-target efficiency as compared to unmodified single guide RNA.²

As an optimal cycle, we recommend using DCI as an activator (30-3150-XX) and a 15 minute coupling time. Following coupling, cap using Unicap (10-4410-XX) with a regular coupling time and then oxidize using 0.5 M CSO for 3 minutes. Alternatively, a 33 minute coupling time using 0.45 M tetrazole, oxidation using low-water iodine (40-4032-XX) followed by capping with 6.5% DMAP as Cap B will give acceptable results. For deprotection, pre-treat the synthesis column with 1.5% DBU in anhydrous acetonitrile for 60 minutes at room temperature to remove 1,1-dimethyl-2-cyanoethyl protecting groups. Rinse the column with acetonitrile, dry under argon and complete the deprotection with 40% aqueous methylamine for 2 hours at room temperature.


Item	Catalog No.	Pack	Price(\$)
2'-OMe-A-PACE Phosphoramidite	10-3150-02	0.25g	110.00
	10-3150-05	0.5g	220.00
	10-3150-10	1.0g	440.00
2'-OMe-Ac-C-PACE Phosphoramidite	10-3151-02	0.25g	110.00
	10-3151-05	0.5g	220.00
	10-3151-10	1.0g	440.00
2'-OMe-G-PACE Phosphoramidite	10-3152-02	0.25g	110.00
	10-3152-05	0.5g	220.00
	10-3152-10	1.0g	440.00
2'-OMe-U-PACE Phosphoramidite	10-3153-02	0.25g	110.00
	10-3153-05	0.5g	220.00
	10-3153-10	1.0g	440.00




2'-OMe-A-PACE



2'-OMe-Ac-C-PACE



2'-OMe-G-PACE



2'-OMe-U-PACE

REFERENCES

1. A. Hendel, *et al.*, *Nat Biotechnol*, 2015, **33**, 985-989.
2. D.E. Ryan, *et al.*, *Nucleic Acids Res*, 2018, **46**, 792-803.

INTELLECTUAL PROPERTY

These products are covered by patents, US 6,693,187 and 7,067,641, and patents pending owned by Metasense Technologies. Purchase of all or any of these products includes a limited license to use the products solely for the manufacture of oligonucleotides for research use only. This license specifically excludes the use of the product or oligonucleotides containing the product for: (a) therapeutic or diagnostic applications (including kits, pools, libraries and other products or services that incorporate oligonucleotides containing the product), (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), or (c) resale (including sale of kits, pools, libraries and other products or services that incorporate the product or oligonucleotides containing the product). If such activities have commercial application, a separate license is required from Metasense Technologies. Neither the product nor any product created through its use may be used in human clinical trials.

A simple agreement must be signed before end-users and custom oligo services may purchase these products for use as defined above.
<http://www.glenresearch.com/Reference/PACE.pdf>

SEE ALSO

DNA PACE on page 37
 DCI on page 30
 UniCap on page 32
 0.5M CSO on page 32

This page is intentionally left blank.

GLEN-PAK™ PURIFICATION

Glen-Pak™ DNA and RNA cartridges have advantages over Poly-Pak cartridges in that a single loading of the diluted crude deprotection solution is all that is necessary. Also, the range of purification has been extended to 100+ using DMT-ON oligos. In addition, Glen-Pak cartridges allow purification of virtually the complete range of dyes and modifiers.

The Glen-Pak DNA Cartridge 3g is a large cartridge capable of purifying 10-20 µmole oligonucleotide syntheses using the standard DMT-ON procedure and Glen-Pak DNA 30mg 96-Well Plates are for parallel purification of up to 50 nmole scale syntheses. The Glen-Pak DNA 3mg 384-Well Plate is designed for use with 384-well plate compatible vacuum manifold systems and can purify up to a 20 nmole scale synthesis. Each well contains 3mg of Glen-Pak DNA resin, which binds about 15 nmoles of full length 40-mer DMT-ON oligo.

Scale suggestions for the Glen-Pak DNA product line are shown below:

Glen-Pak DNA Product	Catalog Number	Synthesis Scale Compatibility
Glen-Pak DNA 50mg Purification Cartridge	60-5000-96	10 nmole – 200 nmole
Glen-Pak DNA Purification Cartridge	60-5100-XX and 60-5200-XX	10 nmole – 1.0 µmole
Glen-Pak DNA Cartridge 3G	60-5300-01	5 µmole – 20 µmole
Glen-Pak DNA 30mg 96-Well Plate	60-5400-01	10 nmole – 50 nmole
Glen-Pak DNA 3mg 384-Well Plate	60-5500-XX	up to 20 nmole

A User Guide to *Glen-Pak™ Purification* describes in detail the process and several applications for DNA and RNA purification. This booklet is available online at: http://www.glenresearch.com/Technical/GlenPak_UserGuide.pdf.

Item	Catalog No.	Pack	Price (\$)
<i>DNA Purification Cartridges</i>			
Glen-Pak™ 50mg DNA Purification Cartridge (For use in vacuum manifolds and high-throughput devices)	60-5000-96	Pack of 96	415.00
Glen-Pak™ DNA Purification Cartridge (For use in vacuum manifolds and high-throughput devices)	60-5100-10	Pack of 10	80.00
	60-5100-30	Pack of 30	200.00
	60-5100-96	Pack of 96	475.00
Glen-Pak™ DNA Purification Cartridge (For use with disposable syringes)	60-5200-01	each	8.00
	60-5200-10	Pack of 10	80.00
Glen-Pak™ DNA Cartridge 3g	60-5300-01	Pack of 1	150.00
Glen-Pak™ DNA 30mg 96-Well Plate	60-5400-01	Pack of 1	475.00
Glen-Pak™ DNA 3mg 384-Well Plate	60-5500-01	Pack of 1	675.00
	60-5500-10	Pack of 10	6750.00



SEE ALSO

Poly-Pak Reagents on page 149

GLEN-PAK™ PURIFICATION (CONT.)

Item	Catalog No.	Pack	Price (\$)
<i>RNA Purification Cartridges</i>			
Glen-Pak™ RNA Purification Cartridge (For use in vacuum manifolds and high-throughput devices)	60-6100-10	Pack of 10	95.00
	60-6100-30	Pack of 30	225.00
	60-6100-96	Pack of 96	575.00
Glen-Pak™ RNA Purification Cartridge (For use with disposable syringes)	60-6200-01	each	9.50
	60-6200-10	Pack of 10	95.00
<i>Reagents</i>			
RNA Quenching Buffer	60-4120-82	250mL	80.00
	60-4120-80	1L	200.00
<i>Racks and Seals</i>			
Adapter Rack (For use with 96 well manifolds)	60-0010-01	each	20.00
Seal for Adapter Rack (For use on 96 well adapter rack)	60-0020-01	each	30.00

POLY-PAK™ PURIFICATION

The use of Poly-Pak™ packings in cartridges or barrels overcomes several disadvantages usually associated with reverse phase (RP) cartridges. The packing is stable in the pH range 1-13, thus the ammonium hydroxide solution, diluted with water, is loaded directly onto the packing. Also, after elution of failure sequences, the trityl group is removed and washed from the support-bound oligonucleotide. The fully deprotected product can then be eluted and isolated by lyophilization. Poly-Pak™ Cartridges may also be used for desalting normal or labeled oligonucleotides. The original Poly-Pak cartridge and barrel are designed for 0.2 μmole syntheses or less. Poly-Pak II cartridges and barrels are designed for use with 1 μmole syntheses. A booklet, *User Guide To Poly-Pak™ Cartridge Purification*, describes in detail the process and several applications. This booklet is available online at: <http://www.glenresearch.com/Technical/PolyPakBooklet.pdf>.

Item	Catalog No.	Pack	Price (\$)
<i>Packing, Cartridges and Barrels</i>			
Poly-Pak™ Packing	60-1000-05	5g	70.00
	60-1000-25	25g	350.00
Poly-Pak™ Cartridge	60-1100-01	each	8.00
	60-1100-10	Pack of 10	80.00
Poly-Pak™ II Cartridge	60-3100-01	each	12.00
	60-3100-10	Pack of 10	120.00
<i>Reagents</i>			
2.0M Triethylamine Acetate (TEAA) HPLC Grade	60-4110-52	200mL	60.00
	60-4110-57	450mL	120.00
	60-4110-60	960mL	200.00
	60-4110-62	2L	400.00
2% Aqueous Trifluoroacetic Acid	60-4040-57	450mL	36.00



Poly-Pak Cartridge Used Manually

PHYSICAL DATA

The physical data table contains information which is unique to each monomer phosphoramidite. The molecular weight (MW) is the formula weight of the fully-protected monomer phosphoramidite. The MW is used to calculate the volume of solvent required to dilute 0.25g of the monomer to give a final 0.1M concentration. This figure is also shown in the table. The unit molecular weight (Unit FW) is the formula weight of each monomer once inserted into an oligonucleotide with all protecting groups removed. To obtain the molecular weight of a specific oligonucleotide, the following formula is used:

$$\text{Oligonucleotide MW} = \text{Sum of Unit FW} - 61.96$$

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution 10.1M)
10-0001	dA-5'-CE Phosphoramidite	857.95	313.21	0.25g/2.91mL
10-0101	dC-5'-CE Phosphoramidite	833.93	289.18	0.25g/3.00mL
10-0301	dT-5'-CE Phosphoramidite	744.83	304.2	0.25g/3.36mL
10-1000	dA-CE Phosphoramidite	857.95	313.21	0.25g/2.91mL
10-1001	7-Deaza-dA-CE Phosphoramidite	856.96	312.22	0.25g/2.92mL
10-1003	N6-Me-dA-CE Phosphoramidite	767.86	327.24	0.25g/3.26mL
10-1004	3'-dA-CE Phosphoramidite	857.95	313.21	0.25g/2.91mL
10-1006	Etheno-dA-CE Phosphoramidite	777.86	337.23	0.25g/3.21mL
10-1007	8-Br-dA-CE Phosphoramidite	887.81	392.11	0.25g/2.82mL
10-1008	8-oxo-dA-CE Phosphoramidite	873.95	329.21	0.25g/2.86mL
10-1010	dC-CE Phosphoramidite	833.93	289.18	0.25g/3.00mL
10-1014	pdC-CE Phosphoramidite	907.1	327.23	0.25g/2.76mL
10-1015	Ac-dC-CE Phosphoramidite	771.85	289.18	0.25g/3.24mL
10-1016	TMP-F-dU-CE Phosphoramidite	866.97	307.18	0.25g/2.88mL
10-1017	Pyrrolo-dC-CE Phosphoramidite	767.85	327.23	0.25g/3.26mL
10-1018	5-Me-dC Brancher Phosphoramidite	942.1	402.36	0.25g/2.65mL
10-1019	Amino-Modifier C6 dC	1049.14	457.42	0.25g/2.38mL
10-1020	dG-CE Phosphoramidite	839.92	329.21	0.25g/2.98mL
10-1021	7-deaza-dG-CE Phosphoramidite	823.93	328.22	0.25g/3.03mL
10-1027	8-Br-dG-CE Phosphoramidite	903.9	408.1	0.25g/2.77mL
10-1028	8-oxo-dG-CE Phosphoramidite	855.93	345.21	0.25g/2.92mL
10-1029	dmf-dG-CE Phosphoramidite	824.92	329.21	0.25g/3.03mL
10-1030	dT-CE Phosphoramidite	744.83	304.2	0.25g/3.36mL
10-1031	5'-OMe-dT-CE Phosphoramidite	456.48	318.22	0.25g/5.48mL
10-1032	O4-Me-dT-CE Phosphoramidite	758.85	318.22	0.25g/3.29mL
10-1034	4-Thio-dT-CE Phosphoramidite	813.95	320.26	0.25g/3.07mL
10-1035	Carboxy-dT	814.88	360.22	0.25g/3.07mL
10-1036	2-Thio-dT-CE Phosphoramidite	879.02	320.26	0.25g/2.84mL
10-1037	Amino-Modifier C2 dT	938.94	402.3	0.25g/2.66mL
10-1038	Biotin-dT	1285.55	684.7	0.25g/1.94mL
10-1039	Amino-Modifier C6 dT	995.05	458.41	0.25g/2.51mL
10-1040	dl-CE Phosphoramidite	754.79	314.19	0.25g/3.31mL
10-1041	2'-DeoxyNebularine-CE Phosphoramidite (Purine)	738.82	298.19	0.25g/3.38mL
10-1042	O6-Phenyl-dl-CE Phosphoramidite	830.92	Varies	0.25g/3.01mL
10-1044	5-Nitroindole-CE Phosphoramidite	780.86	340.23	0.25g/3.20mL
10-1046	2-Aminopurine-CE Phosphoramidite	809.01	313.21	0.25g/3.09mL
10-1047	dP-CE Phosphoramidite	771.85	330.23	0.25g/3.24mL
10-1048	dK-CE Phosphoramidite	853.96	358.25	0.25g/2.93mL
10-1050	dU-CE Phosphoramidite	730.8	290.17	0.25g/3.42mL
10-1051	O4-Triazolyl-dU-CE Phosphoramidite	781.84	varies	0.25g/3.20mL
10-1052	4-Thio-dU-CE Phosphoramidite	799.93	306.23	0.25g/3.13mL
10-1053	5-OH-dU-CE Phosphoramidite	788.83	306.17	0.25g/3.17mL
10-1054	pdU-CE Phosphoramidite	768.85	328.22	0.25g/3.25mL

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution 10.1M)
10-1055	2'-deoxypseudoU-CE Phosphoramidite	730.8	290.17	0.25g/3.42mL
10-1056	Fluorescein-dT Phosphoramidite	1425.57	815.71	0.25g/1.75mL
10-1057	TAMRA-dT	1311.48	870.85	0.25g/1.91mL
10-1058	Dabcyl-dT	1150.32	709.7	0.25g/2.17mL
10-1059	EDTA-C2-dT-CE Phosphoramidite	1201.32	676.53	0.25g/2.08mL
10-1060	5-Me-dC-CE Phosphoramidite	847.9	303.21	0.25g/2.95mL
10-1061	5-Me-2'-deoxyZebularine-CE Phosphoramidite	728.82	288.19	0.25g/3.43mL
10-1062	5-Hydroxymethyl-dC-CE Phosphoramidite	917	319.21	0.25g/2.73mL
10-1063	5-OH-dC-CE Phosphoramidite	954.03	305.18	0.25g/2.62mL
10-1064	3'-dC-CE Phosphoramidite	833.92	289.18	0.25g/3.00mL
10-1065	dmf-5-Me-isodC-CE Phosphoramidite	798.91	303.21	0.25g/3.13mL
10-1066	5-Carboxy-dC-CE Phosphoramidite	905.97	333.19	0.25g/2.76mL
10-1068	N4-Et-dC-CE Phosphoramidite	757.87	317.42	0.25g/3.30mL
10-1070	O6-Me-dG-CE Phosphoramidite	853.97	343.24	0.25g/2.93mL
10-1072	6-thio-dG-CE Phosphoramidite	934.97	345.26	0.25g/2.67mL
10-1073	7-Deaza-8-aza-dG-CE Phosphoramidite (PPG)	824.91	329.2	0.25g/3.03mL
10-1074	3'-dG-CE Phosphoramidite	824.92	329.21	0.25g/3.03mL
10-1076	7-deaza-dX-CE Phosphoramidite	769.83	329.21	0.25g/3.25mL
10-1078	dmf-isodG-CE Phosphoramidite	1020.13	329.21	0.25g/2.45mL
10-1079	8-Amino-dG-CE Phosphoramidite	895.01	344.22	0.25g/2.79mL
10-1080	5-Br-dC-CE Phosphoramidite	912.82	368.08	0.25g/2.74mL
10-1081	5-l-dC-CE Phosphoramidite	959.83	415.08	0.25g/2.60mL
10-1082	2-F-dl-CE Phosphoramidite	921.96	varies, 2F=332.18	0.25g/2.71mL
10-1083	7-deaza-8-aza-dA-CE Phosphoramidite	808.91	313.2	0.25g/3.09mL
10-1084	3'-dT-CE Phosphoramidite	744.83	304.2	0.25g/3.36mL
10-1085	2-Amino-dA-CE Phosphoramidite	1047.33	328.22	0.25g/2.39mL
10-1086	8-Amino-dA-CE Phosphoramidite	879.01	328.22	0.25g/2.84mL
10-1088	3-deaza-dA-CE Phosphoramidite	856.95	312.22	0.25g/2.92mL
10-1089	Amino-Modifier C6 dA	1068.14	427.4	0.25g/2.34mL
10-1090	5-Br-dU-CE Phosphoramidite	809.69	369.07	0.25g/3.09mL
10-1091	5-l-dU-CE Phosphoramidite	856.69	416.07	0.25g/2.92mL
10-1092	5-F-dU-CE Phosphoramidite	748.79	308.16	0.25g/3.34mL
10-1093	5-Hydroxymethyl-dU-CE Phosphoramidite	802.86	320.19	0.25g/3.11mL
10-1096	Thymidine Glycol CE Phosphoramidite	1007.36	338.21	0.25g/2.48mL
10-1097	AP-dC-CE Phosphoramidite	974.97	438.33	0.25g/2.56mL
10-1098	8,5'-Cyclo-dA CE Phosphoramidite	855.92	311.19	0.25g/2.92mL
10-1100	dA-Me Phosphoramidite	802.91	311.24	0.25g/3.11mL
10-1115	Ac-dC-Me Phosphoramidite	716.81	287.21	0.25g/3.49mL
10-1120	dG-Me Phosphoramidite	784.89	327.24	0.25g/3.19mL
10-1130	dT-Me Phosphoramidite	689.79	302.23	0.25g/3.62mL
10-1140	dA-PACE Phosphoramidite	928.02	354.24	0.25g/2.69mL
10-1150	Ac-dC-PACE Phosphoramidite	841.93	330.21	0.25g/2.97mL
10-1160	dG-PACE Phosphoramidite	910.01	370.24	0.25g/2.75mL
10-1170	dT-PACE Phosphoramidite	814.9	345.22	0.25g/3.07mL
10-1200	dA-H-Phosphonate, TEA Salt	822.9	313.21	0.25g/3.04mL
10-1210	dC-H-Phosphonate, DBU Salt	849.35	289.18	0.25g/2.94mL
10-1220	dG-H-Phosphonate, TEA Salt	804.88	329.21	0.25g/3.11mL
10-1230	dT-H-Phosphonate, TEA Salt	709.78	304.2	0.25g/3.52mL
10-1301	Pac-dA-Me Phosphoramidite	848.93	327.23 (Methyl triester)	0.25g/2.94mL
10-1315	Ac-dC-Me Phosphoramidite	732.81	303.21 (Methyl triester)	0.25g/3.41mL
10-1321	iPr-Pac-dG-Me Phosphoramidite	907.01	343.23 (Methyl triester)	0.25g/2.76mL
10-1330	dT-Me Phosphoramidite	705.79	318.22 (Methyl triester)	0.25g/3.54mL

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution (0.1M)
10-1440	CleanAmp™-Pac-dA-CE Phosphoramidite	1045.25	523.56 (triester)	0.25g/2.39mL
10-1450	CleanAmp™-Ac-dC-CE Phosphoramidite	929.13	499.54 (triester)	0.25g/2.69mL
10-1460	CleanAmp™-Pac-dG-CE Phosphoramidite	1061.25	539.56 (triester)	0.25g/2.36mL
10-1470	CleanAmp™-dT-CE Phosphoramidite	902.11	514.55 (triester)	0.25g/2.77mL
10-1501	1-Me-dA-CE Phosphoramidite	814.31	328.24	0.25g/3.07mL
10-1503	N6-Ac-N6-Me-dA-CE Phosphoramidite	809.89	327.23	0.25g/3.09mL
10-1504	def-dA-CE Phosphoramidite	836.97	313.21	0.25g/2.99mL
10-1510	5-Hydroxymethyl-dC II-CE Phosphoramidite	785.82	319.21	0.25g/3.18mL
10-1511	5-aza-5,6-dihydro-dC-CE Phosphoramidite	787.89	292.18	0.25g/3.17mL
10-1513	N4-Ac-N4-Et-dC-CE Phosphoramidite	799.89	317.24	0.25g/3.13mL
10-1514	5-Formyl-dC-CE Phosphoramidite	915.96	317.19 (formyl) 349.23 (diol)	0.25g/2.73mL
10-1516	tC-CE Phosphoramidite	835.95	395.33	0.25g/2.99mL
10-1517	tC ⁺ -CE Phosphoramidite	819.88	379.26	0.25g/3.05mL
10-1518	tNitro-CE Phosphoramidite	880.94	440.32	0.25g/2.84mL
10-1529	N2-Amino-Modifier C6 dG	965.01	428.38	0.25g/2.59mL
10-1530	5,6-Dihydro-dT-CE Phosphoramidite	746.84	306.21	0.25g/3.35mL
10-1531	N3-Cyanoethyl-dT	797.88	357.26	0.25g/3.13mL
10-1532	5'-Dabsyl-dT-CE Phosphoramidite	729.78	591.53	0.25g/3.43mL
10-1534	N-POM Caged-dT-CE Phosphoramidite	967.99	527.38 (N-POM-dT)	0.25g/2.58mL
10-1535	NHS-Carboxy-dT	897.91	varies, -CO ₂ H=360.22	0.25g/2.78mL
10-1536	Fmoc Amino-Modifier C6 dT	1121.28	458.41(NH ₂)	0.25g/2.23mL
10-1537	dX-CE Phosphoramidite	1069.1	330.19	0.25g/2.34mL
10-1538	S-Bz-Thiol-Modifier C6-dT	1091.26	546.53	0.25g/2.29mL
10-1539	DBCO-dT-CE Phosphoramidite	1214.57	773.77	0.25g/2.06mL
10-1540	C8-Alkyne-dT-CE Phosphoramidite	834.94	394.32	0.25g/2.99mL
10-1541	C8-TIPS-Alkyne-dC-CE Phosphoramidite	1094.4	393.33	0.25g/2.28mL
10-1542	C8-TMS-Alkyne-dC-CE Phosphoramidite	1010.24	393.33	0.25g/2.47mL
10-1543	C8-Alkyne-dC-CE Phosphoramidite	938.06	393.33	0.25g/2.67mL
10-1544	C8-TIPS-Alkyne-dT-CE Phosphoramidite	991.28	394.32	0.25g/2.52mL
10-1545	C8-TMS-Alkyne-dT-CE Phosphoramidite	907.12	394.32	0.25g/2.76mL
10-1550	5,6-Dihydro-dU-CE Phosphoramidite	732.81	292.19	0.25g/3.41mL
10-1554	5-Ethynyl-dU-CE Phosphoramidite	754.81	314.19	0.25g/3.31mL
10-1555	TIPS-5-Ethynyl-dU-CE Phosphoramidite	911.15	314.19	0.25g/2.74mL
10-1560	Ac-5-Me-dC-CE Phosphoramidite	785.86	303.21	0.25g/3.18mL
10-1564	5-Formyl dC III CE Phosphoramidite	950.02	317.19 375.27 (acetal)	0.25g/2.63mL
10-1576	Ferrocene-dT-CE Phosphoramidite	1125.07	684.45	0.25g/2.22mL
10-1585	Pac-2-Amino-dA-CE Phosphoramidite	1042.21	328.22	0.25g/2.40mL
10-1590	Pyrene-dU-CE Phosphoramidite	955.04	514.42	0.25g/2.62mL
10-1591	Perylene-dU-CE Phosphoramidite	1005.1	564.48	0.25g/2.49mL
10-1598	8,5'-Cyclo-dG-CE Phosphoramidite	619.65	327.19	0.25g/4.03mL
10-1601	Pac-dA-CE Phosphoramidite	887.97	313.21	0.25g/2.82mL
10-1621	iPr-Pac-dG-CE Phosphoramidite	946.05	329.21	0.25g/2.64mL
10-1700	dA-Thiophosphoramidite	955.09	345.34(dithioate)	0.25g/1.75mL
10-1710	dC-Thiophosphoramidite	931.07	321.31(dithioate)	0.25g/1.79mL
10-1720	dG-Thiophosphoramidite	937.07	361.34(dithioate)	0.25g/1.78mL
10-1730	dT-Thiophosphoramidite	841.97	336.32(dithioate)	0.25g/1.98mL
10-1900	Chemical Phosphorylation Reagent	656.77	79.98	0.25g/3.81mL
10-1901	Chemical Phosphorylation Reagent II	722.82	79.98	0.25g/3.46mL
10-1902	Solid Chemical Phosphorylation Reagent II	692.79	79.98	0.25g/3.61mL
10-1905	5'-Amino-Modifier 5	577.71	167.1	0.25g/4.33mL

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution (0.1M)
10-1906	5'-Amino-Modifier C6	589.76	179.16	0.25g/4.24mL
10-1907	5'-DMS(O)MT-Amino-Modifier C6	681.34	179.16	0.25g/3.67mL
10-1908	5'-Hexynyl Phosphoramidite	298.36	160.11	0.25g/8.38mL
10-1909	Spacer Phosphoramidite 9	652.77	212.14	0.25g/3.83mL
10-1910	1-Ethynyl-dSpacer CE Phosphoramidite	644.74	204.12	0.25g/3.88mL
10-1912	5'-Amino-Modifier C12	673.92	263.32	0.25g/3.71mL
10-1913	Spacer Phosphoramidite C3	578.69	138.06	0.25g/4.32mL
10-1914	dSpacer CE Phosphoramidite	620.73	180.1	0.25g/4.03mL
10-1915	Pyrrrolidine-CE Phosphoramidite	841.97	178.1	0.25g/2.97mL
10-1916	5'-Amino-Modifier C6-TFA	413.42	179.16	0.25g/6.05mL
10-1917	5'-Amino-Modifier TEG CE-Phosphoramidite	489.47	255.21	0.25g/5.11mL
10-1918	Spacer Phosphoramidite 18	784.93	344.3	0.25g/3.18mL
10-1919	5'-Amino-Modifier-11-CE Phosphoramidite	711.82	271.21	0.25g/3.51mL
10-1920	Symmetric Doubler Phosphoramidite	1095.32	351.31	0.25g/2.28mL
10-1922	Treble Phosphoramidite	1417.72	370.33	0.25g/1.76mL
10-1923	5'-Amino-Modifier C3-TFA	371.34	137.08	0.25g/6.73mL
10-1925	Long Treble Phosphoramidite	1475.78	428.41	0.25g/1.69mL
10-1926	5'-Thiol-Modifier C6	576.78	196.2	0.25g/4.33mL
10-1927	Abasic II Phosphoramidite	750.98	196.1	0.25g/3.33mL
10-1928	Spacer C12 CE Phosphoramidite	704.93	264.3	0.25g/3.55mL
10-1931	5'-I-dT-CE Phosphoramidite	552.35	414.09	0.25g/4.53mL
10-1932	5'-Amino-dT-CE Phosphoramidite	713.81	303.21	0.25g/3.50mL
10-1933	5'-Aldehyde-Modifier C2 Phosphoramidite	480.58	228.14	0.25g/5.20mL
10-1934	5-Formylindole-CE Phosphoramidite	763.86	323.24	0.25g/3.27mL
10-1935	5'-Carboxy-Modifier C10	485.56	varies, -CO ₂ H = 250.23	0.25g/5.15mL
10-1936	Thiol-Modifier C6 S-S	769.05	328.4 (disulfide) 196.2 (thiol)	0.25g/3.25mL
10-1938	5'-Maleimide-Modifier Phosphoramidite	437.47	299.22 (pre-retro-DA) 203.09 (maleimide)	0.25g/5.71mL
10-1939	Spermene Phosphoramidite	1233.17	408.52	0.25g/2.03mL
10-1941	5'-DBCO-TEG Phosphoramidite	708.82	570.57	0.25g/3.53mL
10-1945	5'-Carboxy-Modifier C5	595.11	180.1	0.25g/4.20mL
10-1946	5'-Bromohexyl Phosphoramidite	381.29	243.04 (bromide) 205.15 (azide)	0.25g/6.56mL
10-1947	5'-Amino-Modifier C6-PDA	478.57	179.15	0.25g/5.22mL
10-1948	5'-Amino-Modifier C12-PDA	562.7	263.32	0.25g/4.44mL
10-1949	5'-Amino-Modifier TEG PDA	554.62	255.21	0.25g/4.51mL
10-1952	DesthiobiotinTEG Phosphoramidite	980.19	539.56	0.25g/2.55mL
10-1953	Biotin Phosphoramidite	876.1	435.48	0.25g/2.85mL
10-1955	BiotinTEG Phosphoramidite	1010.24	569.61	0.25g/2.47mL
10-1963	Fluorescein Phosphoramidite	1207.5	598.56	0.25g/2.07mL
10-1964	6-Fluorescein Phosphoramidite	1176.35	566.48	0.25g/2.13mL
10-1973	Acridine Phosphoramidite	891.53	450.86	0.25g/2.80mL
10-1974	5'-GalNAc C3 Phosphoramidite	1206.38	609.61	0.25g/2.07mL
10-1975	Cholesteryl-TEG Phosphoramidite	1196.6	755.97	0.25g/2.09mL
10-1976	5'-Cholesteryl-TEG Phosphoramidite	820.13	682.89	0.25g/3.05mL
10-1977	a-Tocopherol-TEG Phosphoramidite	1139.56	698.91	0.25g/2.19mL
10-1979	Stearyl Phosphoramidite	470.71	332.46	0.25g/5.31mL
10-1981	Asymmetric Doubler (Lev) Phosphoramidite	891.04	352.32	0.25g/2.81mL
10-1982	Psoralen C2 Phosphoramidite	502.55	364.29	0.25g/4.97mL
10-1983	Psoralen C6 Phosphoramidite	558.65	420.4	0.25g/4.48mL
10-1985	DNP-TEG Phosphoramidite	950.00	509.41	0.25g/2.63mL

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution I0.1M)
10-1986	5'-Trimethoxystilbene Cap Phosphoramidite	571.65	433.39	0.25g/4.37mL
10-1987	5'-Pyrene Cap Phosphoramidite	501.6	363.35	0.25g/4.98mL
10-1991	Dithiol Serinol Phosphoramidite	853.08	412.46	0.25g/2.93mL
10-1992	Alkyne-Modifier Serinol Phosphoramidite	758.88	318.26	0.25g/3.29mL
10-1993	Protected Biotin Serinol Phosphoramidite	1051.28	450.45	0.25g/2.38mL
10-1994	6-Fluorescein Serinol Phosphoramidite	1191.3	582.45	0.25g/2.10mL
10-1995	Protected BiotinLC Serinol Phosphoramidite	1298.57	697.74	0.25g/1.93mL
10-1996	COT Serinol Phosphoramidite	822.97	382.35	0.25g/3.04mL
10-1997	Amino-Modifier Serinol Phosphoramidite	887.01	224.15	0.25g/2.82mL
10-1998	DBCO-Serinol Phosphoramidite	909.08	468.45	0.25g/2.75mL
10-2000	Bz-A-LA-CE Phosphoramidite	885.96	341.22	0.25g/2.82mL
10-2011	5-Me-Bz-C-LA-CE Phosphoramidite	875.96	331.22	0.25g/2.85mL
10-2029	dmf-G-LA-CE Phosphoramidite	852.93	357.22	0.25g/2.93mL
10-2030	T-LA-CE Phosphoramidite	772.84	332.20	0.25g/3.23mL
10-3000	Pac-A-CE Phosphoramidite	1018.23	329.21	0.25g/2.46mL
10-3003	Bz-A-CE Phosphoramidite	988.21	329.21	0.25g/2.53mL
10-3004	A-TOM-CE Phosphoramidite	998.24	329.21	0.25g/2.50mL
10-3005	N6-Methyl-A-CE Phosphoramidite	1032.25	343.23	0.25g/2.42mL
10-3011	Zebularine-CE Phosphoramidite	845.05	290.17	0.25g/2.96mL
10-3012	Pyridin-2-one-CE Phosphoramidite	844.06	289.18	0.25g/2.96mL
10-3014	C-TOM-CE Phosphoramidite	974.22	305.18	0.25g/2.57mL
10-3015	Ac-C-CE Phosphoramidite	902.11	305.18	0.25g/2.77mL
10-3017	Pyrrolo-C-TOM-CE Phosphoramidite	970.23	343.27	0.25g/2.58mL
10-3021	iPr-Pac-G-CE Phosphoramidite	1076.31	345.21	0.25g/2.32mL
10-3024	G-TOM-CE Phosphoramidite	1014.24	345.21	0.25g/2.46mL
10-3025	Ac-G-CE Phosphoramidite	941.43	345.21	0.25g/2.66mL
10-3030	U-CE Phosphoramidite	861.06	306.17	0.25g/2.90mL
10-3034	U-TOM-CE Phosphoramidite	933.17	306.17	0.25g/2.68mL
10-3039	Amino-Modifier C6-U Phosphoramidite	1197.41	474.4	0.25g/2.09mL
10-3040	I-CE Phosphoramidite	885.08	330.19	0.25g/2.82mL
10-3050	5-Me-U-CE Phosphoramidite	875.08	320.19	0.25g/2.86mL
10-3052	4-Thio-U-TOM-CE Phosphoramidite	1002.29	322.22	0.25g/2.49mL
10-3055	PseudoUridine-CE Phosphoramidite	861.05	306.17	0.25g/2.90mL
10-3056	1-Methyl-PseudoUridine Phosphoramidite	875.07	320.19	0.25g/2.86mL
10-3064	5-Me-C-TOM-CE Phosphoramidite	988.25	319.21	0.25g/2.53mL
10-3070	2-Aminopurine-TBDMS-CE Phosphoramidite	954.19	329.21	0.25g/2.62mL
10-3072	6-Thio-G-CE Phosphoramidite	1039.31	361.26	0.25g/2.41mL
10-3083	8-Aza-7-deaza-A-CE Phosphoramidite	939.16	329.21	0.25g/2.66mL
10-3085	2,6-Diaminopurine-TOM-CE Phosphoramidite	1113.36	344.22	0.25g/2.25mL
10-3090	Br-U-CE Phosphoramidite	939.96	385.06	0.25g/2.66mL
10-3091	5-I-U-CE Phosphoramidite	986.96	432.07	0.25g/2.53mL
10-3100	2'-OMe-A-CE Phosphoramidite	887.97	343.24	0.25g/2.82mL
10-3110	2'-OMe-C-CE Phosphoramidite	863.95	319.21	0.25g/2.89mL
10-3111	2'-OMe-TMP-5-F-U-CE Phosphoramidite	897.08	337.2	0.25g/2.79mL
10-3115	2'-OMe-Ac-C-CE Phosphoramidite	801.88	319.21	0.25g/3.12mL
10-3116	2'-OMe-3-deaza-5-aza-C-CE Phosphoramidite	816.91	319.21	0.25g/3.06mL
10-3120	2'-OMe-ibu-G-CE Phosphoramidite	869.97	359.24	0.25g/2.87mL
10-3121	2'-OMe-G-CE Phosphoramidite	854.93	359.24	0.25g/2.92mL
10-3123	2'-OMe-2-Aminopurine-CE Phosphoramidite	839.04	343.24	0.25g/2.98mL
10-3124	2'-OMe-2,6-Diaminopurine-CE Phosphoramidite	924.05	358.25	0.25g/2.71mL
10-3130	2'-OMe-U-CE Phosphoramidite	760.82	320.2	0.25g/3.29mL
10-3131	2'-OMe-5-Me-U-CE Phosphoramidite	774.84	334.22	0.25g/3.23mL

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution I0.1M)
10-3132	2'-OMe-5-F-U-CE Phosphoramidite	778.78	338.19	0.25g/3.21mL
10-3140	2'-OMe-I-CE Phosphoramidite	784.85	344.22	0.25g/3.19mL
10-3150	2'-OMe-A-PACE Phosphoramidite	958.07	385.27	0.25g/2.61mL
10-3151	2'-OMe-Ac-C-PACE Phosphoramidite	871.97	361.25	0.25g/2.87mL
10-3152	2'-OMe-G-PACE Phosphoramidite	940.05	401.27	0.25g/2.66mL
10-3153	2'-OMe-U-PACE Phosphoramidite	830.92	362.23	0.25g/3.01mL
10-3160	2'-OMe-5-Me-C-CE Phosphoramidite	815.9	333.24	0.25g/3.06mL
10-3170	2'-OMe-A-Thiophosphoramidite	985.12	375.36	0.25g/1.69mL
10-3171	2'-OMe-C-Thiophosphoramidite	899.02	351.34	0.25g/1.85mL
10-3172	2'-OMe-G-Thiophosphoramidite	967.1	391.36	0.25g/1.72mL
10-3173	2'-OMe-U-Thiophosphoramidite	857.97	352.32	0.25g/1.94mL
10-3190	2'-OMe-5-Br-U-CE Phosphoramidite	839.72	399.09	0.25g/2.98mL
10-3400	2'-F-A-CE Phosphoramidite	875.93	331.2	0.25g/2.85mL
10-3415	2'-F-Ac-C-CE Phosphoramidite	789.84	307.18	0.25g/3.17mL
10-3420	2'-F-G-CE Phosphoramidite	857.91	347.19	0.25g/2.91mL
10-3430	2'-F-U-CE Phosphoramidite	748.79	308.16	0.25g/3.34mL
10-3501	1-Me-A-CE Phosphoramidite	944.57	344.24	0.25g/2.65mL
10-3517	Ribo-tC ⁺ Phosphoramidite	950.16	395.26	0.25g/2.63mL
10-3601	2'-OMe-Pac-A-CE Phosphoramidite	917.99	343.24	0.25g/2.72mL
10-3621	2'-OMe-iPr-Pac-G-CE Phosphoramidite	976.07	359.24	0.25g/2.56mL
10-3800	2'-FANA-A-CE Phosphoramidite	875.93	331.2	0.25g/2.85mL
10-3815	2'-FANA-Ac-C-CE Phosphoramidite	789.83	307.17	0.25g/3.16mL
10-3820	2'-FANA-G-CE Phosphoramidite	857.91	347.19	0.25g/2.91mL
10-3830	2'-FANA-U-CE Phosphoramidite	748.79	308.16	0.25g/3.34mL
10-3914	rSpacer CE Phosphoramidite	823.09	196.09	0.25g/3.04mL
10-3915	rSpacer TBDMS CE Phosphoramidite	750.99	196.09	0.25g/3.33mL
10-4410	UniCap Phosphoramidite	334.39		0.25g/7.48mL
10-4906	PC Amino-Modifier Phosphoramidite	605.59	371.32	0.25g/4.13mL
10-4913	PC Spacer Phosphoramidite	784.88	344.26	0.25g/3.19mL
10-4920	PC Linker Phosphoramidite	699.78	259.15	0.25g/3.57mL
10-4950	PC Biotin Phosphoramidite	1038.25	597.62	0.25g/2.41mL
10-4960	3-Cyanovinylcarbazole Phosphoramidite (CNVK)	836.95	396.33	0.25g/2.99mL
10-5800	Azobenzene Phosphoramidite	815.94	375.32	0.25g/3.06mL
10-5901	5'-Fluorescein Phosphoramidite	843.95	537.46	0.25g/2.96mL
10-5902	5'-Hexachloro-Fluorescein Phosphoramidite	1050.62	744.13	0.25g/2.38mL
10-5903	5'-Tetrachloro-Fluorescein Phosphoramidite	981.73	675.24	0.25g/2.55mL
10-5905	SIMA (HEX) Phosphoramidite	1065.02	759.54	0.25g/2.35mL
10-5906	5'-Dichloro-dimethoxy-Fluorescein Phosphoramidite II	972.88	666.4	0.25g/2.57mL
10-5912	5'-DabcyI Phosphoramidite	568.69	430.18	0.25g/4.40mL
10-5913	Cyanine 3 Phosphoramidite	953.64	507.59	0.25g/2.62mL
10-5914	Cyanine 3.5 Phosphoramidite	1053.76	607.7	0.25g/2.37mL
10-5915	Cyanine 5 Phosphoramidite	979.68	533.63	0.25g/2.55mL
10-5916	Cyanine 5.5 Phosphoramidite	1171.25	633.74	0.25g/2.13mL
10-5920	Redmond Red [®] Phosphoramidite	971.09	445.34	0.25g/2.57mL
10-5921	Yakima Yellow [®] Phosphoramidite	1023.81	718.33	0.25g/2.44mL
10-5923	5'-AquaPhluor [®] 593 CE Phosphoramidite	1239.17	787.82	0.25g/2.02mL
10-5924	5'-CDPI ₃ MGB [™] Phosphoramidite	1323.42	872.96	0.25g/1.89mL
10-5925	Eclipse [®] Quencher Phosphoramidite	978.5	537.89	0.25g/2.55mL
10-5931	5'-BHQ-1 Phosphoramidite	676.75	538.49	0.25g/3.69mL
10-5932	5'-BHQ-2 Phosphoramidite	678.72	540.47	0.25g/3.68mL
10-5934	5'-BBQ-650 [®] -CE Phosphoramidite	802.9	665.65	0.25g/3.11mL
10-5941	BHQ-1-dT	1401.56	960.93	0.25g/1.78mL

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution (0.1M)
10-5942	BHQ-2-dT	1403.53	962.91	0.25g/1.78mL
10-5944	BBQ-650®-dT-CE Phosphoramidite	1441.57	1000.95	0.25g/1.73mL
10-5945	SIMA (HEX)-dT Phosphoramidite	1646.64	1037.79	0.25g/1.52mL
10-5950	5'-Biotin Phosphoramidite	846.08	405.45	0.25g/2.95mL
10-5961	Methylene Blue II Phosphoramidite	967.67	489.57	0.25g/2.58mL
10-7001	2',3'-ddA-CE Phosphoramidite	574.7	297.21	0.25g/4.35mL
10-7101	2',3'-ddC-CE Phosphoramidite	550.68	273.18	0.25g/4.54mL
10-7201	2',3'-ddG-CE Phosphoramidite	506.54	313.2	0.25g/4.94mL
10-7301	2',3'-ddT-CE Phosphoramidite	426.45	288.19	0.25g/5.86mL
10-9201	dmf-dG-5'-CE Phosphoramidite	824.92	329.21	0.25g/3.03mL
11-1330	Cis-syn Thymine Dimer Phosphoramidite	1024.01	608.39	0.25g/2.44mL
13-1000	AAA Trimer Phosphoramidite	1911.5		0.25g/1.31mL
13-1001	AAC Trimer Phosphoramidite	1887.5		0.25g/1.32mL
13-1011	ACC Trimer Phosphoramidite	1863.5		0.25g/1.34mL
13-1013	ACT Trimer Phosphoramidite	1774.5		0.25g/1.41mL
13-1020	AGA Trimer Phosphoramidite	1893.5		0.25g/1.32mL
13-1031	ATC Trimer Phosphoramidite	1774.5		0.25g/1.41mL
13-1032	ATG Trimer Phosphoramidite	1780.5		0.25g/1.40mL
13-1102	CAG Trimer Phosphoramidite	1869.5		0.25g/1.34mL
13-1103	CAT Trimer Phosphoramidite	1774.5		0.25g/1.41mL
13-1110	CCA Trimer Phosphoramidite	1863.5		0.25g/1.34mL
13-1112	CCG Trimer Phosphoramidite	1845.5		0.25g/1.35mL
13-1122	CGG Trimer Phosphoramidite	1851.5		0.25g/1.35mL
13-1123	CGT Trimer Phosphoramidite	1756.5		0.25g/1.42mL
13-1132	CTG Trimer Phosphoramidite	1756.5		0.25g/1.42mL
13-1200	GAA Trimer Phosphoramidite	1893.5		0.25g/1.32mL
13-1201	GAC Trimer Phosphoramidite	1869.5		0.25g/1.34mL
13-1203	GAT Trimer Phosphoramidite	1780.5		0.25g/1.40mL
13-1210	GCA Trimer Phosphoramidite	1869.5		0.25g/1.34mL
13-1212	GCG Trimer Phosphoramidite	1851.5		0.25g/1.35mL
13-1213	GCT Trimer Phosphoramidite	1756.5		0.25g/1.42mL
13-1223	GGT Trimer Phosphoramidite	1762.5		0.25g/1.42mL
13-1230	GTA Trimer Phosphoramidite	1780.5		0.25g/1.40mL
13-1233	GTT Trimer Phosphoramidite	1667.5		0.25g/1.50mL
13-1301	TAC Trimer Phosphoramidite	1774.5		0.25g/1.41mL
13-1313	TCT Trimer Phosphoramidite	1661.4		0.25g/1.50mL
13-1321	TGC Trimer Phosphoramidite	1756.5		0.25g/1.42mL
13-1322	TGG Trimer Phosphoramidite	1762.5		0.25g/1.42mL
13-1331	TTC Trimer Phosphoramidite	1661.4		0.25g/1.50mL
13-1333	TTT Trimer Phosphoramidite	1572.4		0.25g/1.59mL
20-0002	dA-5'-CPG		313.21	
20-0102	dC-5'-CPG		289.18	
20-0202	dG-5'-CPG		329.21	
20-0302	dT-5'-CPG		304.2	
20-2000	dA-CPG 500		313.21	
20-2001	dA-CPG 1000		313.21	
20-2002	dA-CPG 2000		313.21	
20-2004	3'-dA-CPG		313.21	
20-2010	dC-CPG 500		289.18	
20-2011	dC-CPG 1000		289.18	
20-2012	dC-CPG 2000		289.18	
20-2013	Ac-dC-CPG 500		289.18	
20-2015	Ac-dC-CPG 1000		289.18	

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution (0.1M)
20-2017	2',3'-ddC-CPG		273.19	
20-2019	3'-Amino-Modifier C6 dC CPG		457.42	
20-2020	dG-CPG 500		329.21	
20-2021	dG-CPG 1000		329.21	
20-2022	dG-CPG 2000		329.21	
20-2029	dmf-dG-CPG		329.21	
20-2030	dT-CPG 500		304.2	
20-2031	dT-CPG 1000		304.2	
20-2032	dT-CPG 2000		304.2	
20-2040	dI-CPG 500		314.19	
20-2041	dI-CPG 1000		314.19	
20-2050	dU-CPG 500		290.17	
20-2051	dU-CPG 1000		290.17	
20-2056	3'-Fluorescein-dT CPG		815.71	
20-2064	3'-dC-CPG		289.18	
20-2074	3'-dG-CPG		329.21	
20-2084	3'-dT-CPG		304.2	
20-2090	5-Br-dU-CPG		369.07	
20-2101-61	dA-CPG 1000		313.21	
20-2101-62	dA-CPG 1000		313.21	
20-2101-65	dA-CPG 1000		313.21	
20-2115-61	Ac-dC-CPG 1000		289.18	
20-2115-62	Ac-dC-CPG 1000		289.18	
20-2115-65	Ac-dC-CPG 1000		289.18	
20-2129-61	dmf-dG-CPG		329.21	
20-2129-62	dmf-dG-CPG		329.21	
20-2129-65	dmf-dG-CPG		329.21	
20-2131-61	dT-CPG 1000		304.2	
20-2131-62	dT-CPG 1000		304.2	
20-2131-65	dT-CPG 1000		304.2	
20-2601	Pac-dA-CPG		313.21	
20-2621	iPr-Pac-dG-CPG		329.21	
20-2900	3'-Phosphate CPG		79.98	
20-2902	3'-Glyceryl CPG		154.06	
20-2903	3'-CPR II CPG		79.98	
20-2913	3'-Spacer C3 CPG		138.06	
20-2933	3'-Thiol-Modifier C3 S-S CPG	154.12 (thiol), 244.27 (disulfide)		
20-2938	3'-Thiol-Modifier 6 S-S CPG	198.18 (thiol), 332.37 (disulfide)		
20-2952	DesthiobiotinTEG-CPG		539.56	
20-2954	3'-PT-Amino-Modifier C3 CPG		137.07	
20-2955	3'-BiotinTEG CPG		569.61	
20-2956	3'-PT-Amino-Modifier C6 CPG		179.15	
20-2958	3'-Amino-Modifier C7 CPG 1000		209.18	
20-2961	3'-(6-FAM) CPG		569.46	
20-2963	3'-Fluorescein CPG		598.56	
20-2964	3'-(6-Fluorescein) CPG		566.48	
20-2973	3'-Acridine CPG		450.86	
20-2974	GalNAc C3 CPG		609.61	
20-2975	3'-Cholesteryl-TEG CPG		755.97	
20-2980	3'-Uaq Cap CPG		539.39	
20-2981	3'-Amino-dT CPG		303.21	
20-2982	3'-Propargyl-5-Me-dC CPG		341.26	
20-2991	3'-Dithiol Serinol CPG		412.46	

PHYSICAL DATA

Cat. No.	Item	Phosphoramidite MW	Unit FW	Dilution (0.1M)
20-2992	3'-Alkyne-Modifier Serinol CPG		334.26	
20-2993	3'-Protected Biotin Serinol CPG		450.45	
20-2994	3'-6-Fluorescein Serinol CPG		584.47	
20-2995	3'-Protected BiotinLC Serinol CPG		697.74	
20-2997	3'-Amino-Modifier Serinol CPG		224.15	
20-3300	Pac-A-RNA-CPG		329.21	
20-3303	Bz-A-RNA-CPG		329.21	
20-3304	Ac-A-RNA-CPG		329.21	
20-3315	Ac-C-RNA-CPG		305.18	
20-3321	iPr-Pac-G-RNA-CPG		345.21	
20-3324	Ac-G-RNA-CPG		345.21	
20-3330	U-RNA-CPG		306.17	
20-3600	2'-OMe-A-RNA-CPG		343.24	
20-3610	2'-OMe-C-RNA-CPG		319.21	
20-3615	2'-OMe-Ac-C-RNA-CPG		319.21	
20-3621	2'-OMe-G-RNA-CPG		359.24	
20-3630	2'-OMe-U-RNA-CPG		320.2	
20-4040	Puromycin-CPG		533.48	
20-5910	3'-TAMRA CPG		623.6	
20-5911	3'-Dabsyl CPG		498.49	
20-5912	3'-Dabcyl CPG		462.44	
20-5913	Cyanine 3 CPG		507.59	
20-5915	Cyanine 5 CPG		533.63	
20-5920	Redmond Red® CPG		445.34	
20-5921	Yakima Yellow® CPG		718.33	
20-5923	AquaPhluor® 593 CPG		900.93	
20-5924	CDPI ₃ MGB™ CPG		831.87	
20-5925	Eclipse® Quencher CPG		537.89	
20-5931	3'-BHQ-1 CPG		554.49	
20-5932	3'-BHQ-2 CPG		556.47	
20-5933	3'-BHQ-3 CPG		597.63	
20-5934	BBQ-650® CPG		667.63	
21-2000	dA-Q-CPG 500		313.21	
21-2010	dC-Q-CPG 500		289.18	
21-2013	Ac-dC-Q-CPG 500		305.18	
21-2029	dmf-dG-Q-CPG 500		329.21	
21-2030	dT-Q-CPG 500		304.2	
25-2000	dA-High Load-CPG		313.21	
25-2010	dC-High Load-CPG		289.18	
25-2020	dG-High Load CPG		329.21	
25-2030	dT-High Load-CPG		304.2	
25-2900	3'-Phosphate CPG (High Load)		79.98	
26-2600	dA PS		313.21	
26-2610	dC PS		289.18	
26-2629	dmf-dG PS		329.21	
26-2630	dT-PS		304.2	
26-2900	3'-Phosphate PS		79.98	
26-2955	3'-BiotinTEG PS		569.61	
26-2956	3'-PT-Amino-Modifier C6 PS		179.15	
26-2961	3'-(6-FAM) PS		569.46	
26-5910	3'-TAMRA PS		623.6	
26-5912	3'-Dabcyl PS		462.44	

PHYSICAL DATA

Cat. No.	Item	MW	Unit FW
50-1904	Azidobutyrate NHS Ester	226.19	113.12
50-1905	Alkyne-NHS Ester	225.2	110.11
50-1941	DBCO-sulfo-NHS Ester	532.5	316.37
50-2000	BiotinTEG Azide	444.55	
50-2001	DesthiobiotinTEG Azide	414.5	
50-2002	Dipivaloyl 6-FAM-TEG Azide	744.79	
50-2003	6-FAM-TEG Azide	576.55	
50-2004	Coumarin Azide	203.15	
50-2005	6-HEX Azide	665.09	
50-2006	6-TET Azide	596.2	
50-2007	TEMPO Azide	197.26	
50-2008	TEMPO-TEG Azide	373.47	
50-2009	Psoralen Azide	283.28	
50-2010	Disulfo-Cyanine 7 Azide	829.08	
50-5910	TAMRA NHS Ester	527.53	413.45

Index

Symbols

2'-5' Linkages

- 3'-dA-CE Phosphoramidite 65
- 3'-dA-CPG 65
- 3'-dC-CE Phosphoramidite 65
- 3'-dC-CPG 65
- 3'-dG-CE Phosphoramidite 65
- 3'-dG-CPG 66
- 3'-dT-CE Phosphoramidite 65
- 3'-dT-CPG 66

2'-5' Linked Oligonucleotides 64, 66

5' -> 3' SYNTHESIS

- 5'-CE Phosphoramidites 34

A

A

- 1-Me-A-CE Phosphoramidite 135
- 1-Me-dA-CE Phosphoramidite 66
- 2',3'-ddA 56
- 2'-F-A-ANA CE Phosphoramidite 144
- 2'-F-A-CE Phosphoramidite 143
- 2'-OMe-A-CE Phosphoramidite 137
- 2'-OMe-A-PACE Phosphoramidite 145
- 2'-OMe-A-RNA 139
- 2'-OMe-Pac-A-CE Phosphoramidite 138
- 3'-dA-CE Phosphoramidite 65
- 3'-dA-CPG 55, 66
- 7-Deaza-dA-CE Phosphoramidite 57
- 8-Amino-dA-CE Phosphoramidite 59
- 8-Br-dA-CE Phosphoramidite 60
- 8-Oxo-dA-CE Phosphoramidite 61
- Ac-A-RNA-CPG 126
- Amino-Modifier C6 dA 77
- A-TOM-CE Phosphoramidite 126
- Bz-A-CE Phosphoramidite 128
- Bz-A-LNA-CE Phosphoramidite 41
- Bz-A-RNA-CPG 125, 129
- dA-CE Phosphoramidite 8, 12, 15, 16, 18, 20
- dA-H-Phosphonate 39
- dA-Me Phosphoramidite 36
- dA-Me Phosphoramidite 38
- dA-PACE Phosphoramidite 37
- dA-Thiophosphoramidite 40, 142
- def-dA-CE Phosphoramidite 22
- dma-dA-CE Phosphoramidite 22
- N6-Ac-N6-Me-dA-CE Phosphoramidite 47, 66
- N6-Me-A-CE Phosphoramidite 135
- N6-Me-dA-CE Phosphoramidite 47, 66
- Pac-A-CE Phosphoramidite 129
- Pac-A-RNA-CPG 129
- Pac-dA-CE Phosphoramidite 23

A-2-Amino

- 2-Amino-A-TOM-CE Phosphoramidite 131
- 2-Amino-dA-CE Phosphoramidite 47

- 2'-OMe-2-Amino-A-CE Phosphoramidite 140
- Pac-2-Amino-dA-CE Phosphoramidite 47

Abasic Site 62, 84

- Abasic II Phosphoramidite 62
- dSpacer Phosphoramidite 62
- Pyrrolidine-CE Phosphoramidite 63

Acridine Labelling

- 3'-Acridine CPG 114
- Acridine Phosphoramidite 114

Activator (Powder)

- 4,5-Dicyanoimidazole 30
- 5-Benzylthio-1H-tetrazole 30
- 5-Ethylthio-1H-tetrazole 26, 30, 70, 110, 112
- Saccharin 1-Methylimidazole 30

Adamantane Carbonyl Chloride 39

Affinity Chromatography 151

ÄKTA oligopilot 18, 19

Aldehyde Modifier

- 5'-Aldehyde-Modifier C2 Phosphoramidite 83
- 5-Formyl-dC-CE Phosphoramidite 50
- Formylindole CE Phosphoramidite 83

Alternative Solvents and Reagents 30

Amino-dA

- 8-Amino-dA-CE Phosphoramidite 59

Amino-dG

- 8-Amino-dG-CE Phosphoramidite 62

Amino-Modifiers

- 3'-Amino-Modifier C6 dC CPG 81
- 3'-Amino-Modifier C6 dT CPG 81
- 3'-Amino-Modifier Serinol CPG 79
- 3'-PT-Amino-Modifier C3 CPG 79
- 3'-PT-Amino-Modifier C6 CPG 79
- 3'-PT-Amino-Modifier C6 PS 79
- 5'-Amino-dT-CE Phosphoramidite 55
- 5'-Amino-Modifier 5 74
- 5'-Amino-Modifier C3-TFA 74, 75
- 5'-Amino-Modifier C6 74, 75
- 5'-Amino-Modifier C6-PDA 75
- 5'-Amino-Modifier C6-TFA 74, 75
- 5'-Amino-Modifier C12 74
- 5'-Amino-Modifier C12-PDA 75
- 5'-Amino-Modifier TEG 74
- 5'-Amino-Modifier-TEG-PDA 75
- 5'-DMS(O)MT-Amino-Modifier C6 74
- Amino-Modifier C2 dT 77
- Amino-Modifier C6 dA 77
- Amino-Modifier C6 dC 77
- Amino-Modifier C6 dT 77
- Amino-Modifier C6-U Phosphoramidite 132
- Amino-Modifier Serinol Phosphoramidite 78
- Fmoc-Amino-Modifier C6 dT 78
- N2-Amino-Modifier C6 dG 77
- PC Amino-Modifier Phosphoramidite 76, 86, 95

AminoOxy-Modifier

- 5'-AminoOxy-Modifier 11 76

Aminopurine

- 2'-OMe-2-Aminopurine-CE Phosphoramidite 140

Anthraquinone

- 3'-Uaq Cap CPG 49

Applied Biosystems Instruments

- AB 3900 1000Å CPG Columns 10
- AB 3900 Polystyrene Columns 10
- AB 3900 Polystyrene Modifier Columns 11
- CE Phosphoramidites 8
- Solvents/Reagents 8
- Supports and Columns 9

Aptamer Development 73

AquaPhluor® 593

- 5'-AquaPhluor® 593 Phosphoramidite 108
- AquaPhluor® 593 CPG 109

Aza-dC

- 5-Aza-5,6-dihydro-dC-CE Phosphoramidite 71

Azidobutyrate NHS Ester 89

Azobenzene

- Azobenzene Phosphoramidite 122

B

Benzylthio-1H-tetrazole 30

Biocompatible Chemical Ligation 64

Biotin Labelling

- 3'-BiotinTEG CPG 101
- 3'-BiotinTEG PS 101
- 3'-Protected BiotinLC Serinol CPG 96, 101
- 3'-Protected Biotin Serinol CPG 95, 101
- 5'-Biotin Phosphoramidite 100
- Biotin-dT 100
- Biotin Phosphoramidite 99
- BiotinTEG Azide 92
- BiotinTEG Phosphoramidite 99
- DesthiobiotinTEG Azide 92
- DesthiobiotinTEG-CPG 101
- DesthiobiotinTEG Phosphoramidite 100
- PC Biotin Phosphoramidite 86, 100
- Protected BiotinLC Serinol Phosphoramidite 95, 99
- Protected Biotin Serinol Phosphoramidite 94, 99

BlackBerry® Quencher

- 3'-BBQ-650® CPG 112
- 5'-BBQ-650® Phosphoramidite 112
- BBQ-650®-dT 112

Black Hole Quencher™ Dyes

- 3'-BHQ-1 CPG 111
- 3'-BHQ-2 CPG 111
- 3'-BHQ-3 CPG 111, 112
- 5'-BHQ-1 Phosphoramidite 110
- 5'-BHQ-2 Phosphoramidite 70, 110, 112
- BHQ-1-dT 110, 112
- BHQ-2-dT 110

Brancher Phosphoramidite

- dC Brancher Phosphoramidite 85

Br-dA

- 8-Br-dA-CE Phosphoramidite 60

Br-dC

- 5-Br-dC-CE Phosphoramidite 60

Br-dG

- 8-Br-dG-CE Phosphoramidite 60

Br-dU

- 5-Br-dU-CE Phosphoramidite 60
- 5-Br-dU-CPG 60

Bromoethyl Phosphoramidite 89

Br-U

- 2'-OMe-5-Br-U-CE Phosphoramidite 141
- 5-Br-U-CE Phosphoramidite 133

C

C

- 2',3'-ddC 56
- 2',3'-ddC-CPG 56
- 2'-F-Bz-C-ANA CE Phosphoramidite 144
- 2'-OMe-5-Me-C-CE Phosphoramidite 141
- 2'-OMe-Ac-C-CE Phosphoramidite 137, 138
- 2'-OMe-Ac-C-PACE Phosphoramidite 145
- 2'-OMe-Ac-C-RNA 139
- 2'-OMe-C-CE Phosphoramidite 137
- 2'-OMe-C-RNA 139
- 3'-Amino-Modifier C6 dC CPG 81
- 3'-dC-CE Phosphoramidite 65
- 3'-dC-CPG 55, 65
- 5-Br-dC-CE Phosphoramidite 60
- 5-Carboxy-dC-CE Phosphoramidite 50
- 5-Formyl-dC-CE Phosphoramidite 50
- 5-Hydroxymethyl-dC-CE Phosphoramidite 50
- 5-Hydroxymethyl-dC II-CE Phosphoramidite 50
- 5-I-dC-CE Phosphoramidite 60
- 5-Me-Bz-C-LNA-CE Phosphoramidite 41
- 5-Me-dC-CE Phosphoramidite 46
- 5-OH-dC-CE Phosphoramidite 61
- Ac-C-CE Phosphoramidite 128
- Ac-C-RNA-CPG 125, 127, 130
- Ac-dC-CE Phosphoramidite 8, 12, 15, 16, 18, 20, 23
- Ac-dC-Me Phosphoramidite 36
- Ac-dC-PACE Phosphoramidite 37
- Amino-Modifier C6 dC 77
- AP-dC 68
- AP-dC-CE Phosphoramidite 46
- C8-Alkyne-dC-CE Phosphoramidite 87, 90
- C8-TIPS-Alkyne-dC-CE Phosphoramidite 87
- C8-TMS-Alkyne-dC-CE Phosphoramidite 87
- C-TOM-CE Phosphoramidite 126
- dC Brancher Phosphoramidite 85
- dC-CE Phosphoramidite 8, 12, 15, 16, 18, 20
- dC-H-Phosphonate 39
- dC-Me Phosphoramidite 38
- dC-Thiophosphoramidite 40, 142
- N4-Et-dC-CE Phosphoramidite 47
- pdC-CE Phosphoramidite 46
- Pyrrolo-dCTP 68

tC-CE Phosphoramidite 70
tC°-CE Phosphoramidite 70

Camphorsulfonyloxaziridine (CSO) 32

Cap CPG
3'-Uaq Cap CPG 49, 64

Cap Phosphoramidite
5'-Pyrene Cap Phosphoramidite 49
5'-Trimethoxystilbene Cap Phosphoramidite 49

Capping Reagent
UniCap Phosphoramidite 32

Carboxy-dC
5-Carboxy-dC-CE Phosphoramidite 50

Carboxy-Modifiers
5'-Carboxy-Modifier C5 76
5'-Carboxy-Modifier C10 76
Carboxy-dT 77

Chain Terminators 56

Chelates
EDTA-C2-dT-CE Phosphoramidite 118

Chemical Phosphorylation 82

Cholesterol Labelling
3'-Cholesteryl-TEG CPG 115
5'-Cholesteryl-TEG Phosphoramidite 115
Cholesteryl-TEG Phosphoramidite 75, 115

CleanAmp™ Technology
CleanAmp™ Primers 54

Click Chemistry 87
1,2,3-triazoles 87
1-Ethynyl-dSpacer CE Phosphoramidite 90
3'-Alkyne-Modifier Serinol CPG 80, 89, 97
3'-Propargyl-5-Me-dC CPG 64
5'-Bromohexyl Phosphoramidite 89
5-Ethynyl-dU-CE Phosphoramidite 88
5'-Hexynyl Phosphoramidite 89
5'-I-dT-CE Phosphoramidite 89
Alkyne-Modifier Serinol Phosphoramidite 89, 95
Alkyne-NHS Ester 89
Azides 89
Azidobutyrate NHS Ester 89
baseclick Oligo-Click-M-Biotin 90
baseclick Oligo-Click-M-Fluorescein 90
baseclick Oligo-Click-M-Reload 90
baseclick Oligo-Click-M-TAMRA 90
C8-Alkyne-dC-CE Phosphoramidite 87, 90
C8-Alkyne-dT-CE Phosphoramidite 87
C8-TIPS-Alkyne-dC-CE Phosphoramidite 87, 88
C8-TIPS-Alkyne-dT-CE Phosphoramidite 88
C8-TMS-Alkyne-dC-CE Phosphoramidite 87, 88
C8-TMS-Alkyne-dT-CE Phosphoramidite 88
Click DNA and RNA Ligation 64
Copper-free Click Chemistry 89
THPTA Ligand 88
TIPS-5-Ethynyl-dU-CE Phosphoramidite 88

Convertible 2-dG
2-F-dl-CE Phosphoramidite 67

Convertible dA
O6-Phenyl-dl-CE Phosphoramidite 67

Convertible dU
O4-Triazolyl-dU-CE Phosphoramidite 67

Convertible F-dC
TMP-F-dU-CE Phosphoramidite 67

Convertible Nucleosides 67

Copper-free Click Chemistry 90
5'-DBCO-TEG Phosphoramidite 91
DBCO-dT-CE Phosphoramidite 91
DBCO-sulfo-NHS Ester 91

Cross-linking 22, 58, 60, 93, 118, 123

Custom Doping 51

Cyanine Labelling
Cyanine 3.5 Phosphoramidite 106
Cyanine 3 CPG 107
Cyanine 3 Phosphoramidite 106
Cyanine 5.5 Phosphoramidite 106
Cyanine 5 CPG 107
Cyanine 5 Phosphoramidite 106
Disulfo-Cyanine 7 Azide 93, 107

Cyanovinylcarbazole
3-Cyanovinylcarbazole Phosphoramidite 123
CNVK 123

Cyclo-dA
5',8-Cyclo-dA CE Phosphoramidite 63

Cyclo-dG
5',8-Cyclo-dG CE Phosphoramidite 63

Cyclooctatetraene
COT Serinol Phosphoramidite 97

Cytosine Arabanoside 71

D

Dabcyl Labelling
3'-Dabcyl CPG 98
3'-Dabcyl PS 98
3'-Dabsyl CPG 70, 98, 110, 112
5'-Dabcyl Phosphoramidite 98
Dabcyl-dT 98

DBCO
5'-DBCO-TEG Phosphoramidite 91
DBCO-dT-CE Phosphoramidite 91
DBCO-Serinol Phosphoramidite 91
DBCO-sulfo-NHS Ester 91

DCI (4,5-Dicyanoimidazole) 30

Deaza-5-aza-C
2'-OMe-3-deaza-5-aza-C-CE Phosphoramidite 141

Deaza-8-aza-A
7-deaza-8-Aza-A-CE Phosphoramidite 134
7-Deaza-8-aza-dA-CE Phosphoramidite 57

Deaza-8-aza-G
7-deaza-8-Aza-dG-CE Phosphoramidite 57

Deaza-A
3-Deaza-dA-CE Phosphoramidite 57
7-Deaza-dA-CE Phosphoramidite 57

Deaza-G
7-Deaza-dG-CE Phosphoramidite 57

Deaza-X
7-Deaza-dX-CE Phosphoramidite 50 57

Dendrimers
Asymmetric Doubler (LEV) Phosphoramidite 85
Long Trebler Phosphoramidite 85
Symmetric Doubler Phosphoramidite 85
Trebler Phosphoramidite 85

Depurination Resistant CE Phosphoramidites 22

Desthiobiotin
DesthiobiotinTEG-CPG 101
DesthiobiotinTEG Phosphoramidite 100

Deuterated Nucleosides 60

Diaminopurine 47, 131, 140

Dicyanoimidazole 30

Dideoxynucleoside, 2',3'- 55

Dideoxynucleosides
2',3'-ddA-CE Phosphoramidite 56
2',3'-ddC-CE Phosphoramidite 56
2',3'-ddC-CPG 56
2',3'-ddG-CE Phosphoramidite 56
2',3'-ddT-CE Phosphoramidite 56

Dihydro-dT
5,6-Dihydro-dT-CE Phosphoramidite 61

Dihydro-dU
5,6-Dihydro-dU-CE Phosphoramidite 61

Distributors 6

Dithiol
3'-Dithiol Serinol CPG 80
Dithiol Serinol Phosphoramidite 76

DNA Damage/Repair 62–63, 63

DNA Methylation
5-Carboxy-dC-CE Phosphoramidite 50
5-Formyl-dC-CE Phosphoramidite 50
5-Me-dC-CE Phosphoramidite 46

DNA Methyltransferases 71

DNP Labelling
DNP-TEG Phosphoramidite 114

Doubler Phosphoramidite
Symmetric 85

Dr. Oligo Synthesizers
CE Phosphoramidites 20
Solvents and Reagents 20
Supports and Columns 21

Duplex Stabilization 46, 47, 48, 49

E

Eclipse® Quencher
Eclipse® Quencher CPG 109
Eclipse® Quencher Phosphoramidite 108

EDTA-dT
EDTA-C2-dT-CE Phosphoramidite 118

EdU
5-Ethynyl-dU-CE Phosphoramidite 88
TIPS-5-Ethynyl-dU-CE Phosphoramidite 88

ELITechGroup Dyes and Quencher 108

Epigenetics 50, 135

Et-dC-CE Phosphoramidite 47

Etheno-A
Etheno-dA-CE Phosphoramidite 68

Ethylthiotetrazole 30

Excimers 120

Expedite™ Instruments
CE Phosphoramidites 12
Solvents and Reagents 12
Supports and Columns 13

F

FAM
3'-(6-FAM) CPG 104
3'-(6-FAM) PS 104
6-FAM 102
6-FAM-TEG Azide 92
Dipivaloyl 6-FAM-TEG Azide 92

F-ANA Monomers
2'-F-A-ANA CE Phosphoramidite 144
2'-F-Ac-C-ANA CE Phosphoramidite 144
2'-F-Bz-C-ANA CE Phosphoramidite 144
2'-F-G-ANA CE Phosphoramidite 144
2'-F-U-ANA CE Phosphoramidite 144

F-C
2'-F-Ac-C-CE Phosphoramidite 49
2'-OMe-5-F-C Precursor 141
F-dC Precursor 67

F-dl
2-F-dl-CE Phosphoramidite 67

F-dU
5-F-dU-CE Phosphoramidite 60
5-Fluoro-dU 60

Ferrocene Labelling
Ferrocene-dT-CE Phosphoramidite 119

Fluorescein Labelling
3'-(6-FAM) CPG 104
3'-(6-FAM) PS 104
3'-(6-Fluorescein) CPG 104, 107
3'-6-Fluorescein Serinol CPG 96, 104
3'-Fluorescein CPG 104

3'-Fluorescein-dT CPG 104
 5'-Dichloro-dimethoxy-Fluorescein 102
 5'-Fluorescein Phosphoramidite 102
 5'-Hexachloro-Fluorescein 102
 5'-Tetrachloro-Fluorescein 102
 6-Fluorescein Phosphoramidite 103
 6-Fluorescein Serinol Phosphoramidite 94, 103
 Dichloro-diphenyl-fluorescein 105
 Fluorescein-dT Phosphoramidite 103
 SIMA (HEX) 105

Fluorescent Nucleosides 68

2-Aminopurine-CE Phosphoramidite 58
 5-Me-2'-deoxyZebularine-CE Phosphoramidite 71
 AP-dC-CE Phosphoramidite 46
 Etheno-dA-CE Phosphoramidite 68
 Perylene-dU-CE Phosphoramidite 69
 Pyrrolo-C-TOM-CE Phosphoramidite 132
 Pyrrolo-dC-CE Phosphoramidite 68
 Pyrrolo-dCTP 68
 tC-CE Phosphoramidite 70
 tC^o-CE Phosphoramidite 70

Formyl-dC

5-Formyl-dC-CE Phosphoramidite 50
 5-Formyl-dC III-CE Phosphoramidite 50

Formylindole CE Phosphoramidite 83

Free Radicals 61

F-RNA Monomers

2'-F-Ac-C-CE Phosphoramidite 143
 2'-F-A-CE Phosphoramidite 143
 2'-F-G-CE Phosphoramidite 143
 2'-F-U-CE Phosphoramidite 143

F-U

2'-OMe-5-F-U-CE Phosphoramidite 141
 2'-OMe-TMP-5-F-U-CE Phosphoramidite 141

G

G

2',3'-ddG 56
 2'-F-G-ANA CE Phosphoramidite 144
 2'-F-G-CE Phosphoramidite 143
 2'-OMe-G-CE Phosphoramidite 137
 2'-OMe-G-PACE Phosphoramidite 145
 2'-OMe-G-RNA 139
 2'-OMe-iPr-Pac-G-CE Phosphoramidite 138
 3'-dG-CE Phosphoramidite 65
 3'-dG-CPG 55, 66
 6-Thio-dG-CE Phosphoramidite 58
 6-Thio-G-CE Phosphoramidite 133
 7-Deaza-8-aza-dG-CE Phosphoramidite 57
 7-Deaza-dG-CE Phosphoramidite 57
 8-Amino-dG-CE Phosphoramidite 62
 8-Br-dG-CE Phosphoramidite 60
 8-Oxo-dG-CE Phosphoramidite 61
 Ac-G-CE Phosphoramidite 128
 Ac-G-RNA-CPG 127, 130
 dG-CE Phosphoramidite 8, 12, 15, 16, 18, 20

dG-H-Phosphonate 39
 dG-Me Phosphoramidite 36
 dG-Me Phosphoramidite 38
 dG-PACE Phosphoramidite 37
 dG-Thiophosphoramidite 40, 142
 dmf-dG-5'-CE Phosphoramidite 34
 dmf-dG-CE Phosphoramidite 8, 12, 15, 16, 18, 20
 dmf-G-LNA-CE Phosphoramidite 41
 G-TOM-CE Phosphoramidite 126
 iPr-Pac-dG-CE Phosphoramidite 23
 iPr-Pac-G-RNA-CPG 130
 N2-Amino-Modifier C6 dG 77
 O6-Me-dG-CE Phosphoramidite 66

GalNAc

5'-GalNAc C3 Phosphoramidite 116
 GalNAc C3 CPG 116

G-Clamp 46, 68

GE Healthcare Life Sciences Instruments

CE Phosphoramidite 18
 Solvents and Reagents 19

Glen Gel-Pak™ Purification

Glen Gel-Pak™ 150

Glen-Pak™ Purification

Adapter Rack 148
 Glen-Pak™ DNA Purification Cartridge 147
 Glen-Pak™ RNA Purification Cartridge 148
 RNA Quenching Buffer 148
 Seal for Adapter Rack 148

Glen UnySupport™

Glen UnySupport CPG 24, 25
 Glen UnySupport FC CPG 25
 Glen UnySupport PS 24, 25

Glyceryl CPG 80

Gold

Conjugation to gold surfaces 94

G-Quadruplex 72

H

Halogenated Nucleosides

2'-OMe-5-F-U-CE Phosphoramidite 141
 2'-OMe-TMP-5-F-U-CE Phosphoramidite 141
 5-Br-dC-CE Phosphoramidite 60
 5-Br-dU-CE Phosphoramidite 60
 5-Br-dU-CPG 60
 5-Br-U-CE Phosphoramidite 133
 5-F-dU-CE Phosphoramidite 60
 5-I-dC-CE Phosphoramidite 60
 5-I-dU-CE Phosphoramidite 60
 5-I-U-CE Phosphoramidite 133
 8-Br-dA-CE Phosphoramidite 60
 8-Br-dG-CE Phosphoramidite 60

HEX 102

6-HEX Azide 92, 93

Hexynyl Phosphoramidite 89

High Load CPG 29

H-Phosphonate Chemistry

Monomers 39
 Reagents for ABI synthesizers 39

Hydrogen Bonding 57

Hydroxy-C

5-OH-dC-CE Phosphoramidite 61

Hydroxymethyl-dC

5-Hydroxymethyl-dC-CE Phosphoramidite 50
 5-Hydroxymethyl-dC II-CE Phosphoramidite 50

Hydroxymethyl-dU

5-Hydroxymethyl-dU-CE Phosphoramidite 61

Hydroxy-U

5-OH-dU-CE Phosphoramidite 61

I

I

2-F-dI-CE Phosphoramidite 67
 2'-OMe-I-CE Phosphoramidite 141, 142
 dI-CE Phosphoramidite 51
 dI-CPG 51
 I-CE Phosphoramidite 133
 O6-Phenyl-dI-CE Phosphoramidite 67

I-dC

5-I-dC-CE Phosphoramidite 60

I-dT

5'-I-dT-CE Phosphoramidite 89

I-dU

5-I-dU-CE Phosphoramidite 60

i-Motif DNA structures 72

Introduction 1, 3, 5, 6, 8, 10, 12, 14, 16, 18, 20

Ionizing Radiation 61

isodC

dmf-5-Me-isodC-CE Phosphoramidite 53

isodG

dmf-isodG-CE Phosphoramidite 53

Isopropyl Phosphite 39

I-U

5-I-U-CE Phosphoramidite 133

J

JOE

5'-Dichloro-dimethoxy-Fluorescein Phosphoramidite II 102

K

K

dK-CE Phosphoramidite 52

L

Labelling of MicroRNAs 64

Large Scale Synthesis

N3-Cyanoethyl-dT 71

Locked Analog Phosphoramidites

5-Me-Bz-C-LA-CE Phosphoramidite 41
 Bz-A-LA-CE Phosphoramidite 41
 dmf-G-LA-CE Phosphoramidite 41
 T-LA-CE Phosphoramidite 41

Locked Nucleic Acid (LNA) 41

M

Maleimide-Modifier

5'-Maleimide-Modifier Phosphoramidite 76

Me-C

2'-OMe-5-Me-C-CE Phosphoramidite 141
 3'-Propargyl-5-Me-dC CPG 64
 5-Me-C-TOM-CE Phosphoramidite 131
 5-Me-dC-CE Phosphoramidite 46
 Ac-5-Me-dC-CE Phosphoramidite 46

MerMade Instruments

CE Phosphoramidites 16
 Solvents and Reagents 16
 Supports and Columns 17

Methylated Nucleosides

1-Me-A-CE Phosphoramidite 135
 1-Me-dA-CE Phosphoramidite 66
 1-Me-Pseudouridine Phosphoramidite 135
 N6-Ac-N6-Me-dA-CE Phosphoramidite 47, 66
 N6-Me-A-CE Phosphoramidite 135
 N6-Me-dA-CE Phosphoramidite 47, 66
 O4-Me-dT-CE Phosphoramidite 47, 66
 O6-Me-dG-CE Phosphoramidite 66

Methylene Blue

Methylene Blue II Phosphoramidite 119
 Methylene Blue NHS Ester 119

Methyl Phosphoramidites 36

Me-U

2'-OMe-5-Me-U-CE Phosphoramidite 140
 5-Me-U-CE Phosphoramidite 133

MGB

3'-CDPI3 MGB™ CPG 48
 5'-CDPI3 MGB™ Phosphoramidite 48, 117
 CDPI3 MGB™ CPG 117

MicroRNA 64

Minor 2'-OMe-RNA Phosphoramidites 140–142

Minor Groove 58

Minor Groove Binder (MGB) 117

Mixed Base Combinations 51

Modifiers 74, 75, 78

Mutagenesis 66

N

Nebularine

2'-DeoxyNebularine-CE Phosphoramidite 51

Nitroindole

5-Nitroindole-CE Phosphoramidite 52

Non-canonical Structures 72

O

Oligo-Affinity Support

OAS PS 151

OMe-RNA Synthesis

2'-OMe-RNA Phosphoramidites 137

2'-OMe-RNA Supports 139

Minor 2'-OMe-RNA Phosphoramidites 140, 141

OMe-T 140

Oxo-dA

8-Oxo-dA-CE Phosphoramidite 61

Oxo-dG

8-Oxo-dG-CE Phosphoramidite 61

P

P

dP-CE Phosphoramidite 52, 54

PACE Phosphoramidites

2'-OMe-RNA-PACE Phosphoramidites 145

DNA PACE phosphoramidites 37

PCR/Sequencing Utilities 51

Perylene

Perylene-dU-CE Phosphoramidite 69, 120

Phenothiazine

tC-CE Phosphoramidite 70

Phenoxazine

tC°-CE Phosphoramidite 70

Phosphonocarboxylate Monomers 37

Phosphorylation

3'-CPR II CPG 82

3'-Phosphate CPG 82

3'-Phosphate CPG - High Load 82

3'-Phosphate PS 82

Chemical Phosphorylation Reagent 82

Chemical Phosphorylation Reagent II 82

CPR II 82

Solid CPR II 82

Photoaffinity Labelling 58

Photocleavable Monomers

PC Amino-Modifier Phosphoramidite 76, 86, 95

PC Biotin Phosphoramidite 86, 100

PC Linker Phosphoramidite 86

PC Spacer Phosphoramidite 84, 86

Photo cross-linking 58, 123

Photo-Regulation of DNA Function 70

NPOM-Caged-dT-CE Phosphoramidite 70

Photo-responsive DNA

Azobenzene Phosphoramidite 122

Phthalimide (PT)

3'-PT-Amino-Modifier C3 CPG 79

3'-PT-Amino-Modifier C6 CPG 79

3'-PT-Amino-Modifier C6 PS 79

Poly-Pak™ Purification

Poly-Pak™ Cartridge 149

Poly-Pak™ II Cartridge 149

Poly-Pak™ Packing 148, 149

Reagents 148, 149

Polystyrene Supports

3'-(6-FAM) PS 104

3'-BiotinTEG PS 101

3'-Dabcyl PS 98

3'-Phosphate PS 82

3'-PT-Amino-Modifier C6 PS 79

3'-TAMRA PS 113

Glen UnySupport PS 24

Universal Support III PS 26

PPG 57

Propyne Derivatives

pdC-CE Phosphoramidite 46

pdU-CE Phosphoramidite 46

Protein-DNA Interaction 58

PseudoU

1-Me-Pseudouridine Phosphoramidite 135

2'-deoxypseudoU-CE Phosphoramidite 58

PseudoUridine-CE Phosphoramidite 134

Psoralen Labelling

Psoralen Azide 22, 93, 118

Psoralen C2 Phosphoramidite 118

Psoralen C6 Phosphoramidite 118

Purification

Glen-Pak™ Purification 147, 148

Poly-Pak™ Purification 149

Purine

Purine 52

Puromycin

Puromycin CPG 121

Pyrene

5'-Pyrene Cap Phosphoramidite 49

Pyrene-dU-CE Phosphoramidite 69, 70, 120

Pyridin-2-one-CE Phosphoramidite 134

Pyrrolidine-CE Phosphoramidite 63

Pyrrolo-C

Pyrrolo-C-TOM-CE Phosphoramidite 132

Pyrrolo-dC-CE Phosphoramidite 68

Pyrrolo-dCTP 68

Q

Q-Supports 27, 28

Quenched Autoligation (QUAL) Probes 121

5'-Dabsyl-dT-CE Phosphoramidite 121

R

Redmond Red®

Redmond Red® CPG 109

Redmond Red® Phosphoramidite 108

Repair Enzyme 61

Reverse Synthesis 34

Rhodamine 113

RNA Supports

for 3' DNA Modification 125

RNA Synthesis

Minor RNA Phosphoramidites 130, 135

RNA Phosphoramidites 128, 129

RNA Supports 129, 130

RNA Supports for TOM-RNA Synthesis 126, 127

TOM-Protected Minor RNA Phosphoramidites 127, 131, 132, 133

TOM-Protected RNA Phosphoramidites 126

S

Saccharin 1-Methylimidazole 30

SBC Oligos 48

Sequence Modifiers 78

Serinol Backbone

3'-6-Fluorescein Serinol CPG 96, 104

3'-Alkyne-Modifier Serinol CPG 80, 89, 97

3'-Amino-Modifier Serinol CPG 79, 96

3'-Dithiol Serinol CPG 80, 97

3'-Protected BiotinLC Serinol CPG 96, 101

3'-Protected Biotin Serinol CPG 95, 101

6-Fluorescein Serinol Phosphoramidite 94, 103

Alkyne-Modifier Serinol Phosphoramidite 89, 96

Amino-Modifier Serinol Phosphoramidite 78, 95

COT Serinol Phosphoramidite 97

Dithiol Serinol Phosphoramidite 76, 95

Protected BiotinLC Serinol Phosphoramidite 95, 100

Protected Biotin Serinol Phosphoramidite 94, 100

SIMA

SIMA (HEX)-dT Phosphoramidite 105

SIMA (HEX) Phosphoramidite 105

SMI 30

Spacer Modifiers

1-Ethynyl-dSpacer CE Phosphoramidite 90

3'-Spacer C3 CPG 84

dSpacer CE Phosphoramidite 84

PC Spacer Phosphoramidite 84, 86

rSpacer CE Phosphoramidite 132

rSpacer TBDMS CE Phosphoramidite 134

Spacer C12 CE Phosphoramidite 84

Spacer Phosphoramidite 9 84

Spacer Phosphoramidite 18 84

Spacer Phosphoramidite C3 84

Spermine Phosphoramidite 48

Spin Labels

TEMPO Azide 93

TEMPO-TEG Azide 93

Stearyl Labelling 115

5'-Stearyl Phosphoramidite 115

Sterling

Introduction 7

Structural Studies 57

Structure/Activity Relationship 57

Sulfurizing Reagent 33

Sulfurizing Reagent II 33

T

T

2',3'-ddT 56

2-Thio-dT-CE Phosphoramidite 58

3'-Amino-dT CPG 64

3'-Amino-Modifier C6 dT CPG 81

3'-dT-CE Phosphoramidite 65

3'-dT-CPG 55, 66

3'-Fluorescein-dT CPG 104

4-Thio-dT-CE Phosphoramidite 58

5,6-Dihydro-dT-CE Phosphoramidite 61

5'-Amino-dT-CE Phosphoramidite 55

5'-Dabsyl-dT 121

5'-I-dT-CE Phosphoramidite 89

5'-OMe-dT-CE Phosphoramidite 55

Amino-Modifier C2 dT 77

Amino-Modifier C6 dT 77

C8-Alkyne-dT-CE Phosphoramidite 87

C8-TIPS-Alkyne-dT-CE Phosphoramidite 88

C8-TMS-Alkyne-dT-CE Phosphoramidite 88

DBCO-dT-CE Phosphoramidite 91

dT-CE Phosphoramidite 8, 12, 15, 16, 18, 20

dT-H-Phosphonate 39

dT-Me Phosphonamide 36

dT-Me Phosphoramidite 38

dT-PACE Phosphoramidite 37

dT-Thiophosphoramidite 40

EDTA-C2-dT-CE Phosphoramidite 118

Ferrocene-dT-CE Phosphoramidite 119

Fluorescein-dT 103

N3-Cyanoethyl-dT 71

NPOM-Caged-dT 70

O4-Me-dT-CE Phosphoramidite 47, 66

S-Bz-Thiol-Modifier C6-dT 78

TAMRA-dT 113

Thymidine Glycol CE Phosphoramidite 62

T-LNA-CE Phosphoramidite 41

TAMRA Labelling

3'-TAMRA CPG 113

3'-TAMRA PS 113

TAMRA-dT 113
TAMRA NHS Ester 113

tC

tC-CE Phosphoramidite 70

tCnitro

tCnitro-CE Phosphoramidite 70

tCo

Ribo-tC^o-CE Phosphoramidite 136
tC^o-CE Phosphoramidite 70

TEMPO

TEMPO Azide 93
TEMPO-TEG Azide 93

Termination, 3'

2',3'-ddA 56
2',3'-ddC 56
2',3'-ddC-CPG 56
2',3'-ddG 56
2',3'-ddT 56
3'-3' linkage 35, 55
3'-dA-CPG 55
3'-dC-CPG 55
3'-dG-CPG 55
3'-dT-CPG 55
3'-Spacer C3 CPG 84

Termination, 5'

5'-OMe-dT-CE Phosphoramidite 55

Terminus Modifiers 74**TET 102**

6-TET Azide 92, 93

Thio-dT

2-Thio-dT-CE Phosphoramidite 58
4-Thio-dT-CE Phosphoramidite 58

Thio-dU

4-Thio-dU-CE Phosphoramidite 58

Thio-G

6-Thio-dG-CE Phosphoramidite 58
6-Thio-G-CE Phosphoramidite 133

Thiol-Modifiers

3'-Dithiol Serinol CPG 80
3'-Thiol-Modifier 6 S-S CPG 80
5'-Thiol-Modifier C6 76
Dithiol Serinol Phosphoramidite 76
S-Bz-Thiol-Modifier C6-dT 78
Thiol-Modifier C6 S-S 76

Thiophosphoramidites

2'-OMe-RNA Thiophosphoramidites 142
DNA Thiophosphoramidites 40

Thio-U

4-Thio-U-TOM-CE Phosphoramidite 131

Thymidine Glycol

Thymidine Glycol CE Phosphoramidite 62

Thymine Dimer

Cis-syn Thymine Dimer Phosphoramidite 63

Tm Modulation 53**Tocopherol**

a-Tocopherol-TEG Phosphoramidite 115

TOM-Protecting-Group

Ac-A-RNA-CPG 126
Ac-C-RNA-CPG 127
Ac-G-RNA-CPG 127
A-TOM-CE Phosphoramidite 126
C-TOM-CE Phosphoramidite 126
G-TOM-CE Phosphoramidite 126
U-RNA-CPG 127
U-TOM-CE Phosphoramidite 126

Trebler Phosphoramidite

Trebler 85

Trimer phosphoramidites 42**Trimethoxystilbene**

5'-Trimethoxystilbene Cap Phosphoramidite 49

Triphosphate Nucleotides

Pyrrolo-dCTP 68

Triplex 57**Triplex-forming oligonucleotides 72****U****U**

2'-F-U-ANA CE Phosphoramidite 144
2'-OMe-5-Br-U-CE Phosphoramidite 141
2'-OMe-5-F-U-CE Phosphoramidite 141
2'-OMe-5-Me-U-CE Phosphoramidite 140
2'-OMe-TMP-5-F-U-CE Phosphoramidite 141
2'-OMe-U-CE Phosphoramidite 137
2'-OMe-U-PACE Phosphoramidite 145
2'-OMe-U-RNA 139
3'-Uaq Cap CPG 49, 64
4-Thio-dU-CE Phosphoramidite 58
5,6-Dihydro-dU-CE Phosphoramidite 61
5-Br-dU-CE Phosphoramidite 60
5-Br-dU-CPG 60
5-Ethynyl-dU-CE Phosphoramidite 88
5-F-dU-CE Phosphoramidite 60
5-Hydroxymethyl-dU-CE Phosphoramidite 61
5-I-dU-CE Phosphoramidite 60
5-I-U-CE Phosphoramidite 133
5-OH-dU-CE Phosphoramidite 61
Amino-Modifier C6-U Phosphoramidite 132
Br-U-CE Phosphoramidite 133
dU-CE Phosphoramidite 51
dU-CPG 500 51
dU-CPG 1000 51
O4-Triazolyl-dU-CE Phosphoramidite 67
pdU-CE Phosphoramidite 46
Perylene-dU-CE Phosphoramidite 69
Pyrene-dU-CE Phosphoramidite 69, 70
TIPS-5-Ethynyl-dU-CE Phosphoramidite 88
TMP-F-dU-CE Phosphoramidite 67

U-CE Phosphoramidite 128, 129
U-RNA-CPG 125, 127, 130
U-TOM-CE Phosphoramidite 126

UltraMILD Deprotection

2'-OMe-Ac-C-CE Phosphoramidite 138
2'-OMe-iPr-Pac-G-CE Phosphoramidite 138
2'-OMe-Pac-A-CE Phosphoramidite 138
Ac-C-CE Phosphoramidite 129
Ac-dC-CE Phosphoramidite 23
Cap Mix A 23, 38, 130, 138
iPr-Pac-dG-CE Phosphoramidite 23
iPr-Pac-G-CE Phosphoramidite 129
Pac-A-CE Phosphoramidite 129
Pac-dA-CE Phosphoramidite 23
Potassium Carbonate in Methanol 23, 38, 130, 138

UniCap Phosphoramidite 32**Universal Support III**

Universal Support III PS 26

Unnatural base pairs 48**Unnatural Base Pairs**

5-Me-isodC 53
isodG 53

V**Vitamin E 115****X****X**

2'-dX-CE Phosphoramidite 59
7-deaza-dX-CE Phosphoramidite 57

X-ray crystallography 60**Y****Yakima Yellow®**

Yakima Yellow® CPG 109
Yakima Yellow® Phosphoramidite 108

Z**Zebularine**

5-Me-2'-deoxyZebularine-CE Phosphoramidite 71
Zebularine-CE Phosphoramidite 134

Zip Nucleic Acid 48

Spermine Phosphoramidite 48

ZNA® 48

Spermine Phosphoramidite 48

ORDERING

Orders may be placed by mail, telephone, fax or email to:

Glen Research Corporation
22825 Davis Drive, Sterling, VA 20164, USA
(800) 327-GLEN (USA)
(703) 437-6191 (International)
(800) 934-2490 (Fax USA)
(703) 435-9774 (Fax)
orders@glenres.com (Internet email)
<http://www.glenres.com> (World Wide Web Home Page)

There is no minimum order requirement. Our office hours are 9:30 a.m. to 5:30 p.m. EST, Monday through Friday. During this time, Glen Research representatives will be available to take your order and answer any questions you may have regarding our product line or DNA synthesis in general. Messages may be left on voice mail at any other time. Orders placed by 3:00 p.m. EST are normally shipped the same day.

DISCOUNTS

Glen Research offers a 10% discount on catalog products to educational institutions. In addition, discounts based on volume of usage are available, so please request a specific quotation. We are also happy to quote on bulk quantities or custom packaging of our products and on custom synthesis of unique materials.

TERMS AND CONDITIONS OF SALE

Delivery will be F.O.B. Glen Research, Sterling, VA. Freight charges will be prepaid by Glen Research and will be included on the invoice sent to the Customer.

Glen Research will issue an invoice to the Customer upon shipment of Products which shall be due and payable in full within thirty (30) days of the date of the invoice. Amounts not paid when due may incur an interest charge of one and one-half (1-1/2) percent per month until fully paid. Customer shall be responsible for all taxes on Products sold to it, including state and local sales and use taxes. We reserve the right to require payment in advance before shipping any materials ordered.

Glen Research warrants that each product conforms to the specifications for such product. If the Customer notifies Glen Research within thirty (30) days of its receipt of a product that the product does not conform to the specifications, Glen Research will, at its option, replace the product or return the purchase price paid by Customer. No replacement or refund will be made if the Customer does not notify Glen Research of a non-conforming product within said thirty (30) day period.

Glen Research makes no other express or implied warranty with respect to the products. Glen Research disclaims all implied warranties, including the warranties of merchantability and fitness for a particular purpose. Glen Research's liability for breach of warranty is limited to refund of purchase price paid. In no event shall Glen Research be liable for lost profits, loss of goodwill or any incidental, special or consequential damages even if advised of the possibility of same.

It is the responsibility of the Buyer to determine the suitability of any of the listed products for any specific purpose. Responsibility for accidents arising from the handling and use of the listed products rests solely with the Buyer. All of the products should be handled only by trained personnel who are aware of the potential hazards and have ready access to suitable safety equipment. All products supplied are primarily for research and development purposes and Glen Research products are not intended for use in foods or drugs, or for any commercial purpose.

PATENTS

As a research-oriented company, we realize the desirability of patents to cover original research and it is our policy to avoid infringing any approved patents. Accordingly, it is possible that some of our products may have to be withdrawn or adjusted in price as patents are approved and issued.

Oligo synthesis success. The first time and every time.

US Headquarters

Glen Research, LLC

22825 Davis Drive, Suite 100

Sterling, VA 20164

Phone: 703-437-6191

Fax: 703-435-9774

glenresearch.com

© 2018 Glen Research. All rights reserved. For research use only.
Not intended for animal or human therapeutic or diagnostic use.

